# 1,25-(OH)<sub>2</sub>-Vitamin D ELISA Kit

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Hersteller

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# 1,25 (OH)<sub>2</sub> Vitamin D ELISA Kit

***For the in vitro determination of 1,25-(OH)<sub>2</sub> vitamin D in plasma and serum***

*EU: IVD / CE*

*US: Research Use Only. Not for use in diagnostic procedures.*


Valid from 02.09..2014

K 2112

48

## 1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of 1,25-dihydroxy vitamin D in serum and plasma. It is for *in vitro* diagnostic use only.

## 2. INTRODUCTION

Vitamin D is either produced in the skin (under the influence of UV light) or taken up from nourishment. The storage type of vitamin D, namely 25-hydroxy vitamin D, is formed in the liver. The hormone 1,25-dihydroxy vitamin D (D hormone) is formed in a second hydroxylation step in the kidney. The responsible enzyme, the kidney 1 $\alpha$ -hydroxylase, is subjected to a rigid control through hormones (especially parathyroid hormone) and its activity is influenced by the serum concentrations of calcium and phosphate.

The serum concentration of 1,25-dihydroxy vitamin D normally re-adjusts itself to the demands of metabolism. Deviations from the normal range of 1,25-dihydroxy vitamin D must therefore always be interpreted in the context of the remaining parameters of the calcium metabolism. The serum concentration of 1,25-dihydroxy vitamin D decreases only in seldom cases of vitamin D deficiency. For the diagnosis of vitamin D deficiency the precursor metabolite, 25-hydroxyvitamin D, should be measured.

The reason for a non-physiological deficiency of 1,25-dihydroxy vitamin D can be found in metabolic disturbances, caused either by genetic defects of the enzyme 1 $\alpha$ -hydroxylase (rare) or kidney malfunctions (more common). Even a slightly impaired kidney function can lead to a decrease of the 1,25-dihydroxy vitamin D concentration.

Since 1,25-dihydroxy vitamin D has important functions in calcium metabolism as well as supplementing secretion of parathyroid hormone from the parathyroid glands, increasing kidney malfunctioning leads to development of renal osteopathy, which is characterized by osteomalacia and osteitis fibrosa.

Treatment of renal osteopathy consists of the administration of 1,25-dihydroxy vitamin D (calcitriol) or the prohormone 1 $\alpha$ -hydroxy vitamin D. In renal tubules malfunctions decreased or relatively low levels of 1,25-dihydroxy vitamin D (e.g. diabetes insipidus, Fanconi-syndrom) are found. A

non-physiological over-production of 1,25-dihydroxy vitamin D arises in granulomatosis (e.g. sarcoidosis), where extra-renal synthesis of 1,25-dihydroxy vitamin D occurs. This can lead to hypercalcaemia. Also in idiopathic hypercalciuria a relatively high level of 1,25-dihydroxy vitamin D is found. Increased concentrations of 1,25-dihydroxy vitamin D can be seen in case of non-functional vitamin D receptors (rare), during calcium deficient nutrition, as well as a result from overproduction of parathyroid hormone (primary hyperthyroidism).

### Indications

- Defect of kidney functions  
Chronic kidney failure  
Haemodialysis following kidney transplantation
- Renal osteopathy
- Osteomalacia from various types of vitamin D metabolism disturbances
- Kidney tubules function disturbances (diabetes insipidus, Fanconi-Syndrom)
- Monitoring of therapy with active vitamin D metabolites
- Ideopathic hypercalciuria
- Hypercalcaemia

**3 MATERIAL SUPPLIED**

Catalogue No	Content	Kit Components	Quantity
K 2112MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 2112WP	WASHBUF	ELISA wash buffer concentrate 10x	100 ml
K 2112E	ETHANOL	Ethanol, ready to use	1.5 ml
K 2112TP	TRIS-HCL	Tris-HCL buffer, ready to use	30 ml
K 2112A1	AB	Detection antibody, anti 1,25-(OH) <sub>2</sub> vitamin D, ready to use	25 ml
K 2112ST	STD	Standard incl. NSB, ready to use (for range see specification or label)	7 x 2.5 ml
K 2112KO	CTRL	Controls, ready to use (for range see specification)	2 x 2.5 ml
K 2112K	CONJ	Conjugate, polyclonal peroxidase- labeled antibody, ready to use	24 ml
K 2112TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	2 x 15 ml
K 2112AC	STOP	ELISA stop solution, ready to use	15 ml
K 2112FOL	FOL	Foil to cover the microtiter plate	2 x 1

#### 4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- 48 Chromabond columns (Catalog No.: Se2112)
- 48 Silica Cartridges (solid phase extraction cartridges, Catalog No.: Sb2221)
- Diisopropylether (p.A.) 99.0 %
- Isopropanol (p.A.) 99.9 %
- n-Hexan (p.A.) 98.3 %
- Methanol (p.A.) 99.9 %
- 75 x 12 mm glass tubes (no plastic)
- Extraction rack (Catalog No.: K2221sv), containing three plastic stands and one nitrogen distributor
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Vacuum centrifuge or nitrogen distributor
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)


## 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label. **Prepare just the appropriate amount necessary for the assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be redissolved at room temperature or at 37°C using a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for no longer than one month.**
- Care must be taken when **preparing the solvent mixtures** because of the volume contraction. 1. Measure first the solvents in separate measuring cylinders and then mix in a storage bottle. 2. Use glass storage bottles and glass measuring cylinders; plastic vessels are not suitable.
- **Note: Microtiter strips:** Once the vacuum-sealed-aluminium-bag has been opened, all unused strips must be covered with the foil supplied and put back into the aluminium bag. Seal the aluminium bag and store at 2-8°C.
- Care must be taken not to assemble wells of different microtiter plates for analysis, even if they are of the same batch. Furthermore, once the vacuum-sealed-aluminium-bag has been opened, the microtiter plates are exposed to different conditions compared with the microtiter plates stored in the vacuum-sealed- aluminium-bags.
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

## 6. SAMPLE PREPARATION

### Serum/plasma samples

1. Fresh collected blood should be centrifuged within one hour. Vitamin D is an inert substance. The serum samples can be stored at room temperature. However, serum storage at 2-8°C is recommended when the analysis is performed within 72 h after collection. Otherwise, the serum samples must be stored at -20°C until analyzed. Avoid repeated freeze-thaw cycles.
2. Serum samples can be shipped at 2-8°C (for example with Coolpacks) and remain stable for up to 3 days. If the serum samples are not analyzed within this time period, they must be frozen and stored at -20°C.
3. Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.
4. We recommend to apply 1000 µl sample per cartridge. If the sample volume is less than 1000 µl, load an appropriate amount of Tris-HCl buffer in the column, and then the sample (minimum 500 µl) for a total volume of 1000 µl. To do so, first apply 500 µl of Tris-HCl to the column, allow to be completely absorbed and then pipette 500 µl of sample. To calculate the actual concentration, each result should be multiplied with the respective dilution factor.

*Extraction of the samples*

- The extraction unit consists of three parts, which were put on top of each other.
- The upper part is used for the Chromabond columns (extraction rack I), the lower part for the silica cartridges (extraction rack II).
- During sample application and the entire washing procedure, the whole unit should be put into a container big enough to collect the extraction solvents (extraction rack I and extraction rack II). After the first extraction step (ether) remove the extraction rack I with the Chromabond columns.
- It is recommended to place the glass-tubes (extraction rack III) directly under the cartridges (extraction rack II) for the last elution step. The tubes can then be used directly for the next step of the assay.

## 7. ASSAY PROCEDURE

### *Principle of the test*

The assay utilizes of a competitive Enzyme-Immuno-Assay (EIA) technique with a selected monoclonal antibody recognizing 1,25-dihydroxy vitamin D.

Standards, NSB (non-specific binding), controls and patient samples which are assayed for 1,25-dihydroxy vitamin D are incubated after the extraction step with the detection antibody. The pre-incubated solution is then transferred to the microplate coated with 1,25-dihydroxy vitamin D. During this incubation step, 1,25-dihydroxy vitamin D in the sample and a fixed amount of 1,25-dihydroxy vitamin D bound to the microtiter well compete for the binding of the detection antibodies. Then a peroxidase-conjugated anti-mouse antibody is added into each microplate well and a complex of 1,25-dihydroxy vitamin D - detection antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 1,25-dihydroxy vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 1,25-dihydroxy vitamin D in the samples is determined from this curve.

## Extraction

1. Apply **1000 µl** of standards, NSB, control and sample (plasma or serum) on the Chromabond columns and incubate for 10 minutes. For sample volumes less than 1000 µl wet the cartridges with Tris-HCl buffer, e.g. pipette 500 µl Tris-HCl buffer in the cartridge and 500 µl sample (results in **factor 2** for the sample evaluation).
2. Extract vitamin D from the chromabond columns with **4 x 1 ml** diisopropylether (3 min for each elution). The eluate should drip from the chromabond column directly on an untreated and dry silica cartridge. **After the extraction the chromabond columns should be removed (extraction rack I).**  
**Ensure during steps 3 and 4 that the columns should not be allowed to remain dry for longer than 5 min.**
3. Wash the silica cartridges (extraction rack II) with **5 x 2 ml Isopropanol/Hexane (4/96 v/v).**
4. Wash the silica cartridges (extraction rack II) with **3 x 2 ml Isopropanol/Hexane (6/94 v/v).**  
**After step 4, place the extraction rack unit II with the silica cartridges on the extraction rack unit III with the glass tubes.**  
**Note:** The glass tubes (extraction rack unit III) should be placed directly under the silica cartridges.
5. Elute 1,25-dihydroxy vitamin D from the silica cartridges with **2 x 2 ml Isopropanol/Hexane (25/75 v/v).**
6. Evaporate the eluate under a nitrogen stream at 37°C or in a vacuum centrifuge.
7. Before starting the pre-incubations, let the glass tubes cool down to room temperature.

## Pre-incubation

1. Add **20 µl** of ethanol into each glass tube. Immediately after adding ethanol gently vortex app. 1 sec each tube to avoid any possible evaporation.
2. Add **450 µl** antibody solution into each glass tube. Mix thoroughly for at least 10 sec.
3. Cover glass tubes with a plastic film and incubate exactly for **1 hour** at room temperature.

### *Test procedure*

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well
2. Mark the positions of STD/NSB/SAMPLE/CTRL (Standards/non-specific binding/Sample/Control) on a protocol sheet
3. Take microtiter strips out of the kit. Unused strips must be stored covered at 2-8°C. Strips are stable until the expiry date stated on the label (see P.5. Preparation and storage of reagents)
4. Add <b>200 µl of STD/NSB/SAMPLE/CTRL</b> in duplicate into respective well. All these solutions are viscous; pipette slowly and carefully. We recommend to rinse the inside of the pipette tip with the pre-incubate before pipetting it. The transfer of the samples to the microtiter plate should not take longer than 20 minutes.
5. Cover the plate tightly and incubate for 18-22 hours <b>at 6 - 10°C*</b>
6. Discard the contents of each well. Wash <b>5 times by dispensing 250 µl</b> of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
7. Add <b>200 µl CONJ</b> (conjugate) into each well
8. Cover the plate tightly and incubate for exactly <b>1 hour</b> at room temperature (18-26°C) on a horizontal mixer
9. Discard the contents of each well. Wash <b>5 times by dispensing 250 µl</b> of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
10. Add <b>200 µl of SUB</b> (substrate) into each well. <b>ATTENTION:</b> Mind the time needed for adding the substrate. It should not take longer than 2 min. Adding the stop solution in step 12 should

span the same time to ensure an identical incubation time with substrate for all samples.

11. Incubate for **15 - 25 min** at room temperature (18-26°C) in the dark\*\*

12. Add **50 µl of STOP** (stop solution) into each well, mix thoroughly. Keep the same pipetting order and time as in step 10.

13. Determine **absorption** immediately with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm as a reference

\*As with any competitive immunoassay, consistent incubation times and temperature are essential for accurate plate-to-plate comparisons. Fluctuations in overnight incubation can lead to increased inter-assay CV's. It is therefore recommended to use always the same incubation time, i.e. 20 hours.

\*\*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

### 1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

### 2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with

the used program, a control of the paired values should be done manually.

## 9. LIMITATIONS

Samples with 1,25-dihydroxy vitamin D levels greater than the highest standard value should be re-assayed. Apply instead of 1 ml sample a volume of 500 µl or 750 µl on the pre-buffered Chromabond column. Recalculate the results with the dilution factor, e.g.:

250 µl TRIS-HCL + 750 µl sample	results in factor 1.3
500 µl TRIS-HCL + 500 µl sample	results in factor 2.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### **Normal range (plasma or serum):**

Suckling babes and children up to 12:	30 - 180 pg/ml
Teenager and young adults (age 13-21):	30 - 100 pg/ml
Healthy adults (age 21-60):	17 - 65 pg/ml
Persons older than 61:	10 - 50 pg/ml
Pregnant women (8-42 week):	higher values than adults

The normal range is independent of the season.

We recommend each laboratory to establish its own norm concentration range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

Intra-Assay (n=20)		
Probe	1,25 (OH) <sub>2</sub> Vitamin D [pg/ml]	CV [%]
1	55.3	6.69

Inter-Assay (n=20)		
Probe	1,25 (OH) <sub>2</sub> Vitamin D [pg/ml]	CV [%]
1	39	9

### *Sensitivity*

The sensitivity was set as  $B_0 - 2SD$ . The zero-standard was measured 20 times.

Sample	1,25 (OH) <sub>2</sub> Vitamin D mean value [OD]	Standard variation (SD)	Detection limit [pg/ml]
1	1.201	0.025	4.8

### *Cross reactivity*

1,25-(OH) <sub>2</sub> Vit D <sub>3</sub>	100 %
1,25-(OH) <sub>2</sub> Vit D <sub>2</sub>	41 %
Vit D <sub>2</sub> & D <sub>3</sub>	< 0.01 %
24,25-(OH) <sub>2</sub> Vit D <sub>3</sub>	< 0.1 %
25-OH Vit D <sub>2</sub>	< 0.1 %
25-OH Vit D <sub>3</sub>	< 0.01 %
Alfacalcidol	< 0.003 %

### Linearity

Different volumes of two patient samples were analyzed. The results are shown below:

n= 2

Sample	Dilution	Expected [pg/ml]	Measured [pg/ml]
A	undiluted	21.2	21.2
	250 µl Buffer + 750 µl Sample	15.9	15.9
	500 µl Buffer + 500 µl Sample	10.6	11.2
	750 µl Buffer + 250 µl Sample	5.3	8.5
B	undiluted	29.6	29.6
	250 µl Buffer + 750 µl Sample	22.2	27.3
	500 µl Buffer + 500 µl Sample	14.8	16.9
	750 µl Buffer + 250 µl Sample	7.4	7.2

## 12. REGENERATION OF THE SILICA CARTRIDGES

The silica cartridges (lower columns) must be regenerated directly after sample extraction as follows:

- **2 x 2 ml** methanol
  - **2 x 2 ml** n-hexan
  - **Dry** the columns in the hood
- The silica cartridges must be dry before the next use.
  - Up to 5 regeneration cycles are possible.
  - Use equally often regenerated columns in the same run.
  - Store dried columns in plastic bags with drying agent at room temperature.

### 13. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

### 14. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

## 16. REFERENCES

1. Wildermuth S. et al., 1993; Clinica chimica Acta, 220, 61
2. Schilling M. et al., 1987; Clinical Chemistry, 33, 187
3. Armbruster F.P. et al., 1990; Ärztl. Lab., 36, 75
4. Durham B.W. et al., 1995; Ann. Clin. Biochem., 32, 77
5. Hollis B.W., 1996; Calcif. Tissue Int., 58, 4
6. Hollis B.W., 1995; Clinical Chemistry, 41, 1313
7. Iqbal S.J. et al., 1996; Clinical Chemistry, 42, 112
8. Withold W. et al., 1995; Eur. J. Clin. Chem. Clin. Biochem., 33, 15

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number