25(OH)-Vitamin D Xpress ELISA Kit

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15. LITERATUR

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone-

rone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
25(OH)-Vitamin D Xpress ELISA Kit

*For the determination of 25(OH) vitamin D
in human serum*

Valid from 2015-03-31

REF K 2107

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	19
2. CLINICAL RELEVANCE	19
3. MATERIAL SUPPLIED	20
4. MATERIAL REQUIRED BUT NOT SUPPLIED	21
5. PREPARATION AND STORAGE OF REAGENTS	21
6. SPECIMEN COLLECTION AND PREPARATION	22
7. ASSAY PROCEDURE	22
<i>Principle of the test</i>	22
<i>Test procedure</i>	23
8. RESULTS	24
9. LIMITATIONS	25
10. QUALITY CONTROL	25
<i>Reference ranges for 25(OH) vitamin D₃</i>	25
11. PERFORMANCE CHARACTERISTICS	26
<i>Precision and reproducibility</i>	26
<i>Analytical Sensitivity</i>	28
<i>Dilution recovery</i>	28
<i>Correlation data</i>	30
12. PRECAUTIONS	30
13. TECHNICAL HINTS	31
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	31
15. REFERENCES	32

1. INTENDED USE

The Immundiagnostik ELISA is intended for the quantitative determination of the 25-OH-Vitamin D in serum and fresh EDTA plasma (see chapter 6 "Specimen collection and preparation"). For *in vitro* diagnostic use only.

2. CLINICAL RELEVANCE

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. There are two different forms of vitamin D, named D₃ and D₂, which are very similar in structure. The D₂ is a synthetic product, which is predominantly absorbed by fortified food.

Physiological vitamin D₃ levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. In the liver, the vitamin is hydroxylated to 25-hydroxyvitamin D (25(OH)-vitamin D), the major circulating metabolite of vitamin D. Although 1,25-(OH)₂ vitamin D portrays the biological active form of vitamin D, which is synthesized in the kidney, it is widely accepted that the measurement of circulating 25(OH)-vitamin D provides better information with respect to patients vitamin D status and allows its use in diagnose hypovitaminosis.

Preanalytical stability of 25(OH)-vitamin D₃ in human blood or serum at room temperature: solid as a rock (Wielders and Wijnberg, 2009).

The concentration of 25(OH)-vitamin D decreases with age and a deficiency is common among elderly persons.

Clinical applications of 25(OH)-vitamin D measurements are the diagnosis and therapy control of postmenopausal osteoporosis, rickets, osteomalacia, renal osteodystrophy, pregnancy, neonatal hypocalcemia and hyperparathyroidism. In addition, a prevalence of subclinical vitamin D deficiency has been discussed in different European countries.

Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia, hypercalcuria and nephrocalcinosis in susceptible infants.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 2107MTP	PLATE	Microtiter plate frame with strips with 25-hydroxyvitamin D antigen linked to the inner surface of the polystyrene wells	12 x 8 wells
K 2107WP	WASHBUF	10x ELISA wash buffer concentrate phosphate buffered saline containing Triton X-100 and ProClin 300	2 x 100 ml
K 2107RS	RECSOL	Reconstitution solution phosphate buffered saline for dissociating 25-hydroxyvitamin D from binding protein	2 x 20 ml
K 2107RR	RELREAG	Releasing reagent lyophilised phosphate buffered saline containing a vitamin D binding protein inhibitor	2 x 1 vial
K 2107AK	AB	Anti 25(OH)-vitamin D antibody, ready to use phosphate buffered saline containing monoclonal mouse antibodies, stabilisers and preservative	1 x 10 ml
K 2107CAL	CAL	Calibrator, ready to use buffered human serum containing 25-hydroxyvitamin D and 0.09% sodium azide	1 vial à 300 µl
K 2107KO	CTRL A CTRL B	Controls, ready for use see specification for range buffered human serum containing 25-hydroxyvitamin D and 0.09% sodium azide	300 µl 300 µl
K 2107K	CONJ	Conjugate, peroxidase labeled, ready to use phosphate buffered saline containing polyclonal anti mouse-HRP, antibodies stabilisers and preservative	1 x 20 ml
K 2107TMB	SUB	TMB substrate aqueous formulation of tetramethylbvenzidine and hydrogen peroxide	1 x 20 ml
K 2107SD	SAMDIL	Sample dilution buffer phosphate buffered saline containing stabilisers and preservative	1 x 15 ml
K 2107AC	STOP	ELISA stop solution, ready to use sulfuric acid, 0,4 M	1 x 20 ml
K 2107FOL	FOL	Foil to cover the microtiter plate	2 x 1

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Freezer (-20°C)
- Calibrated precision pipettors and 10–1000 µl tips
- Multi-channel pipets or repeater pipets
- Vortex
- Microplate thermoshaker at 37°C with shaking function (available via Immundiagnostik upon request)
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- Allow all reagents to come to **room temperature** before use in the assay.
- The test kit is designed for 96 single determinations. A reduction of the sample or buffer volumes results in erroneous values.
- To run assay more than once, ensure that reagents are stored at conditions stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8°C for two weeks**.
- Reconstitute **RELREAG** (releasing reagent) in **16 ml RECSOL** (reconstitution solution), mix gently by careful swinging (do not vortex). **Unused volume of RELREAG can be stored at -20°C upon further use**. It is then stable until the date of expiry (see label). Frozen RELREAG can be defrosted and used only once. Let frozen reagent come to ambient temperature prior to use.
- Bring **AB** (antibody) at room temperature at least one hour before use.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

- **Note: Microtiter strips:** Once the vacuum-sealed aluminium bag has been opened, all unused strips must be covered with the foil supplied and put back into the aluminium bag. Seal the bag and store it at 2–8 °C. Strips handled in such a way can be used within 4 weeks.

6. SPECIMEN COLLECTION AND PREPARATION

1. a) **Freshly collected blood** should be centrifuged within one hour. Vitamin D is an inert substance. However, serum storage at 2–8 °C **is recommended** when the analysis is performed within 24 h after collection. Otherwise, the serum samples must be stored at -20 °C until analyzed. Avoid repeated freeze-thaw cycles.
b) **Freshly collected EDTA plasma** should be measured within 8 days. It should be stored at 2–8 °C until the measurement. Do not freeze plasma samples.
2. Serum samples and EDTA plasma samples can be shipped at 2–8 °C (for example with coolpacks) and remain stable for up to 3 days.
3. Serum is the preferred sample matrix; **whole blood is not suitable**.
4. Indicated incubation times and temperatures must be strictly observed.
5. Mix samples well before use.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive ELISA technique with a selected monoclonal antibody recognizing 25(OH)-vitamin D. For a reliable determination of 25(OH)-vitamin D, it is necessary to release it from the 25(OH)-vitamin D-VDBP-complex.

The calibrator, controls and patient samples which are assayed for 25(OH)-vitamin D are incubated with the releasing reagent. The pre-incubated solutions are then transferred to the microplate coated with 25(OH)-vitamin D, and an anti-25(OH)-vitamin D antibody is added. During this incubation step, 25(OH)-vitamin D in the sample and a fixed amount of 25(OH)-vitamin D bound to the microtiter well compete for the binding of the antibody. Then a peroxidase-conjugated antibody is added into each microplate well. A complex of 25(OH)-vitamin D – anti-25(OH)-vitamin D antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the 25(OH) vitamin D concentration of the sample. Samples are quantified by referring their optical density to a lot-dependent master

calibration curve and the use of the calibrator that is run with each test.

Test procedure

1.	Prior to use in the assay allow all reagents and samples to come to room temperature (18–26 °C). For this purpose, open the kit, take out the needed individual components and mix gently avoiding foam formation.
2.	Mark the positions of CAL (calibrator, in duplicate)/SAMPLE/CTRLs (control A and B) on a protocol sheet.
3.	Label V-tubes (e. g. 1.5 ml Eppendorf-tubes).
4.	Pipet 10 µl of CAL (calibrator)/ SAMPLE/CTRLs (control A and B) respectively, into the corresponding tube. Instead of dilution in polypropylene reaction tubes, you can also dilute the samples in Deep Well® DSX dilution tubes (available as single tubes or 8 well strips). This offers the additional advantage that you can transfer the diluted samples with a multi-channel pipet directly to the microtiter stripes used for the test.
5.	Reconstitute RELREAG (releasing reagent) (see chapter 5., page 19).
6.	Add 300 µl of RELREAG (releasing reagent) into each tube, vortex shortly.
7.	Take microtiter strips out of the microtiter module. Unused strips must be covered with the enclosed foil, stored at 2–8 °C and used within 4 weeks.
8.	Transfer 100 µl of CAL (calibrator, in duplicate)/ SAMPLE/CTRL (control A and B) from the V-tubes to respective well.
9.	Incubate in the thermoshaker for 75 minutes at 37 °C while shaking at 900 rpm. Note: do not seal the microtiter strips for this step.
10.	Add 100 µl of SAMDIL (sample dilution buffer) into each well.
11.	Add 50 µl of AB (anti 25(OH)-vitamin D antibody) into each well.
12.	Incubate microtiter plate in the thermoshaker for 30 minutes at 37 °C while shaking at 900 rpm. Note: do not seal the microtiter strips for this step.

13.	Aspirate and wash the wells 5 x with 350 µl of diluted wash buffer. The use of an 8-channel pipet is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash. For TECAN and Dynex instruments a programming protocol can be requested from Immundiagnostik AG.
14.	Add 150 µl of CONJ (conjugate) into each well.
15.	Incubate microtiter plate in the thermoshaker for 30 minutes at 37°C while shaking at 900 rpm. Note: do not seal the microtiter strips for this step.
16.	Aspirate and wash the wells 5 x with 350 µl of diluted wash buffer. The use of an 8-channel pipet is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash.
17.	Add 150 µl of SUB (substrate) into each well.
18.	Incubate for 20 minutes at room temperature (18–26°C) in the dark. Note: do not shake the microtiter strips.
19.	Add 150 µl of STOP (stop solution) into each well.
20.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm as a reference. If no reference wavelength is available, read only at 450 nm.

8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A,B,C and D. In both cases the optical density of the calibrator (CAL) is essential. Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Note: The parameters A, B, C, D and CAL must be entered exactly into the evaluation program. Even small deviations of the values can cause significant interference with the evaluation of the results.

The plausibility of the pairs of values should always be examined after the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually by the operator.

9. LIMITATIONS

Whole blood is not suitable as a sample.

10. QUALITY CONTROL

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference ranges for 25(OH) vitamin D₃

Information from ASBMR 2011

Deficiency (seriously deficient)	< 20 ng/ml	or	< 50 nmol/l
Insufficiency (deficient)	20–29 ng/ml	or	50–74 nmol/l
Sufficiency (adequately supplied)	> 30 ng/ml	or	> 75 nmol/l

Society of Osteologists SACHSEN E. V.

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1 ng/ml = 2.5 nmol/l

1 nmol/l = 0.4 ng/ml

Reference intervals for 25(OH)-vitamin D₃ (ng/ml), males and females*

Age	n	2.5%	97.5%
0 to < 3 months	131	5	42
3 to < 6 months	135	9	60
6 months to < 1 year	147	18	58
1 to < 3 years	394	15	54
3 to < 10 years	619	14	46
10 to < 13 years	286	11	50

13 to < 15 years	275	10	44
15 to < 18 years	390	8	45
> 18 years	421	8	56

*All seasons

(Soldin et al., 2009)

Note

The vitamin D production in the skin is high variable and depends on the season and daily time, degree of latitude, age, sun protection etc. The normal ranges depend on the method used (e.g. vitamin D release from the vitamin D binding protein, VDBP) and serve only as orientation.

We recommend each laboratory to establish an own reference concentration range based on their local population.

Literature references

The following literature references for the vitamin D reference values can be found in the references on page 28: Grant et al., Soldin et al., Visser et al., Wicherts et al.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Two reference samples were measured ten times by using three different kit lots.

Reference sample	P1			P2		
Kit lot	130125	130201	130418	130125	130201	130418
Test nr.	Result in nmol/l					
1	46.1	41.7	37.0	123.8	97.4	103.5
2	41.6	39.4	45.8	110.8	108.0	106.1
3	42.5	47.4	43.1	117.9	114.7	118.9
4	46.2	49.1	40.1	117.9	111.3	101.5
5	41.1	48.3	38.21	121.0	118.5	105.1
6	43.9	44.0	39.1	110.9	102.4	99.3
7	40.9	40.7	42.2	104.8	117.7	90.7
8	49.4	45.7	42.7	114.8	113.7	107.2
9	44.2	51.9	42.2	111.1	123.6	120.7

Reference sample	P1			P2		
10	49.1	50.8	42.6	113.9	111.5	105.4
Mean	44.6	45.9	41.3	114.7	111.9	105.8
SD	3.1	4.3	2.7	5.6	7.7	8.8
Intra-assay CV in %	6.9	9.4	6.4	4.9	6.9	8.3
Mean	43.9			110.8		
SD	3.6			8.2		
Inter-assay CV in %	8.8			7.4		

P1: Repeatability (intra-assay) CV is for lot 1: 6.9%, lot 2: 9.4%, lot 3: 6.4%. Inter-batch (inter-assay) CV is 8.8%.

P2: Repeatability (intra-assay) CV is for lot 1: 4.9%, lot 2: 6.9%, lot 3: 8.3%. Inter-batch (inter-assay) CV is 7.4%.

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to 25(OH)-vitamin D₃. The specificity is calculated in percent, based on the cross-reactivity of these compounds with the anti-25(OH)-vitamin D₃ antibody compared to the 25(OH)-vitamin D₃ antigen.

- 25(OH) vitamin D₃ 100.0%
- 25(OH) vitamin D₂ 67%
- 24,25(OH) vitamin D₃ ≥ 100.0%
- Vitamin D₂ (ergocalciferol) 0.3%
- Vitamin D₃ (cholecalciferol) 0.8%

Analytical Sensitivity

Limit of detection, LoD


5.6 ng/ml 14 nmol/l

Linear range

6.40–100 ng/ml 16–250 nmol/l

Two samples with known concentrations were diluted to 11 decreasing levels of concentration within the expected range of measurement. A linear regression analysis was performed.

theoretical nmol/l	measured nmol/l	recovery in %
310	307	99
286	279	97
263	246	94
239	226	94
219	207	95
192	182	95
157	146	93
142	133	94
127	121	95
68	64	93
36	35	98
21	22	107
15	15	101

$r = 0,9991$

The linear range of the assay was found to be 16–250 nmol/l, respectively 6.4–100 ng/ml

Dilution recovery

Samples with concentrations higher than 320 nmol/L must be diluted as follows:

To dilute a sample please pipet a lower volume of the prediluted sample (10 µl sample + 300 µl RELREAG) to the well.

The lowered volume of prediluted sample (10 µl sample + 300 µl RELREAG) together with the fixed volumes of all other reagents leads to the following dilution factors that need to be applied for calculating the original sample concentration.

Sample	Volume of prediluted Sample (10 µl+300 µl RelReag)	target concentration [nmol/l]	observed concentration [nmol/l]	recovery [in %]	dilution factor
A	100 µl	126.00	126.00		1
	90 µl	118.13	116.5	98.6	1.0666
	80 µl	109.57	101.18	92.3	1.15
	70 µl	100.24	94.4	94.2	1.257
	60 µl	90.00	85.5	95.0	1.4
	50 µl	78.75	76.19	96.7	1.6
	40 µl	66.32	72.7	109.6	1.9
B	100 µl	206.03	206.03		1
	90 µl	193.17	191.7	99.2	1.0666
	80 µl	179.16	188.02	104.9	1.15
	70 µl	163.91	183.66	112.1	1.257
	60 µl	147.16	164.04	111.5	1.4
	50 µl	128.77	146.76	114.0	1.6
	40 µl	108.44	122.01	112.5	1.9

When using lower volume of prediluted sample while using the original volumes of all other reagents we get an effective dilution within the reaction well compared to the original procedure. The original procedure uses 100 µl of prediluted sample (x) + 100 µl SD + 50 µl AB = 250 µl total in the reaction well with 150 µl total of reagents. Using the formula $(x + 150) / x$ we get the fraction of the sample volume used per total volume in the reaction well. The formula is:

$$\text{Dilution factor} = [(x + 150) / x] / [(100 + 150) / 100]$$

For example the dilution factor for the use of x = 50 µl prediluted sample is:



$$(50 + 150) / 50 = 4$$

$$(100 + 150) / 100 = 2.5$$

$$\Rightarrow 4 / 2.5 = \text{factor } 1.6$$

Correlation data

Excellent correlation with 25(OH)-Vitamin D direct ELISA und LC-MS/MS results.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- The quality control guidelines should be observed.
- Control samples should be analyzed with each run.

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- **Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions than sealed ones.**
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by