25(OH)-Vitamin D *dried blood* ELISA Kit

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15. LITERATUR

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone-

rone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
25(OH) vitamin D *dried blood* ELISA Kit

For the determination of 25(OH) vitamin D in dried blood

Valid from 2015-06-29

REF K 2108DBS

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. SUMMARY AND EXPLANATION OF THE TEST	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPERATION AND STORAGE OF REAGENTS	19
6. SPECIMEN COLLECTION AND PREPARATION	20
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Test procedure</i>	21
8. RESULTS	22
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Reference ranges for 25(OH) vitamin D₃</i>	23
11. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Specificity</i>	24
<i>Analytical Sensitivity</i>	25
12. PRECAUTIONS	25
13. TECHNICAL HINTS	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	26

1. INTENDED USE

This Immundiagnostik ELISA is intended for the quantitative determination of the 25-OH-vitamin D in dried blood. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. There are two different forms of vitamin D, named D₃ and D₂, which are very similar in structure. The D₂ is a synthetic product, which is predominantly absorbed by fortified food.

Physiological vitamin D₃ levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. In the liver, the vitamin is hydroxylated to 25-hydroxyvitamin D (25(OH)-vitamin D), the major circulating metabolite of vitamin D. Although 1,25-(OH)₂ vitamin D is the biological active form of vitamin D, which is synthesized in the kidney, it is widely accepted that the measurement of circulating 25(OH)-vitamin D provides better information with respect to patients vitamin D status and allows its use in diagnose hypovitaminosis (1,2).

Preanalytical stability of 25(OH)-vitamin D₃ in human blood or serum at room temperature: solid as a rock (Wielders and Wijnberg, 2009).

The concentration of 25(OH)-vitamin D decreases with age and a deficiency is common among elderly persons.

Clinical applications of 25(OH)-vitamin D measurements are the diagnosis and therapy control of postmenopausal osteoporosis, rickets, osteomalacia, renal osteodystrophy, pregnancy, neonatal hypocalcemia and hyperparathyroidism. In addition, a prevalence of subclinical vitamin D deficiency has been discussed in different European countries.

Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia, hypercalcuria and nephrocalcinosis in susceptible infants.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 2108DBS	PLATE	Microtiter plate frame with strips with 25-hydroxyvitamin D antigen linked to the inner surface of the polystyrene wells	12 x 8 wells
K 2108DBS	WASHBUF	10x ELISA wash buffer concentrate phosphate buffered saline containing Triton X-100 and ProClin 300	2 x 100 ml
K 2108DBS	AB	Anti 25(OH)-vitamin D antibody, ready to use phosphate buffered saline containing monoclonal mouse antibodies, stabilisers and preservative	18 ml
K 2108DBS	STD	Standards, ready to use buffered human serum containing 25-hydroxyvitamin D and 0.09% sodium azide	6 vials
K 2108DBS	CTRL 1 CTRL 2	Controls, ready for use see specification for range	each 1 vial
K 2108DBS	CONJ	Conjugate, peroxidase labeled, ready to use phosphate buffered saline containing anti mouse-HRP, stabilisers and preservative	24 ml
K 2108DBS	SUB	TMB substrate aqueous formulation of tetramethylbenzidine and hydrogen peroxide	2 x 15 ml
K 2108DBS	STOP	ELISA stop solution, ready to use sulfuric acid	1 x 15 ml
K 2108DBS	ACTSOL	Activation solution Salt solution	10 ml
K 2108DBS	ASYBUF	Assay buffer phosphate buffered saline, stabilisers and preservative	15 ml
K 2108DBS	Ethanol	Ethanol	2 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10-1000 µl tips

- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 g, 2–8 °C
- Standard laboratory vials, cups, etc. (single-use products) made of polypropylene
- Glass tubes, 12 × 75 mm
- Microtiter plate reader (450 nm, reference wave length 620 or 690 nm)
- Refrigerator
- Water bath, 37 °C
- Acetonitril

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPERATION AND STORAGE OF REAGENTS

- Allow all reagents to come to room temperature before use in the assay.
- The test kit is designed for 96 single determinations. A reduction of the sample or buffer volumes results in erroneous values.
- To run assay more than once, ensure that reagents are stored at conditions stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- Bring **AB** (antibody) and **ASYBUF** (assay buffer) at room temperature (15–30 °C) at least one hour before use.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2–8 °C**.
- **Note: Microtiter strips:** Once the vacuum-sealed aluminum bag has been opened, all unused strips must be covered with the foil supplied and put back into the aluminium bag. Close the aluminium bag and store it at 2–8 °C. Strips handled in such a way can be used within 4 weeks.

6. SPECIMEN COLLECTION AND PREPARATION

We recommend DrySpot-ID (catalogue no. DZ9025ID) as carrier for dried blood samples.

1.	Prepare the amount of polypropylene reaction tubes needed.
2.	Cut one of each control (CTRL 1 und 2) diagonally, put it into an reaction tube and add each 50 µl ACTSOL (activation solution).
3.	Put dried blood samples (SAMPLE) in reaction tubes and add each 50 µl ACTSOL (activation solution).
4.	Put 25 µl STD (standard) in reaction tubes.
5.	Vortex SAMPLE and CTRL shortly.
6.	Add 300 µl acetonitril to STD, CTRL, SAMPLE , mix and incubate for 30 min in a water bath at 37°C . Then centrifuge for 10 min at 2–8°C and 3000 g .
7.	Pipet 200 µl supernatant in a glass tube (12 × 75 mm) and evaporate .
8.	Add 10 µl EtOH (ethanol), mix shortly.
9.	Add 120 µl ASYBUF (assay buffer), mix well .

The samples are now ready to be used in the test.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive ELISA technique with a monoclonal antibody recognizing 25(OH)-vitamin D. For a reliable determination of 25(OH)-vitamin D, it is necessary to release it from the 25(OH)-vitamin D-VDBP-complex.

Standards, controls and dried blood samples which are assayed for 25(OH)-vitamin D are extracted with acetonitril and transferred to the microplate coated with 25(OH)-vitamin D. Then, an anti-25(OH)-vitamin D antibody is added. During an incubation step, 25(OH)-vitamin D in the sample and a fixed amount of 25(OH)-vitamin D bound to the microtiter well compete for the binding of the antibody. Then a peroxidase-conjugated antibody is added into each microplate well. A complex of 25(OH)-vitamin D – anti-25(OH)-vitamin D antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow.

The intensity of the yellow color is inversely proportional to the concentration of 25(OH)-vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 25(OH)-vitamin D in the samples is determined from this curve.

Test procedure

Prior to use in the assay allow all reagents and samples to come to room temperature. For this purpose, open the kit, take out the needed individual components and mix gently avoiding foam formation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Mark the positions of STD (standards)/CTRL (control) and SAMPLE on a protocol sheet.
2.	Take microtiter strips out of the microtiter module. Unused strips must be covered with a foil, stored at 2–8 °C and used within 4 weeks
3.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
4.	Add 50 µl resuspended STD/CTRL/SAMPLE to all wells.
5.	Add 110 µl AB (anti 25(OH) vitamin D antibody) into each well
6.	Cover the plate again tightly with a foil and incubate for 45 min at room temperature* in the dark .
7.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add 200 µl CONJ (conjugate) into each well
9.	Cover the plate again tightly with the enclosed foil and incubate for 45 min at room temperature* in the dark .
10.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.

11.	Add 200 µl SUB (substrate) into each well
12.	Incubate for 15–20 minutes at room temperature* in the dark
13.	Add 50 µl STOP (stop solution) into each well
14.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e.g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

* **The optimal ambient temperature range is 22–25 °C.** All other temperatures result in strong deviations from the optical densities described in the QC data sheet.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend the use of the „4-parameter-algorithm“.

1. 4-parameter-algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Dried blood samples

We recommend DrySpot-ID (catalogue no. DZ9025ID) as carrier for dried blood samples. When using DrySpot-ID, the result of the samples must be **multiplied** by the **factor of 1.5** to obtain the 25(OH) vitamin D concentration. Standards and controls must not be multiplied.

Immundiagnostik AG cannot guarantee accuracy of test results if other dried blood carriers are used.

9. LIMITATIONS

This test is supposed to be used with dried blood samples resulting from 50 µl of capillary blood. Deviations from this value influence the test result.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference ranges for 25(OH) vitamin D₃

Information from ASBMR 2011

Deficiency (seriously deficient)	< 20 ng/ml	or	< 50 nmol/l
Insufficiency (deficient)	20–29 ng/ml	or	50–74 nmol/l
Sufficiency (adequately supplied)	≥ 30 ng/ml	or	> 75 nmol/l

Society of Osteologists SACHSEN E. V.

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1 ng/ml = 2.5 nmol/l

1 nmol/l = 0.4 ng/ml

Reference intervals for 25(OH)-vitamin D₃ (ng/ml), males and females*

Age	n	2.5%	97.5%
0 to < 3 months	131	5	42
3 to < 6 months	135	9	60
6 months to < 1 year	147	18	58
1 to < 3 years	394	15	54
3 to < 10 years	619	14	46
10 to < 13 years	286	11	50
13 to < 15 years	275	10	44

15 to < 18 years	390	8	45
> 18 years	421	8	56

*All seasons

(Soldin et al., 2009)

Note

The vitamin D production in the skin is high variable and depends on the season and daily time, degree of latitude, age, sun protection etc. The normal ranges depend on the method used (e.g. vitamin D release from the vitamin D binding protein, VDBP) and serve only as orientation.

Literature references

The following literature references for the vitamin D reference values can be found in the references on page 28: Grant et al., Soldin et al., Visser et al., Wicherts et al.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)

Sample	25(OH)-vitamin D [ng/ml]	CV [%]
1	48.3	8.8
2	26.9	6.0

Inter-Assay (Day-to-Day variation; n = 20)

Sample	25(OH)-vitamin D [ng/ml]	CV [%]
1	21.3	11.8
2	36.2	7.9

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to 25(OH)-vitamin D₃. The specificity is calculated in percent, based on the cross-reactivity of these compounds with the anti-25(OH)-vitamin D₃ antibody compared to the 25(OH)-vitamin D₃ antigen.

- | | |
|---|----------|
| • 25(OH) vitamin D ₃ | 100,0% |
| • 25(OH) vitamin D ₂ | 67,8% |
| • 24, 25(OH) vitamin D ₃ | ≥ 100,0% |
| • Vitamin D ₂ (ergocalciferol) | 0,3% |

Analytical Sensitivity

The Zero-standard was measured 20 times. The detection limit was set as $B_0 + 2 SD$ and estimated to be 2,7 ng/ml.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions than sealed ones.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.

- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.


15. REFERENCES

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**,

1584–5 (2009).

7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by