25(OH)-Vitamin D *direct day* ELISA Kit

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15. LITERATUR

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxyvitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

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
Manual

25(OH)-Vitamin D *direct day* ELISA Kit

For the determination of 25(OH)-vitamin D in human serum

Valid from 2015-05-26

REF K 2108

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. SUMMARY AND EXPLANATION OF THE TEST	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	19
5. PREPERATION AND STORAGE OF REAGENTS	19
6. SPECIMEN COLLECTION AND PREPARATION	20
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Test procedure</i>	21
8. RESULTS	23
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Reference ranges for 25(OH) vitamin D₃</i>	24
11. PERFORMANCE CHARACTERISTICS	25
<i>Precision and reproducibility</i>	25
<i>Dilution recovery</i>	25
<i>Spiking Recovery</i>	26
<i>Analytical Sensitivity</i>	26
<i>Specificity</i>	27
12. PRECAUTIONS	27
13. TECHNICAL HINTS	27
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	28
15. REFERENCES	28

1. INTENDED USE

This Immundiagnostik ELISA is intended for the quantitative determination of the 25-OH-Vitamin D in serum and fresh plasma. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. There are two different forms of vitamin D, named D₃ and D₂, which are very similar in structure. The D₂ is a synthetic product, which is predominantly absorbed by fortified food.

Physiological vitamin D₃ levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. In the liver, the vitamin is hydroxylated to 25-hydroxyvitamin D (25(OH)-vitamin D), the major circulating metabolite of vitamin D. Although 1,25-(OH)₂ vitamin D is the biological active form of vitamin D, which is synthesized in the kidney, it is widely accepted that the measurement of circulating 25(OH)-vitamin D provides better information with respect to patients vitamin D status and allows its use in diagnose hypovitaminosis (1,2).

Preanalytical stability of 25(OH)-vitamin D₃ in human blood or serum at room temperature: solid as a rock (Wielders and Wijnberg, 2009).

The concentration of 25(OH)-vitamin D decreases with age and a deficiency is common among elderly persons.

Clinical applications of 25(OH)-vitamin D measurements are the diagnosis and therapy control of postmenopausal osteoporosis, rickets, osteomalacia, renal osteodystrophy, pregnancy, neonatal hypocalcemia and hyperparathyroidism. In addition, a prevalence of subclinical vitamin D deficiency has been discussed in different European countries.

Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia, hypercalcuria and nephrocalcinosis in susceptible infants.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 2108MTP	PLATE	Microtiter plate frame with strips with 25-hydroxyvitamin D antigen linked to the inner surface of the polystyrene wells	12 x 8 wells
K 2108WP	WASHBUF	10x ELISA wash buffer concentrate phosphate buffered saline containing Triton X-100 and ProClin 300	2 x 100 ml
K 2108RS	RECSOL	Reconstitution solution phosphate buffered saline for dissociating 25-hydroxyvitamin D from binding protein; red coloured	2 x 20 ml
K 2108RR	RELREAG	Releasing reagent lyophilised phosphate buffered saline containing a vitamin D binding protein inhibitor	2 x 1 vial
K 2108AK	AB	Anti 25(OH)-vitamin D antibody, ready to use phosphate buffered saline containing monoclonal mouse antibodies, stabilisers and preservative	18 ml
K 2108ST	STD	Standards, ready to use buffered human serum containing 25-hydroxyvitamin D and 0.09% sodium azide	6 vials 300 µl each
K 2108KO	CTRL A CTRL B	Controls, ready for use see specification for range buffered human serum containing 25-hydroxyvitamin D and 0.09% sodium azide	300 µl 300 µl
K 2108K	CONJ	Conjugate, peroxidase labeled, ready to use phosphate buffered saline containing anti mouse-HRP, stabilisers and preservative	24 ml
K 2108TMB	SUB	TMB substrate aqueous formulation of tetramethylbenzidine and hydrogen peroxide	2 x 15 ml
K 2108AC	STOP	ELISA stop solution, ready to use sulfuric acid	1 x 15 ml
K 2108FOL	FOL	Foil to cover the microtiter plate	3 x 1

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Deep freezer (-20°C)
- Calibrated precision pipettors and 10-1000 µl tips
- Multi-channel pipets or repeater pipets
- Vortex
- Standard laboratory vials, cups, etc. (single-use products) made of polypropylene
- Microtiter plate reader (450 nm, reference wave length 620 or 690 nm)
- Refrigerator

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles >0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPERATION AND STORAGE OF REAGENTS

- Allow all reagents to come to room temperature before use in the assay.
- The test kit is designed for 96 single determinations. A reduction of the sample or buffer volumes results in erroneous values.
- To run assay more than once, ensure that reagents are stored at conditions stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved in a water bath at 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for 2 weeks**.
- Make sure that **RECSOL** (reconstitution solution) has reached room temperature before use. For example by warming it up for 10 minutes in a waterbath at 37 °C.
- Reconstitute **RELREAG** (releasing reagent) in **16 ml RECSOL** (reconstitution solution), mix gently by careful swinging (do not vortex). 16 ml RELREAG are enough to perform 48 single determinations, for 96 determinations both vials of RELREAG are needed. After reconstitution, RELREAG is stable for 4 weeks at 2–8 °C in the dark. For long term storage, freeze aliquots at -20 °C. It is then stable until the date of expiry (see label). Frozen RELREAG can be defrosted

and used only once. Pre-heat frozen releasing reagent to room temperature before use.

- Bring **AB** (antibody) at room temperature at least one hour before use.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2–8 °C**.
- **Note: Microtiter strips:** Once the vacuum-sealed aluminum bag has been opened, all unused strips must be covered with the foil supplied and put back into the aluminium bag. Close the aluminium bag and store it at 2–8 °C. Strips handled in such a way can be used within 4 weeks.

6. SPECIMEN COLLECTION AND PREPARATION

1. Fresh collected blood should be centrifuged within one hour. Vitamin D is an inert substance. However, serum storage at 2-8°C is recommended when the analysis is performed within 24 h after collection. Otherwise, the serum samples must be stored at -20°C until analyzed. Avoid repeated freeze-thaw cycles.
2. **Serum** samples can be shipped at 4-8 °C (for example with Coolpacks) and remain stable for up to 3 days.
3. Serum is the preferred sample matrix; **whole blood is not suitable**.
4. Indicated **incubation times and temperatures** must be strictly observed. The room temperature should be checked with a thermometer. **Attention: We advise not to use a heatable incubator.**
5. Mix samples well before use.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive ELISA technique with a selected monoclonal antibody recognizing 25(OH)-vitamin D. For a reliable determination of 25(OH)-vitamin D, it is necessary to release it from the 25(OH)-vitamin D-VDBP-complex.

Standards, controls and patient samples which are assayed for 25(OH)-vitamin D are prediluted with the releasing reagent and transferred to the microplate coated with 25(OH)-vitamin D. After an incubation to release the 25(OH) vitamin D, an anti-25(OH)-vitamin D antibody is added. During an incubation step, 25(OH)-vitamin D in the sample and a fixed amount of 25(OH)-vitamin D bound to the microtiter well compete for the binding of the antibody. Then a peroxidase-conjugated antibody

is added into each microplate well. A complex of 25(OH)-vitamin D – anti-25(OH)-vitamin D antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 25(OH)-vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 25(OH)-vitamin D in the samples is determined from this curve.

Test procedure

Prior to use in the assay allow all reagents and samples to come to room temperature. For this purpose, open the kit, take out the needed individual components and mix gently avoiding foam formation.

Sample dilution

1	Prepare the amount of polypropylene reaction tubes needed.
2.	Pipet 20 µl of STD (standard)/ SAMPLE/CTRLs (control A and B) respectively, into the corresponding tube
3.	Reconstitute RELREAG (releasing reagent) (see chapter 5, page 19)
4.	Add 300 µl of RELREAG (releasing reagent) into each tube. The transfer of the diluted samples to the microtiter plate has to be done within 5–10 min.

Instead of dilution in polypropylene reaction tubes, you can also dilute the samples in Deep Well® DSX dilution tubes (available as single tubes or 8 well strips). This offers the additional advantage that you can transfer the diluted samples with a multi-channel pipet directly to the microtiter stripes used for the test.

Test procedure

1.	Mark the positions of STD (standards)/ SAMPLE/CTRLs (control A and B) on a protocol sheet
2.	Take microtiter strips out of the microtiter module. Unused strips must be covered with the enclosed foil, stored at 2–8 °C and used within 4 weeks
3.	Put 20 µl of the samples , diluted as described in the previous steps, speedily into the wells of the microtiter plate, cover the stripes with the enclosed foil and incubate at room temperature (18–28 °C) for 60 min* .

4.	Add 150 µl of AB (anti 25(OH)-vitamin D antibody) into each well
5.	Cover the plate again tightly with the enclosed foil and incubate for 45 min at room temperature.
6.	Aspirate and wash the wells 5 x with 250 µl of diluted wash buffer. Remove remaining wash buffer by hitting the plate against paper towel after the last wash For TECAN and Dynex instruments a programming protocol can be requested from Immundiagnostik AG
7.	Add 200 µl CONJ (conjugate) into each well
8.	Cover the plate again tightly with the enclosed foil and incubate for 45 min at room temperature.
9.	Aspirate and wash the wells 5 x with 250 µl of diluted wash buffer. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
10.	Add 200 µl of SUB (substrate) into each well
11.	Incubate for 10–15 minutes at room temperature in the dark
12.	Add 50 µl of STOP (stop solution) into each well
13.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

* **The optimal ambient temperature range is 22–25°C.** All other temperatures result in strong deviations from the optical densities described in the QC data sheet.

8. RESULTS

We recommend the 4-parameter algorithm for result calculation.

4-parameter algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The usage of the 4-parameter algorithm is strongly recommended. If it is not possible to use the 4-parameter algorithm for result calculation, it is possible to switch to a point-to-point calculation.

Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

9. LIMITATIONS

Samples with 25(OH)-vitamin D concentrations higher than the highest standard should be diluted e.g. 1+1 with standard 1 (= 0 nmol/l) (e.g. 50 µl sample + 50 µl standard 1) and re-assayed.

Whole blood is not suitable as a sample.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference ranges for 25(OH) vitamin D₃

Information from ASBMR 2011

Deficiency (seriously deficient)	< 20 ng/ml	or	< 50 nmol/l
Insufficiency (deficient)	20–29 ng/ml	or	50–74 nmol/l
Sufficiency (adequately supplied)	> 30 ng/ml	or	> 75 nmol/l

Society of Osteologists SACHSEN E. V.

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1 ng/ml = 2.5 nmol/l

1 nmol/l = 0.4 ng/ml

Reference intervals for 25(OH)-vitamin D₃ (ng/ml), males and females*

Age	n	2.5%	97.5%
0 to < 3 months	131	5	42
3 to < 6 months	135	9	60
6 months to < 1 year	147	18	58
1 to < 3 years	394	15	54
3 to < 10 years	619	14	46
10 to < 13 years	286	11	50
13 to < 15 years	275	10	44
15 to < 18 years	390	8	45
> 18 years	421	8	56

*All seasons

(Soldin et al., 2009)

Note

The vitamin D production in the skin is high variable and depends on the season and daily time, degree of latitude, age, sun protection etc. The normal ranges depend on the method used (e.g. vitamin D release from the vitamin D binding protein, VDBP) and serve only as orientation.

Literature references

The following literature references for the vitamin D reference values can be found in the references on page 28: Grant et al., Soldin et al., Visser et al., Wicherts et al.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)

Sample	25(OH)-vitamin D [nmol/l]	CV [%]
1	71.6	5.0
2	20.1	9.6

Inter-Assay (Day-to-Day variation; n = 20)

Sample	25(OH)-vitamin D [nmol/l]	CV [%]
1	33.3	12.1
2	48.5	9.2
3	79.6	8.6
4	109.9	8.8

Dilution recovery

Two samples were diluted with standard 1 (0 nmol/l) and used in the test. The results are shown in the following table.

Sample	Dilution	Observed [nmol/l]	Expected [nmol/l]	Recovery [%]
A	undiluted	80.7	80.7	
	75%	62.6	60.5	103.5
	50%	41.7	40.4	103.2
	25%	18.9	20.2	93.6
B	undiluted	100.8	100.8	
	75%	78.1	75.6	103.3
	50%	52.3	50.4	103.8
	25%	26.7	25.2	106.0

Spiking Recovery

This table shows the recovery rate of 25(OH) vitamin D₃ which was added to 3 different serum samples. Standard 5 (220 nmol/l) and standard 6 (600 nmol/l) were used for this test. All samples were measured in triplicates on three days. The table shows the average target values and obtained values.

Sample	Obtained value [nmol/l]	Target value [nmol/l]	Recovery [%]
A	57.5	56.5	101.8
	68.8	65.1	105.7
	79.1	75.5	104.8
	108.1	103.1	104.8
	161.2	158.4	101.8
B	33.2	32.2	103.1
	40.5	42.1	96.2
	47.5	51.2	92.8
	75.4	80.1	94.1
	129.4	137.9	93.8
C	36.6	35.6	102.8
	43.2	45.3	95.4
	52.8	54.6	96.7
	82.0	83.3	98.4
	142.8	140.7	101.5

Analytical Sensitivity

Limit of blank, LoB	1,08 ng/ml	2,7 nmol/l
Limit of detection, LoD	2,64 ng/ml	6,6 nmol/l
Limit of quantitation, LoQ	6,08 ng/ml	15,2 nmol/l
Measurement range	6–160 ng/ml	15–400 nmol/l

The evaluation was performed according to the CLSI guideline EP-17-A.

NIST standard reference material (SRM) 972 traceable.

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to 25(OH)-vitamin D₃. The specificity is calculated in percent, based on the cross-reactivity of these compounds with the anti-25(OH)-vitamin D₃ antibody compared to the 25(OH)-vitamin D₃ antigen.

- | | |
|---|----------|
| • 25(OH) vitamin D ₃ | 100,0% |
| • 25(OH) vitamin D ₂ | 67,8% |
| • 24, 25(OH) vitamin D ₃ | ≥ 100,0% |
| • Vitamin D ₂ (ergocalciferol) | 0,3% |

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- **Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions than sealed ones.**
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of the enclosed plate sealers during incubation steps is necessary.

- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).

7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxyvitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by