# ADMA ELISA Kit

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## 15. LITERATUR

- Vallance P, Leone A, Calver A, Collier J, Moncada S. Accumulation of an endogenous inhibitor of NO synthesis in chronic renal failure. *Lancet* 1992; **339**: 572 – 575
- Kielstein JT, Böger RH, Bode-Böger SM, et al. Asymmetric dimethylarginine plasma concentrations differ in patients with end-stage renal disease: Relationship to treatment method and atherosclerotic disease. *J. Am. Soc. Nephrol.* 1999; **10**: 594 – 600
- Böger RH, Bode-Böger SM, Szuba A, Tangphao O, Tsao PS, Chan JR, Blaschke TF, Cooke JP. Asymmetric dimethylarginine: a novel risk factor for endothelial dysfunction. Its role in hypercholesterolemia. *Circulation* 1998; **98**: 1842 – 1847
- Stühlinger M, Abbasi F, Chu JW, Lamendola C, McLaughlin TL, Cooke JP, Reaven GM, Tsao PS. Relationship between insulin resistance and an endogenous nitric oxide synthase inhibitor. *J. Am. Med. Assoc.* 2002; **287**: 1420-1426
- Zoccali C, Bode-Böger SM, Mallamaci F, Benedetto FA, Tripepi G, Malatino L, Cataliotti A, Bellanuova I, Fermo I, Frölich JC, Böger RH. Asymmetric dimethylarginine (ADMA): An endogenous inhibitor of nitric oxide synthase predicts mortality in end-stage renal disease (ESRD). *Lancet* 2001; **358**: 2113-2117
- Nijveldt RJ, Teerlink T, Van der Hoven B, Siroen MP, Kuik DJ, Rauwerda JA, van Leeuwen PA. Asymmetrical dimethylarginine (ADMA) in critically ill patients: high plasma ADMA concentration is an independent risk factor of ICU mortality. *Clin. Nutr.* 2003; **22**: 23-30
- Savvidou MD, Hingorani AD, Tsikas D, Frolich JC, Vallance P, Nicolaidis KH. Endothelial dysfunction and raised plasma concentrations of asymmetric dimethylarginine in pregnant women who subsequently develop pre-eclampsia. *Lancet* 2003; **361**: 1511-1517
- Böger RH. The emerging role of asymmetric dimethylarginine as a novel cardiovascular risk factor. *Cardiovasc. Res.* 2003; **59**: 824-833

- Lu TM, Ding YA, Lin SJ, Lee WS, Tai HC. Plasma levels of asymmetrical dimethylarginine and adverse cardiovascular events after percutaneous coronary intervention. *Eur Heart J.* 2003; **24**: 1912-1919.

### Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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Manual

# ADMA ELISA Kit

***For the determination of ADMA in human serum, citrate and EDTA plasma***



*EU: IVD / CE*

*US: Research Use Only. Not for use in diagnostic procedures.*

Valid from 07.07.2014

**K 7828**

**Immundiagnostik AG**

## Table of Contents

<b>1. INTENDED USE</b>	<b>19</b>
<b>2. INTRODUCTION / CLINICAL RELEVANCE</b>	<b>19</b>
<b>3. MATERIAL SUPPLIED</b>	<b>20</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>20</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>21</b>
<b>6. SPECIMEN COLLECTION AND PREPARATION</b>	<b>22</b>
<b>7. ASSAY PROCEDURE</b>	<b>22</b>
<i>Principle of the test</i>	22
<i>Sample preparation procedure</i>	23
<i>Test procedure</i>	24
<b>8. RESULTS</b>	<b>25</b>
<b>9. LIMITATIONS</b>	<b>26</b>
<b>10. QUALITY CONTROL</b>	<b>26</b>
<i>Reference Range</i>	27
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>27</b>
<i>Precision and reproducibility</i>	27
<i>Spiking Recovery</i>	27
<i>Dilution Recovery</i>	28
<i>Analytical Sensitivity</i>	28
<i>Specificity</i>	29
<i>Correlation with HPLC-MS</i>	29
<b>12. PRECAUTIONS</b>	<b>29</b>
<b>13. TECHNICAL HINTS</b>	<b>30</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>30</b>
<b>15. REFERENCES</b>	<b>31</b>

## 1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of asymmetric dimethyl arginine (ADMA) in human serum, citrate and EDTA plasma. For *in vitro* diagnostic use only.

## 2. INTRODUCTION / CLINICAL RELEVANCE

Asymmetric dimethyl arginine (ADMA) is an endogenous inhibitor of NO-synthase. It is formed during proteolysis of methylated proteins and removed by renal excretion or metabolic degradation by the enzyme dimethylargininedimethylaminohydrolase (DDAH). Several cell types, including human endothelial and tubular cells are capable of synthesizing and metabolizing ADMA. Elevated ADMA concentrations in the blood are found in numerous diseases associated with endothelial dysfunction. For example, elevated ADMA levels in blood of dialysis patients correlate significantly with the degree of arteriosclerosis and cardiovascular risk. Furthermore, elevated ADMA levels are found in patients with hypercholesterolemia, hypertension, arteriosclerosis, chronic renal failure and chronic heart failure, and are associated with restrictions in endothelial vasodilatation.

During the last years, the important clinical relevance of the regulation of vascular tone and structure by nitric oxide (NO) has been shown. Moreover, there were reports that human endothelial cells produce ADMA as well as nitric oxide, which points to an endogenous endothelial NO-regulation by ADMA. Therefore it was assumed that hypertension, arteriosclerosis and immunological dysfunction in patients with chronic renal failure are connected to a dysfunction of the L-arginine/NO-metabolism and to ADMA accumulation. The reasons for the deregulation of the L-arginine/NO-metabolism could only partially be elucidated. Certainly, there are multiple factors involved in the L-arginine/NO-metabolism regulation as for example elevation of free superoxide radicals ( $O_2^-$ ), ADMA accumulation and reduced NO-synthase activity.

Prospective clinical studies of the last years demonstrate the increased importance of ADMA as a novel cardiovascular risk factor.

### Indication

- Arteriosclerosis
- Hypertension
- Chronic heart failure
- Coronary artery disease

- Hypercholesterolemia
- Chronic renal failure
- Diabetes mellitus
- Peripheral arterial occlusive disease

### 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 7828MTP	PLATE	Holder with precoated strips	12 x 8 wells
K 7828ST	STD	Standards diluted in reaction buffer, ready to use (0, 0.1, 0.25, 0.5, 1.0, 2.0 $\mu$ M)	6 x 1 vial
K 7828KO1 K 7828KO2	CTRL 1 CTRL 2	Controls diluted in reaction buffer, ready to use	2 x 1 vial
K 7828WP	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 7828AK	AB	ADMA antibody, lyophilized	2 x 1 vial
K 7828K	CONJ	Conjugate (peroxidase-labeled), concentrate	1 x 120 $\mu$ l
K 7828CSP	CONJBUF	Conjugate stabilizing buffer, ready to use	1 x 24 ml
K 7828RP	DERBUF	Reaction buffer, ready to use	1 x 15 ml
K 7828DR	DER	Derivatization reagent	2 x 50 mg
K 7828LM	DMSO	Dimethylsulfoxide (DMSO)	1 x 7 ml
K 7828SL	CODIL	Dilution buffer after derivatization, ready to use	1 x 28 ml
K 7828TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 25 ml
K 7828AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Calibrated precision pipets and 10-1000  $\mu$ l tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker

- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm or 405 nm

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **ELISA wash buffer concentrate** (WASHBUF) with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The buffer concentrate is stable at 2-8°C until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month**.
- **Standards (STD)** and **Controls (CTRL1, CTRL2)** are already diluted in reaction buffer (DERBUF). Store Standards and Controls frozen at **-20°C**, thaw before use in the test, and re-freeze immediately after use. They can be re-frozen up to 3 times
- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (50 mg)** in **3 ml DMSO**. Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use**. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.

- Dissolve the content of one vial of **ADMA antibody (AB)** in **5.6 ml** of **diluted wash buffer**. For this, reconstitute the content of one AB vial with 0.6 ml of diluted wash buffer for 5 minutes. Then transfer quantitatively the obtained AB solution into a separate vial and add 5 ml of diluted wash buffer. The ELISA kit can be separated into two performances by providing two AB vials. Diluted ADMA antibody can be stored at **2-8°C for one month**.
- Dilute the **peroxidase conjugate (CONJ) 1:201** with conjugate stabilizing buffer (CONJBUF) (e.g. **110 µl CONJ + 22 ml CONJBUF**; prepare only the required amount). The undiluted POD conjugate is stable at 2-8°C until the expiry date stated on the label. Diluted POD conjugate can be stored at **2-8°C for 1 week**.
- All other test reagents are ready for use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 6. SPECIMEN COLLECTION AND PREPARATION

### Serum, citrate and EDTA plasma

- Venous fasting blood is suited for this test system. Samples are stable for one week at 2-8°C. For longer storage samples should be frozen at -20°C.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- Samples with visible amounts of precipitates should be centrifuged.
- The serum, citrate and EDTA plasma samples are analyzed **undiluted**.  
If the sample volume is less than 50 µl, we recommend a 1:2 dilution in DERBUF (reaction buffer) (25 µl sample + 25 µl DERBUF). This dilution factor must be considered in data evaluation.
- For sample preparation a derivatization reagent (DER) for derivatization of ADMA is added (details are given in the sample preparation procedure).

## 7. ASSAY PROCEDURE

### *Principle of the test*

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization-reagent for ADMA derivatization. Afterwards, the treated samples and the polyclonal ADMA-antiserum are incubated in wells of a

microtiter plate coated with ADMA-derivative (tracer). During the incubation period, the target ADMA in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The ADMA in the sample displaces the antibodies out of the binding to the tracer. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the ADMA concentration in the sample. During the second incubation step a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-ADMA antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow color is inverse proportional to the ADMA concentration in the sample; this means, high ADMA concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. ADMA present in the patient samples is determined directly from this curve.

### *Sample preparation procedure*

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials) as follows:

1. Bring <b>all reagents and samples to room temperature</b> (15-30°C) and mix well.
2. Add <b>200 µl</b> of ready to use <b>standards (STD)</b> , <b>200 µl</b> of ready to use <b>controls (CTRL)</b> and <b>50 µl of samples (SAMPLE)</b> in the corresponding vials.
3. Add <b>150 µl of reaction buffer (DERBUF) only</b> to the <b>samples (SAMPLE)</b> .
4. Add <b>50 µl</b> of freshly prepared <b>derivatization reagent (DER)</b> into each vial (STD, CTRL, SAMPLE), <b>mix thoroughly</b> and incubate for <b>45 min</b> at <b>room temperature</b> (15-30°C) on a horizontal <b>shaker</b> (180-240 rpm).
5. Afterwards add <b>250 µl of dilution buffer (CODIL)</b> into each vial, mix well and incubate for <b>45 min</b> at <b>room temperature</b> (15-30°C) on a horizontal <b>shaker</b> (180-240 rpm).

**2 x 100 µl of each treated sample (STD, CTRL, SAMPLE)** are used in the ELISA as duplicates.

### *Test procedure*

6.	Mark the positions of standards (STD)/ controls (CTRL)/samples (SAMPLE) in duplicate on a <b>protocol sheet</b> .
7.	Take as many <b>microtiter strips (PLATE)</b> as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
8.	Wash the microtiter strips <b>5x with 250 µl of diluted ELISA wash buffer before use</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
9.	For the analysis in duplicate take <b>2 x 100 µl</b> of the <b>derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE)</b> out of the vials and add into the respective wells.
10.	Add <b>100 µl</b> of diluted <b>ADMA antibody (AB)</b> into each well. Cover the plate tightly.
11.	Incubate <b>overnight (15-20 hours) at 2-8°C</b> .
12.	Aspirate the contents of each well. Wash each well <b>5 x with 250 µl of diluted wash buffer</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
13.	Add <b>200 µl</b> of diluted <b>peroxidase conjugate (CONJ)</b> into each well.
14.	Cover the plate tightly and incubate for <b>1 hour at room temperature (15-30°C)</b> on a horizontal shaker (180-240 rpm).
15.	Aspirate the contents of each well. Wash each well <b>5 x with 250 µl of diluted wash buffer</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
16.	Add <b>200 µl</b> of <b>TMB substrate (SUB)</b> into each well.

17. Incubate for **6-10 min** at **room temperature** (15-30°C) in the dark\*.
18. Add **100 µl** of **stop solution (STOP)** into each well, mix thoroughly.
19. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm (690 nm) as a reference.

\* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

## 8. RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls, and blood samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results**. Exception: Consider the dilution factor for 1:2 pre-diluted samples.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.


### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

The concentration of controls and patient samples can be determined directly from the calibration curve in  $\mu\text{mol/l}$ . In the following, an example of a calibration curve is given; do not use it for the calculation of your results.

### Example of calibration curve

## 9. LIMITATIONS

Samples with concentrations above the measurement range should be diluted in reaction buffer (DERBUF) and re-assayed.

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

*Reference range*

Based on internal studies with plasma samples of apparently healthy persons (n=70) a mean value of 0.45  $\mu\text{mol/l}$  was calculated. The standard deviation was 0.095  $\mu\text{mol/l}$ .

**Plasma mean value  $\pm$  2 x standard deviation:** **0.45  $\pm$  0.19  $\mu\text{mol/l}$**

**Normal range:** **0.26 – 0.64  $\mu\text{mol/l}$**

We recommend each laboratory to establish its own reference range.

**11. PERFORMANCE CHARACTERISTICS***Precision and reproducibility***Intra-assay (n=12)**

Sample	ADMA [ $\mu\text{mol/l}$ ]	CV [%]
1	0.27	7.8
2	0.78	5.3

**Inter-assay (n=6)**

Sample	ADMA [ $\mu\text{mol/l}$ ]	CV [%]
1	0.33	8.8
2	0.79	5.7

*Spiking recovery*

One serum sample was spiked with different ADMA concentrations and measured in this assay. The mean recovery rate for all concentrations was 104 % (n=5).

spike [ $\mu\text{mol/l}$ ]	ADMA expected [ $\mu\text{mol/l}$ ]	ADMA measured [ $\mu\text{mol/l}$ ]	recovery [%]
0		0.419	
0.25	0.669	0.694	104
0.5	0.919	0.948	103

### *Dilution recovery*

One spiked sample was diluted with DERBUF. The mean recovery was 103 % (n=10).

dilution	ADMA expected [ $\mu\text{mol/l}$ ]	ADMA measured [ $\mu\text{mol/l}$ ]	recovery [%]
original		0.468	
1:2	0.234	0.255	109
1:4	0.117	0.114	97

### *Analytical sensitivity*

The zero-standard was measured 6 times. The detection limit was set as  $B_0 - 1 \text{ SD}$  and estimated to be  $0.05 \mu\text{mol/l}$ .

Sample	ADMA mean value [OD]	Standard deviation (SD)	Detection limit [ $\mu\text{mol/l}$ ]
zero-standard	2.78	0.12	0.05


## Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to ADMA. The specificity is calculated in percent in relation to the ADMA binding activity.

L-Arginin	< 0.02 %
SDMA	< 0.5 %
NMMA	< 0.5 %

## Correlation with HPLC-MS

10 samples were measured with this ELISA and HPLC-MS. The correlation was  $r = 0.98$ .

## 12. PRECAUTIONS

- All reagents in the package are for *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.

- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

### **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch, as wells from already opened microtiter plates are exposed to different conditions than sealed ones.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure which is not coordinated with the producer may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

- Vallance P, Leone A, Calver A, Collier J, Moncada S. Accumulation of an endogenous inhibitor of NO synthesis in chronic renal failure. *Lancet* 1992; **339**: 572 – 575
- Kielstein JT, Böger RH, Bode-Böger SM, et al. Asymmetric dimethylarginine plasma concentrations differ in patients with end-stage renal disease: Relationship to treatment method and atherosclerotic disease. *J. Am. Soc. Nephrol.* 1999; **10**: 594 – 600
- Böger RH, Bode-Böger SM, Szuba A, Tangphao O, Tsao PS, Chan JR, Blaschke TF, Cooke JP. Asymmetric dimethylarginine: a novel risk factor for endothelial dysfunction. Its role in hypercholesterolemia. *Circulation* 1998; **98**: 1842 – 1847
- Stühlinger M, Abbasi F, Chu JW, Lamendola C, McLaughlin TL, Cooke JP, Reaven GM, Tsao PS. Relationship between insulin resistance and an endogenous nitric oxide synthase inhibitor. *J. Am. Med. Assoc.* 2002; **287**: 1420-1426
- Zoccali C, Bode-Böger SM, Mallamaci F, Benedetto FA, Tripepi G, Malatino L, Cataliotti A, Bellanuova I, Fermo I, Frölich JC, Böger RH. Asymmetric dimethylarginine (ADMA): An endogenous inhibitor of nitric oxide synthase predicts mortality in end-stage renal disease (ESRD). *Lancet* 2001; **358**: 2113-2117
- Nijveldt RJ, Teerlink T, Van der Hoven B, Siroen MP, Kuik DJ, Rauwerda JA, van Leeuwen PA. Asymmetrical dimethylarginine (ADMA) in critically ill patients: high plasma ADMA concentration is an independent risk factor of ICU mortality. *Clin. Nutr.* 2003; **22**: 23-30
- Savvidou MD, Hingorani AD, Tsikas D, Frolich JC, Vallance P, Nicolaidis KH. Endothelial dysfunction and raised plasma concentrations of asymmetric dimethylarginine in pregnant women who subsequently develop pre-eclampsia. *Lancet* 2003; **361**: 1511-1517
- Böger RH. The emerging role of asymmetric dimethylarginine as a novel cardiovascular risk factor. *Cardiovasc. Res.* 2003; **59**: 824-833
- Lu TM, Ding YA, Lin SJ, Lee WS, Tai HC. Plasma levels of asymmetrical dimethylarginine and adverse cardiovascular events after percutaneous coronary intervention. *Eur Heart J.* 2003; **24**: 1912-1919.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number