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# AOPP Kit

*Zur in vitro Bestimmung von Advanced oxidation protein products (AOPP) in EDTA-Plasma*

# AOPP Kit

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### 14. LITERATUR

1. Deschamps-Latscha B et al. (2005) Advanced oxidation protein products as risk factors for atherosclerotic cardiovascular events in nondiabetic predialysis patients. *Am J Kidney Dis* **45**(1):39-47
2. Deschamps-Latscha B et al. (2004) Early prediction of IgA nephropathy progression: proteinuria and AOPP are strong prognostic markers. *Kidney International* **66**: 1606-1612
3. Nguyen-Khoa T et al. (2001) Oxidative stress and haemodialysis: role of inflammation and duration of dialysis treatment. *Nephrol Dial Transplant* **16**: 335-340
4. Witko-Sarsat V et al. (1996) Advanced oxidation protein products as a novel marker of oxidative stress in uremia. *Kidney International* **49**: 1304-1313

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
# AOPP Kit

*For the in vitro determination of Advanced oxidation protein products (AOPP) in EDTA-plasma*

Gültig ab / Valid from 30.04.2012

**K 7811w**

## 1. INTENDED USE

The described Assay is intended for the quantitative determination of **advanced oxidation protein products (AOPPs)** in EDTA-plasma. It is for research use only.

## 2. INTRODUCTION

Increased oxidative stress has been implicated in a wide range of diseases. In haemodialysis patients, massive generation of reactive oxygen species (ROS) is induced during each dialysis session. Proteins are highly susceptible to oxidative stress damage, which could result in formation of **AOPPs (advanced oxidation protein products)** or AGEs (advanced glycation end products, non-enzymatically glycosylated proteins with irreversible chemical modifications). **AOPP** measurement is proposed to be a reliable marker for the extent of protein oxidative damage in uremic patients. In addition, plasma AOPP determination may be useful for monitoring the effect of treatments with drugs reducing oxidative stress.

### Indications

- Monitoring of oxidative stress, e.g. in hemodialyzed patients
- Monitoring of inflammatory processes
- Prognostic marker for progressive IgA-nephropathy

### 3. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 7811w MTP	PLATE	Holder with strips	12 x 8 wells
K 7811w PV	SAMPLEBUF	Sample dilution buffer	50 ml
K 7811w SV	STDBUF	Standard dilution buffer	15 ml
K 7811w ST	STDKONZ	Standard concentrate, lyophilized	4 vials
K 7811w KO	CTRL	Control, lyophilized	2 vials
K 7811w DEL	DELIP	Delipidation reagent	2.5 ml

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000  $\mu$ l
- Multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 340

## 5. PREPARATION AND STORAGE OF REAGENTS

- **SAMPLEBUF** (sample dilution buffer) and **STDBUF** (standard dilution buffer) are ready-to-use.
- **DELIP** (delipidation reagent) is ready-to-use.
- The lyophilized **STDKONZ** (standard concentrate) and the lyophilized **CTRL** (control) are stable at **2-8°C** until the expiry date stated on the label. **CTRL** (Control) and **STDKONZ** (standard concentrate) must be reconstituted with bidistilled water (volume and concentration, see respective product specification). Allow the vial content to solve for 10 minutes and then mix gently. Reconstituted standard and control are **not stable** and cannot be stored.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 6. SAMPLE PREPARATION

**Centrifugate** fresh collected EDTA-Plasma in 1.5 ml cups at 3 000 g for 30 sec before analysis

**Mix 125 µl centrifugated EDTA-Plasma** with **25 µl DELIP** (delipidation reagent), vortex

**Dilution 1:1.2**

**Incubate for 10 min** at room temperature

Afterwards, **centrifugate** at 3 000 g for 5 min

Mix **100 µl delipidated EDTA-Plasma** with **400 µl SAMPLEBUF** (sample dilution buffer) in 1.5 ml cup, vortex

**Final dilution 1:6**

## 7. ASSAY PROCEDURE

### *Principle of the test*

The assay is based on the spectroscopic analysis of modified proteins at 340 nm. Standards, controls and patient samples assayed for AOPP are placed in each well of a 96-well microtiter plate. The absorbance is read at 340 nm. The chloramine-T (CT) absorbance at 340 nm is linear within the range of 0 to 100  $\mu\text{mol/l}$ , AOPP concentrations are expressed as CT equivalents.

### *Test procedure*

1. Bring all <b>reagents and samples</b> to <b>room temperature</b> (18-26 °C) and mix well
2. Mark the positions of <b>STD/SAMPLE/CTRL</b> (Standards/Sample/Control) in duplicate on a protocol sheet
3. Add <b>200 <math>\mu\text{l}</math> STD</b> (100 $\mu\text{mol/l}$ ; 50 $\mu\text{mol/l}$ ; 25 $\mu\text{mol/l}$ ; 12.5 $\mu\text{mol/l}$ ; 6.25 $\mu\text{mol/l}$ ) in duplicate into respective well
4. Add <b>200 <math>\mu\text{l}</math> STDBUF</b> (Standard dilution buffer) as <b>standard 0 <math>\mu\text{mol/l}</math></b> in duplicate into respective well
5. Add <b>200 <math>\mu\text{l}</math> Sample</b> in duplicate into respective well
6. Add <b>200 <math>\mu\text{l}</math> CTRL</b> in duplicate into respective well
7. Determine directly the <b>absorption</b> of standards and samples at <b>340 nm</b>
8. Determine the <b>AOPP concentration</b> in diluted samples directly from the linear standard curve

## 8. RESULTS

The estimated AOPP value must be multiplied by **6** to obtain the concentration in the patient samples.

## 9. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### **Normal ranges**

We recommend each laboratory to establish its own norm concentration range.

## 10. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay-Variation**

The precision (intra-assay variation) was calculated from 16 determinations on each one of three samples.

Intra-Assay CV n= 16

<b>Sample</b>	<b>AOPP mean value [<math>\mu\text{mol/l}</math>]</b>	<b>Intra-Assay CV [%]</b>
1	28.7	5.6
2	171.9	1.3
3	522.3	3.0

### Inter-Assay-Variation

The total precision (inter-assay variation) was calculated from data on 2 samples measured by different technicians on different days.

Inter-Assay CV n= 12

Sample	AOPP mean value [µmol/l]	Inter-Assay CV [%]
1	84.0	14.3
2	72.5	16.6

## 11. PRECAUTIONS

- For research use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

## 12. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for research use only.
- Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

### 14. REFERENCES

1. Deschamps-Latscha B et al. (2005) Advanced oxidation protein products as risk factors for atherosclerotic cardiovascular events in nondiabetic predialysis patients. *Am J Kidney Dis* **45**(1):39-47
2. Deschamps-Latscha B et al. (2004) Early prediction of IgA nephropathy progression: proteinuria and AOPP are strong prognostic markers. *Kidney International* **66**: 1606-1612
3. Nguyen-Khoa T et al. (2001) Oxidative stress and haemodialysis: role of inflammation and duration of dialysis treatment. *Nephrol Dial Transplant* **16**: 335-340
4. Witko-Sarsat V et al. (1996) Advanced oxidation protein products as a novel marker of oxidative stress in uremia. *Kidney International* **49**: 1304-1313

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number

For research use only