α 1-Antitrypsin ELISA Kit

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15. LITERATUR

1. Amarri, S. et al., 2006. Changes of gut microbiota and immune markers during the complementary feeding period in healthy breast-fed infants. *Journal of pediatric gastroenterology and nutrition*, **42**(5), pp.488–95.
2. Faust, D et al., 2001. Determination of alpha1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line. *Zeitschrift für Gastroenterologie*, **39**(9), pp.769–74.
3. Faust, D. et al., 2002. Regulation of alpha1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells. *Clinical and experimental immunology*, **128**(2), pp.279–84.
4. Hsu, P.-I. et al., 2010. Diagnosis of gastric malignancy using gastric juice alpha1-antitrypsin. *Cancer epidemiology, biomarkers & prevention*, **19**(2), pp.405–11.
5. Lamprecht, M. et al., 2012. Probiotic supplementation affects markers of intestinal barrier, oxidation, and inflammation in trained men; a randomized, double-blinded, placebo-controlled trial. *Journal of the International Society of Sports Nutrition*, **9**(1), p.45.
6. Muss, C., Schütz, B. & Kirkamm, R., 2002. Alpha1-Antitrypsin - ein objektiver Verlaufsparemeter bei entzündlichen Darmerkrankungen. *Ärztzeitschrift für Naturheilverfahren*, 43(4).
7. Oswari, H. et al., 2013. Comparison of stool microbiota compositions, stool al-

- pha1-antitrypsin and calprotectin concentrations, and diarrhoeal morbidity of Indonesian infants fed breast milk or probiotic/prebiotic-supplemented formula. *Journal of paediatrics and child health*, Epub ahead of print.
8. Quint, J.K. et al., 2011. SERPINA1 11478G>A variant, serum α_1 -antitrypsin, exacerbation frequency and FEV1 decline in COPD. *Thorax*, **66**(5), pp.418–24.
 9. Ragab, H.M. et al., 2009. Clinical utility of serum TNF alpha and alpha-1 anti-tryptsin in predicting the stage and progression of lung cancer. *International Journal of Integrative Biology*, **7**(1), pp.45–52.
 10. Roeckel, N. et al., 2009. High frequency of LMAN1 abnormalities in colorectal tumors with microsatellite instability. *Cancer research*, **69**(1), pp.292–9.
 11. Strygler, B. et al., 1990. α_1 -Antitrypsin Excretion in Stool in Normal Subjects and in Patients With Gastrointestinal Disorders. *Gastroenterology*, **99**(5), pp.1380–1387.
 12. Török, E. et al., 2011. Primary human hepatocytes on biodegradable poly(l-lactic acid) matrices: a promising model for improving transplantation efficiency with tissue engineering. *Liver transplantation*, **17**(2), pp.104–14.

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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α 1-Antitrypsin ELISA Kit


For the in vitro determination of α ₁-antitrypsin in stool

EU: IVD / CE

US: Research Use Only. Not for use in diagnostic procedures.

Valid from 04.09.2014

REF **K 6760**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. STORAGE AND PREPARATION OF REAGENTS	19
6. STORAGE AND PREPARATION OF SAMPLES	20
<i>Storage</i>	20
<i>Extraction of the stool samples</i>	20
<i>Dilution of samples</i>	21
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Test procedure</i>	22
8. RESULTS	23
9. LIMITATIONS	23
10. QUALITY CONTROL	24
<i>Reference range</i>	24
11. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Spiking Recovery</i>	25
<i>Dilution recovery</i>	25
<i>Analytical Sensitivity</i>	25
<i>Specificity</i>	26
12. PRECAUTIONS	26
13. TECHNICAL HINTS	26
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	27
15. REFERENCES	27

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of α_1 -antitrypsin in stool. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Intestinal protein loss is a serious consequence of various systemic or local gastro-intestinal diseases (e.g. allergies, chronic inflammation, malignancies). These pathologies damage the mucosal integrity and/or cause lymphostasis, thereby leading to an increased transfer of plasma proteins into the bowel lumen. Subsequently, hypoproteinemia accompanied with edema may develop. This condition is diagnosed by exclusion of other sources of protein loss and by proof of an elevated α_1 -antitrypsin concentration in stool.

In serum, α_1 -antitrypsin represents the majority of serine protease inhibitors and protects tissues from protease damages during inflammation. The protein is synthesized primarily in the liver but also to a small extent in intestinal macrophages, monocytes, and intestinal epithelial cells. Since α_1 -antitrypsin is relatively resistant against enzymatic digestion, the secreted amount in stool reflects the internal concentration of the protein. An elevated α_1 -antitrypsin stool concentration is therefore a widely recognized marker for intestinal protein loss and for an increased mucosal permeability.

In clinical routine, the α_1 -antitrypsin clearance (ratio of the α_1 -antitrypsin ELISA values of stool and serum samples) has been established along with the sole determination of the 24h α_1 -antitrypsin secretion in stool. Thus the group of J. S. Fordtran reports that the sole determination of the α_1 -antitrypsin concentration in stool yielded false positive or false negative results in 21 % of the patients compared to the α_1 -antitrypsin clearance measurement (Strygler et al. 1990).

The analytical quality of Immundiagnostik's α_1 -antitrypsin ELISA surpasses by far the conventional radial immunodiffusion (RID) technique in the determination of serum, stool and tissue culture supernatants. In direct comparison, the concentrations measured with the ELISA were approximately 30% above the corresponding RID levels. Cell culture supernatants of an intestinal cell line as well as fecal samples of lymphostasis patients yielded negative results with RID. Our ELISA could detect α_1 -antitrypsin in all of these samples, in some of them even in very high concentrations.

These results clearly prove that the α_1 -antitrypsin ELISA is far more sensitive than the conventional method and that it recognizes not only hepatic but also enteral α_1 -antitrypsin. The discrepancy of both methods and hence the superiority of the ELISA to RID is especially striking in the analysis of extremely high enteral protein losses. The combination of two specific antibodies in our α_1 -antitrypsin ELISA widely exclu-

des the possibility of false negative results thereby enabling a reliable diagnostics of enteral protein loss.

Indications

- Suspected enteric protein loss
- Crohn's disease
- Necrotic enterocolitis
- Chronic mesenterial ischemia
- Viral, bacterial, allergic, or autoimmune-induced gastrointestinal inflammation

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6760MTP	PLATE	Holder with precoated strips	12 x 8 wells
K 6760WB	WASHBUF	ELISA wash concentrate 10x	2 x 100 ml
K 6760K	CONJ	Conjugate, (goat-anti α_1 -antitrypsin, peroxidase-labelled)	200 μ l
K 6760CAL	CAL	Calibrator, ready to use	1 vial
K 6760KO1	CTRL	Control, ready to use	1 vial
K 6760KO2	CTRL	Control, ready to use	1 vial
K 6760TMB	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 6760AC	STOP	ELISA stop solution, ready to use	1 x 15 ml
K 6760EP	IDK Extract®	Extraction buffer concentrate <i>IDK Extract® 2.5x</i>	2 x 100 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 μ l tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex

- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)
 - * Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 μm) with an electrical conductivity of 0.055 $\mu\text{S}/\text{cm}$ at 25 °C ($\geq 18.2 \text{ M}\Omega\text{cm}$).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month**.
- The **extraction buffer concentrate** *IDK Extract*[®] must be diluted with ultra pure water **1:2.5** before use (100 ml concentrate + 150 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37 °C in a water bath. The **extraction buffer concentrate** *IDK Extract*[®] is stable at **2–8 °C** until the expiry date stated on the label. Diluted **buffer solution** (extraction buffer) can be stored in a closed flask at **2–8 °C for three months**.
- Use **100 μl** of **wash buffer** (diluted WASHBUF) as **BLANK**. Pipette into the respective well.
- The **calibrator** (STD) and **controls** (CTRL) are stable at **2–8 °C** until the expiry date stated on the label.
- The **conjugate concentrate** (CONJ) must be diluted **1:101 in wash buffer** (100 μl CONJ + 10 ml wash buffer). The concentrate is stable at **2–8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Storage

The **sample stability** is as follows:

Raw stool: over night at room temperature (15–30°C), 3 days at 2–8°C or at least 4 weeks at -20°C

Stool extracts: 9 days at room temperature, 2–8°C or -20°C, maximum 3 freeze-thaw cycles

Extraction of the stool samples

Diluted extraction buffer *IDK Extract*[®] is used as a sample extraction buffer. We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instructions for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

SAS with 1.5 ml extraction buffer:

Applied amount of stool: 15 mg

Buffer Volume: 1.5 ml

Dilution Factor: 1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty sample tube** with **1.5 ml** of ready to use *IDK Extract*[®] extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert the orange dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.

- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~ 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

Dilution I: 1:100

Dilution of samples

The supernatant of the sample preparation procedure (dilution I) is further diluted **1:250 in wash buffer**. For example:

20 μ l supernatant (dilution I) + **980 μ l** wash buffer, mix well = **dilution II (1:50)**

200 μ l dilution II + **800 μ l** wash buffer, mix well = **dilution III (1:5)**

This results in a **final dilution of 1:25 000**.

For analysis, pipet **100 μ l** of **dilution III** per well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the “sandwich” technique with two selected polyclonal antibodies that bind to human α -1-antitrypsin.

Calibrator, controls and prediluted samples which are assayed for human α -1-antitrypsin are added into the wells of a micro plate coated with a high affine polyclonal anti-human α -1-antitrypsin antibody. During the first incubation step, α -1-antitrypsin is bound by the immobilized antibody. Then a peroxidase-conjugated polyclonal anti-human α -1-antitrypsin antibody is added into each microtiter well and a “sandwich” of capture antibody - human α -1-antitrypsin – peroxidase-conjugate is formed. Tetramethyl–benzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of α -1-antitrypsin. Samples are quantified by referring their optical density to a Lot-dependant master calibration curve and the use of a Calibrator that is run with each test.

Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

Mark the positions of CAL (calibrator), BLANK (blank), SAMPLE (sample) CTRL (controls) on a **protocol sheet**.

We recommend to carry out the tests in duplicate.

1.	Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
2.	Add 100 µl of CAL (calibrator), BLANK (blank), SAMPLE (sample) and CTRL (controls) into respective well.
3.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.
4.	Discard the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
5.	Add 100 µl CONJ (conjugate) into each well.
6.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.
7.	Discard the contents of each well. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
8.	Add 100 µl of SUB (substrate) into each well.
9.	Incubate for 10–20 minutes at room temperature (15–30°C) in the dark*.
10.	Add 100 µl of STOP (stop solution) into each well, mix thoroughly.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A,B,C and D. In both cases the optical density of the calibrator (CAL) is essential. Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Caution: Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Stool samples

Multiply the obtained results by the dilution factor of 25 000 to get the real concentration.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the calibration curve \times sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB \times sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on Immundiagnostik studies of stool samples of apparently healthy persons (n = 76), a **cut off value of < 26,8 mg/dl** was estimated.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)

The precision (intra-assay variation) of the Immundiagnostik α_1 -antitrypsin ELISA test was calculated from 20 replicate determinations on each one of two samples.

Sample	α_1 -antitrypsin [mg/dl]	CV [%]
1	15.2	4.5
2	42.2	13.1

Inter-Assay (n = 20)

The total precision (inter-assay variation) of the Immundiagnostik α_1 -antitrypsin ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Sample	α_1 -antitrypsin [mg/dl]	CV [%]
1	16.15	9.8
2	54.46	14.8

Spiking Recovery

Two samples were spiked with different α_1 -antitrypsin standard amounts and measured with the assay.

Sample	Unspiked Sample [$\mu\text{g/l}$]	Spike [$\mu\text{g/l}$]	α_1 -antitrypsin expected [$\mu\text{g/l}$]	α_1 -antitrypsin measured [$\mu\text{g/l}$]
A	6.3	1.65	7.95	7.2
		5.0	11.3	11.4
		15.0	21.3	20.9
		45.0	51.3	46.7
B	5.6	1.65	7.3	7.0
		5.0	10.6	10.9
		15.0	20.6	17.5
		45.0	50.6	46.7

Dilution recovery

Two patient serum samples were diluted with wash buffer. The results are shown below ($n = 2$):

Sample	Dilution	α_1 -antitrypsin expected [mg/dl]	α_1 -antitrypsin measured [mg/dl]
A	1:12500	48.0	48.88
	1:25000	24.5	23.25
	1:50000	12.3	11.9
	1:100000	6.1	6.0
B	1:12500	158.4	158.4
	1:25000	79.3	99.0
	1:50000	39.6	33.0
	1:100000	19.8	22.1

Analytical Sensitivity

The Zero-standard was measured 20 times. The detection limit was set as $B_0 + 2 \text{ SD}$ and estimated to be 1.8 mg/dl.

Specificity

No cross reactivity with other plasma proteins in stool was observed.

No cross reactivity with alpha-1-antitrypsin in mouse serum was observed.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Amarri, S. et al., 2006. Changes of gut microbiota and immune markers during the complementary feeding period in healthy breast-fed infants. *Journal of pediatric gastroenterology and nutrition*, **42**(5), pp.488–95.
2. Faust, D et al., 2001. Determination of alpha1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line. *Zeitschrift für Gastroenterologie*, **39**(9), pp.769–74.
3. Faust, D. et al., 2002. Regulation of alpha1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells. *Clinical and experimental immunology*, **128**(2), pp.279–84.
4. Hsu, P.-I. et al., 2010. Diagnosis of gastric malignancy using gastric juice alpha1-antitrypsin. *Cancer epidemiology, biomarkers & prevention*, **19**(2), pp.405–11.
5. Lamprecht, M. et al., 2012. Probiotic supplementation affects markers of intestinal barrier, oxidation, and inflammation in trained men; a randomized, double-blinded, placebo-controlled trial. *Journal of the International Society of Sports Nutrition*, **9**(1), p.45.
6. Muss, C., Schütz, B. & Kirkamm, R., 2002. Alpha1-Antitrypsin - ein objektiver Verlaufsparemeter bei entzündlichen Darmerkrankungen. *Ärztzeitschrift für Naturheilverfahren*, 43(4).
7. Oswari, H. et al., 2013. Comparison of stool microbiota compositions, stool alpha1-antitrypsin and calprotectin concentrations, and diarrhoeal morbidity of Indonesian infants fed breast milk or probiotic/prebiotic-supplemented formula.

Journal of paediatrics and child health, Epub ahead of print.

8. Quint, J.K. et al., 2011. SERPINA1 11478G>A variant, serum α_1 -antitrypsin, exacerbation frequency and FEV1 decline in COPD. *Thorax*, **66**(5), pp.418–24.
9. Ragab, H.M. et al., 2009. Clinical utility of serum TNF alpha and alpha-1 anti-tryptsin in predicting the stage and progression of lung cancer. *International Journal of Integrative Biology*, **7**(1), pp.45–52.
10. Roedel, N. et al., 2009. High frequency of LMAN1 abnormalities in colorectal tumors with microsatellite instability. *Cancer research*, **69**(1), pp.292–9.
11. Strygler, B. et al., 1990. α_1 -Antitrypsin Excretion in Stool in Normal Subjects and in Patients With Gastrointestinal Disorders. *Gastroenterology*, **99**(5), pp.1380–1387.
12. Török, E. et al., 2011. Primary human hepatocytes on biodegradable poly(l-lactic acid) matrices: a promising model for improving transplantation efficiency with tissue engineering. *Liver transplantation*, **17**(2), pp.104–14

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number