Adiponectin (total) ELISA Kit

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15. LITERATUR

1. Blüher (Blueher), M., et al. (2007) Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care* 30(2): p. 280-5
2. Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelities in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
3. Fernandez-Real, J.M., et al. (2004) Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care* 27(3): p. 739-45
4. Halperin F, Beckman JA, Patti ME, Trijillo ME, Garvin M, Craeger MA, Scherer PE and Goldfine AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154
5. Higashiura, K., et al. (2004) Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)* 61(6): p. 753-9
6. Koenig W, Khuseynova N, Baumert J, Meisinger C and Löwel H (2006) Serum Concentrations of Adiponectin and Risk of Type 2 Diabetes Mellitus and Coronary Heart Disease in Apparently Healthy Middle-Aged Men: Results From the 18-Year Follow-Up of a Large Cohort From Southern Germany. *JACC* 48:1369-1377

7. Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. *Diabetes* 55:249-259
8. Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. *Clin Chem* 50:1511-1525
9. Nakamura, Y., et al. (2004) Implications of plasma concentrations of adiponectin in patients with coronary artery disease. *Heart* 90(5): p. 528-33
10. Pajvani, U.B., et al. (2003) Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. *J Biol Chem* 278(11): p. 9073-85
11. Shibata, R., et al. (2004) Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. *J Biol Chem* 279(27): p. 28670-4

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
Adiponectin (total) ELISA Kit

*For the in vitro determination of human adiponectin
in serum and plasma*

Valid from 2014-11-27

REF K 6250

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

| | |
|---|-----------|
| 1. INTENDED USE | 15 |
| 2. INTRODUCTION | 15 |
| 3. MATERIAL SUPPLIED | 16 |
| 4. MATERIAL REQUIRED BUT NOT SUPPLIED | 16 |
| 5. STORAGE AND PREPARATION OF REAGENTS | 17 |
| 6. STORAGE AND PREPARATION OF SAMPLES | 17 |
| 7. ASSAY PROCEDURE | 18 |
| <i>Principle of the test</i> | 18 |
| <i>Test procedure</i> | 18 |
| 8. RESULTS | 19 |
| 9. LIMITATIONS | 20 |
| 10. QUALITY CONTROL | 20 |
| <i>Reference range</i> | 20 |
| 11. PERFORMANCE CHARACTERISTICS | 21 |
| <i>Precision and reproducibility</i> | 21 |
| <i>Spiking Recovery</i> | 21 |
| <i>Dilution recovery</i> | 22 |
| <i>Analytical Sensitivity</i> | 22 |
| 12. PRECAUTIONS | 22 |
| 13. TECHNICAL HINTS | 23 |
| 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE | 23 |
| 15. REFERENCES | 24 |

1. INTENDED USE

The described ELISA is intended for the quantitative determination of human adiponectin in serum and plasma. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Adiponectin consists of 244 amino acids and has an approximate molecular weight of 30 kDa. The adiponectin concentration in blood is 5–30 µg/ml that accounts for about 0.01% of the total plasma proteins. Adiponectin has four domains: a signal peptide, a variable domain, a collagen-like N-terminal domain and a globular C-terminal domain. Adiponectin exists in different oligomer forms *in vivo*. Beside the trimer and dimer also high molecular multimers are present. Adiponectin is mainly expressed by adipocytes, but also muscle cells and hepatocytes can synthesize it. The only, until now known, natural inductor of the adiponectin synthesis is IGF-I. Two different adiponectin receptors are known. Both of them are ubiquitarily expressed, though their distribution in the tissues varies: Adiponectin Receptor 1 (AdipoR1) is synthesized in muscle and AdipoR2 in liver tissue. The adiponectin significance for the human organism is not completely clear until now. Studies show, that adiponectin correlates negatively with BMI and could participate at the energy metabolism through the regulation of fatty acid oxidation. Adiponectin influences further physiological processes like angiogenesis. It is also associated with glucose and lipid metabolism. Adiponectin levels are associated with insulin resistance and accordingly linked with Type II Diabetes. Furthermore, it is involved in inflammatory processes and therewith of importance for development of arteriosclerosis and coronary artery diseases. Blüher et al. (2007) showed a correlation of total serum adiponectin to insulin sensitivity and the ability of total serum adiponectin levels to predict insulin resistance and impaired glucose tolerance. Low adiponectin concentrations result in inhibition of fatty acid oxidation and are associated with insulin resistance and metabolic syndrome as well as arteriosclerosis.

Indications

- Energy metabolism and body weight regulation
- Metabolic syndrome
- Type 2 diabetes
- Coronary artery disease
- Atherosclerosis

3. MATERIAL SUPPLIED

| Cat. No. | Label | Kit components | Quantity |
|----------|-----------|--|--------------|
| K 6250 | PLATE | Precoated microtiter plate | 12 x 8 wells |
| K 6250 | WASHBUF | ELISA wash buffer concentrate 10x | 2 x 100 ml |
| K 6250 | SAMPLEBUF | Sample dilution buffer | 2 x 100 ml |
| K 6250 | STD | Adiponectin standards, lyophilized (0; 1.4; 5.5; 22; 88 ng/ml) | 2 x 5 vials |
| K 6250 | CTRL | Control, lyophilized (see specification for range) | 2 x 1 vial |
| K 6250 | CTRL | Control, lyophilized (see specification for range) | 2 x 1 vial |
| K 6250 | CONJ | Conjugate, ready to use | 15 ml |
| K 6250 | SUB | TMB substrate (Tetramethylbenzidine), ready to use | 15 ml |
| K 6250 | STOP | ELISA stop solution, ready to use | 15 ml |

For reorders of single components, use the catalogue number followed by the label as purchase order number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8°C for one month**.
- The **lyophilized STD (standards) and CTRL (controls)** are stable at **2–8°C** until the expiry date stated on the label. Reconstitution details are given in the **data sheet**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

6. STORAGE AND PREPARATION OF SAMPLES

Serum and Plasma

Samples are diluted **1:1000 in SAMPLEBUF** (sample dilution buffer) before use. E. g.:

Dilution I:

25 µl sample + 975 µl SAMPLEBUF = 1:40

Dilution II:

40 µl of dilution I + 960 µl SAMPLEBUF = 1:25

For analysis, pipette 100 µl of dilution II per well.

Samples can be stored for **one year at -20°C** and up to **three times frozen and thawed**.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the two-site “sandwich” technique with one monoclonal and one polyclonal antibody that bind to human Adiponectin.

Standards, controls and diluted patient samples which are assayed for human Adiponectin are added to wells of microplate coated with a high affine monoclonal anti-human Adiponectin antibody. During the first incubation step, Adiponectin in the samples is bound by the immobilized antibody. Then a peroxidase labeled conjugate is added to each well and the following complex is formed: capture antibody - human Adiponectin – Peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the Adiponectin concentration of sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standard. Adiponectin present in the patient samples, is determined directly from this curve.

Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

| | |
|----|---|
| 1. | Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper. |
| 2. | Add 100 µl of STD (standard) SAMPLE (sample) CTRL (controls) into respective well. |
| 3. | Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker. |
| 4. | Discard the contents of each well. Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper. |
| 5. | Add 100 µl CONJ (conjugate) into each well. |

| | |
|-----|---|
| 6. | Cover the plate tightly and incubate for 1 hour at room temperature (15–30 °C) on a horizontal shaker. |
| 7. | Discard the contents of each well. Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper. |
| 8. | Add 100 µl of SUB (substrate) into each well. |
| 9. | Incubate for 10–20 min.* at room temperature (15–30 °C) in the dark. |
| 10. | Add 100 µl of STOP (stop solution) into each well, mix. |
| 11. | Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference. |

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform.

For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Serum and plasma samples

To obtain the final Adiponectin concentration of the samples, multiply the estimated value by the dilution factor of **1000**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

Analytical sensitivity × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

A significant correlation between adiponectin serum values and age as well as gender of the probands is reported, in turn the correlation between the respective BMI seems to be less significant.

Serum and Plasma (n = 80): 8.79 µg/ml

Based on Immundiagnostik studies of apparently healthy individuals (n= 80) a mean value of 8.79 µg/ml stool was determined.

| Company | Adiponectin mean value* [µg/mL] | | |
|------------|------------------------------------|----------------------------|-----------------|
| | Normal glucose tolerance | Impaired glucose tolerance | Type 2 diabetes |
| LINCO | 8.95 ± 0.55 | 3.38 ± 0.26 | 3.48 ± 0.42 |
| Mediagnost | 8.81 ± 3.43 | 3.51 ± 1.47 | 3.82 ± 2.16 |

*Values from Reference 1.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 30)

| Sample | Adiponectin [ng/ml] | VK [%] |
|--------|---------------------|--------|
| 1 | 5.25 | 3.2 |
| 2 | 5.92 | 3.4 |

Inter-Assay (n = 12)

| Sample | Adiponectin [ng/ml] | VK [%] |
|--------|---------------------|--------|
| 1 | 6.9 | 6.3 |
| 2 | 8.4 | 6.1 |

Spiking Recovery

| Sample | Unspiked sample [ng/ml] | Spike [ng/ml] | Adiponectin expected [ng/ml] | Adiponectin measured [ng/ml] |
|--------|-------------------------|---------------|------------------------------|------------------------------|
| A | 3.8 | 2 | 5.8 | 6.1 |
| | 3.8 | 4 | 7.8 | 7.6 |
| | 3.8 | 6 | 9.8 | 9.2 |

| Sample | Unspiked sample [ng/ml] | Spike [ng/ml] | Adiponectin expected [ng/ml] | Adiponectin measured [ng/ml] |
|--------|-------------------------|---------------|------------------------------|------------------------------|
| B | 2.1 | 5 | 7.1 | 8.1 |
| | 2.1 | 7 | 9.1 | 9.6 |
| | 2.1 | 9 | 11.1 | 10.8 |

Dilution recovery

| Sample | Dilution | Adiponectin expected [ng/ml] | Adiponectin measured [ng/ml] |
|--------|----------|------------------------------|------------------------------|
| A | 1 : 600 | 7.65 | 7.65 |
| | 1 : 1200 | 3.83 | 4.14 |
| | 1 : 2400 | 1.91 | 2.33 |
| | 1 : 4800 | 0.96 | 1.23 |
| B | 1 : 300 | 7.60 | 7.60 |
| | 1 : 600 | 3.80 | 3.87 |
| | 1 : 1200 | 1.90 | 2.11 |
| | 1 : 2400 | 0.95 | 1.17 |

Analytical Sensitivity

The Zero-standard was measured 20 times. The detection limit was set as $B_0 + 2 \text{ SD}$ and estimated to be 0.338 ng/ml.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Blüher (Blueher), M., et al. (2007) Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care* 30(2): p. 280-5
2. Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelties in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
3. Fernandez-Real, J.M., et al. (2004) Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care* 27(3): p. 739-45
4. Halperin F, Beckman JA, Patti ME, Trijillo ME, Garvin M, Craeger MA, Scherer PE and Goldfine AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154
5. Higashiura, K., et al. (2004) Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)* 61(6): p. 753-9
6. Koenig W, Khuseynova N, Baumert J, Meisinger C and Löwel H (2006) Serum Concentrations of Adiponectin and Risk of Type 2 Diabetes Mellitus and Coronary Heart Disease in Apparently Healthy Middle-Aged Men: Results From the 18-Year Follow-Up of a Large Cohort From Southern Germany. *JACC* 48:1369-1377
7. Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. *Diabetes* 55:249-259
8. Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. *Clin Chem* 50:1511-1525
9. Nakamura, Y., et al. (2004) Implications of plasma concentrations of adiponectin in patients with coronary artery disease. *Heart* 90(5): p. 528-33
10. Pajvani, U.B., et al. (2003) Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. *J Biol Chem* 278(11): p. 9073-85
11. Shibata, R., et al. (2004) Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. *J Biol Chem* 279(27): p. 28670-4

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by