Albumin ELISA Kit

Super sensitiv

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
Albumin ELISA Kit

Super sensitive

For the in vitro determination of Albumin in urine and stool

Valid from 2015-01-14

REF K 6330

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	15
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. STORAGE AND PREPARATION OF REAGENTS	16
6. STORAGE AND PREPARATION OF SAMPLES	17
<i>Extraction of the stool samples</i>	17
<i>Dilution of samples</i>	18
7. ASSAY PROCEDURE	18
<i>Principle of the test</i>	18
<i>Test procedure</i>	19
8. RESULTS	20
9. LIMITATIONS	21
10. QUALITY CONTROL	21
<i>Reference range</i>	22
11. PERFORMANCE CHARACTERISTICS	22
<i>Precision and reproducibility</i>	22
<i>Analytical Sensitivity</i>	22
<i>Specificity</i>	22
12. PRECAUTIONS	23
13. TECHNICAL HINTS	23
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
15. REFERENCES	24

1. INTENDED USE

This Immundiagnostik assay is an enzyme immunoassay intended for the quantitative determination of albumin in urine and stool. For *in vitro* diagnostic.

2. INTRODUCTION

Albumin is the major protein in human plasma (40-60%). It is synthesized in the liver depending of the protein uptake. Albumin in fecal samples refers to an inflammatory reaction combined with intestinal bleeding. Elevated levels of albumin and hemoglobin in stool are found not only by colorectal carcinomas, but also with polyps and during chronic inflammatory diseases (Morbus Crohn or Colitis Ulcerosa).

Indications

- Detection of source of bleeding in the lower gastrointestinal tract
- Detection of colorectal carcinoma
- Investigation of high risk patients
- Morbus Crohn, Colitis Ulcerosa

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6330	PLATE	Precoated strips	96 wells
K 6330	WASHBUF	ELISA wash buffer concentrate (10x)	1 x 100 ml
K 6330	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 100 ml
K 6330	CONJBUF	Conjugate dilution buffer, ready to use	1 x 15 ml
K 6330	IDK Extract®	Extraction buffer concentrate <i>IDK Extract</i> ®, 2.5x	1 x 100 ml
K 6330	CONJ	Conjugate, (rabbit-anti-Albumin peroxidase-labeled)	1 x 50 µl
K 6330	STD	Calibrators, lyophilized (0, 12.5, 50, 200, 800 µg/l)	4 x 5 vials
K 6330	CTRL	Control, lyophilized	4 x 1 vial
K 6330	CTRL	Control, lyophilized	4 x 1 vial
K 6330	SUB	TMB-Substrat (Tetramethylbenzidin), gebrauchsfertig	1 x 15 ml
K 6330	STOP	ELISA-Stopplösung, gebrauchsfertig	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month**.
- The **extraction buffer concentrate** *IDK Extract*[®] must be diluted with ultra pure water **1:2.5** before use (100 ml concentrate + 150 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37 °C in a water bath. The **extraction buffer concentrate** *IDK Extract*[®] is stable at **2–8 °C** until the expiry date stated on the label. Diluted **buffer solution** (extraction buffer) can be stored in a closed flask at **2–8 °C for three months**.
- The **lyophilized standards** (STD) and **controls** (CTRL) are stable at **2–8 °C** until the expiry date stated on the label. Before use, the standards and controls

must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted standards and controls are not stable and cannot be stored.**

- The **conjugate concentrate** (CONJ) must be diluted **1:401 in conjugate dilution buffer** (25 µl CONJ + 10 ml CONJBUF). The concentrate is stable at **2–8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Urine

Adjust urine to a pH 6 to 8 with 1 N NaOH. Samples can be stored for two weeks at 2–8 °C or at -20 °C for longtime storage.

For the assay, dilute samples **1:200** with **SAMPLEBUF** (sample dilution buffer), e.g. 10 µl sample + 1990 µl SAMPLEBUF.

Extraction of the stool samples

Diluted extraction buffer *IDK Extract*[®] is used as a sample extraction buffer. We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instructions for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

SAS with 1.5 ml extraction buffer:

Applied amount of stool:	15 mg
Buffer Volume:	1.5 ml
Dilution Factor:	1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, ino-

cultulation loop or similar device.

- b) Fill the **empty sample tube** with **1.5 ml** of ready to use *IDK Extract*[®] extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert the orange dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~ 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

Dilution I: 1:100

Dilution of samples

The supernatant of the sample preparation procedure (dilution I) is diluted **1:2.5 in sample dilution buffer** (SAMPLEBUF). For example:

- **200 µl** supernatant (dilution I) + **300 µl** sample dilution buffer = **1: 2.5 (dilution II)** .

This results in a final dilution of **1:250**.

For analysis, pipet **100 µl** of **dilution II** per well.

7. ASSAY PROCEDURE

Principle of the test

This Enzyme-Linked Immunosorbent Assay (ELISA) is a two step assay for the ultra sensitive determination of human albumin in stool and urine. A polyclonal rabbit antibody specific for human albumin is immobilized on a microtiter plate and a second anti-albumin antibody is conjugated to peroxidase.

In a first incubation step, the albumin in the samples is bound to the immobilized anti-albumin antibodies. A washing step is carried out to remove all unbound substances. In a second incubation step, a peroxidase-labeled anti-albumin antibody is added. After another washing step, to remove all unbound substances, a peroxidase-substrate, tetramethylbenzidine, is added. The enzymatic reaction is terminated by an acidic stop solution, whereby the color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of albumin in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated using the results obtained from the calibrators. Albumin in the patient samples is determined directly from this curve.

Test procedure

Bring all reagents and samples to room temperature (15–30°C) and mix well.

Mark the positions of STD (Standard) SAMPLE (Sample) CTRL (Controls) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (wash-buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
2.	Add 100 µl of STD (standard) SAMPLE (sample) CTRL (controls) into respective well.
3.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.
4.	Discard the contents of each well. Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
5.	Add 100 µl diluted CONJ (conjugate) into each well.
6.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.

7.	Discard the contents of each well. Wash each well 5 x by dispensing 250 µl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
8.	Add 100 µl of SUB (substrate) into each well.
9.	Incubate for 10–20 min* at room temperature (15–30 °C) in the dark.
10.	Add 100 µl of STOP (stop solution) into each well, mix.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform.

For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the

duplicate values should be evaluated manually.

Urine

The obtained albumin levels of urine samples have to be multiplied with the dilution factor of **200**.

Stool

The obtained albumin levels of stool samples have to be multiplied with the dilution factor of **250**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used to get the real concentration.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this higher dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

Analytical sensitivity × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on Immundiagnostik studies of stool samples of apparently healthy persons (n = 76), the following reference range was estimated.

Albumin in stool: < 9.2 mg/l

Albumin in urine: < 20 mg/l*

* Test information, Labor Limbach, Heidelberg, <http://www.labor-limbach.de>

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)

The precision (intra-assay variation) was calculated from 20 replicate determinations on one sample.

Sample	Albumin [mg/l]	CV [%]
1	50	5

Inter-Assay (n = 20)

The total precision (inter-assay variation) of was calculated from data on one sample obtained in 20 different assays.

Sample	Albumin [mg/l]	CV [%]
1	50	8

Analytical Sensitivity

The Zero-standard was measured 20 times. The detection limit was set as $B_0 + 2 \text{ SD}$ and estimated to be 12.5 ng/ml.

Specificity

No cross reactivity to other plasma proteins was observed.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Trasch and Bloch: 1993; Clin. Lab. 39, 479-484
2. John et al.: 1994; Clin. Lab. 40, 77-81

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by