# CRP ELISA Kit

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## 15. LITERATUR

1. Doumatey, Ayo P, Jie Zhou, Adebowale Adeyemo, and Charles Rotimi. 2014. "High Sensitivity C-Reactive Protein (Hs-CRP) Remains Highly Stable in Long-Term Archived Human Serum." *Clinical Biochemistry* **47** (4-5) (March): 315–8. doi:10.1016/j.clinbiochem.2013.12.014.
2. Koenig, Wolfgang, Hannelore Löwel, Jens Baumert, and Christa Meisinger. 2004. "C-Reactive Protein Modulates Risk Prediction Based on the Framingham Score: Implications for Future Risk Assessment: Results from a Large Cohort Study in Southern Germany." *Circulation* **109** (11) (March 23): 1349–53. doi:10.1161/01.CIR.0000120707.98922.E3.
3. Pearson, T. a. 2003. "American Heart Association Guide for Improving Cardiovascular Health at the Community Level: A Statement for Public Health Practitioners, Healthcare Providers, and Health Policy Makers From the American Heart Association Expert Panel on Population and Pre." *Circulation* **107** (4) (February 4): 645–651. doi:10.1161/01.CIR.0000054482.38437.13.
4. Ridker, P M, C H Hennekens, J E Buring, and N Rifai. 2000. "C-Reactive Protein and Other Markers of Inflammation in the Prediction of Cardiovascular Disease in Women." *The New England Journal of Medicine* **342** (12) (March 23): 836–43. doi:10.1056/NEJM200003233421202.
5. Salzer, Jonatan, Göran Hallmans, Maria Nyström, Hans Stenlund, Göran Wadell, and Peter Sundström. 2013. "Vitamin A and Systemic Inflammation as Protective

Factors in Multiple Sclerosis." *Multiple Sclerosis (Houndmills, Basingstoke, England)*  
**19** (8) (July 18): 1046–51. doi:10.1177/1352458512472752.

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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# CRP ELISA Kit


*For the in vitro determination of C-reactive protein in  
serum, plasma, stool and urine*

*EU: IVD / CE*

*US: Research Use Only. Not for use in diagnostic procedures.*

Valid from 26.06.2014

**REF** **K 9710S**

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>17</b>
<b>2. INTRODUCTION</b>	<b>17</b>
<b>3. MATERIAL SUPPLIED</b>	<b>17</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>18</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>19</b>
<b>7. ASSAY PROCEDURE</b>	<b>21</b>
<i>Principle of the test</i>	21
<i>Test procedure</i>	21
<b>8. RESULTS</b>	<b>22</b>
<b>9. LIMITATIONS</b>	<b>23</b>
<b>10. QUALITY CONTROL</b>	<b>24</b>
<i>Reference range</i>	24
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>24</b>
<i>Precision and reproducibility</i>	24
<i>Spiking Recovery</i>	25
<i>Dilution recovery</i>	25
<i>Analytical Sensitivity</i>	26
<i>Specificity</i>	26
<b>12. PRECAUTIONS</b>	<b>26</b>
<b>13. TECHNICAL HINTS</b>	<b>27</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>27</b>
<b>15. REFERENCES</b>	<b>28</b>

## 1. INTENDED USE

The Immundiagnostik Assay is a sandwich Enzyme Immuno Assay intended for the quantitative determination of C-reactive Protein in plasma, serum, stool and urine. It is for *in-vitro* diagnostic use only.

## 2. INTRODUCTION

C-reactive Protein (CRP) is mainly formed in hepatocytes. The synthesis rate of CRP is influenced by the cytokines involved in the inflammatory processes. The biological half-life time is estimated to be 13-16 hours. The serum CRP concentration reflects very sensitive acute fever, pneumonia and myocardial infraction.

Recent studies describe an association between inflammatory reactions and cardiovascular diseases like arteriosclerosis or latent and chronic persistent infections. As a marker for inflammation, CRP high-sensitive can be used to predict the risk of myocardial infarction and stroke.

The CRP determination in urine using the high sensitive ELISA method allows - together with  $\alpha_2$ -Macroglobulin - an early screening diagnose after kidney transplantations. The CRP kit provides an easy to use assay for monitoring anti-rejection therapies. The ELISA-assay was used for years to test hundreds of patients. Its predictive diagnostic value was compared with the gold standard (kidney biopsy).

### Indications

- Prognosis factor for myocardial infarction or stroke
- Inflammatory processes

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9710s MTP	PLATE	Micro titer plate with precoated strips	12 x 8 wells
K 9710s WP	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K 9710s K	CONJ	Conjugate, (rabbit-anti-CRP-antibody, peroxidase-labeled)	1 x 150 $\mu$ l
K 9710s ST	STD	Standards*, ready to use, (0; 1.9; 5.6; 16.7; 50; 150 ng/ml)	6 x 1 ml
K 9710s KO1	CTRL	Control, ready to use	1 x 1 ml
K 9710s KO2	CTRL	Control, ready to use	1 x 1 ml
K 9710s PV	SAMPLEBUF	Sample buffer, ready to use	2 x 100 ml

Cat. No.	Label	Kit components	Quantity
K9710s TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K9710s AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

\*The CRP calibrators were standardized against WHO standard 470.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month.**

- The **conjugate concentrate** (CONJ) must be diluted **1:101 in wash buffer** (100 µl CONJ + 10 ml wash buffer). The concentrate is stable at **2–8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### Serum, plasma

#### *Collection and storage*

**Collection of serum:** Collect sufficient blood (at least 1 ml) by venipuncture into a tube or a plastic syringe, avoid hemolysis, allow to stand for 15 min, centrifuge for 15 min at 1,000 x g and 4 °C and collect the serum.

**Collection of plasma:** Collect sufficient blood (at least 1 ml) by venipuncture into an EDTA venipuncture tube or a plastic syringe, allow to stand for 15 min, centrifuge for 15 min at 1,000 x g and 4 °C and separate the plasma from the cells.

**Storage:** Serum samples can be stored at -80 °C for 11 years.\*

\*A. P. Doumatey et al. 2014.

#### *Sample dilution*

Serum and plasma samples have to be diluted **1:100 or 1:500** before performing the assay.

Add **10 µl serum /plasma to 990 µl SAMPLEBUF** (sample buffer), mix well. **(1:100)**

Use the dilution factor (100 or 500) to calculate the CRP concentration read off the calibration curve.

Patient's samples with **elevated CRP-concentrations** must be diluted **1:4000 – 1:8000**. Samples of other patient collectives must be diluted according to the expected CRP-concentration. The corresponding dilution factor must be used for calculation of the CRP-concentration.

### Urine

Urine samples must be diluted **1:5 with SAMPLEBUF** (sample buffer).

### Stool

#### *Storage*

The samples should be refrigerated and can be stored at **2-8°C for 2 days**. If the test

cannot be performed within this period, the specimen should be stored at  $-20^{\circ}\text{C}$  or colder.

### ***Extraction***

**Diluted WASHBUF** (washbuffer) is used as a sample extraction buffer. We recommend the following sample preparation:

### **Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)**

#### ***Stool sample tube – Instructions for use***

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

#### ***SAS with 0,75 ml extraction buffer:***

Applied amount of stool: 15 mg

Buffer Volume: 0,75 ml

Dilution Factor: 1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty sample tube** with **0,75 ml** of ready to use **washbuffer** before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert the orange dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for  $\sim 10$  minutes improves the result.
- e) Allow sample to stand for  $\sim 10$  minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus

the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

**Dilution:** 1:50

The extract can be stored **1 month at -20°C**. Before analyzing it in the test, the suspension should be centrifuged for **2 min at 13,000 rpm** in a 1.5 ml reaction tube.

**100 µl** per well of this supernatant is used in the assay.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This Enzyme Immuno Assay is a sandwich assay for the determination of CRP in serum, plasma, urine and stool samples. The wells of the microtiter plate are coated with polyclonal antibodies directed against C-reactive Protein. In a first incubation step, the CRP in the samples is bound to the coated polyclonal rabbit antibodies (in excess). To remove all unbound substances, a washing step is carried out.

In a second incubation step, a peroxidase-labeled antibody (polyclonal, rabbit-anti-CRP) is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine. An acidic stopping solution is then added. The color converts to yellow. The intensity of the yellow color is directly proportional to the concentration of CRP in the sample. A dose response curve of the absorbance (at 450 nm) unit vs. concentration is generated. CRP, present in the patient samples, is determined directly from this calibration curve.

The combination of two specific antibodies in the CRP ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Mark the positions of STD (Standard) SAMPLE (Sample) CTRL (Controls) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Wash each well <b>5 x</b> by dispensing <b>250 µl of diluted WASHBUF</b> (wash-buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
----	---

2.	Add <b>100 µl of STD</b> (Standard) <b>SAMPLE</b> (Sample) <b>CTRL</b> (Controls) into respective well.
3.	Cover the plate tightly and incubate for <b>1 hour</b> at room temperature (15-30°C) on a horizontal shaker.
4.	Discard the contents of each well. Wash each well <b>5 x</b> by dispensing <b>250 µl of diluted WASHBUF</b> (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
5.	Add <b>100 µl diluted CONJ</b> (conjugate) into each well.
6.	Cover the plate tightly and incubate for <b>1 hour</b> at room temperature (15-30°C) on a horizontal shaker.
7.	Discard the contents of each well. Wash each well <b>5 x</b> by dispensing <b>250 µl of diluted WASHBUF</b> (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
8.	Add <b>100 µl of SUB</b> (substrate) into each well.
9.	Incubate for <b>10 - 20 min.*</b> at room temperature (15-30°C) in the dark.
10.	Add <b>100 µl of STOP</b> (stop solution) into each well, mix.
11.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero

standard must be specified with a value less than 1 (e.g. 0.001).

## 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

## 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

## Serum, plasma

The value read from the calibration curve must be multiplied by **100 or 500** respectively to get the CRP concentration in serum/plasma samples.

If samples were diluted **1:4000** or **1:8000**, the estimated values must be multiplied by 4000 or 8000 respectively.

## Urine

The measured CRP concentration must be multiplied by factor **5** to get the actual concentration of the samples.

## Stool

The measured CRP concentration must be multiplied by factor **50** to get the actual concentration of the samples.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

## 9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

The lower limit of the measurement range can be calculated as:

*detection limit × sample dilution factor to be used*

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

CRP-Concentration, **serum/plasma\***:

< 1 mg/l low CHD-Risk

1-3 mg/l medium CHD-Risk

> 3 mg/l high CHD-Risk

\*Pearson et al., 2003.

CRP-Concentration, **stool**: < 56 ng/ml

CRP-Concentration, **urine**: < 6 ng/ml

If the CRP-concentration is found to be higher than 3 mg/l, a second determination should be made within 2 to 3 weeks. If the CRP-concentration is again high, and other reasons are excluded (acute infection, chronic-inflammatory diseases), the obtained CRP-concentration can be used for risk stratification in coronary heart disease (CHD) patients. If the CHD risk is high, the lifestyle should be changed together with medical treatment. These normal ranges should be used as a guideline only.

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay (n = 20)**

The precision (intra-assay variation) of the Immundiagnostik CRP ELISA test was calculated from 20 replicate determinations on each one of two samples.

Sample	CRP [ng/ml]	CV [%]
1	23.3	6
2	99.4	5.5

**Inter-Assay (n = 15)**

The total precision (inter-assay variation) of the Immundiagnostik CRP ELISA test was calculated from data on 2 samples obtained in 15 different assays by three technicians on two different lots of reagents over a period of three months.

Sample	CRP [ng/ml]	CV [%]
1	22,1	11,6
2	90,4	13,8

**Spiking Recovery**

Two samples with known CRP concentration were spiked with four different amounts of CRP and measured.

Sample	Unspiked Sample [ng/ml]	Spike [ng/ml]	CRP expected [ng/ml]	CRP measured [ng/ml]
1	9.8	37.5	47.3	44.5
		18.8	28.6	27.3
		9.4	19.2	18.2
		4.7	14.5	14.3
2	9.3	37.5	46.8	48.2
		18.8	28.1	26.3
		9.4	18.7	18.0
		4.7	14.0	13.7

**Dilution recovery**

Two patient samples were diluted and analyzed. The results are shown below (n = 2):

Sample	Dilution	CRP expected [ng/ml]	CRP measured [ng/ml]	Recovery [%]
A	1:100	2.90	2.88	99.3
	1:200	1.45	1.55	106.8
	1:400	0.73	0.83	113.7
	1:800	0.36	0.39	108.3
	1:1600	0.18	0.18	100.0

Sample	Dilution	CRP expected [ng/ml]	CRP measured [ng/ml]	Recovery [%]
B	1:200	10.80	10.80	100.0
	1:400	5.40	5.80	107.4
	1:800	2.70	2.90	107.4
	1:1600	1.35	1.61	119.3
	1:3200	0.68	0.83	122.1
	1:6400	0.33	0.35	106.1

### Analytical Sensitivity

The Zero-standard was measured 18 times. The detection limit was set as  $B_0 + 2 \text{ SD}$  and estimated to be 0,921 ng/ml.

### Specificity

No cross reactivity to other serum proteins was observed.

Alpha-1-Antitrypsin	0 %
Lysozym	0 %
Albumin	0 %
Other acute phase proteins	0 %

No cross reactivity with CRP in mouse serum was observed.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any

spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

### **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Doumatey, Ayo P, Jie Zhou, Adebowale Adeyemo, and Charles Rotimi. 2014. "High Sensitivity C-Reactive Protein (Hs-CRP) Remains Highly Stable in Long-Term Archived Human Serum." *Clinical Biochemistry* **47** (4-5) (March): 315–8. doi:10.1016/j.clinbiochem.2013.12.014.
2. Koenig, Wolfgang, Hannelore Löwel, Jens Baumert, and Christa Meisinger. 2004. "C-Reactive Protein Modulates Risk Prediction Based on the Framingham Score: Implications for Future Risk Assessment: Results from a Large Cohort Study in Southern Germany." *Circulation* **109** (11) (March 23): 1349–53. doi:10.1161/01.CIR.0000120707.98922.E3.
3. Pearson, T. a. 2003. "American Heart Association Guide for Improving Cardiovascular Health at the Community Level: A Statement for Public Health Practitioners, Healthcare Providers, and Health Policy Makers From the American Heart Association Expert Panel on Population and Pre." *Circulation* **107** (4) (February 4): 645–651. doi:10.1161/01.CIR.0000054482.38437.13.
4. Ridker, P M, C H Hennekens, J E Buring, and N Rifai. 2000. "C-Reactive Protein and Other Markers of Inflammation in the Prediction of Cardiovascular Disease in Women." *The New England Journal of Medicine* **342** (12) (March 23): 836–43. doi:10.1056/NEJM200003233421202.
5. Salzer, Jonatan, Göran Hallmans, Maria Nyström, Hans Stenlund, Göran Wadell, and Peter Sundström. 2013. "Vitamin A and Systemic Inflammation as Protective Factors in Multiple Sclerosis." *Multiple Sclerosis (Houndmills, Basingstoke, England)* **19** (8) (July 18): 1046–51. doi:10.1177/1352458512472752.

### Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number