Casomorphin ELISA Kit

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1. INTENDED USE	16
2. INTRODUCTION	16
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	17
5. MATERIAL REQUIRED BUT NOT SUPPLIED	18
6. PREPARATION AND STORAGE OF REAGENTS	18
7. PRECAUTIONS	19
8. SAMPLE AND TEST PREPARATION	19
9. ASSAY PROCEDURE	20
<i>TEST PROCEDURE</i>	20
10. EVALUATION OF RESULTS	22
<i>EXPECTED VALUES</i>	23
11. PERFORMANCE CHARACTERISTICS	24
<i>CROSS REACTIVITY</i>	24
<i>PRECISION AND REPRODUCIBILITY</i>	24
<i>SENSITIVITY</i>	24
<i>LINEARITY</i>	25
12. LIMITATIONS	25
13. REFERENCES	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26

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13. LITERATUR

- Dohan FC. (1973) Coeliac disease and schizophrenia. *Br Med J.* Jul 7;3(5870):51-2
- Haileselassie SS, Lee BH, Gibbs BF. (1999) Purification and identification of potentially bioactive peptides from enzyme-modified cheese. *J Dairy Sci.* Aug; 82(8):1612-7
- Hole K, Lingjaerde O, Mørkrid L, Bøler JB, Saelid G, Diderichsen J, Ruud E, Reichelt KL. Attention deficit disorders: a study of peptide-containing urinary complexes. (1988) *J Dev Behav Pediatr.* Aug; 9(4):205-12
- Jinsmaa Y, Yoshikawa M. (1999) Enzymatic release of neocasomorphin and beta-casomorphin from bovine beta-casein. *Peptides.* 20(8):957-62
- Kampa M, Loukas S, Hatzoglou A, Martin P, Martin PM, Castanas E. (1996) Identification of a novel opioid peptide (Tyr-Val-Pro-Phe-Pro) derived from human alpha S1 casein (alpha S1-casomorphin, and alpha S1-casomorphin amide). *Biochem J.* Nov 1;319 (Pt 3):903-8
- Kaplan BJ, McNicol J, Conte RA, Moghadam HK (1989) Dietary replacement in preschool-aged hyperactive boys. *Pediatrics.* Jan; 83(1):7-17
- Nygaard E, Reichelt KL, Fagan JF. (2001) The relation between the psychological functioning of children with Down syndrome and their urine peptide levels and levels of serum antibodies to food proteins. *Downs Syndr Res Pract.* Jul; 6(3):139-45
- Reichelt KL, Knivsberg AM. Can the pathophysiology of autism be explained by the nature of the discovered urine peptides? (2003) *Nutr Neurosci.* Feb; 6(1):19-28
- Sun Z, Cade R. (2003) Findings in normal rats following administration of gliadorphin-7 (GD-7). *Peptides.* Feb; 24(2):321-3
- Teschemacher H, Koch G, Brantl V. (1997) Milk protein-derived opioid receptor ligands. *Biopolymers.* 43(2):99-117
- Trompette A, Claustre J, Caillon F, Jourdan G, Chayvialle JA, Plaisancié P (2003) Milk bioactive peptides and beta-casomorphins induce mucus release in rat jejunum. *J Nutr.* Nov; 133(11):3499-503
- White JF. (2003) Intestinal pathophysiology in autism. *Exp Biol Med (Maywood).* Jun; 228(6):639-49
- Zioudrou C, Streaty RA, Klee WA. Opioid peptides derived from food proteins. The exorphins. *J Biol Chem.* 1979 Apr 10;254(7):2446-9

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Manual



Casomorphin ELISA Kit

For the determination of Casomorphin in urine

Valid from 11.07.2013

K 7010

1. INTENDED USE

This ELISA Kit is intended for the determination of Casomorphin in urine. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Casomorphin is a 7 amino acids peptide derived from the milk protein casein. Casein-derived peptides bind to opioid receptors in the brain and exhibit morphine-like effects, for example like heroin. These compounds have been shown to react with areas of the brain, which are involved in speech and auditory integration.

Urine samples from people with autism and schizophrenia contained high amounts of casomorphin. It is suspected that this peptide may also be elevated in other disorders such as chronic fatigue, fibromyalgia, and depression. Symptom remission has been observed after exclusion of dairy products from the diet.

Indication

- Autism
- Schizophrenia

3. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. Sample dilution buffer is used for sample preparation. Afterwards, the diluted samples and a polyclonal Casomorphin-antiserum are incubated in microtiter plate wells coated with a Casomorphin-derivative (tracer). During the incubation, the target Casomorphin in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. Casomorphin in the sample displaces the antibodies bound to the tracer. Therefore, the concentration of the tracer-bound antibodies is inverse proportional to the Casomorphin concentration in the sample. During the second incubation step, a peroxidase-conjugated antibody, which binds to the polyclonal anti-Casomorphin antibodies, is added to each microtiter well. After washing the unbound components, the peroxidase substrate tetramethylbenzidine (TMB) is added. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow color is inverse proportional to the Casomorphin concentration in the sample; this means high Casomorphin

concentration in the sample reduces the concentration of the tracer-bound antibodies and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. Casomorphin present in the patient samples is determined directly from this curve.

The ELISA results are normalized to the creatinine concentration of the urine sample. For this reason, a parallel determination of the creatinine concentration is required.

4. MATERIAL SUPPLIED

Catalog No	Content	Kit Components	Quantity
K7010MTP	PLATE	One holder with precoated strips	12 x 8 wells
K7010ST	STD	Standards, lyophilized	6 x 1 vial
K7010KO1 K7010KO2	CTRL 1 CTRL 2	Controls, lyophilized	2 x 1 vial
K7010WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K7010CSP	2.ABDIL	Conjugate stabilizing buffer, ready to use	24 ml
K7010AK	AB	Casomorphin antibody, lyophilized	2 x 1 vial
K7010K	2.AB	POD antibody (concentrate 200x)	120 µl
K7010DB	DIL	Dilution buffer, ready to use	50 ml
K7010VR	ABBUF	Antibody dilution buffer (lyophilized)	2 x 1 vial
K7010TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	25 ml
K7010AC	STOP	ELISA stop solution, ready to use	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm or 405nm

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 μ m) with an electrical conductivity of < 0.055 μ S/cm at 25°C (\geq 18.2 M Ω cm).

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μ l** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The WASHBUF is stable at 2-8°C until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month.**
- **Standards (STD)** and **controls (CTRL)** are lyophilized and must be reconstituted in **1 ml of dilution buffer (DIL)** per vial. Put the vials on a horizontal shaker for 5 min. Store the reconstituted STD and CTRL frozen at **-20°C**, re-freeze immediately after use. They can be re-frozen up to 2 times.
- The **antibody dilution buffer (ABBUF)** must be reconstituted in **5.6 ml of diluted wash buffer** per vial (total volume 2 x 5.6 ml = 11.2 ml).

- The **Casomorphin antibody (AB)** is lyophilized and must be reconstituted in **5.5 ml of antibody dilution buffer (ABBUF)** per vial. Put the vial on a horizontal shaker for 5 min. The ELISA kit can be separated into two performances by providing two AB vials. When more than one vial is to be used, combine the reconstituted solutions in a separate vial and mix prior to use. Diluted casomorphin antibody can be stored at **-20°C for one month**.
- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL) (e.g. **110 µl 2.AB + 22 ml 2.ABDIL**, prepare only the required amount). Undiluted POD antibody is stable at 2-8°C until the expiry date stated on the label. Diluted POD antibody can be stored at **2-8°C for 1 week**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at 2-8°C.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SAMPLE AND TEST PREPARATION

- Urine is suited for this test system (recommendation: early morning urine).
- Samples should be sent cooled; but they are stable for 24 h at room temperature.
- Samples must be **diluted** at least **1:10 in DIL** (dilution buffer) prior to analyses (see test procedure).
- Samples are stable for two days at 2-8°C. For longer storage samples should be frozen at -20°C.

9. ASSAY PROCEDURE

Procedural notes

- The assay should always be performed according to the enclosed manual.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The ELISA results are normalized to the creatinine concentration of the urine sample. For this reason, a parallel determination of the creatinine concentration is required.

Test procedure

Dilute samples with **dilution buffer (DIL)** by **factor 1:10**, e.g. **50 µl sample + 450 µl DIL** (dilution buffer).

If the measurement range is exceeded, dilute samples **1:100**, e.g. **10µl sample + 990 µl DIL**.

1. Bring all reagents and samples to room temperature (15-30°C)
2. Mark the positions of standards (STD)/ controls (CTRL)/ samples (SAMPLE) in duplicate on a protocol sheet.
3. Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
4. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5. For the analysis in duplicate, pipette 2 x 100 µl of reconstituted standards (STD)/ controls (CTRL)/ diluted samples (SAMPLE) into the respective wells of the microtiter plate (PLATE).

6.	Add 100 µl reconstituted casomorphin antibody (AB) into each well, cover the plate tightly.
7.	Incubate overnight (15-20 hours) at 2-8°C
8.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
9.	Add 200 µl of diluted POD antibody (2.AB) into each well.
10.	Cover the plate and incubate for 1 hour at room temperature (15-30°C) on a horizontal shaker (180-240 rpm)
11.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
12.	Add 200 µl of TMB substrate (SUB) into each well.
13.	Incubate for 5-15 min at room temperature (15-30°C) in the dark*.
14.	Add 100 µl of stop solution (STOP) into each well, mix thoroughly.
15.	Determine absorption immediately with an ELISA reader at 450 nm . If the highest extinction of the standards (STD) is above the range of the photometer, absorption must be measured immediately at 405 nm and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

10. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls and samples and with correct sample treatment) standards, controls and samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results.**

For samples with 1:100 dilution, the results have to be multiplied by 10.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

Plausibility of the measured pairs of values should be examined before automatically evaluating the results. If this option is not available within the used program, the pairs of values should be controlled manually.

The ELISA results are normalized to the creatinine concentration of the urine sample.

$$\text{Concentration}_{\text{Sample}} [\text{ng}/\mu\text{mol Creatinine}] = \frac{\text{Casomorphin Concentration}_{\text{Sample}} [\text{ng/ml}]}{\text{Creatinine Concentration}_{\text{Sample}} [\text{mmol/l}]}$$


Controls

Control samples should be analyzed with each run. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve in ng/ml if a dilution of 1:10 has been used.

In the following an example of a calibration curve is given. Do not use it for the calculation of your results

Example of calibration curve

Expected values

Based on internal studies of apparently healthy persons (n=223) a normal range of

0 - 0.8 ng casomorphin / μ mol creatinine

was estimated.

We recommend each laboratory to develop its own normal range. The values mentioned above are indicative only and can deviate from other published data.

11. PERFORMANCE CHARACTERISTICS

Cross reactivity

Gliadorphin: No cross reactivity was observed with gliadorphin at a concentration up to 1000 ng/ml in urine.

Casein: No cross reactivity was observed with casein at a concentration up to 10 µg/ml in urine.

Precision and reproducibility

Intra-Assay (n=18)		
Sample	Casomorphin [ng/ml]	Coefficient of variation (CV) [%]
1	3.5	13.8
2	13	11.1

Inter-Assay (n=8)		
Sample	Casomorphin [ng/ml I]	Coefficient of variation (CV) [%]
1	3.7	14.4
2	12.5	10.3

Sensitivity

The detection limit was set as $B_0 - 1SD$. The zero-standard was measured 12 times.

Sample	Casomorphin mean value [OD]	Standard Deviation (SD)	Detection limit [ng/ml]
Zero standard	2.7	0.07	0.5

Linearity

The linearity of the ELISA was determined by the dilution of a urine sample spiked with 24 ng/ml casomorphin. The mean linearity was 127% (n=6).

Dilution	Measured [ng/ml]	Expected [ng/ml]	Recovery [%]
1:2	14.3	12	119
1:4	8.1	6	135

12. LIMITATIONS

Casomorphin can only be determined in urine samples.

13. REFERENCES

- Dohan FC. (1973) Coeliac disease and schizophrenia. *Br Med J.* Jul 7;3(5870):51-2
- Haileselassie SS, Lee BH, Gibbs BF. (1999) Purification and identification of potentially bioactive peptides from enzyme-modified cheese. *J Dairy Sci.* Aug; 82(8):1612-7
- Hole K, Lingjaerde O, Mørkrid L, Bøler JB, Saelid G, Diderichsen J, Ruud E, Reichelt KL. Attention deficit disorders: a study of peptide-containing urinary complexes. (1988) *J Dev Behav Pediatr.* Aug; 9(4):205-12
- Jinsmaa Y, Yoshikawa M. (1999) Enzymatic release of neocasomorphin and beta-casomorphin from bovine beta-casein. *Peptides.* 20(8):957-62
- Kampa M, Loukas S, Hatzoglou A, Martin P, Martin PM, Castanas E. (1996) Identification of a novel opioid peptide (Tyr-Val-Pro-Phe-Pro) derived from human alpha S1 casein (alpha S1-casomorphin, and alpha S1-casomorphin amide). *Biochem J.* Nov 1;319 (Pt 3):903-8
- Kaplan BJ, McNicol J, Conte RA, Moghadam HK (1989) Dietary replacement in preschool-aged hyperactive boys. *Pediatrics.* Jan; 83(1):7-17
- Nygaard E, Reichelt KL, Fagan JF. (2001) The relation between the psychological functioning of children with Down syndrome and their urine peptide levels and levels of serum antibodies to food proteins. *Downs Syndr Res Pract.* Jul; 6(3):139-45

Reichelt KL, Knivsberg AM. Can the pathophysiology of autism be explained by the nature of the discovered urine peptides? (2003) *Nutr Neurosci.* Feb; 6(1):19-28

Sun Z, Cade R. (2003) Findings in normal rats following administration of gliadorphin-7 (GD-7). *Peptides.* Feb; 24(2):321-3

Teschemacher H, Koch G, Brantl V. (1997) Milk protein-derived opioid receptor ligands. *Biopolymers.* 43(2):99-117

Trompette A, Claustre J, Caillon F, Jourdan G, Chayvialle JA, Plaisancié P (2003) Milk bioactive peptides and beta-casomorphins induce mucus release in rat jejunum. *J Nutr.* Nov; 133(11):3499-503

White JF. (2003) Intestinal pathophysiology in autism. *Exp Biol Med (Maywood).* Jun; 228(6):639-49

Zioudrou C, Streaty RA, Klee WA. Opioid peptides derived from food proteins. The exorphins. *J Biol Chem.* 1979 Apr 10;254(7):2446-9

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the test package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the results of the test. Immundiagnostik AG can, therefore, not be held reliable for any damage resulting from this.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number