---

CERTEST  
*Clostridium*  
*perfringens*

---

ONE STEP *Clostridium perfringens*  
CARD TEST

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## ENGLISH

### INTENDED USE

CerTest *Clostridium perfringens* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Clostridium perfringens* in stool samples.

CerTest *Clostridium perfringens* card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium perfringens* infection.

### INTRODUCTION

*Clostridium perfringens* (*C. perfringens*) is a Gram-positive spore bearing anaerobic bacterium. This bacteria is distributed in nature, especially in soil, and in the intestines of many animals' species including the humans.

*C. perfringens* can develop a variety of exotoxins (11 histotoxins and one enterotoxin) and depending on the production of the four major toxins, the species are divided into five types (A, B, C, D and E) of which only the type A (food poisoning) and type C (enteritidis necrotic) affecting humans mainly exogenous infections.

Food poisoning can be caused by *C. perfringens* enterotoxin (CPE) produced by *C. perfringens* spores in the small intestine, which can germinate in foods such as meat and poultry (contaminated pre-cooked food which has been preserved or re-heated improperly). The main symptoms of *C. perfringens* are nausea, diarrhoea and abdominal pain. Symptoms appear after 8 to 24 hours after ingestion and lasts about 24 hours.

### TEST PRINCIPLE

CerTest *Clostridium perfringens* is based on the principle of a qualitative immunochromatographic assay for the determination of *Clostridium perfringens* in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *C. perfringens* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*C. perfringens*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*C. perfringens* monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*C. perfringens* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *C. perfringens* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*C. perfringens* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.


- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Clostridium perfringens* test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

## SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 24 hours prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool sample as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

## MATERIALS

### MATERIALS PROVIDED

- CerTest *Clostridium perfringens* card tests
- Instructions for use
- Stool collection tubes with diluent

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

## TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Clostridium perfringens* card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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## INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>Clostridium perfringens</i>	Interpretation of the results
1.	- GREEN	There is no <i>Clostridium perfringens</i> presence. No infection caused by <i>Clostridium perfringens</i> .
2.	+ GREEN-RED	There is <i>Clostridium perfringens</i> presence. Infection caused by <i>Clostridium perfringens</i> .
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

**NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

**INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens presents in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.

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- CerTest *Clostridium perfringens* should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *Clostridium perfringens* in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a *Clostridium perfringens* determination should be carried out on a sample from an enrichment culture or using another technique (for example PCR and microscopy).

## EXPECTED VALUES

The incidence of *C. perfringens* gastrointestinal illnesses in the United States has been estimated by the Centers for Disease Control and prevention to be 248,000 cases per year, leading to 41 hospitalizations and 7 deaths each year, with 100% of these cases being due to food-borne transmission of the pathogen and with an estimated cost of \$200 per case.

In another report from 1993 to 1997, *C. perfringens* accounted for 2.1% of the outbreaks and 3.2% of the cases of food-borne. Most cases of *C. perfringens* food poisoning are mild and are not reported. In 1994, the total cost of illnesses due to *C. perfringens* was estimated at \$123 million in the United States.

The estimated large number of illnesses due to *C. perfringens* clearly stresses the importance of cooling foods quickly after cooking, with proper refrigeration during shelf-life storage.

The effects is stronger symptoms in older persons residing in nursing homes and in hospitalized patients with compromised immune system or have other disease, because they are more susceptible to bacterial colonization.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using a culture to detect *Clostridium perfringens* infection. The results were as follows:

An evaluation was performed comparing two immunochromatographic tests to detect *Clostridium perfringens* infection (CerTest *Clostridium perfringens*, CerTest vs. *BIO K 195*, *BIO X* Diagnostics). The results were as follows:

		IC test: <i>BIO K 195</i>		
		+	-	Total
IC test: <i>CerTest Clostridium perfringens</i>	+	27	6	33
	-	0	37	37
Total		27	43	70

	Sensitivity	Specificity	PPV	NPV
<i>CerTest Clostridium perfringens</i> vs <i>BIO K 195</i>	>99%	86%	82%	>99%

The results showed a high sensitivity and specificity to detect *Clostridium perfringens* using *CerTest Clostridium perfringens*.

### Cross reactivity

An evaluation was performed to determine the cross reactivity of *CerTest Clostridium perfringens*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter coli</i>	<i>Salmonella enteritidis</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Salmonella paratyphi</i>	<i>Shigella flexneri</i>
<i>E. coli</i> O157:H7	<i>Salmonella typhi</i>	<i>Shigella sonnei</i>
<i>H. pylori</i>	<i>Salmonella typhimurium</i>	<i>Staphylococcus aureus</i>
<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>	<i>Yersinia enterocolitica</i>

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1. WREN, M.W.D, et al. "Laboratory diagnosis of *Clostridium* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". *British Journal of Biomedical Science*, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of *Clostridium* associated diseases". *Indican J. Med. Res.* 131, April 2010, p 487-499.
3. VIJAY, K. JUNEJA, et al. "Pathogens and Toxins in Foods Challenges and Interventions". *ASM Press, Washington, DC*, 2010, p 53-66.
4. MIRANDA, C., ROJO, MD., "Clostridium perfringens PIEL Y TEJIDOS BLANDOS". *Control Calidad SEIMC*, pp 1-10.
5. DOMÍNGUEZ-BERJÓN MF, ET AL. "Brote de toxiinfección alimentaria por *Clostridium perfringens* en un comedor escolar". *Med Clin (Barc)* 2003;121(2):58-60.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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