# IDK<sup>®</sup> DAO ELISA

*Zur in-vitro-Bestimmung von Diaminoxidase (DAO)  
in Serum*

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## 15. LITERATUR

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). *Agents and Actions* **23**:361-65.
2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. *Scand J Clin Lab Invest* **24**:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. *Inflammation Research* **48**:169-70.

4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. *Allergy Proceedings* **15**:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. *Inflammation Research* **47**:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): *The Atopy Syndrome in the Third Millennium. Curr Probl Dermatol*, Basel, Karger, **28**:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches *Clin Exp Allergy* **23**: 982-85.
8. Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel *Allergologie* **9**: 426-30.
9. Jarisch Reinhart, Histamin-Intoleranz. **1. Auflage** (1999), *Thieme-Verlag*, Stuttgart

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
# IDK<sup>®</sup> DAO ELISA

*For the in vitro determination of DAO in serum*

Valid from 2015-07-22

**REF** K 8500

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com) [www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>17</b>
<b>2. INTRODUCTION</b>	<b>17</b>
<b>3. MATERIAL SUPPLIED</b>	<b>17</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>19</b>
<b>6. SAMPLE COLLECTION, STORAGE AND PREPARATION</b>	<b>20</b>
<i>Sample stability</i>	20
<i>Sample preparation</i>	20
<b>7. ASSAY PROCEDURE</b>	<b>20</b>
<i>Principle of the test</i>	20
<i>Test procedure</i>	21
<b>8. RESULTS</b>	<b>22</b>
<b>9. LIMITATIONS</b>	<b>23</b>
<b>10. QUALITY CONTROL</b>	<b>23</b>
<i>Reference range</i>	23
<i>Heparin treatment</i>	24
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>26</b>
<i>Precision and reproducibility</i>	26
<i>Sample dilution</i>	26
<i>Analytical Sensitivity</i>	27
<b>12. PRECAUTIONS</b>	<b>27</b>
<b>13. TECHNICAL HINTS</b>	<b>27</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>28</b>
<b>15. REFERENCES</b>	<b>28</b>

## 1. INTENDED USE

This ELISA is intended for the quantitative determination of diamine oxidase (DAO) in serum. For *in vitro* diagnostic only.

## 2. INTRODUCTION

Diamine oxidase (DAO) is a body's own enzyme that metabolizes histamine. Although DAO is found practically in the whole body, the most important site of its action is the intestine. The enzymatic activity of DAO determines the histamine degradation speed. In the case of DAO deficiency or inhibition, incorporated or endogenous histamine cannot be degraded quickly enough, and the symptoms of histamine intolerance are presented. Millions of people suffer from gastrointestinal problems, migraine, irritations of nasal mucosa and other allergy-like symptoms after consumption of certain nutrients. Too much histamine in the body can be the reason for this wide range of symptoms.

The determination of DAO serum concentration (K 8500) combined with the determination of DAO activity (K 8220 DAO REA) is a suitable marker for the differential diagnosis of histamine intolerance and associated symptoms.

Our IDK® DAO ELISA kit is intended for determination of the diamine oxidase (DAO) concentration in serum.

### Indications

- Frequent headaches or migraine
- Snuffles after consumption of histamine-containing nutrients
- Tissue oedema
- Eyelid turgor
- Skin redness
- Limb aches
- Gastrointestinal discomfort
- Monitoring of a histamine free diet

## 3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 8500	PLATE	Holder with precoated strips	12 x 8 wells
K 8500	WASHBUF	ELISA wash buffer concentrate, 5 x	4 x 100 ml
K 8500	STD	Standards, lyophilized (see specification for concentration range)	4 x 5 vials

Cat. No	Content	Kit Components	Quantity
K 8500	CTRL	Control, lyophilized (see specification for concentration range)	4 x 1 vial
K 8500	CTRL	Control, lyophilized (see specification for concentration range)	4 x 1 vial
K 8500	AB	Detection antibody, (biotinylated), concentrate	1 x 200 µl
K 8500	CONJ	Conjugate (Streptavidin, peroxidase-labeled), concentrate	1 x 200 µl
K 8500	ABBUF	Dilution buffer for AB und CONJ, ready to use	1 x 50 ml
K 8500	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 50 ml
K 8500	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 8500	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Standard laboratory reaction vessels (1.5 ml)
- Standard laboratory reaction vessel (15 ml)
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Vortex mixer
- Ultra pure water\*
- Shaking incubator at 37 °C
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:5** before use (200 ml WASHBUF + 800 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer (1:5 diluted WASHBUF)** can be stored in a closed flask at **2–8 °C for one month.**
- **Please note:**
- **This WASHBUF is intended only for use in the IDK® DAO ELISA. Crystals in the WASHBUF must be completely dissolved before dilution.**
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at 2–8 °C until the expiry date stated on the label. The **STD** (standards) and **CTRL** (controls) must be reconstituted with **500 µl SAMPLEBUF** (sample dilution buffer). Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards are **not stable.**
- **Preparation of the conjugate and the detection antibody:** The **conjugate concentrate (CONJ)** and the **detection antibody concentrate (AB)** must be diluted **1:101** in dilution buffer (100 µl CONJ + 10 ml ABBUF) (100 µl AB + 10 ml ABBUF). The CONJ and the AB are stable at **2–8 °C** until expiry date stated on the label. **Conjugate (1:101 diluted CONJ) and detection antibody (1:101 diluted AB) are not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C.**

## 6. SAMPLE COLLECTION, STORAGE AND PREPARATION

### *Sample stability*

Storage at room temperature (15–30°C):	4 days
Storage at 2–8°C:	9 days
Storage at -20°C:	6 months

### *Sample preparation*

Pipet **50 µl** of sample in a 1.5 ml reaction vial, add **200 µl SAMPLEBUF** (sample dilution buffer) and mix well (corresponds to **1:5 dilution**).

## 7. ASSAY PROCEDURE

### *Principle of the test*

The assay utilizes the sandwich technique with two polyclonal antibodies against recombinant DAO.

Standards, controls and diluted samples which are assayed for DAO are added into the wells of a micro plate coated with polyclonal rabbit anti-DAO antibody. During the first incubation step, DAO is bound by the immobilized primary antibody. Then a biotinylated polyclonal anti-DAO antibody is added into each microtiter well. In the next step, the streptavidin peroxidase conjugate is added and a sandwich of 1st antibody – DAO - biotinylated antibody – streptavidin peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of DAO. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standard.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of STD /SAMPLE/CTRL (standards/sample/controls) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips in the carefully closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash each well <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
2.	Add <b>100 µl</b> of <b>STD</b> (standards)/ <b>SAMPLE/CTRL</b> (controls) into the respective wells.
3.	Cover the plate tightly and incubate for <b>2 hours</b> at <b>37 °C</b> on a horizontal mixer.
4.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5.	Add <b>100 µl of detection antibody</b> into each well, mix gently.
6.	Cover the plate tightly and incubate for <b>1 hour</b> at <b>37 °C</b> on a horizontal mixer.
7.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add <b>100 µl</b> of <b>conjugate</b> into each well.
9.	Cover the plate tightly and incubate for <b>1 hour</b> at <b>37 °C</b> on a horizontal mixer.
10.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
11.	Add <b>100 µl</b> of <b>SUB</b> (substrate) into each well.
12.	Incubate for <b>10–20 minutes*</b> at <b>room temperature</b> (15–30 °C) in the dark.

13.	Add <b>100 µl</b> of <b>STOP</b> (stop solution) into each well, shake well.
14.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### Samples

The obtained DAO concentration must be multiplied by the dilution factor of **5**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used to get the real concentration.

## 9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

The lower limit of the measurement range can be calculated as:

*LoB × sample dilution factor to be used*

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*


< 3 U/ml:	high incidence for HIT (Histamine intolerance)
3 - 10 U/ml:	HIT probable
> 10 U/ml:	low incidence for HIT

Conversion factor: 1 U/ml = 1.25 ng/ml

We recommend each laboratory to establish its own reference range.

### Comparison “Pregnant women” and “Healthy controls”

For the clinical evaluation of this assay we have analyzed samples from pregnant women and apparently healthy controls. The IDK® DAO ELISA detects, as required and expected, higher values in pregnant women than in healthy controls.

### Heparin treatment

Furthermore, DAO levels in healthy study participants increased sharply within 30 minutes of heparin administration. It has been documented in scientific literature that DAO levels rise after heparin administration.

### Treatment outcome before and after heparin administration

IDK® DAO ELISA [U/ml]

	BEFORE administration	30 min AFTER	60 min AFTER
Patient 1	76	219	-
Patient 2	55	152	-
Patient 3	2.5	277	622
Patient 4	18.9	621	555

Since the IDK® DAO ELISA determines DAO concentration while the conventional histamine intolerance tests with putrescine or histamine as substrates determine DAO activity, the correlation coefficient must not necessarily be  $r > 0.8$ . This can be explained by the fact that the activity does not depend on the number of molecules alone, but also on cofactors such as vitamin C, vitamin B6, copper or manganese ions in vitro and in vivo. For the diagnosis of histamine intolerance via DAO activity test we therefore recommend to determine the above mentioned cofactors as well. The problem may not be a low DAO level, but a cofactor deficiency.

The symptoms of histamine intolerance can be caused by low DAO activity because

the above-mentioned cofactors are not sufficiently available. By quantitating the cofactors it can be determined which one needs to be supplemented.

### Medication effects

In addition, histamine intolerance symptoms may be due to low DAO activity caused by medication such as:

Muscle relaxants	Pancuronium, alcuronium, D-tubocurarine
Narcotics	Thiopental
Analgetics	Morphine, pethidine, nonsteroidal anti-inflammatory drugs, acetylsalicylic acid, metamizole
Local anesthetics	Prilocaine
Antihypotonics	Dobutamine
Antihypertensive drugs	Verapamil, alprenolol, dihydralazine
Antiarrhythmics	Propafenone
Diuretics	Amiloride
Drugs influencing gut motility	Metoclopramide
Antibiotics	Cefuroxime, cefotiam, isoniazid, pentamidin, clavulanic acid, choroquine
Mucolytics	Acetylcysteine, ambroxol
Broncholytics	Aminophylline
H <sub>2</sub> -receptor antagonists	Cimetidine
Cytostatics	Cyclophosphamide
Antidepressants	Amitriptyline

If you are taking such medication, you may want to discuss with your physician alternative medication in order to relieve your symptoms.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay (n = 16)**

The precision (intra-assay variation) was calculated from 16 replicate determinations on each of two samples.

Sample	DAO [U/ml]	CV [%]
1	5.0	1.42
2	5.9	1.72

#### **Inter-Assay (n = 8)**

The total precision (inter-assay variation) was calculated from data on 2 samples obtained in 8 different assays by different technicians on three different days.

Probe	DAO [U/ml]	CV [%]
1	23.27	7.9
2	13.36	10.7

### *Sample dilution*

Two patient samples were diluted and assayed. The results are shown below:

Linearity n= 2

Sample	Dilution	Expected [U/ml]	Measured [U/ml]
A	1:5	15.65	15.65
	1:10	7.83	7.68
	1:20	3.92	3.86
	1:40	1.96	1.86
B	1:5	5.58	5.58
	1:10	2.79	2.89
	1:20	1.39	1.45
	1:40	0.7	0.7

### *Analytical Sensitivity*

Limit of blank, LoB	0.079 U/ml
Limit of detection, LoD	0.120 U/ml
Limit of quantitation, LoQ	0.120 U/ml

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20% CV.

## **12. PRECAUTIONS**

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). *Agents and Actions* **23**:361-65.
2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. *Scand J Clin Lab Invest* **24**:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. *Inflammation Research* **48**:169-70.
4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. *Allergy Proceedings* **15**:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. *Inflammation Research* **47**:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): *The Atopy Syndrome in the Third Millennium. Curr Probl Dermatol*, Basel, Karger, **28**:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches *Clin*

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Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by