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Epidermal Growth Factor-Receptor ELISA Kit

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

For the in vitro determination of Epidermal Growth Factor-Receptor in plasma and tissue extracts

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TABLE OF CONTENTS	PAGE
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
TEST PROCEDURE	18
10. RESULTS	19
11. LIMITATIONS	19
12. QUALITY CONTROL	20
13. PERFORMANCE CHARACTERISTICS	20
PRECISION AND REPRODUCIBILITY	20
SENSITIVITY	21
CROSS REACTIVITY	21
SAMPLE DILUTION	21
14. REFERENCES	21
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	21

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In-Vitro-Diagnostikum

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Hersteller

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

Epidermal Growth Factor- Receptor ELISA Kit

*For the in vitro determination of Epidermal Growth Factor-
Receptor in plasma and tissue extracts*

Gültig ab / Valid from 13.02.2008

K 7110

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of EGF-R in plasma and tissue extracts. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Polyclonal growth factors e.g. the epidermal growth factor (EGF) have been recognized as important determinants in the regulation of cellular proliferation and differentiation as shown in cell culture. EGF mediates its effects through receptors.

The EGF-receptor (EGF-R) is a glycoprotein of 170.000 Dalton. It consists of an extra cellular domain for EGF binding, a transmembrane domain for signal transduction and an intracellular domain with tyrosine-kinase activity. Binding of EGF leads to a Ca^{2+} release by activation of Calmodulin-dependent enzymes.

Overexpression or amplification of EGF-R in breast cancer is associated with autonomous tumor growth and poor clinical prognosis.

The most widely used method of receptor identification has been a radio ligand - assay with ^{125}I -labeled EGF. The advantage of the radio ligand assay is the measurement of native receptor. On the other hand, this test system is combined with a lot of disadvantages (e.g. great amount of tumor material, radioactivity, poor reproducibility, time-consuming calculation).

In the offered test system EGF-R is determined by specific antibodies. This immunoassay shows good reproducibility, small sample volume, and is easy to calculate. The sensitivity for h EGF-R in ELISA - technique is 0.4 fMol/ml.

3. PRINCIPLE OF THE TEST

This Enzyme-Linked-Immunosorbent Assay (ELISA) allows the quantitative determination of human Epidermal Growth Factor-Receptor (h EGF-R). The h EGF-Rs in the samples are bound to monoclonal mouse antibodies against h EGF-R immobilized on the surface of the microtiter plates. After a washing step, to remove all foreign substances, the quantification of bound h EGF-R is carried out by adding a rabbit anti human EGF-R antibody. Detection of the bound rabbit antibody is performed by a peroxidase labeled goat anti rabbit antibody (PO-antibody). The amount of converted substrate is directly proportional to the amount of bound h EGF-R and can be determined photometrically at 450 nm.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 7110MTP	PLATE	One holder with strips	96
K 7110WP	WASHBUF	ELISA wash buffer concentrate 10x	100 ml
K 7110AP	ABBUF	Antibody dilution buffer, ready-to-use	25 ml
K 7110A1	1.AB	1 st Antibody (mouse anti hEGF-R)	130 µl
K 7110A2	2.AB	2 nd Antibody (rabbit anti hEGF-R)	80 µl
K 7110K	CONJ	Conjugate, (goat anti rabbit, peroxidase-labeled), ready to use	22 ml
K 7110BeP	COATBUF	Coating buffer, ready-to-use	25 ml
K 7110PV	SAMPLEBUF	Sample dilution buffer, ready-to-use	50 ml
K 7110BP	BLOCKBUF	Blocking buffer, ready-to-use	25 ml
K 7110ST	STD	Calibrator (0; 3; 6; 12; 24 fmol/ml)	2 x 5 vial
K 7110KO	CTRL	Control, lyophilized	2 x 1 vial
K 7110TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	2 x 15 ml
K 7110AC	STOP	ELISA Stop solution, ready-to-use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized water
- Precision pipettors calibrated to deliver 10-1000 µl
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **WASHBUF** (ELISA wash buffer concentrate) must be diluted with aqua dest. **1:10** before use (50 ml WASHBUF + 450 ml aqua dest.). Crystals could occur due to high salt concentration. The crystals must be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The **WASHBUF** is stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **STD** (calibrators) and **CTRL** (control) must be reconstituted with **500 µl SAMPLEBUF** (sample dilution buffer). Allow the vial to stand for minimum 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted Calibrators and Control can not be stored.
- The **1. AB** (1st Antibody) must be diluted **1:200** in **COATBUF** (coating buffer) (100 µl 1.ABi + 20 ml COATBUF). The **1. AB** is stable at -20 °C until expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **2. AB** (2nd Antibody) must be diluted **1:500** in **ABBUF** (antibody dilution buffer) (40 µl 2. AB + 20 ml ABBUF). The **2. AB** is stable at -20 °C until expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop solution contains diluted sulfuric acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Plasma

The **plasma** has to be diluted between 1/2 and 1/10 in dilution buffer to be in the range of standard curve and then measured directly.

Tissue

Tissue material has to be homogenised (e.g. mechanically with a dismembrator). Then, an ultra centrifugation step (1h, 100.000 x g) follows. The supernatant of the centrifugation is the cytosolic fraction (for oestrogen and progesterone receptor measurement) and the pellet includes membrane of the outer and nucleus walls.

For measurement of EGF-R by ELISA the following matrices could be used:

1. The **membrane pellet** has to be solubilized in PBS (approximate 1 ml per mg protein 0.14M NaCl, 2.6 mM KCl, 8 mM Na₂HPO₄, 1.4 mM KH₂PO₄, pH 7.4) 1% Tween 20 with an ultra turrax for 2 min. After incubation for 30 min. at room temperature with constant agitation the samples have to be centrifuged with 10 – 15000 x g for 10 min. Protein in the supernatant should be determined according to Lowry et al. or with BCA Protein Assay [Pierce]. Protein should be diluted to 50 – 100 µg protein / ml in PBS 1% Tween 20 and 2 x 200 µl/sample should be used for the ELISA.
2. The **cytosolic fraction** has to be diluted to 50-200 µg protein/ml and then be measured directly. Protein determination, see above.
3. To prepare the **total lysate** the tissue should be first homogenized as described above without centrifugation step. For further preparation see point 1.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

Carry out the tests in duplicate.

1. Pipette **200 µl** of prediluted **1. AB** (1st antibody) solution into each well.
2. Incubate for **1 hour at room temperature (RT)**.
3. Decant the contents of the plate.
4. Pipette **250 µl** of **BLOCKBUF** (blocking buffer) into each well.
5. Incubate for **1 hour at room temperature** (in order to block remaining binding sites).
6. Decant the contents of the plate and wash the cavities **5 x with 250 µl** of ELISA wash buffer.
7. Pipette **200 µl** of **STD** (standards), **CTRL** (control) or **samples** into the wells of the microtitre plate.
8. Incubate (**16 – 18 h**) **over night at room temperature**.
9. Decant the contents of the plate and wash the cavities **5 x with 250 µl** of ELISA wash buffer.
10. Add **200 µl** of **2. AB** (2nd antibody)
11. Incubate for **1 hour at RT**
12. Decant the contents of the plate and wash the cavities **5x with 250 µl** of ELISA wash buffer.
13. Add **200 µl** of **CONJ** (conjugate, POD-antibody).
14. Incubate for **1 hour at RT**.
15. Decant the contents of the plate and wash the cavities **5x with 250 µl** of ELISA wash buffer.
16. Add **200 µl** of **SUB** (TMB substrate) solution
17. Incubate for **15 –20 minutes at room temperature**, shaking lightly, until sufficient coloring has been achieved.
18. Add **50 µl** of **STOP** (stop solution).
19. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the „4-Parameter-algorithm“.

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Example for calculation

Example for calculation of the EGF-R content in an unknown sample

$$\frac{\text{Result in fMol / ml}}{\text{Protein concentration in mg / ml}} = \text{EGF - Receptor concentration in fMol / mg Protein}$$

Example:

An EGF-R concentration of 4.5 fMol/ml is determined by the calibration curve. The protein concentration of the sample is 50µg/ml.

The final result is:

$$\frac{4.5 \text{ fMol / ml}}{0.050 \text{ mg / ml}} = 90 \text{ fMol EGF - R / mg Protein}$$

11. LIMITATIONS

Samples with EGF-R levels greater than the highest calibrator, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra Assay variation

The precision (intra-assay variation) of the Immundiagnostik EGF-R ELISA test was calculated from 11 replicate determinations on each of one samples.

Intra-Assay CV n = 11

Sample	EGF-R [fmol/ml]	Intra-Assay CV [%]
1	7.8	10.8
2	16.1	5.4

Inter Assay variation

The total precision (inter-assay variation) of the Immundiagnostik EGF-R ELISA test was calculated from data on 2 samples obtained in 11 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 11

Sample	EGF-R [fmol/ml]	Inter-Assay -Assay CV [%]
1	7.8	14.8
2	16.1	9.6

Sensitivity

The detection limit of this ELISA was defined as $B_0 + 2SD$. The value is 0.4 fmol/ml

Cross reactivity

No cross reactions with other cellular proteins have been observed.

Sample dilution

The linearity was calculated from the dilution of EGF-R containing material with dilution buffer. Linearity was achieved in the range of 3-24 fmol/ml.

14. REFERENCES

Tumor Necrosis Factor α Sensitizes Low Epidermal Growth Factor Receptor (EGFR)-expressing Carcinomas for Anti_EGFR Therapy.
Hambeck M. et al. (2001) *Cancer Res* 61:1045-1049

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent; these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.

- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from this.

13.02.2008 23012004_EGF-R.DOC

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number