CERTEST FOB 50+200

ONE STEP human
haemoglobin COMBO
CARD TEST

CERTEST BIOTEC S.L.

CerTest
BIOTEC

For Information Only Purposes

ENGLISH**INTENDED USE**

CerTest FOB 50+200 one step combo card test is a coloured chromatographic immunoassay for the semi-quantitative detection of human haemoglobin (hHb) in stool samples to the following cut-off values: A: 50ng/mL and B: 200ng/mL.

CerTest FOB 50+200 combo card test offers a simple, a highly sensitive and non-invasive screening assay to detect gastrointestinal bleeding (faecal occult blood).

INTRODUCTION

Colorectal cancer is the second leading cause of illness and death in the Western world. The screening using faecal occult blood tests is based on the bleeding detection, as important symptom of colonic neoplasm, either at an early-stage cancer or large adenomatous polyps. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

Haemoglobin (hHb) is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases.

TEST PRINCIPLE

CerTest FOB 50+200 is based on the principle of a semi-quantitative immunochromatographic assay for the determination of human haemoglobin in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human haemoglobin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human haemoglobin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human haemoglobin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human haemoglobin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is human haemoglobin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human haemoglobin monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is human haemoglobin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human haemoglobin monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human haemoglobin antibodies present on the membrane of strip A (test line) and the anti-human haemoglobin antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no human haemoglobin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-human haemoglobin antibodies present on the membrane (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date marked on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest FOB 50+200 test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In that case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (1), to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 15µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (2). The stool collection vial with diluted sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

Precautions: Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples.

MATERIALS

MATERIALS PROVIDED

- CerTest FOB 50+200 combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest FOB 50+200 combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (3) and dispense 4 drops in the circular window marked with the letter A (4) and 4 drops, using the same tube, in the circular window marked with the letter B (5). Avoid adding solid particles with the liquid.
4. **Read the result at 10 minutes.** Do not read the test result after more than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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
Add 4 drops in B-circular window.


(5)


INTERPRETATION OF THE RESULTS (please refer to the illustration below)

A: Green → Negative haemoglobin 50
 B: Green → Negative haemoglobin 200


A: Green-Red → Positive haemoglobin 50
 B: Green-Red → Positive haemoglobin 200


A: Green-Red → Positive haemoglobin 50
 B: Green → Negative haemoglobin 200

INVALID RESULTS

A: Green → Negative haemoglobin 50
 B: Green-Red → Positive haemoglobin 200

| | A (Haemoglobin 50 (concentration 50ng/mL)) | B (Haemoglobin 200 (concentration 200ng/mL)) | Interpretation of the results |
|-----------|---------------------------------------------------------|-----------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. | - GREEN | - GREEN | There is no human haemoglobin presence. Haemoglobin marker is not present in patient sample (<50ng/mL), which might mean no faecal occult blood and no gastrointestinal bleeding. |
| 2. | + GREEN-RED | + GREEN-RED | High haemoglobin concentration (200ng/mL) is present in patient sample, which might mean some gastrointestinal bleeding problem (e.g. caused by colorectal cancer, bacterial infections, inflammatory bowel disease, etc.). Additional confirmatory diagnostic (endoscopy) should be followed up to determine the exact cause and source of the blood. |
| 3. | + GREEN-RED | - GREEN | Haemoglobin marker is present in patient sample (50ng/mL and <200ng/mL), which might mean small amounts of blood loss near the detection limit of the combo test. If the problem still persists, physician should evaluate all clinical and laboratory findings. |
| 4. | - GREEN | + GREEN-RED | Invalid result. If B is positive, A has to be also positive. |
| 5. | Any other result | Any other result | Invalid result either A or B, we recommend repeating the assay using the same sample with another test. |

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the result lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay using a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

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NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of human haemoglobin present in the specimen.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening of the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high haemoglobin concentration to faint when the haemoglobin concentration is close to the cut-off value of the test.
4. CerTest FOB 50+200 should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained
5. Positive results determine the presence of human haemoglobin in faecal samples; nevertheless, it can be due to several causes, besides colorectal bleeding, such as haemorrhoid, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
6. A negative result is not meaningful because of it is possible the haemoglobin concentration in the stool sample is lower than the cut-off values. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
7. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine, or if they have strained during bowel movement.

EXPECTED VALUES

Some studies established 200ng haemoglobin/mL as optimal cut-off value in faecal occult blood test used in colorectal cancer screening programs.

Colorectal cancer affects both men and women of all racial and ethnic groups, being most often found in people aged 50 years or older. Regarding men, colorectal cancer is the third most common cancer after prostate and lung cancers. Regarding women, colorectal cancer is the third most common cancer after breast and lung cancers.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off values of CerTest FOB 50+200 are 50ng/mL (5.1µg hHb/g faeces) for strip A and 200ng/mL (20µg hHb/g faeces) for strip B.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest FOB 50+200, CerTest) and another commercial available test (ACTIM FECAL BLOOD, Medix Biochemica). The results were as follows:

| | | IC Test: ACTIM FECAL BLOOD | | |
|-----------------------------|-------|----------------------------|-----|-------|
| | | + | - | Total |
| IC test: CerTest FOB 50+200 | + | 6 | 0 | 6 |
| | - | 0 | 121 | 121 |
| | Total | 6 | 121 | 127 |

| | Sensitivity | Specificity | PPV | NPV |
|-----------------------------------------|-------------|-------------|------|------|
| CerTest FOB 50+200 vs ACTIM FECAL BLOOD | >99% | >99% | >99% | >99% |

The results showed a high sensitivity and specificity to detect human haemoglobin using CerTest FOB 50+200.

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Cross-Reactivity

An evaluation was performed to determine the cross reactivity of CerTest FOB 50+200; no cross-reactivity against other faecal markers occasionally present in faeces:

| | |
|--------------------|----------------------------|
| Human lactoferrin | Bovine and pig haemoglobin |
| Human transferrin | Bovine and pig transferrin |
| Human calprotectin | |

No special diet is recommended prior to testing. There are not interferences with any foods (Vitamin C, broccoli, carrots...).

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
1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.
4. ITOH M., et al. "Estimation of the optimal cut-off point in a new immunological faecal occult blood test in a corporate colorectal cancer screening programme". J. Med. Screen 1996; 3(2):66-71

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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