

CE

IVD

For Professional Use Only

AmpliSens[®] HPV HCR screen-EPh
PCR kit
Instruction Manual

AmpliSens[®]



Ecoli s.r.o., Studenohorska 12
841 03 Bratislava 47
Slovak Republic
Tel.: +421 2 6478 9336
Fax: +421 2 6478 9040



Federal Budget Institute of
Science "Central Research
Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia

TABLE OF CONTENTS

1. INTENDED USE	3
2. PRINCIPLE OF PCR DETECTION.....	3
3. CONTENT	4
4. ADDITIONAL REQUIREMENTS	5
5. GENERAL PRECAUTIONS.....	5
6. SAMPLING AND HANDLING	6
7. WORKING CONDITIONS.....	7
8. PROTOCOL.....	7
9. DATA ANALYSIS.....	9
10. TROUBLESHOOTING.....	9
11. TRANSPORTATION.....	10
12. STABILITY AND STORAGE	10
13. SPECIFICATIONS.....	10
14. REFERENCES	10
15. QUALITY CONTROL.....	11
16. KEY TO SYMBOLS USED	12

1. INTENDED USE

AmpliSens® HPV HCR screen-EPh PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of *human papillomavirus (HPV)* of high carcinogenic risk (HCR) types 16, 18, 31, 33, 35, 39, 45, 52, 53, 56, 58, 59, 66, 70 DNA in the clinical material (urogenital swabs) by using electrophoretic detection of the amplified products in agarose gel.



The results of PCR analysis are taken into account in complex diagnostics of disease.

2. PRINCIPLE OF PCR DETECTION

HPV HCR detection by the polymerase chain reaction (PCR) is based on the amplification of pathogen DNA specific region using special *HPV* HCR primers. After PCR the amplified product is detected in agarose gel. **AmpliSens® HPV HCR screen-EPh** PCR kit is able to detect three main phylogenetic groups of *HPV*: A7, A6, and A9 that include following types: 16, 18, 31, 33, 35, 39, 45, 52, 53, 56, 58, 59, 66, 70. These types are known for high transformation ability and they are responsible for more than 97% of severe cervical dysplasia and cervical cancer. **AmpliSens® HPV HCR screen-EPh** PCR kit contains the internal endogenous control (fragment of β -globine gene). **AmpliSens® HPV HCR screen-EPh** PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase by using a chemically modified polymerase (TaqF). Chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The test is based on simultaneous amplification (multiplex PCR) of E1-E2 gene fragments of three *HPV* groups and a fragment of β -globin gene which is used as an internal endogenous control. PCR test detecting fourteen HPV types of is run in a single tube. All amplified products have a similar length.

Detection of clinically significant virus quantity by using AmpliSens® HPV HCR screen-EPh PCR kit.

Nowadays epidemiologic studies show that most routine screening tests for dysplastic changes of cervix, vagina, and vulva as well as for evaluation of the risk of their development require detection of *clinically valuable* quantity of high-risk HPV. Believed, that detection of virus in quantity not exceeding certain threshold value is clinically insignificant because 100% of such cases associate with spontaneous recovery. On the contrary, high virus load indicate dysplasia or risk of its development. However, in case of treatment monitoring the diagnosis of even low virus quantity can be an early marker of a relapse. Currently, the level of clinically significant virus quantity estimates at 10^5 GE of HCR *HPV* per cervical swab when standard obtaining of clinical material is provided. Clinical trails with model samples have showed that

clinically significant virus quantity can be defined if following steps are taken:

- urogenital swabs are collected according to the standard procedure (placed in 0.6 ml of transport medium).
- DNA is extracted (DNA-sorb-A or DNA-sorb-AM were used).
- DNA samples are diluted 100x with TE-buffer.
- PCR test is performed.

Clinical trials on samples collected from both healthy patients and patients suffering from severe dysplasia and cervical cancer demonstrated increase of dysplasia detection specificity by 22.9% (from 61.7% without dilution to 84.6% with dilution) whereas the high level of severe dysplasia and cervical cancer diagnosis was maintained (98.9%). Note that level of clinically significant virus quantity wasn't validated in men.

Therefore, AmpliSens® HPV HCR screen-EPh PCR kit can run two formats of high-risk HPV detection:

- **HPV HCR presence** (a sample is tested after DNA extraction).
- **HPV HCR clinically significant quantity** (a sample is tested after DNA extraction and dilution with TE-buffer). Note that standard sampling is necessary.

3. CONTENT

AmpliSens® HPV HCR screen-EPh PCR kit is produced in 1 form:

AmpliSens® HPV HCR screen-EPh PCR kit variant 100 F **REF** V31-100F-CE.

AmpliSens® HPV HCR screen-EPh PCR kit variant 100 F includes:

<i>Reagent</i>	<i>Description</i>	<i>variant 100 F</i>	
		<i>Volume (ml)</i>	<i>Quantity</i>
PCR-mix -1 HPV HCR screen	colorless clear liquid	0.3	2 tubes
2.5x PCR-buffer blue	blue clear liquid	1.15	1 tube
Polymerase (TaqF)	colorless clear liquid	0.06	1 tube
Mineral oil for PCR	colorless viscous liquid	4.0	1 dropper bottle
Positive Control DNA HPV types 31, 39, 56 and DNA human (C+_{HPV 31, 39, 56})	colorless clear liquid	0.2	1 tube
Positive Control DNA human (C+h)	colorless clear liquid	0.2	1 tube
TE-buffer	colorless clear liquid	5.0	4 tubes
Negative Control (C-)*	colorless clear liquid	1.2	1 tube

* must be used in the extraction procedure as Negative Control of Extraction (see DNA-sorb-AM, **REF** K1-12-100-CE, or DNA-sorb-B, **REF** K1-2-100-CE or DNA-sorb-C, **REF**

K1-6-50-CE protocols).

AmpliSens® HPV HCR screen-EPh PCR kit variant 100 F is intended for 110 reactions, including controls.

4. ADDITIONAL REQUIREMENTS

- DNA extraction kit.
- Agarose gel detection kit.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol barriers (up to 200 µl).
- Vortex mixer.
- Desktop centrifuge with rotor for 2 ml reaction tubes.
- PCR box.
- Personal thermocyclers (for example, Terzik (DNA-Technology), GeneAmp PCR System 2400, GeneAmp PCR System 2700 (Applied Biosystems), T-personal (Biometra), PTC-100 (MJ Research)).
- Disposable polypropylene microtubes for PCR with 0.2- 0.5-ml capacity (for example, Axygen, USA).
- Refrigerator for 2-8 °C.
- Deep-freezer for ≤ -16 °C.
- Waste bin for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store and handle amplicons away from all other reagents.
- Thaw all components thoroughly at room temperature before starting detection.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats, protect eyes while samples and reagents handling. Thoroughly wash hands afterward.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all samples and unused reagents in compliance with local authorities requirements.
- Samples should be considered potentially infectious and handled in a biological cabinet in

accordance with appropriate biosafety practices.

- Clean and disinfect all sample or reagent spills using a disinfectant such as 0.5 % sodium hypochlorite, or other suitable disinfectant.
- Avoid contact with the skin, eyes and mucosa. If skin, eyes and mucosa contact immediately flush with water, seek medical attention.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.



Some components of this kit contain Sodium Azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING



Obtaining samples of biological materials for PCR-analysis, transportation, and storage are described in manufacturer's handbook [1]. It is recommended to read this handbook before starting work

AmpliSens[®] HPV HCR screen-EPh PCR kit is intended for analysis of DNA extracted with DNA extraction kits from:

- *Urogenital swabs*

Female: epithelium samples should be obtained as for cytological examination:

Method I - use the sampling kit which includes one/two cervical cytobrushes and the 2-ml tube with 0.5 ml of transport medium TSM.

Endocervical epithelial scrape taken with the first cytobrush and/or exocervical epithelial scrape taken with the second cytobrush should be placed into the tube with a transport medium. Break the effective part of the cytobrush with the sample at the score mark and leave it in the tube.

Method II - use Digene (USA) sampling kit, which contains the cervical cytobrush and 1.0-ml tube with Digene transport medium.

Endocervical epithelial scrape obtained with cytobrush should be placed into the tube with Digene transport media.

Method III - use the sampling kit, which contains combined gynecological probe for simultaneous obtaining of epithelial cells from endo-/exocervix and 5-ml tube with 2.0 ml of transport medium TSM.

Place endocervical and exocervical epithelial scrapes into the tube with transport media. Break the effective part of the cytobrush with the sample at the score mark and leave it in the

tube.

Method IV - use CytoScreen (Italy) or PreservCyt (USA) sampling kits which contains the combined gynecological probe for simultaneous obtaining of epithelium from endo-/exocervix and a vial with transport-fixation medium.

Place endocervical and exocervical epithelial scrapes into the tube with transport-fixation medium. Break the effective part of the cytobrush with the sample at the score mark and leave it in the vial.

Male: Obtain the urethral epithelial scrape by a universal probe and place it into the 2.0-ml tube with 0.5 ml of transport medium TSM.



Only one freeze-thaw cycle of clinical material is allowed.

7. WORKING CONDITIONS

AmpliSens® HPV HCR screen-EPh PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. DNA Extraction

It's recommended to use the following nucleic acid extraction kits:

- DNA-sorb-AM, **REF** K1-12-100-CE.
- DNA-sorb-B, **REF** K1-2-100-CE.
- DNA-sorb-C, **REF** K1-6-50-CE (is used for mucosa biopsy samples).



Please carry out the DNA extraction according to the instructions provided by the manufacturer.

8.2. Preparing the PCR

Total reaction volume - **25 µl**, volume of DNA sample - **10 µl**.

8.2.1 Preparing tubes for PCR

1. Prepare reaction mix for N reactions:

5*(N+1) µl of PCR-mix-1 HPV HCR screen

10*(N+1) µl of 2.5x PCR-buffer blue

0.5*(N+1) µl of polymerase (TaqF)

Calculate the reaction mixture volume taking into account three controls (one negative and two positive) and one extra reaction.



The prepared reaction mixture can be stored for up to 2 hours.

2. Vortex the tube with the reaction mixture. Add 15 µl of the reaction mixture into PCR tubes.

3. Add above 1 drop of **mineral oil for PCR** (about 25 µl). Omit if a thermocycler with constant-temperature cover is used.
4. Using tips with aerosol barrier add **10 µl of DNA samples** obtained from clinical or control samples at the DNA extraction stage, under oil or directly on the level of oil.(in case of *clinically significant quantity of HPV HCR format* use 100x dilution of obtained DNA samples).
5. Carry out the control amplification reactions:
 - NCA - Add 10 µl of **TE-buffer** to the tube for Negative Control of Amplification (NCA).
 - C+_h - Add 10 µl of **Positive Control DNA human** to the tube for Positive Control of human DNA Amplification.
 - C+_{HPV31, 39, 56} - Add 10 µl of **Positive Control DNA HPV types 31, 39, 56 and human DNA** to the tube for Positive Control of HPV Amplification.

8.2.2 Amplification

Run the following program on the thermocycler (see table 1). When the temperature reaches 95 °C (pause regimen), insert tubes to cells of amplifier and press button to continue.

It is recommended to sediment drops from walls of tubes by short vortex (1-3 s) before their insertion in thermocycler.

Table 1

Programming thermocyclers for DNA amplification of DNA of HPV HCR types 16, 18, 31, 33, 35, 39, 45, 52, 53, 56, 58, 59, 66, 70

Step	Thermocyclers with active temperature adjustment:						Thermocyclers with block temperature adjustment:		
	Terzik (DNA-Technology)			GeneAmp PCR System 2400 (Perkin Elmer), GeneAmp PCR System 2700 (Applied Biosystems)			T-personal (Biometra), PTC-100 (MJ Research)		
	Temperature	Time	Cycles	Temperature	Time	Cycles	Temperature	Time	Cycles
1	95 °C	900 s	1	95 °C	900 s	1	95 °C	900s	1
2	95 °C	10 s	42	95 °C	15 s	42	95 °C	30 s	42
	63 °C	20 s		63 °C	30 s		65 °C	40 s	
	72 °C	20 s		72 °C	30 s		72 °C	40 s	
3	72 °C	60 s	1	72 °C	60 s	1	71 °C	60 s	1

The run takes approximately 2 h 30 min to complete in a thermocycler with block temperature adjustment or 1 h 50 min in a thermocycler with active temperature adjustment.

After the reaction is completed, the PCR tubes must be collected and sent to the room for PCR products analysis.

Analysis of amplification products is performed by separation of DNA fragments in agarose gel.

The amplified samples can be stored for 16 h at room temperature, for 1 week at 2–8 °C (make sure that the samples are warmed up to room temperature before running electrophoresis).

9. DATA ANALYSIS

It's recommended that the following detection agarose kit is used:

- EPh variant 200, **REF** K5-200-CE.

Analysis of results is based on the presence or absence of specific bands of amplified DNA in agarose gel (1.7%).

PCR-mix -1 HPV HCR screen includes primers for amplification of 16, 18, 31, 33, 35, 39, 45, 52, 53, 56, 58, 59, 66, 70 types DNA fragments as well as the fragment of human genome DNA (β -globin gene).

The length of specific amplified DNA fragments is:

- Internal Control (fragment of β -globin gene) - 723 bp
- HPV types 16, 18, 31, 33, 35, 39, 45, 52, 53, 56, 58, 59, 66, 70 - from 267 to 325 bp



Use a protective mask or a glass filter when looking at the gel or taking photos

9.1. Results interpretation

Table 2

Results for controls

Control	Step for control	Specific bands in the agarose gel		Interpretation
		267-325 bp	723 bp	
C-	DNA extraction	No	No	OK
NCA	Amplification	No	No	OK
C+h	Amplification	No	Yes	OK
C+ <i>HPV</i> _{31, 39, 56}	Amplification	Yes	Yes	OK

- The sample is considered positive for *HPV* HCR DNA if the band at the level from 267 to 325 bp is present in agarose gel regardless of the band of Internal Control (723 bp).
 - The sample is considered negative for *HPV* HCR DNA if the only 723-bp band is present.
- Besides the specific products the fuzzy bands of primer dimers may appear in the lanes below the 100-bp level.

10. TROUBLESHOOTING

Analysis results are not obtained as per the following examples:

- If the results of the controls do not match with the listed above (Table 2), then the appropriate step of the test should be repeated.
- If the Internal Control band (723 bp) is not observed in the lane of the clinical sample, it can indicate that the insufficient quantity of clinical material was taken or mistakes in clinical processing, DNA extraction, or PCR conducting were made.
- If nonspecific bands are seen at different levels in the lanes, this may be caused by the lack of "hot start" or incorrect temperature profile of the thermocycler.

- If the specific bands appear in the lanes corresponding to the negative controls (NCA, C–) it indicates contamination of the reagents or samples. In this case results of the analysis for all samples are considered invalid. The test must be repeated and measures to detect and eliminate the source of contamination should be taken.

If you have any further questions or if encounter problems, please contact our Authorized representative in the European Community.

11. TRANSPORTATION

AmpliSens® HPV HCR screen-EPh PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the **AmpliSens® HPV HCR screen-EPh** PCR kit are to be stored at 2-8 °C (except for polymerase TaqF) when not in use. All components of the PCR kit are to be stable until the labeled expiration date. The shelf life of reagents before and after the first use is the same, unless otherwise stated.



Polymerase (TaqF) is to be stored at temperature from minus 24 to minus 16 °C when not in use.

13. SPECIFICATIONS

13.1. Sensitivity

Analytical Sensitivity of **AmpliSens® HPV HCR screen-EPh** PCR kit is no less than 5×10^3 genome equivalents per 1 ml of sample (GE/ml) (types 16, 18, 31, 35, 39, 45, 52, 56, 59, 66, 70) and 2.5×10^4 GE/ml (types 33, 58).



Claimed analytical features of **AmpliSens® HPV HCR screen-EPh** PCR kit are guaranteed only when additional reagent kits of, DNA-sorb-AM, DNA-sorb-B, or DNA-sorb-C and EPh, are used.

13.2. Specificity

The analytical specificity of **AmpliSens® HPV HCR screen-EPh** PCR kit is ensured by selection of specific primers and strict reaction conditions.

The clinical specificity of **AmpliSens® HPV HCR screen-EPh** PCR kit was confirmed in laboratory clinical trials.













14. REFERENCES

1. Manual “Sampling, transportation and storage of clinical material for PCR diagnostics”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology”, Moscow, 2008.

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of **AmpliSens® HPV HCR screen-EPh** PCR kit is tested against predetermined specifications to ensure consistent product quality.

16. KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Sufficient for
	<i>In vitro</i> diagnostic medical device		Expiration Date
	Version		Consult instructions for use
	Temperature limitation	C+_{HPV 31, 39, 56}	Positive Control DNA <i>HPV</i> types 31, 39, 56 and DNA human
	Manufacturer	C+h	Positive Control DNA human
	Date of manufacture	NCA	Negative control of amplification
	Authorised representative in the European Community	C-	Negative control of extraction

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
06.07.10	Content	«C+ _{HPV31, 39, 56} » was added for Positive Control DNA <i>HPV</i> types 31, 39, 56 and DNA human, «C+ _h » - for Positive Control DNA human
25.12.10 KM	Through the text	Correcting the document template «Cervical or urethral scrapes» was changed to «urogenital swabs»
	Cover page	The phrase “For Professional Use Only” was added
	Intended use	The phrase “The results of PCR analysis are taken into account in complex diagnostics of disease” was added.
	Content	New sections “Working Conditions” and “Transportation” were added The “Explanation of Symbols” section was renamed to “Key to Symbols Used”
	Stability and Storage	The information about the shelf life of open reagents was added
	Key to Symbols Used	The explanation of symbols was corrected
22.06.11 VV	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science “Central Research Institute for Epidemiology”
23.07.12 IV	Content	AmpliSens [®] <i>HPV</i> HCR screen-EPh PCR kit variant 50 F was deleted
	Footer	Catalogue number REF V31-50F-CE was deleted