Hämoglobin ELISA Kit

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Aus: Hoepffner N et al. (2006) Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. *Aliment Pharmacol Ther* 23 (1):145-54

Aus: Hoepffner N et al. (2006) Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. *Aliment Pharmacol Ther* 23 (1):145-54

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
Manual

Hemoglobin ELISA Kit

For the in vitro determination of hemoglobin in stool

Valid from 2015-03-23

REF K 7816D

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	18
6. STORAGE AND PREPARATION OF SAMPLES	19
<i>Sample stability and storage</i>	19
<i>Extraction of the stool samples</i>	19
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Test procedure</i>	20
8. RESULTS	22
9. LIMITATIONS	22
10. QUALITY CONTROL	23
<i>Reference range</i>	23
11. PERFORMANCE CHARACTERISTICS	23
<i>Precision and reproducibility</i>	23
<i>Analytical Sensitivity</i>	23
<i>Dilution recovery</i>	24
12. PRECAUTIONS	25
13. TECHNICAL HINTS	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	26
<i>General literature</i>	26
<i>Literature using Immundiagnostik hemoglobin ELISA</i>	27

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of human **hemoglobin** in stool. For *in-vitro* diagnostic use only.

2. INTRODUCTION

In contrast to other commercially available Hemoglobin quick tests, this **hemoglobin ELISA** does not require previous adherence to a diet (no raw meat etc.) and recognizes human **hemoglobin** in 100-fold lower concentrations. This avoids false-negative results. Because of the choice of antibodies, false-positive results are almost excluded.

Recent data show, that the clinical specificity and sensitivity can be increased by using Hemoglobin ELISA.

Indications

- Occult blood in stool
- Crohn´s disease; ulcerative colitis
- Suspicion of colon carcinoma
- Polyps in the colorectum

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 7816D	PLATE	Holder with pre-coated strips	12 x 8
K 7816D	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 7816D	IDK Extract®	Extraction buffer concentrate <i>IDK Extract</i> ®, 2.5 x	1 x 100 ml
K 7816D	SAMPLEBUF	Dilution buffer, ready to use	2 x 15 ml
K 7816D	CONJ	Conjugate, (mouse-anti humanHb, peroxidase-labeled), ready to use	1 x 15 ml
K 7816D	STD	Standards, lyophilized (50; 10; 3.3; 0.67; 0 µg/g)	4 x 5 vials
K 7816D	CTRL	Controls, lyophilized	4 x 2 vials
K 7816D	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 7816D	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- **Preparation of the extraction buffer:** The **extraction buffer concentrate IDK Extract®** must be diluted with ultra pure water **1:2.5** before use (100 ml *IDK Extract®* + 150 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37 °C in a water bath. The **IDK Extract®** is stable at **2–8 °C** until the expiry date stated on the label. Extraction buffer (1:2.5 diluted *IDK Extract®*) can be stored in a closed flask at **2–8 °C for three months**.

- The lyophilized **STD** (standards) and **CTRL** (controls) must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards and controls are not stable.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Sample stability and storage

Raw stool can be stored for **1 month** at **-20 °C**.

Stool extract is stable at room temperature (15-30 °C) for one day, at 2–8 °C as well as at -20 °C for five days. Avoid more than three freeze-thaw cycles.

Extraction of the stool samples

We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instructions for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

SAS with 1.5 ml extraction buffer:

Applied amount of stool: 15 mg

Buffer Volume: 1.5 ml

Dilution Factor: 1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube** with **1.5 ml** of ready to use *IDK Extract*[®] extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into

sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.

- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~ 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

The sample suspension is now ready for use.

The sample can also be used in a pipetting automat. Place the sample in the sample rack according to instrument instructions.

7. ASSAY PROCEDURE

Principle of the test

This ELISA is used for quantitative determination of hemoglobin in stool. The hemoglobin in the sample is bound to anti-hemoglobin antibodies (in excess), which are immobilized on the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step an anti-hemoglobin peroxidase labeled antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine. An acidic solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of hemoglobin in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the values obtained from standard. Hemoglobin present in the patient samples is determined directly from this curve.

Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Take as many microtiter strips as needed from kit. Store unused strips covered at

2–8° C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash the pre-coated microtiter plate 5 x with 250 µl wash buffer before use . After the final washing step, the inverted microtiter plate should be tapped on absorbent paper.
2.	Add in each well 50 µl of SAMPLEBUF (sample dilution buffer)
3.	Then 50 µl of CTRL/STD/SAMPLE (control/standard/sample) into respective well (see template given in the test kit).
4.	Incubate for 1 hour at room temperature (15–30 °C).
5.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
6.	Add 100 µl conjugate in each well.
7.	Incubate for 1 hour at room temperature (15–30 °C).
8.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
9.	Add 100 µl SUB (TMB substrate).
10.	Incubate for 10–20 minutes at room temperature (15–30 °C)*
11.	Add 100 µl STOP (stop solution) and mix shortly.
12.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the „4 parameter algorithm“.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Stool

Since the sample dilution is already considered in the calibration curve, the dilution factor is **1**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this higher dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

Analytical sensitivity × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Hemoglobin (stool)* $< 2 \mu\text{g/ml} \triangleq < 2 \mu\text{g/g}$
(1 g stool \triangleq 1 ml)

* H.G. Bischoff et al., 1990. Fecal Occult Blood: A Sensitive and Specific Method for Detection by Immunological Determination of Human Albumin and Hemoglobin. *Ärztl. Lab.*, **36**. pp. 101–112.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 42)

Sample	Hemoglobin [$\mu\text{g/g}$]	CV [%]
1	3.99	4.75
2	7	2.53

Inter-Assay (n = 63)

Sample	Hemoglobin [$\mu\text{g/g}$]	CV [%]
1	0.863	8.16
2	3.501	4.76

Analytical Sensitivity

The sensitivity was set as $B_0 + 2 \text{SD}$. The zero-standard was measured 20 times. A detection limit of $0.042 \mu\text{g/ml}$ was obtained.

Dilution recovery

Two patient samples were diluted and analyzed. The results are shown below (n = 2):

Sample	Dilution	Hämoglobin expected [µg/g]	Hämoglobin measured [µg/g]
A	-	13.83	13.83
	1:2	6.92	6.74
	1:4	3.46	3.52
	1:8	1.73	1.61
B	-	17.74	17.74
	1:2	8.87	8.78
	1:4	4.44	4.47
	1:8	2.22	1.98

Sensitivities, specificities, positive predictive value (PPV) and negative predictive value (NPV) of three fecal occult blood tests (FOBTs) in different diagnostic groups.

Parameters	Guajak-based test	Bedside-IFOBT	ELISA-IFOBT
Sensitivity (%)			
Adenomas	5.56 (0.14-27.29)	18.9 (3.58-41.42)	22.2* (6.41-47.64)
Carcinomas	37.0 (24.29-51.26)	74.0* (60.35-85.04)	77.7* (64.40-87.96)
Carcinomas + adenomas	29.1 (19.05-41.07)	59.7* (47.50-71.12)	63.8* (51.71-74.88)
Colitis ulcerosa	48.2 (29.45-67.47)	93.1 (77.23-99.15)	96.5 (82.24-99.9)
Morbus Crohn	42.1 (26.31-59.18)	94.7* (82.25-99.36)	94.7* (82.25-99.36)
Specifity (%)			
Carcinomas + adenomas	90.2 (84.64-94.32)	94.5 (89.84-97.46)	96.3* (92.21-98.65)
Colitis ulcerosa	90.9 (70.84-98.88)	95.4* (77.16-99.88)	95.4 (77.16-99.88)
Morbus Crohn	89.3 (79.36-95.63)	93.9* (85.20-98.32)	95.4* (87.29-99.05)
PPW (%)			
Adenomas	5.9 (0.15-28.69)	25.0 (5.49-57.19)	40.0 (12.16-73.76)
Carcinomas	55.6 (38.10-72.06)	81.6 (67.98-91.24)	87.5 (74.75-95.27)
Carcinomas + adenomas	56.7 (39.49-72.90)	82.6 (69.67-91.77)	88.4 (76.56-95.65)
Colitis ulcerosa	87.5 (61.65-98.45)	96.4 (81.65-99.91)	96.6 (82.24-99.91)
Morbus Crohn	69.6 (47.08-86.79)	90.0 (76.34-97.21)	92.3 (79.13-98.38)

Parameters	Guajak-based test	Bedside-IFOBT	ELISA-IFOBT
NPW (%)			
Adenomas	89.7 (84.02-93.88)	91.2 (85.86-94.98)	91.9 (86.72-95.48)
Carcinomas	81.3 (74.89-86.70)	91.7 (86.49-95.40)	92.9 (87.99-96.30)
Carcinomas + adenomas	74.4 (67.72-80.28)	84.2 (78.16-89.18)	85.9 (79.99-90.56)
Colitis ulcerosa	57.1 (39.35-73.68)	91.3 (71.96-98.93)	95.4 (77.16-99.88)
Morbus Crohn	72.8 (61.81-82.13)	96.9 (89.16-99.62)	96.7 (89.32-9.63)

* $P < 0.05$ compared with Guaiac-based test. Values in parentheses indicate 95% CI. **IFOBT**, immunochemical fecal occult blood tests; **PPV**, positive predictive value; **NPV**, negative predictive value; **ELISA**, enzyme-linked immunosorbent assay; **CI**, confidence interval.

From: Hoepffner N et al. (2006) Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. *Aliment Pharmacol Ther* 23 (1):145-54

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

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
Literature using Immundiagnostik hemoglobin ELISA

4. Trojan, J. et al., 2002. A new immunological test strip device for the rapid, qualitative detection of faecal occult blood. *Zeitschrift für Gastroenterologie / German Journal of Gastroenterology*, **40**, pp.921–924.
5. Hoepffner, N. et al., 2006. Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. *Alimentary pharmacology & therapeutics*, **23**(1), pp.145–54.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by