# Histamin ELISA Kit

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
Manual

# Histamine ELISA Kit

*For the in vitro determination of histamine in stool*

Gültig ab / Valid from 2015-03-03

**Immundiagnostik AG**, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: ++49 6251 70190-0

Fax: ++ 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

<b>1. INTENDED USE</b>	<b>17</b>
<b>2. MATERIAL SUPPLIED</b>	<b>17</b>
<b>3. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>17</b>
<b>4. PREPARATION AND STORAGE OF REAGENTS</b>	<b>18</b>
<b>5. PREPARATION AND STORAGE OF SAMPLES</b>	<b>19</b>
<i>Extraction of the stool sample</i>	19
<i>Dilution of samples</i>	20
<i>Storage of samples</i>	20
<b>6. ASSAY PROCEDURE</b>	<b>21</b>
<i>Principle of the test</i>	21
<i>Sample preparation procedure</i>	21
<i>Test procedure</i>	22
<b>7. RESULTS</b>	<b>23</b>
<b>8. LIMITATIONS</b>	<b>24</b>
<b>9. QUALITY CONTROL</b>	<b>25</b>
<i>Reference range</i>	25
<b>10. PERFORMANCE CHARACTERISTICS</b>	<b>25</b>
<i>Precision and reproducibility</i>	25
<i>Spiking recovery</i>	25
<i>Dilution recovery</i>	26
<i>Analytical sensitivity</i>	26
<i>Specificity</i>	27
<b>11. PRECAUTIONS</b>	<b>27</b>
<b>12. TECHNICAL HINTS</b>	<b>27</b>
<b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>28</b>

## 1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of histamine in stool. For *in vitro* diagnostic use only.

## 2. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 8213MTP	PLATE	Holder with precoated strips	12 x 8 wells
K 8213KP	COPLATE	Plate for Derivatization	12 x 8 wells
K 8213ST	STD	Standards, ready to use (0; 1; 3; 10; 30; 120 ng/ml)	6 x 2 ml
K 8213KO1 K 8213KO2	CTRL 1 CTRL 2	Controls, ready to use (see specification for range)	2 x 2 ml
K 8213WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 8213EP	EXBUF	Extraction buffer, ready to use	2 x 75 ml
K 8213AB	AB	Histamine antibody, peroxidase-labeled, lyophilized	6 x 1 vial
K 8213VR	ABBUF	Antibody dilution buffer, ready to use	1 x 11 ml
K 8213RP	REABUF	Reaction buffer, ready to use	1 x 50 ml
K 8213DR	DER	Derivatization reagent, lyophilized	4 x 1 vial
K 8213LM	DMSO	Dimethylsulfoxide (DMSO)	1 x 10 ml
K 8213TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 8213AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

## 3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*

- Calibrated precision pipets and 10-1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 *g*
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 6)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 4. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **ELISA wash buffer concentrate** (WASHBUF) with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The buffer concentrate is stable at 2-8°C until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month**.
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20°C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use. They can be re-frozen up to 3 times.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Store the **derivatization reagent (DER)** at **-20°C**. Bring to room temperature before opening. Dissolve the content of one vial in **DMSO** as stated on the label. Allow the content of the vial to dissolve for 10 min and mix thoroughly

with a vortex-mixer. DER must be **prepared immediately before use**. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. The ELISA kit can be separated into 4 performances by providing 4 DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.

- Store the peroxidase conjugated **histamine antibody (AB)** at **-20°C**. Reconstitute with **1.5 ml of antibody dilution buffer (ABBUF)**. Allow the antibody dilution buffer to reach **room temperature** before use. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest after use. The ELISA kit can be separated into several performances by providing 6 AB vials.
- All other test reagents are ready for use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 5. PREPARATION AND STORAGE OF SAMPLES

### *Extraction of the stool sample*

The extraction buffer (EXBUF) is ready to use. We recommend the following sample preparation:

#### **Stool Sample Application System (SAS)** (Cat. No. K 6998SAS)

##### ***Stool sample tube – Instructions for use:***

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer:

##### ***SAS with 1.5 ml buffer:***

Applied amount of stool:	15 mg
Buffer volume (EXBUF):	1.5 ml
Dilution factor:	1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenization using an applicator, inoculation loop or similar device.

- b) Fill the **empty sample tube** with **1.5 ml** of extraction buffer (EXBUF) before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert the yellow dipstick into the sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place the dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped of, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

**Dilution I:            1:100**

### *Dilution of samples*

Dilute the supernatant of the sample extraction (dilution I) **1:2** with extraction buffer (EXBUF). For example:

**100 µl supernatant** (dilution I) + **100 µl EXBUF**, mix well  
= **1:2 (dilution II)**

This results in a final dilution of **1:200**.

To **25 µl of dilution II** (SAMPLE) a derivatization reagent (DER) is added for derivatization of histamine (details are given in the sample preparation procedure).

### *Storage of samples*

**Raw stool** is stable for 2 days at room temperature. For longer storage keep frozen at -20 °C.

**Stool extract** is stable for 2 days at room temperature. For longer storage keep frozen at -20 °C.

## 6. ASSAY PROCEDURE

### *Principle of the test*

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for histamine derivatization. Afterwards, the treated samples and a peroxidase conjugated polyclonal histamine antibody are incubated in wells of a microtiter plate coated with histamine-derivative (tracer). During the incubation period, the target histamine in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the histamine concentration in the sample.

After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the histamine concentration in the sample; this means, high histamine concentration in the sample reduces the concentration of tracer-bound antibody and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Histamine present in the patient samples is determined directly from this curve.

### *Sample preparation procedure*

Derivatization of standards (STD), controls (CTRL) and diluted stool samples (SAMPLE) is carried out in single analysis in the wells of the COPLATE.

Alternatively, reaction vials (e.g. 1.5 ml vials) can be used.

We recommend carrying out 48 derivatizations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1. Bring **all reagents and samples to room temperature** (15-30°C) and mix well.

2. Add **25 µl of standards (STD), 25 µl of controls (CTRL) and 25 µl of samples from dilution II (SAMPLE)** in the wells of the COPLATE or in the corresponding vials.
3. Add **250 µl of reaction buffer (REABUF)** into each well (STD, CTRL, SAMPLE) of the COPLATE or into each vial.
4. Add **50 µl** of freshly prepared **derivatization reagent (DER)** into each well (STD, CTRL, SAMPLE) of the COPLATE and incubate immediately on a **horizontal shaker** (180-240 rpm) **for 30 min at room temperature** (15-30°C).  
*Alternatively:* Add **50 µl** of freshly prepared **derivatization reagent (DER)** into each vial (STD, CTRL, SAMPLE), **mix thoroughly** by repeated inversion or several seconds on a vortex mixer and incubate for **30 min at room temperature** (15-30°C) on a horizontal **shaker** (180-240 rpm).

**2 x 50 µl of each treated sample (STD, CTRL, SAMPLE)** are used in the ELISA as duplicates.

### *Test procedure*

5. Mark the positions of standards (STD)/ controls (CTRL)/samples (SAMPLE) in duplicate on a **protocol sheet**.
6. Take as many microtiter plate strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.  
Please note: **Do not wash the plate**.
7. For the analysis in duplicate take **2 x 50 µl** of the **derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE)** out of the COPLATE or vials and add into the respective wells of the microtiter plate (PLATE).
8. Add **50 µl of histamine antibody (AB)** into each well of the microtiter plate.
9. Cover the plate tightly with foil and incubate for **1 hour at room temperature** (15-30°C) on a horizontal **shaker** (180-240 rpm).

10. Discard the contents of each well and wash **5 times** with **250 µl of diluted wash buffer**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
11. Add **100 µl** of **TMB substrate (SUB)** into each well.
12. Incubate for **12-18 min** at **room temperature** (15-30°C) in the dark\*.
13. Add **100 µl** of **stop solution (STOP)** into each well, mix thoroughly.
14. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm (690 nm) as a reference.

\*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

## 7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.


### Stool samples

The obtained histamine levels of the stool samples have to be multiplied by the dilution factor of **200** (1 ng/ml = 1 ng/g stool).

In case another dilution factor has been used, multiply the obtained result by the dilution factor used.

In the following, an example of a calibration curve is given; do not use it for the calculation of your results.

### Example of calibration curve

## 8. LIMITATIONS

Samples with an OD lower than the OD of the highest standard should be further diluted with extraction buffer (EXBUF) and re-assayed. Please consider this dilution factor when calculating the results.

## 9. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

### *Reference range*

Based on internal studies with stool samples of apparently healthy persons (n=36) a **median of 455 ng/ml** was calculated.

We recommend each laboratory to establish its own reference range.

## 10. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra assay (n=16)**

sample	histamine [ng/ml]	CV [%]
1	355	8.4
2	1234	7.3

#### **Inter assay (n=6)**

sample	histamine [ng/ml]	CV [%]
1	225	10.3
2	1012	5.7
3	2369	6.8

### *Spiking recovery*

One sample was spiked with different histamine concentrations and measured in this assay. The mean recovery rate was 94.8 % (n=9).

spike [ng/ml]	histamine expected [ng/ml]	histamine measured [ng/ml]	recovery [%]
0		262	
300	562	484	86.1
2400	2662	2755	103.5

### *Dilution recovery*

One spiked sample was diluted and measured in this assay. The mean recovery rate was 87.6 % (n=9).

dilution	histamine expected [ng/ml]	histamine measured [ng/ml]	recovery [%]
original		901	
1:2	450	374	83,1
1:4	225	207	92,1

### *Analytical sensitivity*

The zero-standard was measured 75 times. The detection limit was set as  $B_0 - 2 \text{ SD}$  and estimated to be 0,11 ng/ml. Considering the dilution factor, the detection limit is calculated to be 22 ng/ml (22 ng/g stool).

sample	mean value [OD]	2 x standard deviation (SD)	detection limit histamine (stool) [ng/ml]
zero-standard	1.991	0.19	$0.11 \times 200 = 22$

## Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to histamine. The specificity is calculated in percent in relation to the histamine binding activity.

3-methyl-histamine	< 0,1 %
tyramine	<0,001%
L-phenylalanine	<0,0002%
L-histidine	<0,0002%
L-tyrosine	<0,0002%
tryptamine	<0,0002%
5-hydroxy-Indole-acetic acid	<0,0002%
serotonine	<0,0002%

## 11. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

## 12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure which is not coordinated with the producer may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

#### **Used symbols:**

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by