IDK[®] Calprotectin ELISA Kit

***Zur in-vitro-Bestimmung von Calprotectin (MRP 8/14,
S100A8/A9) in Urin***

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
Manual

IDK[®] Calprotectin ELISA Kit

***For the in vitro determination of calprotectin (MRP 8/14,
S100A8/A9) in urine***

Gültig ab / Valid from 2015-03-11

REF **K 6928**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	15
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. STORAGE AND PREPARATION OF REAGENTS	16
6. SPECIMEN COLLECTION AND PREPARATION	17
7. ASSAY PROCEDURE	17
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
8. RESULTS	19
9. LIMITATIONS	20
10. QUALITY CONTROL	20
<i>Reference range</i>	20
11. PERFORMANCE CHARACTERISTICS	20
<i>Precision and reproducibility</i>	20
<i>Analytical Sensitivity</i>	21
<i>Spiking Recovery</i>	21
<i>Dilution recovery</i>	22
<i>Specificity</i>	22
12. PRECAUTIONS	22
13. TECHNICAL HINTS	23
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
15. REFERENCES	24
<i>General literature</i>	24
<i>Literature using Immundiagnostik Calprotectin ELISA</i>	24

1. INTENDED USE

The described calprotectin ELISA is intended for the quantitative determination of calprotectin (MRP (8/14, S100A8/A9) in urine. For *in vitro* diagnostic use only.

2. INTRODUCTION

Alternative names of calprotectin:

- MRP8/14, L1, (p8,14), p34, S100 A8/A9

Alternative names of the two proteins forming the heterocomplex calprotectin:

- S100A8, calgranulin A, MRP8 (migration inhibition factor-related protein-8), CP-10 (in mouse)
- S100A9, calgranulin B, MRP14 (migration inhibition factor-related protein-14)

Calprotectin is a heterodimeric complex, which consists of the calcium-binding proteins S100A8 and S100A9. It is expressed in a tissue-specific manner mainly in cells of the myeloid lineage, including monocytes, neutrophils and early differentiation states of macrophages. However, expression of calprotectin is also inducible in mature macrophages, osteoclasts, keratinocytes, fibroblasts and microvascular endothelial cells (Foell, Roth, 2004).

Calprotectin is involved in innate immunity, leukocyte adhesion, and endothelial transmigration. In addition, it also emerges as important mediator of diverse processes within chronic inflammation, with recent attention being focused on its involvement in inflammation-associated carcinogenesis (Gebhardt et al., 2006).

Indeed, strong up-regulation of the S100A8/A9 proteins has been observed in many tumors, including gastric, colon, pancreatic, bladder, ovarian, thyroid, breast and skin cancers, with both tumor cells and infiltrating immune cells in the tumor stroma being a source of calprotectin (Srikrishna, 2012).

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6928	PLATE	Holder with precoated strips	12 x 8 wells
K 6928	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6928	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 100 ml
K 6928	STD	Calprotectin standards, lyophilized (0; 13; 52; 210; 840 ng/ml)	2 x 5 vials

Cat. No.	Label	Kit components	Quantity
K 6928	CTRL 1	Control, lyophilized (see specification for range)	2 x 1 vial
K 6928	CTRL 2	Control, lyophilized (see specification for range)	2 x 1 vial
K 6928	CONJ	Conjugate, ready to use	15 ml
K 6928	SUB	TMB substrate (Tetramethylbenzidine), ready to use	15 ml
K 6928	STOP	ELISA stop solution, ready to use	15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Multi-channel pipets or repeater pipets
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≤ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions.

The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.

- The **lyophilized STD** (standards) and **CTRL** (controls) are stable at **2–8 °C** until the expiry date stated on the label. Before use, the **STD** (standards) and **CTRL** (controls) must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards and controls can be stored at 2–8 °C for four weeks**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. SPECIMEN COLLECTION AND PREPARATION

Urine samples should be diluted **1:10 with sample dilution buffer** (SAMPLEBUF) before performing the assay.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the two-site sandwich technique with two selected monoclonal antibodies that bind to human calprotectin.

Standards, controls and diluted patient samples which are assayed for human calprotectin are added to wells of microplate coated with a high affine monoclonal anti-human calprotectin antibody. During the first incubation step, calprotectin in the samples is bound by the immobilized antibody. Then a peroxidase labeled conjugate is added to each well and the following complex is formed: capture antibody – human calprotectin – peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the calprotectin concentration of sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the values obtained from standard. Calprotectin present in the patient samples is determined directly from this curve.

Test procedure

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Bring all reagents and samples to room temperature (15–30 °C) and mix well.
2.	Mark the positions of STD / SAMPLE / CTRL (standard / samples / controls) on a protocol sheet.
3.	Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label.
4.	Add 100 µl of STD / SAMPLE / CTRL into the respective wells.
5.	Cover plate tightly and incubate for 30 minutes at room temperature (15–30 °C).
6.	Aspirate the contents of each well. Wash each well 5 x with 250 µl of wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
7.	Add 100 µl CONJ (conjugate) into each well.
8.	Cover plate tightly and incubate for 30 minutes at room temperature (15–30 °C).
9.	Aspirate the contents of each well. Wash each well 5 x with 250 µl of wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
10.	Add 100 µl of SUB (substrate) into each well.
11.	Incubate for 10–20 minutes at room temperature (15–30 °C) in the dark*.
12.	Add 100 µl of STOP (stop solution) into each well, mix thoroughly.

- | | |
|-----|---|
| 13. | Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference. |
|-----|---|

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Urine

For calculation of calprotectin concentration in urine, the result must be multiplied by the dilution factor of **10**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

- The median value of calprotectin in urine in healthy adults is about 51 ng/ml (interquartile range 20.5–86.6 ng/ml).
- A cut-off of 140 ng/ml has been shown to differentiate best between patients with urothelial bladder carcinoma and healthy subjects (Ebbing et al).

We recommend each laboratory to establish its own reference concentration range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 27)

Sample	Calprotectin [ng/ml]	CV [%]
1	416.1	6.499
2	123.9	5.257

Inter-Assay (n = 12)

Sample	Calprotectin [ng/ml]	CV [%]
1	2289.7	10.7
2	202.6	12.0

Analytical Sensitivity

The LoB (limit of blank) was evaluated according to the guideline CLSI EP17-A2 and resulted in 0,987 ng/ml.

The LoD (limit of detection) was evaluated according to the guideline CLSI EP17-A2 and resulted in 1,132 ng/ml.

Spiking Recovery

Two samples were spiked with different calprotectin concentrations and measured using this assay (n = 2).

Sample	Unspiked Sample [ng/ml]	Spike [ng/ml]	Calprotectin expected [ng/ml]	Calprotectin measured [ng/ml]
A	18.5	25	43.5	47.1
	18.5	12.5	31	30.1
	18.5	6.25	24.7	23.2
	18.5	3.125	21.6	20.2
B	43.6	25	68.6	73.1
	43.6	12.5	56.1	56.3
	43.6	6.25	49.8	49.6
	43.6	3.125	46.7	45.8

Dilution recovery

Two samples were diluted and analyzed. The results are shown below (n = 2):

Sample	Dilution	Calprotectin expected [ng/ml]	Calprotectin measured [ng/ml]
A	1:10		1544.1
	1:20	772	733.5
	1:40	386	354.4
	1:80	193	175.0
	1:160	96	75.4
B	1:10		4334.5
	1:20	2167	2180.9
	1:40	1083	1048.9
	1:80	541	489.1
	1:160	270	223.7

Specificity

No cross-reactivity was observed to the following proteins:

- Lysozyme 0%
- PMN-Elastase 0%
- Myeloperoxidase 0%
- Lactoferrin 0%

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be

handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

General literature

1. Foell, Dirk, and Johannes Roth. 2004. "Proinflammatory S100 Proteins in Arthritis and Autoimmune Disease." *Arthritis and Rheumatism* **50** (12) (December): 3762–71. doi:10.1002/art.20631.
2. Gebhardt, Christoffer, Julia Németh, Peter Angel, and Jochen Hess. 2006. "S100A8 and S100A9 in Inflammation and Cancer." *Biochemical Pharmacology* **72** (11) (November 30): 1622–31. doi:10.1016/j.bcp.2006.05.017.
3. Srikrishna, Geetha. 2012. "S100A8 and S100A9: New Insights into Their Roles in Malignancy." *Journal of Innate Immunity* **4** (1) (January): 31–40. doi:10.1159/000330095.

Literature using Immundiagnostik Calprotectin ELISA

4. Ebbing, Jan, Susanne Mathia, Felix S Seibert, Nikolaos Pagonas, Frederic Bauer, Barbara Erber, Karsten Günzel, et al. 2013. "Urinary Calprotectin: A New Diagnostic Marker in Urothelial Carcinoma of the Bladder." *World Journal of Urology* [Epub ahead of print]. doi:10.1007/s00345-013-1227-8.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by