IDK[®] Calprotectin ELISA Kit

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15. LITERATUR

Allgemeine Literatur

1. Stríz, I. & Trebichavský, I. Calprotectin - a pleiotropic molecule in acute and chronic inflammation. *Physiological research / Academia Scientiarum Bohemoslovaca* **53**, 245–53 (2004).
2. Yui, S., Nakatani, Y. & Mikami, M. Calprotectin (S100A8/S100A9), an inflammatory protein complex from neutrophils with a broad apoptosis-inducing activity. *Biological & pharmaceutical bulletin* **26**, 753–60 (2003).

3. Odink, K. et al. Two calcium-binding proteins in infiltrate macrophages of rheumatoid arthritis. *Nature* **330**, 80–2 (1987).
4. Wilkinson, M. M. et al. Expression pattern of two related cystic fibrosis-associated calcium-binding proteins in normal and abnormal tissues. *Journal of cell science* **91** (Pt 2), 221–30 (1988).
5. Müller, F., Frøland, S. S., Aukrust, P. & Fagerhol, M. K. Elevated serum calprotectin levels in HIV-infected patients: the calprotectin response during ZDV treatment is associated with clinical events. *Journal of acquired immune deficiency syndromes* **7**, 931–9 (1994).
6. Hermani, A. et al. Calcium-binding proteins S100A8 and S100A9 as novel diagnostic markers in human prostate cancer. *Clinical cancer research : an official journal of the American Association for Cancer Research* **11**, 5146–52 (2005).

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7. Heller, F., Frischmann, S., Grünbaum, M., Zidek, W. & Westhoff, T. H. Urinary calprotectin and the distinction between prerenal and intrinsic acute kidney injury. *Clinical journal of the American Society of Nephrology : CJASN* **6**, 2347–55 (2011).
8. Lönnkvist, M. H., Theodorsson, E., Holst, M., Ljung, T. & Hellström, P. M. Blood chemistry markers for evaluation of inflammatory activity in Crohn's disease during infliximab therapy. *Scandinavian journal of gastroenterology* **46**, 420–7 (2011).
9. Börekçi, B., Aksoy, H., Öztürk, N. & Kadanali, S. Correlation between Calprotectin and Oxidized LDL in Preeclampsia. *Turkish Journal of Medical Sciences* **39**, 191–195 (2009).

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
Manual

IDK[®] Calprotectin ELISA Kit

***For the in vitro determination of calprotectin (MRP 8/14,
S100A8/A9) in serum and plasma***

Gültig ab / Valid from 2015-03-06

REF **K 6935**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. CLINICAL RELEVANCE	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	19
5. PREPARATION AND STORAGE OF REAGENTS	19
6. SAMPLE COLLECTION AND PREPARATION	20
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Test procedure</i>	21
8. RESULTS	22
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Reference range</i>	23
11. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Specificity</i>	24
<i>Analytical Sensitivity</i>	24
<i>Spiking Recovery</i>	25
<i>Dilution recovery</i>	25
12. PRECAUTIONS	26
13. TECHNICAL HINTS	26
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	27
<i>General literature</i>	27
<i>Literature using Immundiagnostik Calprotectin ELISA</i>	27

1. INTENDED USE

The described **IDK®** Calprotectin ELISA is intended for the quantitative determination of calprotectin (MRP (8/14, S100A8/A9) in serum and plasma. For *in vitro* diagnostic use only.

2. CLINICAL RELEVANCE

Alternative names of calprotectin:

- MRP8/14, L1, (p8,14), p34, S100 A8/A9

Alternative names of the two proteins forming the heterocomplex calprotectin:

- S100A8, Calgranulin A, MRP8 (Migration inhibition factor-related protein-8), CP-10 (in mouse)
- S100A9, Calgranulin B, MRP14 (Migration inhibition factor-related protein-14)

Calprotectin is a calcium-binding protein secreted predominantly by neutrophils and monocytes. The heterocomplex consists of the two proteins, S100A8 (calgranulin A) and S100A9 (calgranulin B), also designated as MRP8 and MRP14, respectively. Expression of S100A8 and S100A9 in epithelial tissues was first described in context with squamous epithelia and with murine and human wound repair. More recently, an association of S100 protein expression with adenocarcinomas in humans has emerged. The genes S100A8 and S100A9 are located in a gene cluster on chromosome 1q21, a region in which several rearrangements that occur during tumor development have been observed.

Elevated MRP8/14 levels have been found in many sites of inflammation and in the extracellular fluid of patients with many types of inflammatory conditions. The concentration of MRP8/14 in blood is increased in patients with rheumatoid arthritis, cystic fibrosis, multiple sclerosis, and HIV infections, while elevated MRP8/14 levels have been detected in stool of patients with Crohn's disease and colorectal cancer [1-5]. Extracellular MRP8/14 has antimicrobial, antigrowth and apoptotic effects. It suppresses the growth of some fungi and bacteria [1,2]. It also suppresses the proliferation of several different types of cells including: macrophages, lymphocytes, hematopoietic progenitors, and tumor cell lines. MRP8/14 can also induce apoptosis of some tumor cell lines [1,2].

Hermani et al. (2005) [6] reported recently that enhanced expression of S100A8 and S100A9 is an early event in prostate tumor genesis and may contribute to development and progression or extension of prostate carcinomas. Furthermore, they tested the value of S100A9 as a serum marker for prostate cancer comparing the serum concentrations of S100A9 in cancer patients with healthy controls or patients with benign prostatic hyperplasia (BPH). Significantly increased S100A9 serum levels in

prostate cancer were found in prostate cancer patients compared to patients with BPH, the latter exhibiting values similar to that obtained for healthy individuals.

Pathological significance and clinical application

The diagnostic value and advantage of MRP8/14 over other disease markers is that they are preformed and released immediately upon activation of the respective cell population. Other markers may be generated in downstream events or need to be synthesized de novo in the liver. Various conditions have shown significant correlation of MRP8/14 (or MRP8, MRP14) levels with disease activity:

- Concentrations of MRP8/14 in serum, and particularly in synovial fluid, correlate strongly with disease activity in rheumatoid arthritis.
- Plasma MRP8/14 levels are very early, specific and sensitive prediction markers for acute rejection in kidney allograft transplantation.
- Serum MRP8/14 concentration is a prognostic marker of recurrent infection and survival in alcoholic liver cirrhosis.
- MRP8/14 is useful for evaluating the extent of periodontal inflammation.
- In cerebral malaria, MRP 8/14 expression correlates with microglial activation in brain.
- MRP8/14 is present in urinary stones and in dental calculus.
- S100A9 in serum may serve as a useful marker for discrimination between prostate cancer and benign prostatic hyperplasia (BPH).

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6935	PLATE	Holder with precoated strips	12 x 8 wells
K 6935	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6935	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 100 ml
K 6935	STD	Calprotectin standards, lyophilized (0; 3,9; 15,6; 62,5; 250 ng/ml)	2 x 5 vials
K 6935	CTRL 1	Control, lyophilized (see specification for range)	2 x 1 vial
K 6935	CTRL 2	Control, lyophilized (see specification for range)	2 x 1 vial
K 6935	CONJ	Conjugate, ready to use	15 ml
K 6935	SUB	TMB substrate (tetramethylbenzidine), ready to use	15 ml

Cat. No.	Label	Kit components	Quantity
K 6935	STOP	ELISA stop solution, ready to use	15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month.**
- The **lyophilized STD** (standards) and **CTRL** (controls) are stable at **2–8 °C** until the expiry date stated on the label. Before use, the **STD** (standards) and **CTRL** (controls) must be reconstituted with **500 µl of ultra pure water.** Allow the

vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards and controls can be stored at 2–8 °C for four weeks.**

- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. SAMPLE COLLECTION AND PREPARATION

Preanalytic handling

Significant differences in the calprotectin levels can be observed due to different sample preparation procedures, e. g. up to 10-fold higher serum levels compared to the plasma calprotectin concentrations. The reasons are as follows:

Granulocytes are activated during serum clotting and release granulocyte-activating markers. The time between serum collecting and analysis as well as repeated freeze-thaw cycles don't cause a calprotectin concentration shift.

On the contrary, in the case of plasma samples, varying the time between sampling and analysis or the number of freeze-thaw cycles will cause variation in the observed calprotectin levels. Therefore, the preanalytical conditions of plasma samples should be held constant. This is a general requirement independent of the used test-system.

Immundiagnostik recommends the use of serum samples for calprotectin determinations.

Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.

Serum samples

Serum samples should be diluted **1:100 with sample dilution buffer** before performing the assay, e. g.

50 µl sample + 450 µl SAMPLEBUF = **dilution I (1:10)**

50 µl dilution I + 450 µl SAMPLEBUF = **dilution II (1:10)**

Final dilution 1:100

Plasma samples

EDTA plasma samples should be diluted **1:30 with sample dilution buffer** before performing the assay (e. g. 20 µl sample + 580 µl SAMPLEBUF).

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the two-site sandwich technique with two selected monoclonal antibodies that bind to human Calprotectin.

Standards, controls and diluted patient samples which are assayed for human calprotectin are added to wells of microplate coated with a high affine monoclonal anti-human calprotectin antibody. During the first incubation step, calprotectin in the samples is bound by the immobilized antibody. Then a peroxidase labeled conjugate is added to each well and the following complex is formed: capture antibody - human calprotectin – Peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the calprotectin concentration of sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standard. The calprotectin concentration of the patient samples is determined directly from this curve.

Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Add 100 µl of STD / SAMPLE / CTRL (standard / sample / controls) into the respective wells.
2.	Cover plate tightly and incubate for 30 min at room temperature .
3.	Aspirate the contents of each well. Wash 5 times by dispensing 250 µl wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
4.	Add 100 µl CONJ (conjugate) into each well.
5.	Cover plate tightly and incubate for 30 min at room temperature .

6.	Aspirate the contents of each well. Wash 5 times by dispensing 250 µl wash buffer into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
7.	Add 100 µl of SUB (substrate) into each well.
8.	Incubate for 10–20 minutes at room temperature (15–30 °C) in the dark* .
9.	Add 100 µl of STOP (stop solution) into each well, mix thoroughly.
10.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the „4 parameter algorithm“.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Serum

For calculation of calprotectin concentration in serum, the result must be multiplied by the **dilution factor of 100**.

EDTA plasma

For calculation of calprotectin concentration in plasma, the result must be multiplied by the **dilution factor of 30**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Calprotectin in serum of healthy persons: < 3 µg/ml (< 3000 ng/ml)

We recommend each laboratory to establish its own reference concentration range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 80)

Sample	Calprotectin [ng/ml]	VK [%]
1	447.7	3.7
2	784.6	7.9

Inter-Assay (n = 10)

Sample	Calprotectin [ng/ml]	VK [%]
1	485.28	5.5
2	816.11	7.0

Specificity

No cross-reactivity with MPR 8/14 in mouse serum was observed.

No cross-reactivity was observed to the following plasma proteins:

- Lysozyme 0%
- PMN-Elastase 0%
- Myeloperoxidase 0%
- Lactoferrin 0%

Analytical Sensitivity

Limit of blank, LoB 0,52 ng/ml

Limit of detection, LoD 0,78 ng/ml

Limit of quantitation, LoQ 0,78 ng/ml

The evaluation was performed according to the CLSI-Guideline:EP-17-A.

Spiking Recovery

Two samples were spiked with different calprotectin concentrations and measured using this assay (n = 2).

Sample	Sample unspiked [ng/ml]	Spike [ng/ml]	Calprotectin expected [ng/ml]	Calprotectin measured [ng/ml]
1	4.29	8.36	12.65	11.39
	4.25	13.19	17.44	16.41
	4.21	16.83	21.04	19.32
	4.13	24.12	28.25	27.21
2	7.35	8.36	15.71	15.42
	7.28	13.19	20.47	20.27
	7.21	16.83	24.04	23.98
	7.08	24.12	31.20	30.27

Dilution recovery

Two serum samples were diluted and analyzed. The results are shown below (n = 2):

Sample	Dilution	Calprotectin expected [ng/ml]	Calprotectin measured [ng/ml]
A	1:50		9.34
	1:100	4.67	4.53
	1:200	2.34	2.27
	1:400	1.17	1.13
B	1:50		18.17
	1:100	9.09	8.9
	1:200	4.54	4.25
	1:400	2.27	2.11

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not to assemble wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Quality control guidelines should be followed.
- IDK® is a trademark of Immundiagnostik AG.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

General literature

1. Stríz, I. & Trebichavský, I. Calprotectin - a pleiotropic molecule in acute and chronic inflammation. *Physiological research / Academia Scientiarum Bohemoslovaca* **53**, 245–53 (2004).
2. Yui, S., Nakatani, Y. & Mikami, M. Calprotectin (S100A8/S100A9), an inflammatory protein complex from neutrophils with a broad apoptosis-inducing activity. *Biological & pharmaceutical bulletin* **26**, 753–60 (2003).
3. Odink, K. et al. Two calcium-binding proteins in infiltrate macrophages of rheumatoid arthritis. *Nature* **330**, 80–2 (1987).
4. Wilkinson, M. M. et al. Expression pattern of two related cystic fibrosis-associated calcium-binding proteins in normal and abnormal tissues. *Journal of cell science* **91** (Pt 2), 221–30 (1988).
5. Müller, F., Frøland, S. S., Aukrust, P. & Fagerhol, M. K. Elevated serum calprotectin levels in HIV-infected patients: the calprotectin response during ZDV treatment is associated with clinical events. *Journal of acquired immune deficiency syndromes* **7**, 931–9 (1994).
6. Hermani, A. et al. Calcium-binding proteins S100A8 and S100A9 as novel diagnostic markers in human prostate cancer. *Clinical cancer research : an official journal of the American Association for Cancer Research* **11**, 5146–52 (2005).

Literature using Immundiagnostik Calprotectin ELISA

7. Heller, F., Frischmann, S., Grünbaum, M., Zidek, W. & Westhoff, T. H. Urinary calprotectin and the distinction between prerenal and intrinsic acute kidney injury. *Clinical journal of the American Society of Nephrology : CJASN* **6**, 2347–55 (2011).
8. Lönnkvist, M. H., Theodorsson, E., Holst, M., Ljung, T. & Hellström, P. M. Blood

chemistry markers for evaluation of inflammatory activity in Crohn's disease during infliximab therapy. *Scandinavian journal of gastroenterology* **46**, 420–7 (2011).

9. Börekçi, B., Aksoy, H., Öztürk, N. & Kadanali, S. Correlation between Calprotectin and Oxidized LDL in Preeclampsia. *Turkish Journal of Medical Sciences* **39**, 191–195 (2009).

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by