# **IDKmonitor<sup>®</sup> Adalimumab free ADA ELISA**

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van Doorn, J. D. Bos, and G.-J. Wolbink, "Extent and clinical consequences of anti-body formation against adalimumab in patients with plaque psoriasis.," *Archives of dermatology*, vol. **146**, no. 2, pp. 127–32, Feb. 2010.

3. P. A. van Schouwenburg, T. Rispens, and G. J. Wolbink, "Immunogenicity of anti-TNF biologic therapies for rheumatoid arthritis.," *Nature reviews. Rheumatology*, vol. **9**, no. 3, pp. 164–72, Mar. 2013.
4. I. Ordás, D. R. Mould, B. G. Feagan, and W. J. Sandborn, "Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms.," *Clinical pharmacology and therapeutics*, vol. **91**, no. 4, pp. 635–46, Apr. 2012.
5. N. K. Bender, C. E. Heilig, B. Dröll, J. Wohlgemuth, F. P. Armbruster, and B. Heilig, "Immunogenicity, efficacy and adverse events of adalimumab in RA patients.," *Rheumatology international*, vol. **27**, no. 3, pp. 269–74, Jan. 2007

### Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
# **IDKmonitor<sup>®</sup> Adalimumab free ADA ELISA**

***For the in vitro determination of free human antibodies against  
adalimumab (e. g. HUMIRA<sup>®</sup>) in EDTA plasma and serum***

Valid from 2015-08-26

**K 9652**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

<b>1. INTENDED USE</b>	<b>13</b>
<b>2. INTRODUCTION</b>	<b>13</b>
<b>3. MATERIAL SUPPLIED</b>	<b>13</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>14</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>14</b>
<b>6. SPECIMEN COLLECTION AND PREPARATION</b>	<b>15</b>
<b>7. ASSAY PROCEDURE</b>	<b>15</b>
<i>Principle of the test</i>	15
<i>Test procedure</i>	16
<b>8. RESULTS</b>	<b>17</b>
<b>9. LIMITATIONS</b>	<b>17</b>
<b>10. QUALITY CONTROL</b>	<b>18</b>
<i>Reference range</i>	18
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>18</b>
<i>Precision and reproducibility</i>	18
<b>12. PRECAUTIONS</b>	<b>18</b>
<b>13. TECHNICAL HINTS</b>	<b>19</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>19</b>
<b>15. REFERENCES</b>	<b>20</b>

## 1. INTENDED USE

This enzyme-linked immunosorbent assay (ELISA) kit is intended for the determination of free anti-drug antibodies (ADA) against the therapeutic TNF $\alpha$  antibody adalimumab (e.g. HUMIRA®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis, rheumatoid arthritis or psoriasis are often treated with anti-TNF $\alpha$  antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF $\alpha$  therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF $\alpha$  antibody [1-3]. Several factors influence the trough level, among them dosage and frequency of anti-TNF $\alpha$  blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation anti-drug antibodies, ADA). It is thought that ADA functionally neutralize the therapeutic antibodies or induce their rapid elimination [3]. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF $\alpha$  antibody application [1, 5].

The IDKmonitor® Adalimumab free ADA ELISA for the detection of antibodies against Adalimumab (e.g. HUMIRA®) measures free anti-Adalimumab antibodies. In combination with the drug level determination of Adalimumab, the IDKmonitor® Adalimumab free ADA ELISA is an opportunity for the treating physician to monitor and optimize the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9652	PLATE	Holder with strips, precoated with (F(ab) <sub>2</sub> )	12 x 8 wells
K 9652	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K 9652	CONJ	POD antibody, (therapy antibody, peroxidase labelled), concentrate	1 x 200 $\mu$ l
K 9652	CTRL NEG	Negative control, lyophilised	4 x 1 vial
K 9652	CTRL POS	Positive control, lyophilised	4 x 1 vial
K 9652	CTRL CUT-OFF	Cut-off control, lyophilised	4 x 1 vial
K 9652	SAMPLEBUF	Sample buffer, ready to use	1 x 30 ml

Cat. No.	Label	Kit components	Quantity
K 9652	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 9652	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Calibrated precision pipettors and 10–1000 µl tips
- Absorbent paper
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month.**

- The lyophilized **CTRL NEG, CTRL POS and CTRL CUT-OFF** (controls) must be reconstituted with **300 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes at room temperature, and mix thoroughly by gentle inversion to insure complete reconstitution. The undiluted controls are stable at **2–8 °C** until the expiry date stated on the label. **Reconstituted controls cannot be stored.**
- **Preparation of the conjugate:** The **conjugate concentrate (CONJ)** must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The CONJ is stable at **2–8 °C** until expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 6. SPECIMEN COLLECTION AND PREPARATION

### EDTA plasma and serum

EDTA plasma or serum samples must be diluted **1:10** before performing the assay, e.g. **25 µl** sample + **225 µl** SAMPLEBUF (sample buffer), mix well.

For testing in duplicates, pipet **2 x 100 µl** of each prepared sample.

### Sample storage

Freshly collected EDTA plasma or serum can be stored for one day at room temperature (15–30 °C) or for longer storage at -20 °C.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This enzyme immunoassay is a sandwich assay for the determination of free antibodies against adalimumab (e.g. HUMIRA®). In a first incubation step, the free anti-therapeutic antibodies from the sample are bound to the adalimumab F(ab)<sub>2</sub> fragments coated on the plate. To remove all unbound substances, a washing step is carried out. In a further incubation step, peroxidase-labeled adalimumab is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added. The colour converts to yellow. The absorbance of the color compound is determined photometrically. The intensity of the color is directly proportional to the amount of bound ADAs (here: anti-adalimumab antibodies) from the sample. The results are evaluated by a cut-off control.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash the microtiter strips <b>5 x with 250 µl wash buffer before use</b> . After the final washing step, the inverted microtiter strips should be tapped on absorbent paper.
2.	Add <b>100 µl of CTRL NEG, CTRL POS, CTRL CUT-OFF</b> (controls), and <b>diluted samples</b> in the wells of the microtiter plate.
3.	Seal the stripes with foil and incubate <b>over night (16–20 h)</b> , on a horizontal mixer, at 2–8 °C.*
4.	Aspirate the content of the plate and wash each well <b>5 x with 250 µl wash buffer</b> . After the final washing step, the inverted microtiter should be firmly tapped on absorbent paper.
5.	Add <b>100 µl conjugate</b> into each well.
6.	Seal the stripes with foil and incubate for <b>1 hour</b> shaking on a horizontal mixer at room temperature (15–30 °C).
7.	Aspirate the content of the plate and wash each well <b>5 x with 250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add <b>100 µl of TMB substrate</b> solution into each well.
9.	Incubate for <b>10–20 minutes</b> at room temperature in the dark.**
10.	Add <b>100 µl stop solution</b> into each well and mix shortly.

- |     |  |
|-----|--|
| 11. | Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> . If the highest extinction of the standards ( <b>STD</b> ) is above the range of the photometer, absorption must be measured immediately at <b>405 nm</b> and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used. |
|-----|--|

\* The above incubation step at 2–8 °C on a horizontal mixer is recommended by the producer. If there is no possibility to incubate at 2–8 °C, while shaking, we recommend to incubate at 2–8 °C without any shaking.

\*\* The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.

## 8. RESULTS

The analysis of the results is done using the cut-off control. Samples with a higher optical density (OD) as the OD of the cut-off control are positive. Samples with an OD lower than the OD of the cut-off control are negative.

### **Cut-off = 10 AU/ml = OD of cut-off control**

For the calculation of the sample concentrations, linear regression using a linear ordinate and abscissa is recommended.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### **Sample calculation for a positive sample**

average OD of patient's sample	0.735
average OD of cut-off control	0,065
Concentration of the patient's sample	$\frac{0.735 \times 10 \text{ AU/ml}}{0.065} = 113 \text{ AU/ml}$

## 9. LIMITATIONS

The lower limit of the measurement range can be calculated as:

$$LoB \times \text{sample dilution factor to be used}$$

Samples with concentrations lower than the measurement range cannot be clearly quantified.

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Provided control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay (n = 30)**

Sample	Adalimumab [AU/ml]	VK [%]
1	87.5	1.7
2	40.8	0.9

#### **Inter-Assay (n = 12)**

Sample	Adalimumab [AU/ml]	VK [%]
1	160.8	6
2	86.8	6
3	49.4	5
4	23.3	7
5	77.6	7

## 12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not to assemble wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions as sealed ones.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

### 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Quality control guidelines should be followed.
- *IDKmonitor*® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure,

which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. W. Afif, E. V Loftus, W. a Faubion, S. V Kane, D. H. Bruining, K. a Hanson, and W. J. Sandborn, "Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease.," *The American journal of gastroenterology*, vol. **105**, no. 5, pp. 1133–9, May 2010.
2. L. L. A. Lecluse, R. J. B. Driessen, P. I. Spuls, E. M. G. J. de Jong, S. O. Stapel, M. B. A. van Doorn, J. D. Bos, and G.-J. Wolbink, "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis.," *Archives of dermatology*, vol. **146**, no. 2, pp. 127–32, Feb. 2010.
3. P. A. van Schouwenburg, T. Rispens, and G. J. Wolbink, "Immunogenicity of anti-TNF biologic therapies for rheumatoid arthritis.," *Nature reviews. Rheumatology*, vol. **9**, no. 3, pp. 164–72, Mar. 2013.
4. I. Ordás, D. R. Mould, B. G. Feagan, and W. J. Sandborn, "Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms.," *Clinical pharmacology and therapeutics*, vol. **91**, no. 4, pp. 635–46, Apr. 2012.
5. N. K. Bender, C. E. Heilig, B. Dröll, J. Wohlgemuth, F. P. Armbruster, and B. Heilig, "Immunogenicity, efficacy and adverse events of adalimumab in RA patients.," *Rheumatology international*, vol. **27**, no. 3, pp. 269–74, Jan. 2007

### Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by