IDKmonitor[®] Adalimumab total ADA ELISA

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15. LITERATUR

1. Afif W, Loftus EV, Jr., Faubion WA, Kane SV, Bruining DH, Hanson KA, Sandborn WJ (2010). Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *Am J Gastroenterol* **105**(5): 1133-1139.
2. Kopylov U, Mazor Y, Yavzori M, Fudim E, Katz L, Coscas D, Picard O, Chowers Y, Eliakim R, Ben-Horin S (2012). Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-infliximab antibodies. *Inflamm Bowel Dis* **18**(9): 1628-1633.
3. Tak PP (2012). A personalized medicine approach to biological treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology* **51**(4): 600-609.
4. Ordas I, Mould DR, Feagan BG, Sandborn WJ (2012) Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clin Pharmacol Ther* **91**(4): 635-646.
5. Bender NK, Heilig CE, Droll B, Wohlgemuth J, Armbruster FP, Heilig B (2007). Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatol Int* **27**(3): 269-274.

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
IDKmonitor[®] Adalimumab total ADA ELISA

*For the in vitro determination of total human antibodies
against adalimumab (e. g. HUMIRA[®]) in EDTA plasma
and serum*

Valid from 2015-03-05

REF K 9651

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	13
2. INTRODUCTION	13
3. MATERIAL SUPPLIED	13
4. MATERIAL REQUIRED BUT NOT SUPPLIED	14
5. STORAGE AND PREPARATION OF REAGENTS	14
6. STORAGE AND PREPARATION OF SAMPLES	15
<i>Storage of samples</i>	15
<i>Preparation of samples and controls</i>	15
7. ASSAY PROCEDURE	16
<i>Principle of the test</i>	16
<i>Test procedure</i>	16
8. RESULTS	17
9. LIMITATIONS	17
10. QUALITY CONTROL	17
<i>Reference range</i>	18
11. PERFORMANCE CHARACTERISTICS	18
<i>Precision and reproducibility</i>	18
<i>Analytical Sensitivity</i>	18
<i>Specificity</i>	18
12. PRECAUTIONS	18
13. TECHNICAL HINTS	19
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	19
15. REFERENCES	20

1. INTENDED USE

This ELISA is intended for the determination of human antibodies against TNF α blocker adalimumab (e.g. HUMIRA®) in the presence of adalimumab in EDTA plasma and serum. For *in vitro* diagnostic use only.

2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis or rheumatoid arthritis are often treated with anti-TNF α antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF α therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF α antibody [1–3]. Several factors influence the trough level, among them dosage and frequency of anti-TNF α blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation of anti-drug antibodies, ADA). It is thought that ADA functionally neutralize the therapeutic antibodies or induce their rapid elimination [4]. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF α antibody application [1, 5].

The IDKmonitor® Adalimumab total ADA ELISA for the detection of total antibodies against adalimumab (e.g. HUMIRA®) measures free and bound antibodies against adalimumab. This assay allows a reliable determination of ADA even in the presence of adalimumab; therefore it is ideal for therapy monitoring when a measurable adalimumab concentration is expected. In combination with the drug level determination, the IDKmonitor® Adalimumab total ADA ELISA is an opportunity for the treating physician to monitor and optimize the therapy early on.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9651	PLATE	Microtiter plate, coated	12 x 8 wells
K 9651	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 9651	CONJ	Conjugate concentrate	1 x 600 μ l
K 9651	TRACER	Tracer concentrate	1 x 600 μ l
K 9651	CTRL NEG	Negative control, lyophilized	4 x 1 vial
K 9651	CTRL POS	Positive control, lyophilized	4 x 1 vial
K 9651	CTRL CUTOFF	Cutoff control, lyophilized	4 x 1 vial
K 9651	ASYBUF	Assay buffer, ready to use	1 x 30 ml

Cat. No.	Label	Kit components	Quantity
K 9651	ABBUF	Antibody dilution buffer, ready to use	1 x 10 ml
K 9651	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 9651	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label.

Wash buffer (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.

- **A few minutes before use, TRACER** (tracer concentrate) and **CONJ** (conjugate concentrate) must be diluted **1:12** in ABBUF (antibody dilution buffer) (e.g. 3000 µl ABBUF + 300 µl TRACER + 300 µl CONJ). TRACER and CONJ are stable at **2–8 °C** until expiry date stated on the label. **Diluted tracer and conjugate are not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Storage of samples

Undiluted samples can be stored **1 month at -20 °C**.

Diluted samples are **not stable** and **cannot be stored**.

Preparation of samples and controls

1.	<p>Dilute samples 1:10 in ASYBUF (assay buffer) by pipetting 25 µl sample into a reaction tube and adding 225 µl ASYBUF. Mix well. Addition of ASYBUF to all the samples should be performed without pause since this step dissociates the antibody–therapeutic antibody complexes.</p> <p>Reconstitute the controls with 250 µl ASYBUF and vortex. Carry out this step simultaneously with sample dilution in order to ensure equal treatment of controls and samples.</p>
2.	<p>Incubate controls and diluted samples in reaction tubes for 20 min with shaking on a horizontal shaker at room temperature.</p> <p>CAUTION: incubation time begins upon addition of ASYBUF.</p>
3.	<p>Add 60 µl TRACER/CONJ/ABBUF solution to 250 µl control/diluted sample. Vortex and incubate for 1 hour with shaking at room temperature (15–30 °C).</p>

7. ASSAY PROCEDURE

Principle of the test

This ELISA serves for the determination of antibodies against TNF α blocker adalimumab (e.g. HUMIRA®). During sample preparation, the anti-drug antibodies (ADA) are separated from the therapeutic antibody in order to acquire free ADA. By adding the peroxidase conjugate (POD-therapeutic antibody) and the tracer (biotinylated therapeutic antibody), the unmarked therapeutic antibodies are replaced and the marked antibodies can form a complex with the ADA. This complex binds via biotin to the streptavidin coated microtiter plate. It is detected via the peroxidase conjugate with the peroxidase converting the substrate TMB to a blue product. The enzymatic reaction is stopped by adding an acidic solution. The samples convert from blue to yellow. The colour change should be measured in a photometer at 450 nm. The interpretation is made using the cutoff control.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips in the aluminium bag at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash the coated microtiter plate before use 5 x with 250 μl wash buffer . After the final washing step, tap the inverted microtiter plate firmly on absorbent paper.
2.	Pipet 100 μl of preincubated controls/samples into the wells of the microtiter plate.
3.	Cover plate and incubate for 1.5 hours shaking on a horizontal mixer at room temperature (15–30 °C).
4.	Decant the content of the plate and wash each well 5 x with 250 μl wash buffer . After the final washing step, tap the inverted microtiter plate firmly on absorbent paper.
5.	Add 100 μl of TMB substrate solution into each well.

6.	Incubate for 10–20 minutes at room temperature (15–30 °C) until a sufficiently large difference in colour occurs.*
7.	Add 100 µl stop solution into each well and mix shortly by using the shake function of the microplate reader .
8.	Determine the absorption immediately with the microplate reader at 450 nm against 620 nm (or 690 nm) as reference wavelength.

* The intensity of the colour change is temperature sensitive. We recommend to observe the process of the colour change and to stop the reaction upon good differentiation.

8. RESULTS

Cutoff = 10 AU/ml = OD cutoff control

Samples which have a higher average optical density (OD) than the cutoff control are positive.

Immundiagnostik recommends linear regression using a linear ordinate and abscissa to calculate the results.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a manual control of the paired values should be made.

Sample calculation for a positive sample

Average OD of patient's sample	0.735	
Average OD of cutoff control	0.085 = 10 AU/ml	
Concentration patient's sample	$\frac{0.735 \times 10 \text{ AU/ml}}{0.085}$	= 86.47 AU/ml

9. LIMITATIONS

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

Samples with concentrations lower than the measurement range cannot be clearly quantified.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 12)

Probe	AU/ml	VK [%]
1	17.45	7.17
2	48.09	4.83
3	213.44	5.53

Inter-Assay (n = 23)

Probe	AU/ml	VK [%]
1	59.91	3.47
2	290.92	5.92

Analytical Sensitivity

The LoB (limit of blank) was evaluated according to the guideline CLSI EP17-A2 and resulted in 2.765 AU/ml.

Specificity

No cross reactivity to ADA to infliximab was observed.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK*® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure,

which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Afif W, Loftus EV, Jr., Faubion WA, Kane SV, Bruining DH, Hanson KA, Sandborn WJ (2010). Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *Am J Gastroenterol* **105**(5): 1133-1139.
2. Kopylov U, Mazor Y, Yavzori M, Fudim E, Katz L, Coscas D, Picard O, Chowers Y, Eliakim R, Ben-Horin S (2012). Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-infliximab antibodies. *Inflamm Bowel Dis* **18**(9): 1628-1633.
3. Tak PP (2012). A personalized medicine approach to biological treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology* **51**(4): 600-609.
4. Ordas I, Mould DR, Feagan BG, Sandborn WJ (2012) Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clin Pharmacol Ther* **91**(4): 635-646.
5. Bender NK, Heilig CE, Droll B, Wohlgemuth J, Armbruster FP, Heilig B (2007). Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatol Int* **27**(3): 269-274.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by