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15. LITERATUR

1. Ordas I, Mould DR, Feagan BG, Sandborn WJ (2012) Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clin Pharmacol Ther* **91**(4): 635-646.
2. Xu, Z. et al. Population pharmacokinetics of golimumab, an anti-tumor necrosis factor-alpha human monoclonal antibody, in patients with psoriatic arthritis. *J. Clin. Pharmacol.* **49**, 1056-70 (2009).
3. Vincent, F. B. et al. Antidrug antibodies (ADAb) to tumour necrosis factor (TNF)-specific neutralising agents in chronic inflammatory diseases: a real issue, a clinical perspective. *Ann. Rheum. Dis.* **72**, 165-78 (2013).

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
IDKmonitor[®] Golimumab drug level ELISA

***For the in vitro determination of free golimumab
concentration (e. g. SIMPONI[®]) in EDTA plasma and serum***

Valid from 2015-08-21

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	15
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. STORAGE AND PREPARATION OF REAGENTS	16
6. STORAGE AND PREPARATION OF SAMPLES	17
7. ASSAY PROCEDURE	17
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
8. RESULTS	19
9. LIMITATIONS	20
10. QUALITY CONTROL	20
<i>Reference range</i>	20
11. PERFORMANCE CHARACTERISTICS	20
<i>Precision and reproducibility</i>	20
<i>Analytical Sensitivity</i>	21
<i>Spiking Recovery</i>	21
<i>Dilution recovery</i>	21
12. PRECAUTIONS	22
13. TECHNICAL HINTS	22
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
15. REFERENCES	23

1. INTENDED USE

This Immundiagnostik assay is an enzyme immunoassay intended for the quantitative determination of free human therapeutic TNF α antibody golimumab (e.g. SIMPONI®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis, rheumatoid arthritis or psoriasis are often treated with anti-TNF α antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF α therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF α antibody. Several factors influence the trough level, among them dosage and frequency of anti-TNF α blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation anti-drug antibodies, ADA) [1, 2]. It is thought that ADA functionally neutralize the therapeutic antibodies or induce their rapid elimination. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF α antibody application [1, 3].

The Immundiagnostik TNF α blocker monitoring ELISA for the determination of free golimumab (e.g. SIMPONI®) measures reliably the effective drug concentration and is an opportunity for the treating physician to monitor and optimize the therapy early on.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K9656	PLATE	Holder with strips, pre-coated	12 x 8 wells
K9656	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K9656	CONJ	Conjugate, peroxidase-labeled, concentrate, 100x	1 x 200 μ l
K9656	STD	Calibrators, lyophilized (0; 4.15; 8.3; 25; 75; 225 ng/ml)	2 x 6 vials
K9656	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K9656	CTRL	Control, lyophilized (see specification for range)	2 x 1 vial
K9656	SAMPLEBUF	Dilution buffer, ready to use	1 x 100 ml
K9656	SUB	TMB substrate, ready to use	1 x 15 ml

Cat. No.	Label	Kit Components	Quantity
K9656	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Absorbent paper
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month.**
- The **lyophilized standards (STD)** and **controls (CTRL)** are stable at **2–8 °C** until the expiry date stated on the label. Before use, the standards and controls must be reconstituted with **500 µl of ultra pure water.** Allow the vial

content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted standards and controls can be stored at -20°C for one month. Avoid repeated thawing and freezing.**

- The **lyophilized standards** (STD) and **controls** (CTRL) are stable at **2–8°C** until the expiry date stated on the label. Before use, the standards and controls must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted standards and controls can be stored at -20°C for one month. Avoid repeated thawing and freezing.**
- **Preparation of the conjugate:** The **conjugate concentrate (CONJ)** must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The CONJ is stable at **2–8°C** until expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

6. STORAGE AND PREPARATION OF SAMPLES

Storage

Freshly collected EDTA plasma or serum can be stored for one day at room temperature (15–30°C) or for longer storage at -20°C.

EDTA plasma and serum

EDTA plasma or serum samples must be diluted **1:100** before performing the assay, e.g. **10 µl** sample + **990 µl** SAMPLEBUF (dilution buffer), mix well.

For analysis, pipet **100 µl** of each prepared sample per well. For the recommended analysis in duplicate, **2 x 100 µl** are required.

7. ASSAY PROCEDURE

Principle of the test

This enzyme immunoassay is designed to determine the quantitative of free therapeutic TNFα antibodies golimumab (e.g. SIMPONI®) in EDTA plasma or serum samples. In a first incubation step, the free therapeutic TNFα antibodies golimumab (e.g. SIMPONI®) from the sample are bound to the TNFα coated on the plate. To remove all unbound substances, a washing step is carried out. In a further incubation step, peroxidase-labeled antibody is added. Tetramethylbenzidine (TMB) is used as a sub-

strate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of free therapeutic TNF α antibodies golimumab (e. g. SIMPONI®) in the sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standard. The concentrations of free therapeutic TNF α antibodies golimumab (e. g. SIMPONI®) in the samples are determined directly from this curve.

Test procedure

Bring **all reagents and samples to room temperature** (15–30°C) and mix well.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8°C. Strips are stable until expiry date stated on the label.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Add 100 μl of STD (standards), SAMPLE (samples) or CTRL (controls) into the respective wells.
2.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.
3.	Discard the contents of each well and wash 5 times with 250 μl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
4.	Add 100 μl conjugate in each well.
5.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30°C) on a horizontal shaker.
6.	Discard the contents of each well and wash 5 times with 250 μl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
7.	Add 100 μl SUB (substrate) into each well.
8.	Incubate for 10–20 min* at room temperature (15–30°C) in the dark.
9.	Add 100 μl STOP (stop solution) into each well, mix.

10.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.
-----	---

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the „4 parameter algorithm“.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

EDTA-plasma and serum samples

The obtained golimumab levels of EDTA plasma and serum samples have to be multiplied with the dilution factor of **100**.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible. Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay

The precision (intra-assay variation) was calculated from 19 replicate determinations on three samples. (n = 19)

Sample	Golimumab (e.g. SIMPONI®) [µg/ml]	Intra-Assay CV [%]
1	6.2	6.6
2	15.2	6.4
3	7.8	4.1

Inter-Assay

The total precision (inter-assay variation) of was calculated from data on two samples obtained in 11 different assays. (n = 11)

Sample	Golimumab (e. g. SIMPO-NI®) [$\mu\text{g/ml}$]	Inter-Assay CV [%]
1	13.8	5.3
2	6.5	5.0

Analytical Sensitivity

The LoB (limit of blank) was evaluated according to the guideline CLSI EP17-A2 and resulted in 3.56 ng/ml.

Spiking Recovery

Sample	Unspiked sample [$\mu\text{g/ml}$]	Spike [$\mu\text{g/ml}$]	Golimumab expected [$\mu\text{g/ml}$]	Golimumab measured [$\mu\text{g/ml}$]
A	0.33	10.0	10.33	10.85
	0.33	5.0	5.33	5.52
	0.33	2.5	2.83	3.20

Dilution recovery

Sample	Dilution	Golimumab expected [$\mu\text{g/ml}$]	Golimumab measured [$\mu\text{g/ml}$]
A	1:200	14.07	14.07
	1:400	7.04	7.60
	1:800	3.52	4.05
	1:1600	1.76	2.07
	1:3200	0.88	1.01
	1:6400	0.44	0.57

Sample	Dilution	Golimumab expected [$\mu\text{g/ml}$]	Golimumab measured [$\mu\text{g/ml}$]
B	1:200	8.22	8.22
	1:400	4.11	4.06
	1:800	2.06	2.13
	1:1600	1.03	1.09
	1:3200	0.51	0.61
C	1:200	6.44	6.44
	1:400	3.22	3.33
	1:800	1.61	1.59
	1:1600	0.81	0.86
	1:3200	0.40	0.40

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.

- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDKmonitor® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Ordas I, Mould DR, Feagan BG, Sandborn WJ (2012) Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clin Pharmacol Ther* **91**(4): 635-646.
2. Xu, Z. et al. Population pharmacokinetics of golimumab, an anti-tumor necrosis factor-alpha human monoclonal antibody, in patients with psoriatic arthritis. *J. Clin. Pharmacol.* **49**, 1056–70 (2009).
3. Vincent, F. B. et al. Antidrug antibodies (ADAb) to tumour necrosis factor (TNF)-specific neutralising agents in chronic inflammatory diseases: a real issue, a clinical perspective. *Ann. Rheum. Dis.* **72**, 165–78 (2013).

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by