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Immunglobulin E ELISA Kit

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TABLE OF CONTENTS	PAGE
1. INTENDED USE	13
2. SUMMARY AND EXPLANATION OF THE TEST	13
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	14
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	15
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	16
9. ASSAY PROCEDURE	16
PROCEDURAL NOTES	16
TEST PROCEDURE	17
10. RESULTS	19
11. LIMITATIONS	19
12. QUALITY CONTROL	20
EXPECTED VALUES	20
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	20

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Hersteller

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
Immunglobulin E ELISA Kit

For the in vitro determination of IgE in serum and stool

Valid from 02.08.2010

K 6511

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **IgE** in serum and stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Immunoglobulin E (IgE) has been demonstrated in several human secretions and fluids, such as respiratory tract mucus and saliva, urine, tears, human milk as well as in intestinal juices.

The IgE is responsible for various hypersensitivity reactions of the immediate type (Type-I-hypersensitivity). Polyvalent antigens trigger B-cells to produce specific IgE in large amounts. The released IgE antibodies bind to the Fc receptor on the surface of mast cells. The next time the sensitized person has contact to the allergen, it cross links the IgE on the surface of the mast cells and thereby triggers the mast cells to release histamine, which cause the various symptoms.

IgE is involved in a number of allergic phenomena such as erythema, anaphylaxis, homocytotropic histamine release, skin-sensitization and atopic reaction. In patients suffering from allergy, the IgE levels increase as much as three to four times compared to the healthy.

Indication

- Diagnosis of allergy and atopic disease in infants with relapsing, spastic bronchitis, pseudocroup, etc.
- Food allergy and intolerances
- Immunodeficiency syndromes
- Differential diagnosis of eczema, respiratory tract disease, rhinopathy, urticaria and other allergy related reactions

3. PRINCIPLE OF THE TEST

In a first incubation step, the IgE in the samples is bound to polyclonal antibodies (in excess), which are immobilized to the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labeled anti-IgE (POD-Antibody) antibody (polyclonal) is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of IgE in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. IgE, present in the patient samples, is determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6511MTP	PLATE	One holder with precoated strips	12 x 8
K 6511WB	WASHBUF	ELISA wash concentrate 10x	2 x 100 ml
K 6511K	CONJ	Conjugate, (anti-IgE, Peroxidase-labeled)	1 x 150 µl
K 6511ST	STD	Standards, ready to use (0; 0.31; 1.25; 5; 20 kU/l)	5 x 1ml
K 6511KO1	CTRL	Control, ready to use	1 x 1 ml
K 6511KO2	CTRL	Control, ready to use	1 x 1 ml
K 6511PV	SAMPLEBUF	Sample dilution buffer	1 x 50 ml
K 6511TMB	SUB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 6511AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled or deionized water
- Precision pipettors calibrated to deliver 50-100 μ l
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μ l** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (500 ml WASHBUF + 450 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The **STD** (Calibrators) and the **CTRL** (controls) are ready to use and can be stored at **2-8°C** until the expiry date given on the label. Repeated thawing and freezing should be avoided.
- The **CONJ** (conjugate) must be diluted **1: 100** in ELISA wash buffer (100 μ l CONJ + 10 ml ELISA wash buffer). The **undiluted CONJ** (conjugate) is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the expiry date (see label of test package) when stored at 2-8°C.

7. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum

Serum must be diluted **1:50** with sample dilution buffer:

20 µl Serum /Plasma + 980 µl SAMPLEBUF (sample dilution buffer).

Faeces

Transfer about **100 mg** of the sample, accurately weighed (note the exact weight for calculation), into a flask with **5 ml** of **ELISA wash buffer** and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. The resulting supernatant can be stored at -20°C for about 1 month. **100 µl** of the supernatant is used for the assay.

Immundiagnostik recommends the use of the sample tubes from Roche Diagnostics (Cat.-No. 745804) for sample preparation.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well
2. Mark the positions of STD (Standard) SAMPLE (Sample) CTRL (Controls) on a protocol sheet
3. Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label
4. Wash each well 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
5. Add 100 µl of STD (Standard) SAMPLE (Sample) CTRL (Controls) in duplicate into respective well
6. Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer

7. Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
8. Add 100 µl CONJ (conjugate) into each well
9. Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer
10. Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted wash buffer into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper
11. Add 100 µl of SUB (substrate) into each well
12. Incubate for 5 - 10 minutes at room temperature (18-26°C) in the dark*
13. Add 50 µl of STOP (stop solution) into each well, mix thoroughly
14. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

*The intensity of the colour change is temperature sensitive. We recommend to observe the colour change and to stop the reaction upon good differentiation.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Serum

For the calculation of the serum values the results must be multiplied by **50**.

Faeces

For the IgE concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml stool = 1g) = 0.08 ml

dilution step 1: 5ml / 0,08ml = 62.5

Multiply the result by **62.5** to get the real concentration. **The dilution factor depends on the weight of the faeces.**

11. LIMITATIONS

Samples with IgE levels greater than the highest calibrator should be further diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal range

IgE (serum): 40 - 120 kU/l

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The assay should always be performed according the enclosed manual.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number