# Lysozym ELISA Kit

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## 15. LITERATUR

1. Langhorst, Jost, Sigrid Elsenbruch, Twyla Mueller, Andreas Rueffer, Guenther Spahn, Andreas Michalsen, and Gustav J Dobos. 2005. "Comparison of 4 Neutrophil-Derived Proteins in Feces as Indicators of Disease Activity in Ulcerative Colitis." *Inflammatory Bowel Diseases* **11** (12) (December): 1085–91.
2. Pettersson, C, H Karlsson, M Ståhlman, T Larsson, B Fagerberg, M Lindahl, O Wiklund, J Borén, and L Fogelstrand. 2011. "LDL-Associated Apolipoprotein J and Lysozyme Are Associated with Atherogenic Properties of LDL Found in Type 2 Diabetes and the Metabolic Syndrome." *Journal of Internal Medicine* **269** (3) (March): 306–21. doi:10.1111/j.1365-2796.2010.02292.x.
3. Schwaab M., Euteneuer M., Lautermann J., Sudhoff H. 2005. "Lysozym Und Laktoferrin in Adenoiden, Hyperplastischen Und Chronisch Entzündeten Tonsillen - Eine Quantitative Analyse." *Laryngo-Rhino-Otologie* **84**: 660–664. doi:10.1055/s.
4. Volman, Julia J, Ronald P Mensink, Wim a Buurman, and Jogchum Plat. 2011. "In Vivo Effects of Dietary (1→3), (1→4)-B-D-Glucans from Oat on Mucosal Immune Responses in Man and Mice." *Scandinavian Journal of Gastroenterology* **46** (5) (May): 603–10. doi:10.3109/00365521.2010.545830.
5. West, N P, D B Pyne, J M Kyd, G M Renshaw, P A Fricker, and A W Cripps. 2010. "The Effect of Exercise on Innate Mucosal Immunity." *British Journal of Sports Medicine* **44** (4) (March): 227–31. doi:10.1136/bjism.2008.046532.

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

*In-Vitro*-Diagnostikum

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Hersteller

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
**Manual**

# Lysozyme ELISA Kit

*For the in vitro determination of lysozyme in serum, urine  
and liquor*

Valid from 2015-03-03

**REF** K 6903

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>16</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>17</b>
<b>7. ASSAY PROCEDURE</b>	<b>17</b>
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
<b>8. RESULTS</b>	<b>19</b>
<b>9. LIMITATIONS</b>	<b>20</b>
<b>10. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>21</b>
<i>Analytical Sensitivity</i>	21
<i>Precision and reproducibility</i>	21
<i>Spiking Recovery</i>	21
<i>Dilution recovery</i>	22
<i>Specificity</i>	22
<b>12. PRECAUTIONS</b>	<b>22</b>
<b>13. TECHNICAL HINTS</b>	<b>22</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>15. REFERENCES</b>	<b>23</b>

## 1. INTENDED USE

This ELISA is intended for the quantitative determination of Lysozyme in human serum, urine and liquor. It is for in vitro diagnostic use only.

## 2. INTRODUCTION

Lysozyme (muramidase) is a protein with a molecular weight of approx. 15 kDa and belongs to the group of alkaline glycosidases. Lysozyme is produced by granulocytes, monocytes and macrophages. The main source for faecal Lysozyme are the intestinal granulocytes. Lysozyme can be detected in all cells of the inflammatory infiltrate during an acute attack of Crohn's disease. To some extent, Lysozyme is also secreted actively by mononuclear cells into the bowel lumen.

### Indications

- Diagnosis and monitoring of Crohn's Disease
- Early diagnosis of rejection reactions in kidney transplantation cases
- Differential diagnosis and monitoring of leukosis
- Diagnosis and treatment monitoring of urinary tract infections in children
- Differential diagnosis between viral and bacterial meningitis in children
- Identification of sepsis in neonates

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6903	PLATE	one holder with precoated strips	12 x 8 wells
K 6903	WASHBUF	ELISA wash concentrate 10 x	2 x 100 ml
K 6903	CONJBUF	Conjugate dilution buffer	1 x 15 ml
K 6903	CONJ	Conjugate concentrate, (rabbit-anti-lysozyme, Peroxidase-labeled)	1 x 50 µl
K 6903	CAL	Calibrator, ready to use	1 x 1 ml
K 6903	CTRL	Control, ready to use	1 x 1 ml
K 6903	CTRL	Control, ready to use	1 x 1 ml
K 6903	SUB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 6903	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month**.
- Use **100 µl** of **diluted WASHBUF** (wash buffer) as **BLANK**. Pipet into the respective well.
- The **conjugate concentrate** (CONJ) must be diluted **1:1000 in conjugate dilution buffer** (100 µl CONJ + 10 ml CONJBUF). The concentrate is stable at **2–8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### Serum

Serum should be centrifuged within one hour after collection. **Store samples at -20 °C** if not assayed on the same day.

Lysozyme should be determined within **8 hours**, after this time the concentration decreases at **room temperature**. The sample can be stored **6 days at 4 °C**.\*

Lipemic or hemolytic samples may give false results. Samples should be mixed well before assaying. Samples are diluted between **1:500** and **1:1000** with wash buffer.

Use this dilution factor to calculate the Lysozyme concentration.

\*Thomas, L., 1998, Labor und Diagnose, 5th ed. TH-Books, Frankfurt.

### Urine

We recommend a dilution of **1:5** with wash buffer for the urine samples before analysis.

### Liquor

We recommend a dilution of **1:50** in wash buffer for the liquor samples before analysis.

## 7. ASSAY PROCEDURE

### *Principle of the test*

The assay utilizes the sandwich technique with two selected antibodies that recognize human lysozyme.

Calibrator, controls and diluted samples, which are assayed for human lysozyme, are added into the wells of a micro plate coated with a high affine anti-human lysozyme antibody. During the first incubation step, lysozyme is bound by the immobilized antibody. Then a peroxidase-conjugated anti-human lysozyme antibody is added into each microtiter well and a "sandwich" of capture antibody - human lysozyme - peroxidase-conjugate is formed. Tetramethylbenzidine is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the enzymatic reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of lysozyme. Samples are quantified by referring their optical density to a lot-dependent master calibration curve and the use of the calibrator that is run with each test.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Mark the positions of <b>STD/CTRL/SAMPLE</b> (standard/control/sample) on a protocol sheet
2.	Wash the wells <b>5x with 250 µl of diluted WASHBUF</b> (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash
3.	Add <b>100 µl of STD/CTRL/SAMPLE</b> (standard/control/sample) into respective well
4.	Cover the plate tightly and incubate for <b>1 hour</b> at room temperature ( <b>15-30 °C</b> ) on a horizontal mixer
5.	Aspirate and wash the wells <b>5x with 250 µl of diluted WASHBUF</b> (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash
6.	Add <b>100 µl diluted CONJ</b> (Conjugate) into each well
7.	Cover the plate tightly and incubate for <b>1 hour</b> at room temperature ( <b>15-30 °C</b> ) on a horizontal mixer
8.	Aspirate and wash the wells <b>5x with 250 µl of diluted WASHBUF</b> (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash
9.	Add <b>100 µl of SUB</b> (Substrate) into each well
10.	Incubate for <b>10 - 20 minutes*</b> at room temperature ( <b>15-30 °C</b> ) <b>in the dark</b>
11.	Add <b>100 µl of STOP</b> (Stop solution) into each well, mix well

12.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of a sample exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.
-----	---

\* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

## 8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A,B,C and D. In both cases the optical density of the calibrator (CAL) is essential. Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

**Caution:** Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### Serum

The result should be multiplied by the corresponding dilution factor to obtain the serum value.

### Urine

The estimated values are multiplied by a dilution factor of 5.

### Liquor

The estimated values are multiplied by a dilution factor of 50.

## 9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

The lower limit of the measurement range can be calculated as:

*detection limit × sample dilution factor to be used*

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

#### **Serum\***

700 - 2580 ng/ml

The lysozyme concentration depends on its production by monocytes, macrophages, granulocytes as well as kidney parenchyma cells.

Lysozyme is elevated in: myelomonocytic leukosis, sarcoidosis

Lysozym is reduced in: newborn sepsis, panmyelopathy

#### **Urine\***

1,7 - 123 ng/ml

Lysozyme in urine is elevated in: myelomonocytic leukemia, urinary passage infection of children

Lysozyme in urine is reduced in: panmyelopathy

#### **Liquor\***

< 62 ng/ml

\* Labor Dr. Limbach, Heidelberg, <http://www.labor-limbach.de>

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Analytical Sensitivity*

The Zero-standard was measured 21 times. The detection limit was set as  $B_0 + 1.645 \cdot SD$  and estimated to be 0.144 ng/ml.

### *Precision and reproducibility*

#### **Intra-Assay (n = 26)**

Sample	Lysozyme [ng/ml]	CV [%]
A	658.4	2.2
B	397.9	2.8

#### **Inter-Assay (n = 18)**

Sample	Lysozyme [ng/ml]	CV [%]
A	696,5	4,6
B	391,4	6,4

### *Spiking Recovery*

Two samples were spiked with different lysozyme concentrations and measured using this assay (n = 2).

Sample	Unspiked Sample [ng/ml]	Spike [ng/ml]	Lysozyme expected [ng/ml]	Lysozyme measured [ng/ml]
A	0.46	1.8	2.26	2.17
		2.7	3.16	3.21
		5.0	5.46	5.51
B	0.74	1.8	2.54	2.48
		2.7	3.44	3.50
		3.7	4.44	4.65

### Dilution recovery

A patient sample was diluted and analyzed. The results are shown below (n = 1):

Sample	Dilution	Lysozyme expected [ng/ml]	Lysozyme measured [ng/ml]
A	1:500	626.3	626.3
	1:1000	313.2	324.2
	1:2000	156.6	173.0

### Specificity

No cross reactivity to other plasma proteins was observed.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.

- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Langhorst, Jost, Sigrid Elsenbruch, Twyla Mueller, Andreas Rueffer, Guenther Spahn, Andreas Michalsen, and Gustav J Dobos. 2005. "Comparison of 4 Neutrophil-Derived Proteins in Feces as Indicators of Disease Activity in Ulcerative Colitis." *Inflammatory Bowel Diseases* **11** (12) (December): 1085–91.
2. Pettersson, C, H Karlsson, M Ståhlman, T Larsson, B Fagerberg, M Lindahl, O Wiklund, J Borén, and L Fogelstrand. 2011. "LDL-Associated Apolipoprotein J and Lysozyme Are Associated with Atherogenic Properties of LDL Found in Type 2 Diabetes and the Metabolic Syndrome." *Journal of Internal Medicine* **269** (3) (March): 306–21. doi:10.1111/j.1365-2796.2010.02292.x.
3. Schwaab M., Euteneuer M., Lautermann J., Sudhoff H. 2005. "Lysozym Und Laktoferrin in Adenoiden, Hyperplastischen Und Chronisch Entzündeten Tonsillen - Eine Quantitative Analyse." *Laryngo-Rhino-Otologie* **84**: 660–664. doi:10.1055/s.
4. Volman, Julia J, Ronald P Mensink, Wim a Buurman, and Jogchum Plat. 2011. "In Vivo Effects of Dietary (1→3), (1→4)-B-D-Glucans from Oat on Mucosal Immu-

- ne Responses in Man and Mice." *Scandinavian Journal of Gastroenterology* **46** (5) (May): 603–10. doi:10.3109/00365521.2010.545830.
5. West, N P, D B Pyne, J M Kyd, G M Renshaw, P A Fricker, and A W Cripps. 2010. "The Effect of Exercise on Innate Mucosal Immunity." *British Journal of Sports Medicine* **44** (4) (March): 227–31. doi:10.1136/bjism.2008.046532.

**Used symbols:**

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by