Laktoferrin ELISA Kit

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15. LITERATUR

1. Levay, P. F. & Viljoen, M. Lactoferrin: a general review. *Haematologica* 80, 252–67 (1995).
2. Gisbert, J. P., McNicholl, A. G. & Gomollon, F. Questions and answers on the role of fecal lactoferrin as a biological marker in inflammatory bowel disease. *Inflammatory bowel diseases* 15, 1746–54 (2009).

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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Lactoferrin ELISA kit

For the in vitro determination of lactoferrin in stool

Valid from 2015-01-07

REF **K 6870**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	18
6. SPECIMEN COLLECTION AND PREPARATION	19
<i>Sample stability and storage</i>	19
<i>Extraction of the stool samples</i>	19
<i>Dilution of samples</i>	21
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Test procedure</i>	21
8. RESULTS	22
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Reference range</i>	23
11. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Spiking recovery</i>	24
<i>Dilution recovery</i>	25
<i>Analytical Sensitivity</i>	25
<i>Specificity</i>	25
12. PRECAUTIONS	25
13. TECHNICAL HINTS	26
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	27

1. INTENDED USE

This Immundiagnostik assay is an enzyme immunoassay intended for the determination of lactoferrin in stool. For *in vitro* diagnostic use only.

2. INTRODUCTION

Lactoferrin is a 76 kDa iron-binding glycoprotein which is synthesized and stored in the secondary granules of neutrophils. It is also present in several secretory fluids, such as milk, saliva, tears, and nasal secretions.

Lactoferrin can exist in different polymeric forms ranging from monomers to tetramers; it tends to polymerize especially at high concentrations.

The physiological activities of lactoferrin include regulation of iron homeostasis, innate defense against a broad range of microbial infections, anti-inflammatory activity, regulation of cellular growth and differentiation and protection against cancer development and metastasis.

During intestinal inflammation, polymorphonuclear neutrophils infiltrate the mucosa and release lactoferrin by degranulation, which results in an increased excretion of lactoferrin into the feces. Fecal lactoferrin is therefore a marker for neutrophilic intestinal inflammation.

Indications

- Detection of intestinal inflammatory activity
- Monitoring of disease activity in inflammatory bowel disease (IBD)
- Prediction of relapse of IBD
- Assessment of IBD treatment response

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6870	PLATE	Mikrotiterplatte, vorbeschichtet	12 x 8 wells
K 6870	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 6870	IDK Extract®	Extraction buffer concentrate <i>IDK Extract®</i> , 2,5x	1 x 100 ml
K 6870	CONJ	Conjugate (rabbit anti human lactoferrin)	1 x 200 µl

Cat. No.	Label	Kit components	Quantity
K 6870	STD	Standards, lyophilised	4 x 5 vials
K 6870	CTRL	Controls	4 x 2 vials
K 6870	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 6870	STOP	ELISA stop solution, ready to use	1 x 15 ml
K 6870	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 50 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10–1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Centrifuge, 3000 g
- Vortex mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.

- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8°C for one month**.
- The **extraction buffer concentrate** *IDK Extract*® must be diluted with ultra pure water **1:2.5** before use (100 ml EXBUF + 150 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37°C in a water bath. The **extraction buffer concentrate** *IDK Extract*® is stable at **2–8°C** until the expiry date stated on the label. Diluted **buffer solution** (extraction buffer) can be stored in a closed flask at **2–8°C for three months**.
- The **lyophilized standards** (STD) and **controls** (CTRL) are stable at **2–8°C** until the expiry date stated on the label. **Reconstitution details** are given in the **data sheet**. **Reconstituted standards and controls are not stable and cannot be stored**.
- The **conjugate** (CONJ) must be diluted **1:100 in wash buffer** (100 µl CONJ + 9,9 ml wash buffer). The undiluted conjugate is stable at **2–8°C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

6. SPECIMEN COLLECTION AND PREPARATION

Sample stability and storage

Stool extract is stable at room temperature (15–30°C), at 2–8°C as well as at –20°C for nine days. Avoid more than three freeze-thaw cycles.

Extraction of the stool samples

Diluted extraction buffer *IDK Extract*® is used as a sample extraction buffer. We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

SAS with 1.5 ml buffer:

Applied amount of stool:	15 mg
Buffer Volume:	1.5 ml
Dilution factor:	1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty sample tube** with **1.5 ml** of ready-to-use *IDK Extract*[®] extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert the orange dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~ 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

Dilution I: 1:100

Dilution of samples

The supernatant of the sample preparation procedure (dilution I) is diluted **1:10 in sample dilution buffer**. For example:

- **50 µl** supernatant (dilution I) + **450 µl** sample dilution buffer, mix well = **1:10 (dilution II)**. This results in a final dilution of 1:1000.

For analysis, pipet **100 µl** of **dilution II** per well.

7. ASSAY PROCEDURE

Principle of the test

This ELISA is intended for the quantitative determination of lactoferrin in stool. In a first incubation step, the lactoferrin in the samples is bound to polyclonal antibodies, immobilized to the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a peroxidase-labeled conjugate (rabbit anti human lactoferrin) is added which recognizes specifically the bound lactoferrin. After another washing step to remove all unbound substances, the solid phase is incubated with the substrate, tetramethyl-benzidine (TMB), which reacts with the peroxidase. An acidic stop solution is added to stop the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of lactoferrin. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Lactoferrin present in the patient samples is determined directly from this curve.

Test procedure

Before use, let **all reagents and samples** reach **room temperature** (15–30 °C), mix well.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to pipet the standards and controls in duplicate.

1.	Mark positions for STD/SAMPLE/CTRL (standard / sample / controls) in the protocol sheet.
2.	Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until the expiry date stated on the label.

3.	Add 100 µl of STD/CTRL/SAMPLE into the respective wells.
4.	Cover the strips and incubate for 30 min at room temperature (15–30°C).
5.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
6.	Add 100 µl conjugate in each well.
7.	Cover the strips and incubate for 30 min at room temperature (15–30°C).
8.	Aspirate the contents of each well. Wash each well 5 times with 250 µl of diluted wash buffer . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
9.	Add 100 µl SUB (TMB substrate) in each well.
10.	Incubate for 10–20 minutes at room temperature (15–30°C) in the dark.
11.	Add 100 µl STOP (ELISA stop solution) and mix well.
12.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the „4 parameter algorithm“.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Stool samples

Multiply the result by the dilution factor used to obtain the concentration of lactoferrin in the sample.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

We recommend each laboratory to establish its own reference concentration range.

Reference values in stool

1 g stool is equivalent to 1 ml.

Normal value: <7,2 µg/ml

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)

Sample	Lactoferrin [µg/ml]	CV [%]
1	6.6	4.4
2	28.7	3.4

Inter-Assay (n = 8)

Sample	Lactoferrin [µg/ml]	CV [%]
1	3.7	13.0
2	16.6	7.0

Spiking recovery

Two samples were spiked with different lactoferrin concentrations and measured using this assay.

Sample	Unspiked Sample [ng/ml]	Spike [ng/ml]	Lactoferrin expected [ng/ml]	Lactoferrin measured [ng/ml]
A	1.27	5.68	6.95	7.43
	1.27	28.41	29.68	29.70
	1.27	56.82	58.09	54.10
	1.27	113.64	114.90	98.37
B	17	5.68	22.18	20.88
	17	28.41	44.91	43.13
	17	56.82	73.36	63.25
	17	113.64	130.18	102.07

Dilution recovery

Two samples were diluted and analyzed. The results are shown below (n = 2):

Sample	Dilution	Lactoferrin expected [µg/ml]	Lactoferrin measured [µg/ml]
A	1:800	29.7	29.7
	1:1600	14.9	16.7
	1:3200	7.4	9.3
	1:6400	3.7	5.1
	1:12800	1.9	2.3
B	1:1600	20.2	20.2
	1:3200	10.1	10.5
	1:6400	5.1	5.2
	1:12800	2.6	2.6

Analytical Sensitivity

The LoB (limit of blank) was evaluated according to the guideline CLSI EP17-A2 and resulted in 2.7 ng/ml.

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to lactoferrin. The specificity is calculated in per cent, based on the cross-reactivity of these compounds with lactoferrin:

- Calprotectin < 1 %
- PMN elastase < 1 %
- Myeloperoxidase (MPO) < 3 %

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be

handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not to assemble wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions as sealed ones.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Levay, P. F. & Viljoen, M. Lactoferrin: a general review. *Haematologica* 80, 252–67 (1995).
2. Gisbert, J. P., McNicholl, A. G. & Gomollon, F. Questions and answers on the role of fecal lactoferrin as a biological marker in inflammatory bowel disease. *Inflammatory bowel diseases* 15, 1746–54 (2009).

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by