Nitrotyrosin ELISA

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Hersteller

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Nitrotyrosine ELISA


*For the in vitro determination of nitrotyrosine in human
EDTA-plasma, serum and stool*

EU: IVD / CE

US: Research Use Only. Not for use in diagnostic procedures.

Valid from 15.05.2014

REF K 7824

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	16
2. INTRODUCTION / CLINICAL RELEVANCE	16
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	17
5. STORAGE AND PREPARATION OF REAGENTS	18
6. STORAGE AND PREPARATION OF SAMPLES	18
<i>Extraction of the stool samples</i>	19
<i>Dilution of samples</i>	20
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Test procedure</i>	20
8. RESULTS	21
9. LIMITATIONS	22
10. QUALITY CONTROL	23
<i>Reference range</i>	23
11. PERFORMANCE CHARACTERISTICS	23
<i>Precision and reproducibility</i>	23
<i>Spiking Recovery</i>	24
<i>Analytical Sensitivity</i>	24
12. PRECAUTIONS	25
13. TECHNICAL HINTS	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	25
15. REFERENCES	26

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of protein-bound nitrotyrosine in human EDTA-plasma, serum and stool. It is for *in vitro* use only.

2. INTRODUCTION / CLINICAL RELEVANCE

Nitrotyrosine is the nitrated form of the amino acid tyrosine. The accumulation of protein bound nitrotyrosine is associated with cardiovascular diseases that are based on inflammatory processes (e.g., atherosclerosis, myocardial infarction, diabetic vasculopathy, hypertension, or coronary heart diseases). A growing number of studies have also associated the accumulation of nitrotyrosine with neurological diseases (Alzheimer's disease, Parkinson's disease, multiple sclerosis, stroke). With treatment of some of the associated diseases the levels of nitrated tyrosines have been shown to decrease, so nitrotyrosine has been stated to be a marker of nitrosative stress.

During inflammatory processes, large amounts of nitric oxide ($\bullet\text{NO}$) are locally released from L-arginine. This reaction is catalyzed by the enzyme NO-synthase (NOS). Other causes for the increased $\bullet\text{NO}$ production are exposure to chemicals or heavy metals, drugs, nicotine, or physical and psychological stress, as well as extraordinary physical strain with increased oxygen consumption.

In high concentrations, $\bullet\text{NO}$ that is not trapped by mitochondrial superoxide dismutase (MnSOD) reacts with superoxide ($\bullet\text{OO}^-$) to form peroxynitrite (ONOO^-). Peroxynitrite is implicated as a key oxidant species in several pathologies and is known to be cytotoxic (nitrosative stress).

Peroxynitrite is highly reactive and shows a high affinity to aromatic amino acids, e.g., to the phenolic ring of tyrosine. The nitration of tyrosine in general is a natural process within the post-translational protein modification.

Nitrotyrosine is a stable product and might be seen as a correlate of peroxynitrite production, and its accumulation in cells and tissues is a marker of oxidative stress and nitrosative stress, respectively.

Indications

- Cardiovascular diseases
- Neurological diseases
- Thyroid disturbances
- Blockade of biochemical pathways
- Mitochondriopathy

Consequences of nitrosative stress

- Modification of lipids and proteins (for example structural proteins in mitochondria)
- Inhibition of respiratory chain enzymes in the mitochondria
- Glutamate overload
- Disturbances in ion channels
- Calcium overload
- Initiation of apoptosis processes

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 7824MTP	PLATE	Microtiter plate, precoated	12 x 8 wells
K 7824WP	WASHBUF	ELISA wash concentrate, 10x	1 x 100 ml
K 7824AP	ASYBUF	Assay buffer, ready to use	1 x 25 ml
K 7824ST	STD	Standards, lyophilized	4 x 5 vials
K 7824KO1	CTRL	Control, lyophilized	4 vials
K 7824KO2	CTRL	Control, lyophilized	4 vials
K 7824K	CONJ	Conjugate (goat anti-nitrotyrosine, peroxidase-labeled)	1 x 200 µl
K 7824TMB	SUB	TMB substrate (Tetramethylbenzidine),	1 x 15 ml
K 7824AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8°C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8°C for one month**.
- The **lyophilized standards** (STD) and **controls** (CTRL) are stable at **2–8°C** until the expiry date stated on the label. Before use, the standards and controls must be reconstituted with **500 µl of ultra pure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted standards and controls are not stable and cannot be stored.**
- The **conjugate concentrate** (CONJ) must be diluted **1:101 in wash buffer** (100 µl CONJ + 10 ml wash buffer). The concentrate is stable at **2–8°C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- Store **unused microtiter strips** sealed in the aluminium foil bag with desiccating agent at **2–8°C**. Strips are stable until the expiry date stated on the label.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

6. STORAGE AND PREPARATION OF SAMPLES

EDTA-plasma or serum samples

Pipet **50 µl** of fresh EDTA-plasma or serum sample in a 1.5 ml reaction vial, add **200 µl ASYBUF** (assay buffer) and mix well (dilution **1:5**).

Extraction of the stool samples

Diluted wash buffer is used as a sample extraction buffer. We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instructions for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

SAS with 1.5 ml extraction buffer:

Applied amount of stool: 15 mg

Buffer Volume: 0.75 ml

Dilution Factor: 1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty sample tube** with **0.75 ml** of ready-to-use **wash buffer** before using it with the sample. Important: Allow the wash buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert the orange dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~ 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

The sample suspension is not stable.

Dilution I: 1:50

Dilution of samples

Stool samples

The supernatant of the extraction (dilution step I) is diluted **1:20** with **diluted WASH-BUF** (wash buffer). For example:

30 µl supernatant (dilution I) + **570 µl** diluted wash buffer (dilution II)

Final dilution **1 : 1000**

For analysis, **100 µl of dilution II** is pipetted per well.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the “sandwich” technique with two polyclonal antibodies against nitrated proteins.

Standards, controls and diluted samples which are assayed for nitrotyrosine are added into the wells of a microtiter plate coated with polyclonal anti-nitrotyrosine antibody. During the first incubation step, nitrated proteins are bound by the immobilized primary antibody. Then a peroxidase-conjugated polyclonal anti-nitrotyrosine antibody is added into each microtiter well and a “sandwich” of

primary antibody - nitrated protein – peroxidase-conjugate

is formed. Tetramethylbenzidine is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of nitrotyrosine. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standards.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips at 2–8 °C in the aluminium bag with desiccant. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Wash the pre-coated microtiter plate 5 x with 250 µl ELISA wash buffer . After the final washing step remove residual buffer by tapping the plate on absorbent paper.
2.	Add 100 µl of STD/SAMPLE/CTRL (standard/sample/controls) into respective well.
3.	Cover plate or strips with foil tightly and incubate for 2.5 h at room temperature (15 - 30°C) on the horizontal shaker.
4.	Discard the contents of each well. Wash 5 x with 250 µl ELISA wash buffer . After the final washing step remove residual buffer by tapping the plate on absorbent paper.
5.	Add 100 µl conjugate into each well.
6.	Cover plate or strips with foil tightly and incubate for 1 h at room temperature (15 - 30°C) on the horizontal shaker.
7.	Discard the contents of each well. Wash 5 x with 250 µl ELISA wash buffer . After the final washing step remove residual buffer by tapping the plate on absorbent paper.
8.	Add 100 µl of SUB (TMB substrate) into each well.
9.	Incubate for 10-20* min at room temperature (15-30°C) in the dark .
10.	Add 100 µl of STOP (stop solution) into each well, mix thoroughly in a microtiter plate reader.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

EDTA-plasma and serum samples

The obtained nitrotyrosine concentration must be multiplied by the **dilution factor of 5**.

Stool samples

The obtained nitrotyrosine concentration must be multiplied by the **dilution factor of 1000**.

Controls

The nitrotyrosine concentration can be read directly from the calibration curve. The concentration range is given on the enclosed data sheet specification.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on Immundiagnostik studies of serum samples of apparently healthy persons (n = 78), the following values were estimated:

Min: **48 nM**

Max: **1533 nM**

Median: **207 nM**

For 95% of this collective (95 percentile) a nitrotyrosine concentration of 553 nM and less was obtained.

For 10% of this collective a nitrotyrosine concentration < LoB (Limit of Blank) was obtained.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 26)

Each of two serum samples was measured 26 times by one technician within the same assay.

Sample	Nitrotyrosin [nM]	CV [%]
1	1832	2.7
2	543	9.6

Inter-Assay (n = 20)

Each of two serum sample was independently measured by six, respectively seven different technicians on three different days.

Sample	Nitrotyrosin [nM]	CV [%]
1	1876	4.7
2	596	11.0

Spiking Recovery

Two samples were spiked with different nitrotyrosin concentrations and measured using this assay (n = 2).

Sample	Unspiked Sample [nM]	Spike [nM]	Nitrotyrosin expected [nM]	Nitrotyrosin measured [nM]
A	110	300	410	407
		600	710	731
		900	1010	1098
B	119	600	719	666
		900	1019	975
		1200	1319	1247

Analytical Sensitivity

The calculated detection limit (LoB; Limit of Blank) was set as $B_0 + 1.645 \cdot SD$. Standard 1 (blank) was measured 168 times.

LoB = 20.7 nM

This LoB value was estimated based on the concentration of the calibration curve without considering the sample dilution factor.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Peluffo, Gonzalo, and Rafael Radi. 2007. "Biochemistry of Protein Tyrosine Nitration in Cardiovascular Pathology." *Cardiovascular Research* **75** (2) (July 15): 291–302. doi:10.1016/j.cardiores.2007.04.024.
2. Gonsette, R E. 2008. "Neurodegeneration in Multiple Sclerosis: The Role of Oxidative Stress and Excitotoxicity." *Journal of the Neurological Sciences* **274** (1-2) (November 15): 48–53. doi:10.1016/j.jns.2008.06.029.
3. Ischiropoulos, Harry. 2009. "Protein Tyrosine Nitration--an Update." *Archives of Biochemistry and Biophysics* **484** (2) (April 15): 117–21. doi:10.1016/j.abb.2008.10.034.
4. Köse, Fadime Aydın, Meltem Seziş, Fehmi Akçiçek, and Aysun Pabuççuoğlu. 2011. "Oxidative and Nitrosative Stress Markers in Patients on Hemodialysis and Peritoneal Dialysis." *Blood Purification* **32** (3) (January): 202–8. doi:10.1159/000328030.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		