Osteoprotegerin ELISA Kit

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15. LITERATUR

1. Szulc, P. et al., 2011. Cortical bone status is associated with serum osteoprotegerin concentration in men: the STRAMBO study. *The Journal of clinical endocrinology and metabolism*, **96**(7), pp.2216–26.
2. Samelson, E.J. et al., 2008. Increased plasma osteoprotegerin concentrations are associated with indices of bone strength of the hip. *The Journal of clinical endocrinology and metabolism*, **93**(5), pp.1789–95.
3. Terpos, E. et al., 2009. The clinical significance of serum markers of bone turnover in NSCLC patients: surveillance, management and prognostic implications. *Anti-cancer research*, **29**(5), pp.1651–7.
4. Madarász, E. et al., 2009. Osteoprotegerin levels in women with prior gestational diabetes mellitus. *Diabetes care*, **32**(1), p.e5.
5. Kearns, A.E., Khosla, S. & Kostenuik, P.J., 2008. Receptor activator of nuclear factor kappaB ligand and osteoprotegerin regulation of bone remodeling in health and disease. *Endocrine reviews*, **29**(2), pp.155–92.
6. Pepene, C.E. et al., 2011. Circulating osteoprotegerin and soluble receptor activator of nuclear factor κ B ligand in polycystic ovary syndrome: relationships to insulin resistance and endothelial dysfunction. *European journal of endocrinology / European Federation of Endocrine Societies*, **164**(1), pp.61–8.
7. Van Tuyl, L.H.D. et al., 2010. Baseline RANKL:OPG ratio and markers of bone and cartilage degradation predict annual radiological progression over 11 years in rheumatoid arthritis. *Annals of the rheumatic diseases*, **69**(9), pp.1623–8.
8. Anastasilakis, A.D. et al., 2008. Acute changes in serum osteoprotegerin and receptor activator for nuclear factor-kappaB ligand levels in women with established osteoporosis treated with teriparatide. *European journal of endocrinology / European Federation of Endocrine Societies*, **158**(3), pp.411–5.
9. Faienza, M.F. et al., 2009. Osteoclastogenesis in children with 21-hydroxylase deficiency on long-term glucocorticoid therapy: the role of receptor activator of nuclear factor-kappaB ligand/osteoprotegerin imbalance. *The Journal of clinical*

- endocrinology and metabolism*, **94**(7), pp.2269–76.
10. Park, J.S. et al., 2011. Effect of pioglitazone on serum concentrations of osteoprotegerin in patients with type 2 diabetes mellitus. *European journal of endocrinology / European Federation of Endocrine Societies*, **164**(1), pp.69–74.
 11. Nybo, M. et al., 2011. Rosiglitazone decreases plasma levels of osteoprotegerin in a randomized clinical trial with type 2 diabetes patients. *Basic & clinical pharmacology & toxicology*, **109**(6), pp.481–5.
 12. Semb, A.G. et al., 2009. Osteoprotegerin and soluble receptor activator of nuclear factor-kappaB ligand and risk for coronary events: a nested case-control approach in the prospective EPIC-Norfolk population study 1993-2003. *Arteriosclerosis, thrombosis, and vascular biology*, **29**(6), pp.975–80.
 13. Svensson, M. et al., 2012. Osteoprotegerin as a predictor of renal and cardiovascular outcomes in renal transplant recipients: follow-up data from the ALERT study. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*, **27**(6), pp.2571–5.
 14. Stavroulopoulos, A. et al., 2011. Evolution of coronary artery calcification in patients with chronic kidney disease Stages 3 and 4, with and without diabetes. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*, **26**(8), pp.2582–9.
 15. Lieb, W. et al., 2010. Biomarkers of the osteoprotegerin pathway: clinical correlates, subclinical disease, incident cardiovascular disease, and mortality. *Arteriosclerosis, thrombosis, and vascular biology*, **30**(9), pp.1849–54.
 16. Ueland, T. et al., 2011. Osteoprotegerin predicts progression of chronic heart failure: results from CORONA. *Circulation. Heart failure*, **4**(2), pp.145–52.
 17. Jiang, J.-Q. et al., 2011. Serum osteoprotegerin measurement for early diagnosis of chronic kidney disease-mineral and bone disorder. *Nephrology (Carlton, Vic.)*, **16**(6), pp.588–94.
 18. Kudlacek, S. et al., 2003. Serum levels of osteoprotegerin increase with age in a healthy adult population. *Bone*, **32**(6), pp.681–6.

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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Osteoprotegerin ELISA Kit


For the in vitro determination of OPG in serum and plasma

EU: IVD / CE

US: Research Use Only. Not for use in diagnostic procedures.

Gültig ab / Valid from: 07.05.2014

REF K 1011

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	15
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. PREPARATION AND STORAGE OF REAGENTS	16
6. STORAGE AND PREPARATION OF SAMPLES	17
7. ASSAY PROCEDURE	17
<i>Principle of the test</i>	17
<i>Test procedure</i>	17
8. RESULTS	18
9. LIMITATIONS	19
10. QUALITY CONTROL	19
<i>Reference range</i>	19
<i>Conversion factor</i>	19
<i>OPG-Werte from Kudlacek et al. 2003</i>	20
11. PERFORMANCE CHARACTERISTICS	20
<i>Precision and reproducibility</i>	20
12. PRECAUTIONS	21
13. TECHNICAL HINTS	21
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	22
15. REFERENCES	22

1. INTENDED USE

This Immundiagnostik assay is a sandwich ELISA intended for the quantitative determination of osteoprotegerin in serum, EDTA plasma, heparin plasma or citrate plasma. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

Osteoprotegerin (OPG) or osteoclast inhibitory factor (OCIF) is a glycoprotein of the TNF receptor superfamily 11b (gene name TNFRSF11B) <http://www.uniprot.org/uniprot/O00300>. OPG is synthesized as a monomer of 380 amino acids and is assembled as a homodimer within the cell, and then secreted mainly as a disulfide-linked homodimer into the extracellular compartment. OPG is produced by many different tissues and cell types including osteoblasts. OPG is a negative regulator of bone resorption by acting as decoy receptor for RANKL, thus neutralizing its function in osteoclastogenesis. This glycoprotein is also involved in the regulation of vascular calcification.

Indications

- Osteoporosis ^{1,2}
- Diseases with locally incr. resorption activity ³⁻⁶
- Arthritis ^{7,8}
- Therapy monitoring ⁹⁻¹¹
- Cardiovascular Disease ¹²⁻¹⁷

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 1011MTP	PLATE	Goat polyclonal anti OPG antibody, pre-coated microtiter strips in a strip-holder, packed in an aluminium bag with desiccant	12 x 8 wells
K 1011WP	WASHBUF	Wash buffer concentrate 20x, natural cap	1 x 50 ml
K 1011K	CONJ	Conjugate, (streptavidin-HRP), amber bottle, amber cap, ready to use	1 x 22 ml
K 1011ST	STD	Standards 1-6, (0; 1.25; 2.5; 5; 10; 20 pmol/l), white caps, ready to use	6 x 300 µl
K 1011KO	CTRL	Control, yellow cap, ready to use, (exact concentration on the label)	300 µl

Cat. No.	Label	Kit components	Quantity
K 1011AB	AB	Mouse monoclonal anti OPG antibody, biotin labelled, green cap, yellow dye, ready to use	1 x 7 ml
K 1011AP	ASYBUF	Assay Buffer, red cap, ready to use	1 x 25 ml
K 1011TMB	SUB	Substrate (TMB solution), amber bottle blue cap, ready to use	1 x 22 ml
K 1011AC	STOP	Stop solution, white cap, ready to use	1 x 7 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Plate washer is recommended for washing, alternatively multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩcm).

5. PREPARATION AND STORAGE OF REAGENTS

- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:20 in ultra pure water** before use (50 ml concentrate + 950 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

- Collect venous blood samples by using standardized blood collection tubes for serum or plasma. We recommend performing plasma or serum separation by centrifugation as soon as possible (e.g. 20 min at 2000 *g*, preferably at 4 °C (2–8 °C). If this is not possible, store the samples at 4 °C (2–8 °C) prior to centrifugation (up to one day).
- The acquired plasma or serum samples should be measured as soon as possible. For longer storage aliquot samples and store at -20 °C or lower. Samples are at least stable for 4 freeze-thaw cycles. Lipemic or haemolysed samples may give erroneous results.
- Samples should be mixed well before assaying.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the “sandwich” technique with specific antibodies against OPG. Standards, controls and diluted samples which are assayed for OPG are added into the wells of a micro plate coated with polyclonal goat-anti-OPG-antibody. During the first incubation step, OPG is bound by the immobilized primary antibody. Then a biotinylated monoclonal mouse-anti-OPG-antibody, is added into each microtiter well. In the next step, the streptavidin-peroxidase-conjugate is added and a “sandwich” of 1st antibody – OPG - biotinylated antibody – streptavidin-peroxidase-conjugate is formed. Tetramethylbenzidine is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of OPG. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standard.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips in the aluminium bag with desiccant at 2–8 °C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Pipet 150 µl ASYBUF (assay buffer, red cap) into each well. Add further 100 µl ASYBUF in the well for the BLANK .
	Add 20 µl STD / SAMPLE / CTRL (standard / sample / control) into respective wells.
	Add 50 µl AB (biotinylated anti OPG antibody, natural cap) into each well, except BLANK , mix gently.
	Cover tightly and incubate for 4h at room temperature (18 - 26 °C).
	Aspirate and wash wells 5x with 300 µl diluted wash buffer . After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
	Add 200 µl CONJ (conjugate, amber cap) into each well.
	Cover tightly and incubate for 1 hour at room temperature (18 - 26 °C).
	Aspirate and wash wells 5x with 300 µl diluted Wash buffer. After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel.
	Add 200 µl SUB (substrate, blue cap) into each well.
	Incubate for 30 min* at room temperature (18 - 26 °C) in the dark .
	Add 50 µl STOP (stop solution, white cap) into each well.
12.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithm-

mic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted with STD 1 or OPG negative human serum and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on a laboratory intern study of samples of apparently healthy persons (n = 60), a median of 2,7 pmol/ml was estimated.

We recommend each laboratory to establish its own reference range.

Conversion factor

1 pg/ml = 0.05 pmol/l

OPG-Werte from Kudlacek et al. 2003

OPG reference data⁸ (n = 1134, 50.83 ± 51.47 pg/ml, median 36, 2-584;
n = 1134, 2.54 ± 2.57 pmol/L, median 1.8, 0.1-29.2)

Alter	OPG (mean ± SD)		OPG (median, Bereich)		n
	pg/ml	pmol/L	pg/ml	pmol/L	
Female population (n = 687, 51.9 ± 51.3 pg/ml, median 39; 2.60 ± 2.57 pmol/L, median 1.95)					
< 30	44.5 ± 21.2	2.23 ± 1.06	41.7, 87.7	2.09, 4.39	48
31-40	44.1 ± 1.25	2.21 ± 1.69	36, 260	1.8, 13	198
41-50	41.2 ± 25.0	2.06 ± 33.8	36, 204	1.8, 10.2	217
51-60	39.5 ± 22.3	1.98 ± 1.12	35.5, 160	1.78, 8	150
61-70	69.6 ± 62.9	3.48 ± 3.15	60.5, 342	3.03, 17.1	34
71-80	134.0 ± 70.0	6.70 ± 3.50	131, 244	6.55, 12.2	16
> 81	227.0 ± 100.0	11.35 ± 5.00	206, 479	10.3, 23.95	24
Male Population (n = 447, 50 ± 51.7 pg/ml, median 34.8; 2.5 ± 2.59 pmol/L, median 1.74)					
< 30	40.7 ± 26.3	2.04 ± 1.32	32, 96	1.6, 4.8	19
31-40	41.4 ± 30.5	2.07 ± 1.53	33.5, 170	1.68, 8.5	72
41-50	36.3 ± 31.8	1.82 ± 1.59	36.2, 316	1.81, 15.8	116
51-60	36.9 ± 20.2	1.85 ± 1.01	34, 134	1.7, 6.7	176
61-70	41.0 ± 22.7	2.05 ± 1.14	29.8, 160	1.49, 8	32
71-80	119.0 ± 71.0	5.95 ± 3.55	38, 103	1.9, 5.15	12
> 81	226.0 ± 78.0	11.30 ± 3.90	130, 219	6.5, 10.95	20

11. PERFORMANCE CHARACTERISTICS*Precision and reproducibility***Intra-Assay (n = 5)**

Sample	Mean [pmol/l]	CV [%]	Standard deviation
1	3.2	2	0.05
2	10.1	3	0.34

Inter-Assay (n = 12)

Sample	Mean [pmol/l]	CV [%]	Standard deviation
1	3.2	3	0.10
2	9.9	5	0.50

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Szulc, P. et al., 2011. Cortical bone status is associated with serum osteoprotegerin concentration in men: the STRAMBO study. *The Journal of clinical endocrinology and metabolism*, **96**(7), pp.2216–26.
2. Samelson, E.J. et al., 2008. Increased plasma osteoprotegerin concentrations are associated with indices of bone strength of the hip. *The Journal of clinical endocrinology and metabolism*, **93**(5), pp.1789–95.
3. Terpos, E. et al., 2009. The clinical significance of serum markers of bone turnover in NSCLC patients: surveillance, management and prognostic implications. *Anti-cancer research*, **29**(5), pp.1651–7.
4. Madarász, E. et al., 2009. Osteoprotegerin levels in women with prior gestational diabetes mellitus. *Diabetes care*, **32**(1), p.e5.
5. Kearns, A.E., Khosla, S. & Kostenuik, P.J., 2008. Receptor activator of nuclear factor kappaB ligand and osteoprotegerin regulation of bone remodeling in health and disease. *Endocrine reviews*, **29**(2), pp.155–92.
6. Pepene, C.E. et al., 2011. Circulating osteoprotegerin and soluble receptor activator of nuclear factor kB ligand in polycystic ovary syndrome: relationships to insulin resistance and endothelial dysfunction. *European journal of endocrinology / European Federation of Endocrine Societies*, **164**(1), pp.61–8.
7. Van Tuyl, L.H.D. et al., 2010. Baseline RANKL:OPG ratio and markers of bone and cartilage degradation predict annual radiological progression over 11 years in rheumatoid arthritis. *Annals of the rheumatic diseases*, **69**(9), pp.1623–8.

8. Anastasilakis, A.D. et al., 2008. Acute changes in serum osteoprotegerin and receptor activator for nuclear factor-kappaB ligand levels in women with established osteoporosis treated with teriparatide. *European journal of endocrinology / European Federation of Endocrine Societies*, **158**(3), pp.411–5.
9. Faienza, M.F. et al., 2009. Osteoclastogenesis in children with 21-hydroxylase deficiency on long-term glucocorticoid therapy: the role of receptor activator of nuclear factor-kappaB ligand/osteoprotegerin imbalance. *The Journal of clinical endocrinology and metabolism*, **94**(7), pp.2269–76.
10. Park, J.S. et al., 2011. Effect of pioglitazone on serum concentrations of osteoprotegerin in patients with type 2 diabetes mellitus. *European journal of endocrinology / European Federation of Endocrine Societies*, **164**(1), pp.69–74.
11. Nybo, M. et al., 2011. Rosiglitazone decreases plasma levels of osteoprotegerin in a randomized clinical trial with type 2 diabetes patients. *Basic & clinical pharmacology & toxicology*, **109**(6), pp.481–5.
12. Semb, A.G. et al., 2009. Osteoprotegerin and soluble receptor activator of nuclear factor-kappaB ligand and risk for coronary events: a nested case-control approach in the prospective EPIC-Norfolk population study 1993-2003. *Arteriosclerosis, thrombosis, and vascular biology*, **29**(6), pp.975–80.
13. Svensson, M. et al., 2012. Osteoprotegerin as a predictor of renal and cardiovascular outcomes in renal transplant recipients: follow-up data from the ALERT study. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*, **27**(6), pp.2571–5.
14. Stavroulopoulos, A. et al., 2011. Evolution of coronary artery calcification in patients with chronic kidney disease Stages 3 and 4, with and without diabetes. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*, **26**(8), pp.2582–9.
15. Lieb, W. et al., 2010. Biomarkers of the osteoprotegerin pathway: clinical correlates, subclinical disease, incident cardiovascular disease, and mortality. *Arteriosclerosis, thrombosis, and vascular biology*, **30**(9), pp.1849–54.
16. Ueland, T. et al., 2011. Osteoprotegerin predicts progression of chronic heart failure: results from CORONA. *Circulation. Heart failure*, **4**(2), pp.145–52.
17. Jiang, J.-Q. et al., 2011. Serum osteoprotegerin measurement for early diagnosis of chronic kidney disease-mineral and bone disorder. *Nephrology (Carlton, Vic.)*, **16**(6), pp.588–94.
18. Kudlacek, S. et al., 2003. Serum levels of osteoprotegerin increase with age in a healthy adult population. *Bone*, **32**(6), pp.681–6.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number