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Pankreatische Amylase ELISA Kit

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
Pancreatic amylase ELISA Kit

*For the in vitro determination of human pancreatic amylase in
stool*

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Table of content	Seite/Page
	2
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
TEST PROCEDURE	18
10. RESULTS	20
11. LIMITATIONS	20
12. QUALITY CONTROL	20
EXPECTED VALUES	20
13. PERFORMANCE CHARACTERISTICS	21
SENSITIVITY	21
14. REFERENCES	21
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	21

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Hersteller

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
Pancreatic amylase ELISA Kit

For the in vitro determination of human pancreatic amylase in stool

Valid from 21.10.2010

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **pancreatic Amylase** in stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Like lipase and elastase the **pancreatic amylase** is synthesized in the pancreas. In case of a chronic pancreatitis, **pancreatic amylase** is decreased. Some forms of disease (alcoholism, accidental trauma) are associated with a pancreatic insufficiency. In the routine laboratory **pancreatic amylase** has proved itself to be a genuine alternative to Elastase-I in stool diagnosis.

Indication

- Chronic pancreatitis
- Exocrine pancreas insufficiency

3. PRINCIPLE OF THE TEST

In a first incubation step, the pancreatic amylase in the samples is bound to monoclonal mouse antibodies (in excess), which are immobilized to the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, an anti-Pancreatic Amylase (POD-monoclonal antibody) antibody is added. After another washing step, to remove all unbound substances, the substrate, tetramethylbenzidine (TMB) is added. An acidic stop solution is then added to stop the reaction. The colour converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Pancreatic Amylase in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. Pancreatic amylase, present in the patient samples, is determined directly from this curve.

The combination of two specific antibodies in the pancreatic amylase ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 6410MTP	PLATE	One holder with precoated strips	12 x 8
K 6410WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6410K	CONJ	Conjugate (mouse anti-pancreatic amylase, peroxidase-labeled)	1 x 30 µl
K 6410ST	STD	Calibrators, lyophilized (0; 440; 1750; 7000; 28000 mU/l)	1 x 5 vials
K 6410 KO	CTRL	Control, lyophilized	1 x 1 vial
K 6410TMB	SUB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 6410AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

On request we will send you 3 calibrator sets free of charge. Any further calibrator sets will be charged.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0,2 µm) with an electrical conductivity of 0,055µS/cm at 25°C (≤18,2MΩ cm).

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label. **Prepare just the appropriate amount necessary for the assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **ELISA WASHBUF** (wash buffer concentrate) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be dissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for one month.**

The diluted wash buffer is used as an extraction buffer (see Sample preparation).

- The **lyophilized STD** (standards) and **CTRL** (control) are stable **at 2-8°C** until the expiry date stated on the label. The **STD** (standards) and **CTRL** (control) must be reconstituted with **250 µl** of ultra pure water. Allow the vial content to dissolve for **10 minutes** and mix thoroughly by gentle inversion to insure complete reconstitution.
- The **CONJ** (POD antibody) must be diluted **1:1000** in wash buffer (10 µl CONJ + 10 ml wash buffer). The antibody is stable at 2 -4 °C until expiry date given on the label. **Diluted antibody solution is not stable and could not be stored.**
- All other test reagents are ready to use. The test reagents are stable at **2-8°C** until the expiry date stated on the label.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SAMPLE PREPARATION

Extraction of the stool sample

1a. Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

SAS with 0.75 ml Buffer:

Applied amount of stool:	15 mg
Buffer Volume:	0.75 ml
Dilution Factor:	1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw Stool Sample has to be thawed. For remarkably inhomogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube** with **0.75 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off and leave 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.

- e) Allow sample to stand for app. 10 minutes until sediment has settled down. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure, the sediment will not be dispersed again.

1b. Sample preparation kit from Roche Diagnostics, Mannheim, Germany (Cat. No. 10 745 804 322)

Alternatively, other stool sample preparation kits (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany) can be used. In the Roche sample preparation kit, 100 mg of stool sample are suspended in 5 ml of extraction buffer using a vibrator mixer (e.g. Vortex mixer). Centrifugation of the suspension is recommended.

Dilution I (1a. or 1b.)

1:50

Dilution of samples

Stool samples

The supernatant of the centrifugation (dilution step I) is diluted **1:40 in wash buffer** (WASHBUF). For example:

25 µl supernatant of dilution I + 975 µl WASHBUF = 1:40 (**dilution step II**)

For analysis, pipette **100 µl** of the supernatant of **dilution step II** per well.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any variations of the test procedure, without consulting the producer with the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage.
- The assay should always be performed according the enclosed manual.

Test procedure

Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well.

Wash the PLATE (precoated microtiter plate) 5 x with 250 µl ELISA diluted wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.

Carry out the tests in duplicate.

1. Add **100 µl STD** (standards), **CTRL** (control) and **SAMPLE** (prediluted patient samples) in duplicate into the respective well.
2. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** diluted ELISA wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
4. Add **100 µl of diluted CONJ** (peroxidase-labeled anti-pancreatic amylase antibody).
5. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** diluted ELISA wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
7. Add **100 µl of SUB** (TMB substrate).
8. Incubate for **10-20 minutes** at room temperature.
9. Add **50 µl of STOP** (stop solution) and mix shortly.
10. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Faeces

Multiply the result with **2000** to get the real concentration of Pancreatic Amylase in the corresponding sample.

11. LIMITATIONS

Samples with pancreatic amylase levels greater than the highest calibrator should be diluted and re-assayed.

12. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Pancreatic Amylase (stool):	> 1000 U/l
Grey area	1000 – 1500 U/l

13. PERFORMANCE CHARACTERISTICS

Sensitivity

The detection limit was defined as $B_0 + 2SD$.

n=20

Sample	Pancreatic Amylase Mean value [OD]	Standard variation	Detection limit [mU/l]
1	0.065	0.030	237

14. REFERENCES

1. Bishop et al.: 1996; *Pancreas*, 13 (3), 226
2. Katschinsky et al.: 1997; *Pancreas*, 15 (2), 1991

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.

- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

10/28/2010 Amylase__15112004.DOC

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number