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# Intaktes Proinsulin ELISA

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Langenfeld MR, et al. (2004) IRIS II Study: Sensitivity and specificity of intact proinsulin, adiponectin and the proinsulin/adiponectin ratio as markers for insulin resistance. *Diabetes Technology & Therapeutics*. **6**:836-843

Pfützner A, et al. (2005) IRIS II Study: Intact proinsulin is confirmed as a highly specific indicator for insulin resistance in a large cross-sectional study design. *Diabetes Technology & Therapeutics*. **7**:478-486

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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# Intact Proinsulin ELISA

*For the in vitro determination of intact Proinsulin in serum and plasma*


Gültig ab / Valid from 01.03.2011

K 7821

96

## 1. INTENDED USE

The Intact Proinsulin Assay is an immunometric assay for the quantitative measurement of intact proinsulin in human serum and plasma samples. Measurements of proinsulin are used in the diagnosis and treatment of patients with type 2 diabetes. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Proinsulin is a precursor molecule for insulin and is synthesized by the pancreatic  $\beta$ -cells. Under normal circumstances, virtually all proinsulin is cleaved at residues 32-33 and 65-66 to produce insulin during the formation of secretory granules. Some unmodified proinsulin is released into the circulation, though it is believed to have little or no biological activity. Increased concentrations of circulating proinsulin may occur in insulin-resistant syndromes such as type 2 diabetes and in patients with insulinoma. When used in conjunction with a highly specific insulin assay, it may provide useful information on changes in the processing of insulin in such situations.

### Indications

- Diabetes type I and II
- Insulinomas, functional hypoglycaemia, hyperinsulinemia
- Kidney dysfunction, cirrhosis and hyperthyreosis
- Cardiovascular diseases

### 3. MATERIAL SUPPLIED

Cat Nr.	Label	Kit components	quantity
K 7821MTP	PLATE	Microtiter plate, coated	12 x 8 wells
K 7821WP	WASHBUF	ELISA Wash buffer concentrate 10x	1 x 100 ml
K 7821AK	AB	Biotin conjugated Proinsulin Antibody	1 x 12 ml
K 7821PV	SAMPLEBUF	Sample buffer, ready to use	1 x 12 ml
K 7821K	CONJ	HRP-conjugate, ready to use	1 x 12 ml
K 7821ST	STD	Standards, each 1ml, lyophilized	5 vials
K 7821KO	CTRLA	Control A, 1ml, lyophilized	1 vial
K 7821KO	CTRLB	Control B, 1ml, lyophilized	1 vial
K 7821TMB	SUB	TMB Substrate (Tetramethylbenzidin), ready to use	1 x 15 ml
K 7821AC	STOP	ELISA Stop solution, ready to use	1 x 15 ml

The used standards have been calibrated on the WHO reference material IRR 84/611.

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm and 405nm (reference wave length 620 or 690 nm)

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C using a water bath before dilution. The **buffer concentrate is stable at 2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **STD** (standards) and the **CTRL** (control A or B) must be reconstituted with 1 ml aqua bidest. before use. Allow these to stand for 5 minutes, then mix gently to ensure all solids are dissolved. Stability of the reconstituted Standards is four (4) weeks when stored at 2-8°C. The lyophilized STD and CTRL are stable at 2-8°C until the expiry date stated on the label.
- All other test reagents are ready to use. The test reagents are stable until the expiry date stated on the label when stored at **2-8°C**.

## 6. SAMPLE PREPARATION

### *Dilution of samples*

No sample dilution is necessary.

### **Sample matrices**

Serum, Heparin- or EDTA-Plasma can be used in this assay. Do not use severely hemolyzed specimens.

### **Plasma Collection**

Whole blood should be collected into a tube containing EDTA or heparin as anticoagulant, allowed to stand for maximal 20 min and then centrifuged for at least 15 min at 2000 to 3000 rpm. The clear cell-free supernatant should be immediately collected and used in the assay.

### **Serum Collection**

Serum should be collected into a tube without any anticoagulant, allowed to stand for maximal 20 min and then centrifuged for at least 15 min at 2000 to 3000 rpm. The clear cell-free supernatant should be immediately collected and used in the assay.


### Specimen Storage

After collection, the specimens should be capped and may be stored for up to 24 hours at 2-8°C prior to assaying. Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.


We recommend duplicate analyses for each sample.

### Sample Stability

Generally, freshly collected plasma or serum should be used. Samples can be stored for maximal 24 hours at 2-8°C prior to assaying. For a longer storage, the samples should be frozen at -20°C.

## long term stability at -20°C

## 7. ASSAY PROCEDURE

### *Principle of the test*

The Intact Proinsulin Assay is a two-site immunoassay, employing a specific solid phase antibody immobilized on microtiter wells and a soluble antibody labeled with biotin. The sample is incubated in the microtiter well together with a buffer and, after a wash step, the labeled antibody solution is added. After a second incubation and wash step, HRP labeled streptavidin is added. A third incubation and wash is followed by the addition of substrate solution. Following color development "stop reagent" is added and the color intensity measured in a 96-well ELISA reader.

### Test procedure

1.	Prior to use in the assay, allow all <b>reagents and samples</b> to come to <b>room temperature</b> (18-26 °C) and mix well
2.	Mark the <b>positions of STD</b> (Standard), <b>SAMPLE</b> (Sample) and <b>CTRL</b> (Control) on a protocol sheet
3.	<b>Take microtiter strips</b> out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label
4.	Add <b>50 µl of SAMPLEBUF</b> (sample buffer) into each well
5.	Add <b>50 µl of STD</b> (Standard), <b>SAMPLE</b> (Sample) and <b>CTRL</b> (Control) in duplicate into respective well
6.	Cover the plate tightly and incubate for <b>2 hours at room temperature</b> (18-26°C)
7.	Discard the contents of each well. Remove residual fluids by tapping the plate on absorbent paper. Fill Plate with <b>300 µl of diluted wash buffer</b> . Empty plate and remove residual buffer by tapping the plate on absorbent paper. <b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the plate on absorbent paper after each washing step
8.	Add <b>100 µl AB</b> (biotin conjugated antibody) into each well
9.	Cover the plate tightly and incubate for <b>1 hour at room temperature</b> (18-26°C)
10.	Discard the contents of each well. Remove residual fluids by tapping the plate on absorbent paper. Fill Plate with <b>300 µl of diluted wash buffer</b> . Empty plate and remove residual buffer by tapping the plate on absorbent paper. <b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the plate on absorbent paper after each washing step
11.	Add <b>100 µl CONJ</b> (HRP-Conjugate) into each well

12. Cover the plate tightly and incubate for <b>30 minutes at room temperature</b> (18-26°C)
13. Discard the contents of each well. Remove residual fluids by tapping the plate on absorbent paper. Fill Plate with <b>300 µl of diluted wash buffer</b> . Empty plate and remove residual buffer by tapping the plate on absorbent paper. <b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the plate on absorbent paper after each washing step
14. Add <b>100 µl of SUB</b> (substrate) into each well
15. Incubate for <b>30 minutes at room temperature</b> (18-26°C) in the dark
16. Add <b>100 µl of STOP</b> (stop solution) into each well
17. Determine <b>absorption</b> immediately with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm). If no reference wavelength is available, read only at 450 nm

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "cubic-spline-algorithm".

### 1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

### 2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline-algorithm

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.


## Comparison of Serum and Plasma Samples

In general, plasma samples are recommended for determination of intact proinsulin, because the marker is more stable in plasma than in serum. However, measurement of intact proinsulin in serum samples is possible.

The comparison of the results for serum and plasma samples (n = 40) revealed the following regression equation:

$$Y (\text{plasma concentration}) = 1.3 * (\text{serum}) + 0.76$$

Therefore, to ensure comparability with data from the literature (see Point, 9. Expected Results), the above regression equation must always be considered by the evaluation of the results from measurements of intact proinsulin in serum samples.

## 9. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## Expected values

It is strongly recommended that each laboratory determines its own normal and abnormal values. Studies have been performed with the Intact Proinsulin Kit with adult males and females that had been diagnosed as having type 2 diabetes previously and were being treated with oral anti-diabetes drugs (1-3). Samples from patients with type 2 diabetes with oral medication or dietary treatment were collected from 149 sites that participated in the IRIS-II study. In total, 2,146 male and 2,124 female patients with type 2 diabetes without insulin therapy participated in the study. In an additional study 10 groups of 50 patients, each with incremental homeostasis model assessment (HOMA) scores, were randomly chosen out of a 4,265-person cohort in order to investigate intact proinsulin and adiponectin over a wide range of insulin resistance. Another study evaluated 48 patients with type 2 diabetes and on oral anti-diabetic treatment. Twenty women and 28 men, aged 60 ( $\pm$  9 years), were studied by means of an intravenous glucose tolerance test. Determinations of fasting values of intact proinsulin, insulin, resistin, adiponectin, and glucose were performed. The results of these studies showed that a fasting intact proinsulin concentration of  $\geq 10$  pmol/l predicts the presence of insulin resistance in patients with type 2 diabetes mellitus at a very high specificity and high sensitivity. Fasting proinsulin levels in normal subjects were found to be  $< 10$  pmol/l. To increase the specificity of the detection, and based on recent studies, fasting concentrations less than 11 pmol/L should be considered as normal.

In general, the measurement should be repeated within 3-6 months, if the values in patients with Type 2 Diabetes are in the upper normal range (7-11 pmol/L).

Type 2 Diabetes patients with values  $> 11$  pmol/L may have insulin resistance and insulin secretory dysfunction.

Patients with no Type 2 Diabetes and values  $> 11$  pmol/L should be examined for diabetes mellitus or insulinoma as well for cardiovascular risk factors.

(See also different laboratory service providers, e.g. [www.labor-limbach.de](http://www.labor-limbach.de), [www.bioscientia.de](http://www.bioscientia.de), [www.labor28.de](http://www.labor28.de))

## Limitations

- The values obtained from this assay are intended to aid in diagnosis only. As with all serological tests, interpretation of results obtained with this test must be used in conjunction with the patient's clinical symptoms, medical history and other clinical and/or laboratory findings.

- Only if test instructions are rigidly followed will optimum results be achieved.
- Use fresh plasma or specimens frozen and thawed no more than twice. Specimens that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.
- Reproducible results depend on careful pipetting, observation of incubation periods and temperature, as well as thorough mixing of all prepared solutions.
- While rinsing, check that all wells are filled evenly with Washing Solution, and that there are no residues in the wells.
- Instructions for using appropriate luminometers are to be observed. Check that the instrument has the correct measurement protocol installed.

### *Interfering Substances*

Interferences were studied in accordance with CLSI recommendations (CLSI EP7-A2). To study the effect of lipaemia, test pools were prepared by spiking plasma samples with a commercial lipid emulsion (Intralipid Sigma). Test samples for investigating the effect of haemolysis were obtained by osmotic shock. Icteric samples were prepared by spiking plasma samples with commercial bilirubin (Sigma).

No effect of lipaemia was observed at a lipaemic index up to 975. Interference due to haemolysis was not apparent at a haemolysis index up to 467. Bilirubin produced no apparent interference up to an icterus index of 1065.

## **10. PERFORMANCE CHARACTERISTICS**

### *Precision and reproducibility*

#### **Intra-Assay-Variation**

Three plasma pools were analyzed in 5 individual assays, each in duplicate. The following results were obtained:

Intact Proinsulin (pmol/l)	CV%	n
3.38	2.61	5
27.6	4.47	5
57.2	3.57	5

## Recovery

Five plasma samples containing low endogenous intact proinsulin were spiked with recombinant proinsulin at 3 concentration levels. Recoveries are shown as percentages of the expected result for samples within the range of 9 to 22 pmol/l.

Sample	1	2	3	4	5
Spike 5%	102.4	107.5	100.4	98.8	97.6
Spike 10%	105.1	107.1	102.8	101.9	96.1
Spike 15%	104.4	107.5	102.1	101.3	100.4

Mean spiking recovery was 102.4%.

## Linearity

Four patient samples containing elevated proinsulin concentrations were diluted in Sample dilution buffer. The following table shows the measured intact proinsulin concentrations of the undiluted and diluted specimens.

Measured proinsulin (pmol/l)				
Dilution Factor	Sample 1	Sample 2	Sample 3	Sample 4
0	46.1	46.7	22.8	48.6
1:2	24.1	26.6	12.6	27.5
1:4	11.0	12.7	6.3	12.4
1:8	5.4	5.3	3.3	6.0

## Sensitivity

Sensitivity was estimated as two standard deviations from the mean of 20 replicates of a zero standard. Calculated in this way, analytical sensitivity of the Intact Proinsulin Assay is 0.02 pmol/l. The dynamic range of the assay is 0.02-100 pmol/l.

## High Dose Hook Effect

Because of the assay architecture, which employs separate incubations with solid phase and labeled antibodies, no high dose hook effect is experienced.


### Cross reactivity

Cross reactivities of related proteins were investigated at concentrations of 100 pmol/l. Results are expressed as percentages of the reactivity of an identical concentration of intact proinsulin.

Peptide	CR [%]
Intact proinsulin	100
Insulin	0.0
32-33 split proinsulin	5.6
Des 31-32 split proinsulin	1.4
65-66 split proinsulin	37
Des 65-66 split proinsulin	63
C-peptide	0.0

### Comparison with the Gold Standard Chemoluminescence

#### Correlation Proinsulin intact ELISA / LUMI

## 11. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

## 12. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

## 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

## 14. REFERENCES

Pfützner A, Kunt T, Hohberg C, Mondok A, Pahler S, Konrad T, Lübben G, Forst T (2004) Fasting intact proinsulin is a highly specific predictor of insulin resistance in type 2 diabetes. *Diabetes Care* **27**:682-687

Langenfeld MR, et al. (2004) IRIS II Study: Sensitivity and specificity of intact proinsulin, adiponectin and the proinsulin/adiponectin ratio as markers for insulin resistance. *Diabetes Technology & Therapeutics*. **6**:836-843

Pfützner A, et al. (2005) IRIS II Study: Intact proinsulin is confirmed as a highly specific indicator for insulin resistance in a large cross-sectional study design. *Diabetes Technology & Therapeutics*. **7**:478-486

### Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		