---

**CERTEST**  
*RSV+*  
*Adenovirus Resp.*

---

ONE STEP  
*RSV and Adenovirus Resp.*  
COMBO CARD TEST

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## ENGLISH

### INTENDED USE

CerTest *RSV+Adenovirus Resp.* one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Respiratory Syncytial Virus* and *Adenovirus* antigens from nasal swab, nasopharyngeal wash or aspirate specimens.

CerTest *RSV+Adenovirus Resp.* combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Respiratory Syncytial Virus* and/or *Adenovirus* respiratory infection. The clinical presentation of infections caused by the heterogeneous group of the respiratory viruses can be very similar.

### INTRODUCTION

*Respiratory Syncytial Virus (RSV)* is the most common cause of bronchiolitis and pneumonia among infants and children under 1 year of age. Illness begins most frequently with fever, runny nose, cough and sometimes wheezing. Severe lower respiratory tract disease may occur at any age, especially among the elderly or among those with compromised cardiac, pulmonary or immune systems. *RSV* is spread from respiratory secretions through close contact with infected persons or contact with contaminated surfaces or objects.

Adenoviruses most commonly cause respiratory illness; however, depending on the infecting serotype, they may also cause various other illnesses, such as gastroenteritis, conjunctivitis, cystitis and rash illness. Symptoms of respiratory illness caused by *Adenovirus* infection range from the common cold syndrome to pneumonia, croup and bronchitis. Patients with compromised immune systems are especially susceptible to severe complications of *Adenovirus* infection. *Adenovirus* is transmitted by direct contact, fecal-oral transmission and occasionally waterborne transmission. Some types are capable of establishing persistent asymptomatic infections in tonsils, adenoids and intestines of infected hosts and shedding can occur for months or years.

### TEST PRINCIPLE

CerTest *RSV+Adenovirus Resp.* is based on the principle of a qualitative immunochromatographic assay for the determination of *Respiratory Syncytial Virus* and *Adenovirus* from nasal swab, nasopharyngeal wash or aspirate specimens.

**Strip A** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Respiratory Syncytial Virus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*RSV*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

**Strip B** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Adenovirus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Adenovirus*) conjugated to blue polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is *RSV* positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-*RSV* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is *Adenovirus* positive, the antigens of the diluted sample react with the blue-coloured conjugated complex (anti-*Adenovirus* monoclonal antibodies-blue polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*RSV* antibodies present on the membrane of strip A (test line) and the anti-*Adenovirus* antibodies present on the membrane of strip B (test line) capture the coloured conjugate and one red and/or one blue line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *RSV* and *Adenovirus* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red/blue-coloured conjugate complex. The anti-*RSV* and anti-*Adenovirus* antibodies present on the membranes (test lines) will not capture the antigen-red/blue-coloured conjugate complex (not formed), for which the red/blue lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

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### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *RSV+Adenovirus Resp.* combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


### SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

#### Specimen preparation (see illustration):

##### Nasal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.

Specimen collection

##### Nasopharyngeal aspirate method (suction apparatus, sterile suction catheter):

###### For adult:

1. Place the irrigator up to the nose.
2. Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side. Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
3. Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

Specimen collection

###### For children:

1. Use an aspiration bulb or bulb syringe to instil the saline water into one nostril leaning the children head.
2. Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
3. Repeat for the other nostril and transfer the fluid into the same specimen container.

### MATERIALS

#### MATERIALS PROVIDED

- CerTest *RSV+Adenovirus Resp.* combo card tests
- Reagent B (sample diluent)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- CerTest *RSV* Positive Control swab + Instructions for use
- CerTest *Adenovirus Resp.* Positive Control swab + Instructions for use

#### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer
- Vortex



## TEST PROCEDURE

**Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.**

### - Procedure A using nasal swab samples:


1. Add 15 drops (1) Reagent B and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest *RSV+Adenovirus Resp.* combo card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense exactly 4 drops from the testing tube, into the circular window marked with the letter A (3) and 4 drops, using the same tube, into the circular window marked with the letter B (4).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.

### RSV Strip A-procedure


Add 4 drops in A-circular window

(3)

### Adenovirus Strip B-procedure

Add 4 drops in B-circular window

(4)

### Procedure B using nasopharyngeal wash or aspirate samples:

1. Add 6 drops (1) of the nasopharyngeal wash or aspirate samples with a pipette and 9 drops (2) of Reagent B in a testing tube. Mixer with vortex at least 1 minute to homogenize. Best results are obtained when the specimen is vigorously extracted in the solution (3).
2. Remove the CerTest *RSV+Adenovirus Resp.* combo card test from its sealed bag just before using it.
3. Dispense exactly 4 drops from the testing tube, into the circular window marked with the letter A (4) and 4 drops, using the same tube, into the circular window marked with the letter B (5).
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.


### RSV Strip A-procedure

### Adenovirus Strip B-procedure


- (1) Add 6 drops of nasopharyngeal wash/aspirate
- (2) Add 9 drops of Reagent B
- (3) Mix the solution with vortex 1 minute

Add 4 drops in A-circular window

Add 4 drops in B-circular window

(4)

(5)

### INTERPRETATION OF THE RESULTS (please refer to the illustration below)

1. 

A: Green → Negative RSV  
B: Green → Negative Adenovirus
2. 

A: Green/Red → Positive RSV  
B: Green/Blue → Positive Adenovirus
3. 

A: Green/Red → Positive RSV  
B: Green → Negative Adenovirus
4. 

A: Green → Negative RSV  
B: Green/Blue → Positive Adenovirus

### INVALID RESULTS

- 5.

	A (RSV)	B (Adenovirus)	Interpretation of the results
1.	- GREEN	- GREEN	There is no <i>Respiratory Syncytial Virus (RSV)</i> and <i>Adenovirus</i> presence. No infection caused by <i>RSV</i> and <i>Adenovirus</i> .
2.	+ GREEN-RED	+ GREEN-BLUE	There is <i>Respiratory Syncytial Virus (RSV)</i> and <i>Adenovirus</i> presence. Simultaneous infection caused by <i>RSV</i> and <i>Adenovirus</i> .
3.	+ GREEN-RED	- GREEN	There is <i>Respiratory Syncytial Virus</i> presence. Infection caused by <i>RSV</i> .
4.	- GREEN	+ GREEN-BLUE	There is <i>Adenovirus</i> presence. There is respiratory infection caused by <i>Adenovirus</i> .
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

**INVALID:** Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED/BLUE). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

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## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red and/or blue coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

## LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. CerTest RSV+Adenovirus Resp. combo card test should be used only with nasal swab, nasopharyngeal wash or aspirate specimens. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of RSV and/or Adenovirus respiratory infection. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. A negative result is not meaningful because of it is possible the antigens concentration in the nasal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture or PCR.

## EXPECTED VALUES

Respiratory Syncytial Virus (RSV) is the leading cause of emergency visits and hospitalization for acute lower respiratory tract infections (LRTI) in infants and young children worldwide.

In adults, older and healthy children, the symptoms of Respiratory Syncytial Virus are mild and typically mimic the common cold. Self-care measures are usually all that's needed to relieve any discomfort.

Infection with Respiratory Syncytial Virus can be severe in some cases, especially in premature babies and infants with underlying health conditions. RSV can also become serious in older adults, adults with heart and lung diseases, or anyone with a very weak immune system (immunocompromised).

The exact prevalence and incidence of adenoviral infections are unknown, because most cases are seen by general practitioners and optometrists. Adenovirus is a very common infection, estimated to be responsible for between 2% and 5% of all respiratory infections. In winter, infection with type 4 or 7 causes recognisable illness in military recruits, with about 25% requiring hospitalisation for fever and lower respiratory tract disease.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity (Strip A)

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest RSV+Adenovirus Resp., CerTest) and another commercial available immunochromatographic test (BinaxNOW® RSV, Alere). The results were as follows:

IC test: CerTest RSV+Adenovirus Resp.( RSV Strip A)	IC test: BinaxNOW® RSV		
	+	-	Total
+	18	0	18
-	1	10	11
Total	19	10	29

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### Clinical sensitivity and specificity (Strip B)

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest RSV+Adenovirus Resp., CerTest), another commercial available immunochromatographic test (Adenovirus Respi, CorisBioConcept) and a commercial available Immunofluorescence test (PathoDx®Adenovirus, Remel). The results were as follows:

		IFI test: PathoDx®Adenovirus		
		+	-	Total
IC test: RSV+Adenovirus Resp.(Adenovirus Strip B)	+	20	0	20
	-	0	5	5
	Total	20	5	25
		IC test: Adenovirus Respi		
		+	-	Total
IC test: RSV+Adenovirus Resp.(Adenovirus Strip B)	+	20	0	20
	-	0	5	5
	Total	20	5	25

	Sensitivity	Specificity	PPV	NPV
CerTest RSV+Adenovirus Resp. vs PathoDx®Adenovirus and Adenovirus Respi	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Respiratory Syncytial Virus* and *Adenovirus* using CerTest RSV+Adenovirus Resp..

### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest RSV+Adenovirus Resp.; no cross reactivity against organisms that cause other respiratory infections:

- RSV Strip A no cross reaction with:
  - Influenza type A
  - Influenza type B
  - Adenovirus
- Adenovirus Strip B no cross reaction with:
  - Influenza type A
  - Influenza type B
  - Respiratory Syncytial Virus

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- MARCELA ECHAVARRIA, JOSE L. SANCHEZ, et al. "Rapid Detection of Adenovirus in Throat Swab Specimens by PCR during Respiratory Disease Outbreaks among Military Recruits", *Journal of Clinical Microbiology*, Feb. 2003, Vol. 41, No. 2, p. 810–812
- MARILYN J. AUGUST AND ANN L. WARFORD; "Evaluation of a Commercial Monoclonal Antibody for Detection of Adenovirus Antigen" *Journal of Clinical Microbiology*, Nov. 1987, Vol. 25, No. 11; p. 2233-2235
- ANN R. FALSEY and EDWARD E. WALSH. "Respiratory Syncytial Infection in Adults", *Clinical Microbiology Reviews*, July 2000, Vol. 13, No.3, p. 371-384.
- DIANE C. HALSTEAD, SANDRA TODD, and GALE FRITCH; "Evaluation of Five Methods for Respiratory Syncytial Virus Detection" *Journal of Clinical Microbiology*, May 1990, Vol. 28, No. 5; p. 1021-1025.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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