anti-heTG IgA ELISA

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
***For the in vitro determination of anti-epidermal
transglutaminase IgA antibodies in serum and plasma***

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CONTENT	PAGE
1. INTENDED USE	15
2. PRINCIPLE OF THE TEST	15
3. MATERIAL SUPPLIED	16
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. PREPARATION AND STORAGE OF REAGENTS	17
5. PREPARATION AND STORAGE OF REAGENTS	18
6. PRECAUTIONS	19
7. SPECIMEN COLLECTION AND PREPARATION	19
8. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
9. RESULTS	21
10. LIMITATIONS	21
11. QUALITY CONTROL	21
EXPECTED VALUES	21
12. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
13. TECHNICAL HINTS	23
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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

anti-heTG IgA ELISA

For the in vitro determination of anti-epidermal Transglutaminase IgA antibodies in serum and plasma

Valid from 23.02.2012

K 9396

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **IgA anti-epidermal-transglutaminase antibodies** in serum and plasma. It is for *in vitro* diagnostic use only.

2. PRINCIPLE OF THE TEST

This Enzyme Immuno Assay is a sandwich assay for the quantitative determination of anti-epidermal transglutaminase antibodies (IgA anti-heTG). The wells of the microtiter plate are coated with the antigen. In a first incubation step, the anti-heTG antibodies are bound to the coated antigen. To remove all unbound substances, a washing step is carried out.

In a further incubation step, a peroxidase-labeled antibody is added. After another washing step to remove all unbound substances, the solid phase is incubated with the peroxidase substrate, tetramethylbenzidine (TMB). Finally, an acidic stop solution is added to terminate the enzymatic reaction, whereby the colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the IgA anti-epidermal-transglutaminase antibody concentration. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standards. IgA anti-epidermal-transglutaminase antibodies present in the samples are determined directly from this curve.

Indication

- See our proposal for stepwise diagnosis of gluten sensitivity on page 17.


3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 9396MTP	PLATE	One holder with pre-coated strips	8 x 12 wells
K 9396WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9396ST	STD	heTG-Standard, lyophilized	4 x 1 vial
K 9396KO1	CTRL	Control, lyophilized	4 x 1 vial
K 9396KO2	CTRL	Control, lyophilized	4 x 1 vial
K 9396K	CONJ	Conjugate (rabbit-anti-IgA antibody, peroxidase labeled)	1 x 200 µl
K 9396TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9396AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker with 37 °C incubator
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0,2 µm) with an electrical conductivity of 0,055µS/cm at 25°C (≤18,2MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be re-dissolved at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The **CONJ** (conjugate) must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted conjugate is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- The lyophilized **STD** (standard) and **CTRL** (controls) are stable at 2-8 °C until expiry date stated on the label. Volume and concentration see product specification. **Reconstituted standard and controls are not stable and cannot be stored.**

The preparation of the **standard curve** (volumes and concentrations) is described in the product specification.

- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

6. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

7. SPECIMEN COLLECTION AND PREPARATION

Serum/Plasma

Serum and plasma have to be diluted 1:250 in ELISA wash buffer.

Example:

The dilution has to be carried out/taken in two dilution steps.

Dilution step 1: 10 µl serum + 90 µl wash buffer

Dilution step 2: 40 µl dilution 1 + 960 µl wash buffer

8. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the pre-coated microtiter plate **5 x with 250 µl** of wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** of **STD** (standards), **CTRL** (controls) and **diluted samples**.
2. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of diluted **CONJ** (conjugate).
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
7. Add **100 µl** of **SUB** (TMB substrate).
8. Incubate for **10 - 20 minutes** at room temperature.
9. Add **100 µl STOP** (ELISA stop solution) and mix shortly.
10. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

9. RESULTS

Calculation

Since the sample dilution of 1:250 is already considered in the calibration curve, the dilution factor is 1.

If a sample dilution different from 1:250 has been used, a corresponding dilution factor must be considered, e.g., at a dilution of 1000 the results must be multiplied by 4.

10. LIMITATIONS

Samples with IgA anti-heTG levels higher than the highest standard should be further diluted with wash buffer and re-assayed.

11. QUALITY CONTROL

Controls should be analyzed with each run. Results, generated from the analysis of controls, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality controls are outside the acceptable limits.

Expected values

Preliminary normal value

Serum > 22 AU/ml positive

Serum < 16 AU/ml negative

Results between 16 - 22 AU/ml are in the grey range.

We recommend each laboratory to establish its own normal range values.

12. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The intra-assay variation of the assay was calculated from 64 determinations of one sample by one person.

Intra-Assay CV n = 64

Sample	IgA anti-heTG mean value [AU/ml]	Intra-Assay CV [%]
1	83,81	4,1

The inter-assay variation of the assay was calculated from data on 2 samples obtained in 15 or 16 determinations by different persons on two different days.

Inter-Assay CV n = 15

Sample	IgA anti-heTG mean value [AU/ml]	Inter-Assay CV [%]
1 (n=15)	26,18	11,4
2 (n=16)	12,38	15,9

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in-vitro* diagnostics only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number