Thymosin α 1 ELISA KIT

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14. LITERATUR

1. Melatonin is responsible for the nocturnal increase observed in serum and thymus of thymosin alpha1 and thymulin concentrations: observation in rats and humans. Molinero P et al. (2000) J Neuroimmunol 103:180-188
2. Effect of Thymosin α 1 on Hypothalamic Hormone Release. Milenkovic et al. (1992): Neuroendocrinol 56: 674-679
3. Anti-Tumor Effect of Combined Treatment with Thymosin alpha 1 and Interleukin-2 after 5-Fluorouracil in Liver Metastases from Colorectal Cancer in Rats. Rasi et al. (1994) Int J Cancer 57:701-705
4. Antitumor Effect of Thymosin α 1/Interleukin-2 or Thymosin α 1/Interferon α , β Following Cyclophosphamide in Mice Injected with Highly Metastatic Friend Erythroleukemia Cells. Garaci et al. (1993) J Immunotherapy 13:7-17
5. Determination of Thymosin α 1 with enzyme-immunoassay in colorectal cancer patients. Jevromovic et al. (1997) Archive of Oncology 5:193
6. A randomized trial to evaluate the immunorestorative propertise of synthetic Thymosin α 1 in patients with lung cancer. Schulof et al. (1985) J Biol resp Med 4:147-158

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
Thymosin α 1 ELISA KIT

*For the in vitro determination of thymosin α 1
in serum and thymus extract*

Valid from 2014-12-15

REF K 9510

RUO

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	12
2. INTRODUCTION	12
3. MATERIAL SUPPLIED	12
4. MATERIAL REQUIRED BUT NOT SUPPLIED	13
5. STORAGE AND PREPARATION OF REAGENTS	13
6. STORAGE AND PREPARATION OF SAMPLES	14
7. ASSAY PROCEDURE	14
<i>Principle of the test</i>	14
<i>Precubation</i>	14
<i>Test procedure</i>	14
8. RESULTS	16
9. LIMITATIONS	16
10. QUALITY CONTROL	16
<i>Reference range</i>	17
11. PRECAUTIONS	17
12. TECHNICAL HINTS	17
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	18
14. REFERENCES	18

1. INTENDED USE

Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) Kit is intended for the quantitative determination of Thymosin α 1 in serum and thymus extract. For research use only. Not for use in diagnostic procedures.

2. INTRODUCTION

Thymosin α 1 was the first single peptide isolated from thymus fraction 5. It acts on T-helper and NK-cells. Thymosin α 1 has been reported for exert effects on hormone regulating the hypothalamus (1). Thymosin α 1 has been demonstrated to have beneficial effects in animal models of liver and colon carcinoma (2) or leukaemia (3). Its use as a prognostic factor in human studies, e.g. colon carcinoma (4) has been discussed. Thymosin α 1 has been successfully used as component combined chemotherapy in bronchial carcinoma (5).

Indications

- Disorder of immune system
- Control of immune status in association with a chemotherapy
- Disorder of endocrinium
- Quality control of thymus extracts

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K9510	PLATE	One holder with precoated strips	12 x 8 wells
K9510	WASHBUF	ELISA wash buffer concentrate, 10x	1 x 100 ml
K9510	AB	Antibody (rabbit anti Thymosin α 1), ready to use	3 x 3,5 ml
K9510	CONJ	Conjugate (goat anti rabbit, Peroxidase-labeled), ready to use	22 ml
K9510	STDBUF	Standard dilution buffer, ready to use	50 ml
K9510	STD	Calibrator concentrate, lyophilized	3 x 1 vial
K9510	SUB	TMB substrate (Tetramethylbenzidine), ready to use	2 x 15 ml
K9510	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as purchase order number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 μ l tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 *g*
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 μ m) with an electrical conductivity of 0.055 μ S/cm at 25 °C (\geq 18.2 M Ω cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 μ l** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted **1:10 in ultra pure water** before use (100 ml concentrate + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or 37 °C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2–8 °C** until the expiry date stated on the label. **Diluted buffer solution** (wash buffer) can be stored in a closed flask at **2–8 °C for one month**.
- The **AB** (antibody) is ready to use. It can be stored at 2–8 °C up to 4 weeks. Long time storage until the expiry date given on the label has to be at **-20 °C**.
- The **lyophilized standard** (STD) is stable at **2–8 °C** until the expiry date stated on the label. Before use, the standard must be reconstituted with **STDBUF** (standard dilution buffer) (volume, concentration and dilution schema for the calibration curve see product specification). Allow the vial content to solve for 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standard is not stable and can not be stored.**

- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Serum

Serum can be used **without dilution**. Store samples at -20 °C.

Thymusextract

Thymus extracts have varying compositions. Please contact the supplier when using thymus extracts.

7. ASSAY PROCEDURE

Principle of the test

In this enzyme immunoassay (EIA) for the determination of Thymosin α 1 ($T\alpha$ 1) polyclonal rabbit antibodies directed against synthetic $T\alpha$ 1 are used. The test principle is based on a competition between antigen in the sample or standards and the antigen coated on the wells of microplate. A peroxidase-conjugated antibody is used for detection and quantification, and tetramethylbenzidine (TMB) as a peroxidase substrate. The enzymatic reaction is terminated by an acidic stop solution. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from standard. $T\alpha$ 1 present in the patient samples is determined directly from this curve.

Precubation

1.	Add 200 μl of the AB (1st antibody) to 400 μl of STD (calibrator) or sample solutions into one way test tubes.
2.	Incubate for 18 hours at 2-8°C on a shaker.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many **microtiter strips** as needed from kit. Store unused strips in the aluminium foil bag with desiccant at 2–8 °C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Wash each well 5 x by dispensing 250 μl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
2.	Add 200 μl of the preincubated mixture into respective well.
3.	Cover the plate tightly and incubate for 90 min at 2–8°C in the dark on a horizontal shaker.
4.	Discard the contents of each well. Wash each well 5 x by dispensing 250 μl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
5.	Add 200 μl CONJ (conjugate) into each well.
6.	Cover the plate tightly and incubate for 1 hour at room temperature (15-30°C) on a horizontal shaker.
7.	Discard the contents of each well. Wash each well 5 x by dispensing 250 μl of diluted WASHBUF (washbuffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
8.	Add 200 μl of SUB (substrate) into each well.
9.	Incubate for 20–30 min.* at room temperature (15-30°C) in the dark.
10.	Add 50 μl of STOP (stop solution) into each well, mix shortly.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform.

For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Thymus extract

For the calculation of the sample values the results from the microplate reader has to be multiplied with the selected dilution factor.

9. LIMITATIONS

Samples with an OD lower than the OD of the highest standard should be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

We recommend each laboratory to establish its own baseline values.

Baseline values depend on the patient's age and vary between different individuals.

11. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

14. REFERENCES

1. Melatonin is responsible for the nocturnal increase observed in serum and thymus of thymosin alpha1 and thymulin concentrations: observation in rats and humans. Molinero P et al. (2000) *J Neuroimmunol* 103:180-188
2. Effect of Thymosin α 1 on Hypothalamic Hormone Release. Milenkovic et al. (1992): *Neuroendocrinol* 56: 674-679
3. Anti-Tumor Effect of Combined Treatment with Thymosin alpha 1 and Interleukin-2 after 5-Fluorouracil in Liver Metastases from Colorectal Cancer in Rats. Rasi et al. (1994) *Int J Cancer* 57:701-705
4. Antitumor Effect of Thymosin α 1/Interleukin-2 or Thymosin α 1/Interferon α , β Following Cyclophosphamide in Mice Injected with Highly Metastatic Friend Erythroleukemia Cells. Garaci et al. (1993) *J Immunotherapy* 13:7-17
5. Determination of Thymosin α 1 with enzyme-immunoassay in colorectal cancer patients. Jevromovic et al. (1997) *Archive of Oncology* 5:193
6. A randomized trial to evaluate the immunorestorative propertise of synthetic Thymosin α 1 in patients with lung cancer. Schulof et al. (1985) *J Biol resp Med* 4:147-158

Used symbols:

Temperature limitation

Catalogue Number

For research use only

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by