# Thymosin $\beta_4$ EIA Kit

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<b>1. INTENDED USE</b>	<b>12</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>12</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>12</b>
<b>4. MATERIAL SUPPLIED</b>	<b>13</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>13</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>14</b>
<b>7. PRECAUTIONS</b>	<b>15</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>15</b>
<b>9. ASSAY PROCEDURE</b>	<b>15</b>
PROCEDURAL NOTES	15
TEST PROCEDURE	16
<b>10. RESULTS</b>	<b>17</b>
TYPICAL CALIBRATION CURVE	17
<b>11. LIMITATIONS</b>	<b>18</b>
<b>12. QUALITY CONTROL</b>	<b>18</b>
EXPECTED VALUES	18
<b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>18</b>

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Hersteller

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# Thymosin $\beta_4$ EIA Kit


*For the in vitro determination of Thymosin  $\beta_4$  in serum and thymus extract*

*For research use only*

Valid from 16.08.2011

K 9520

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Thymosin  $\beta_4$**  in serum and thymus extract. It is for research use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

The beta-thymosins are a family of highly conserved polar 5 kDa peptides originally thought to be thymus hormones. Further studies demonstrated that beta-thymosins are ubiquitous; they have been found in a variety of tissues and cell lines. The highest concentrations have been detected in spleen, thymus, lung, and peritoneal macrophages.

**Thymosin  $\beta_4$**  binds monomeric actin in a 1:1 complex acting as an actin buffer, preventing polymerization into actin filaments but supplying a pool of actin monomers when needed by the cell. Changes in the expression of **Thymosin  $\beta_4$**  appear to be related to the differentiation of cells. **Thymosin  $\beta_4$**  is detected outside of cells in blood plasma and wound fluid. Several biological effects are attributed to **Thymosin  $\beta_4$** , like induction of metalloproteinases, chemotaxis, angiogenesis and inhibition of inflammation as well as inhibition of bone marrow stem cell proliferation.

## 3. PRINCIPLE OF THE TEST

This enzyme immuno assay (EIA) can be used for the determination of **Thymosin  $\beta_4$  ( $\text{TB}_4$ )** in serum and thymus preparations.

The test principle is based on a competition between antigen in the sample or standards and the antigen coated on the wells of microplate. A peroxidase-conjugated antibody is used for detection and quantification, and tetramethylbenzidine (TMB) as a peroxidase substrate. The enzymatic reaction is terminated by acidic stop solution. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from standard.  **$\text{TB}_4$**  present in the patient samples is determined directly from this curve.

#### 4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 9520MTP	MTP	Microtiter plate, precoated	12 x 8 wells
K 9520WB	WASHBUF	ELISA wash concentrate 10 x	100 ml
K 9520 PV	SAMPLEBUF	Sample dilution buffer	100 ml
K 9520 A	AB	Antibody (Anti-Thymosin $\beta_4$ antibody)	16 ml
K 9520K	CONJ	Conjugate (Peroxidase-labeled), ready to use	22 ml
K 9520ST	STD	Calibrators, ready to use	3 x 5 vials
K 9520 KO1	CTRL	control	3 vials
K 9520 KO2	CTRL	control	3 vials
K 9520TMB	SUB	TMB substrate (Tetramethylbenzidine)	2 x 15 ml
K 9520AC	STOP	ELISA stop solution, ready to use	15 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 5-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm  
(reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100  $\mu$ l** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at **2-8°C** until the expiry date stated on the label. Reconstitution details are given in the data sheet.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

## 7. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Serum

Serum can be used without dilution. Store samples at -20 °C.

### Thymus extract

Thymus extracts have varying compositions. Please contact the supplier when using thymus extracts.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Wash the precoated microtiter plate **5 x with 250  $\mu$ l** ELISA wash buffer.  
Carry out the tests in duplicate.


1. Add **50  $\mu$ l** of **STD** (standards), **CTRL** (controls) or **samples** into each well in duplicate
2. Add **150  $\mu$ l** of **AB** (anti-thymosin  $\beta_4$  antibody solution) into each well
3. Incubate for **1 hour**, shaking on a horizontal mixer in the dark, at room temperature
4. Aspirate and wash the wells **5 x with 250  $\mu$ l** ELISA wash buffer
5. Add **200  $\mu$ l** of **CONJ** (conjugate) into each well
6. Incubate for **1 hour**, shaking on a horizontal mixer in the dark, at room temperature
7. Aspirate and wash the wells **5 x with 250  $\mu$ l** ELISA wash buffer
8. Add **200  $\mu$ l** of **SUB** (substrate solution)
9. Incubate for **10-20 minutes** at room temperature
10. Add **50  $\mu$ l** of **STOP** (stop solution) and mix shortly
11. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

## 10. RESULTS

A calibration curve is constructed from the calibrators. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

<b>Concentration [µg/ml]</b>	10	2.5	0.62	0.16	0
<b>OD mean values</b>	0.320	0.833	1.308	1.894	2.073

The data is for demonstration only and cannot be used for the evaluation of test results.

### Thymus extract

For the calculation of the sample values the results from the microplate reader has to be multiplied with the selected **dilution factor**.

## 11. LIMITATIONS

Samples with levels greater than the highest standard value should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercially control samples as internal quality control.

**Control samples** should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

We recommend each laboratory to establish its own baseline values.

Baseline values depend on the patient's age and vary between different individuals.

## 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number