Tryptophan ELISA Kit

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15. LITERATUR

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *Am J Med Genet B Neuropsychiatr Genet.* 2009 Mar 5; 150B (2): 202-8.

Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. *International Journal of Tryptophan Research* 2009; 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry.* 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophan. *Er Neurol.* 1986 ; 25 (Suppl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr.* 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. J Affect Disord. 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Relation of serum cholesterol, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. J Psychiatry Neurosci. 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; Neurogastroenterol Motil 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; Stress, Volume 11, Issue 3 2008, 198-209.

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

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Hersteller

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Manual

Tryptophan ELISA Kit



*For the in vitro determination of Tryptophan in human
EDTA plasma, serum and urine*

For research use only

Valid from 15.07.2014

K 7730

Immundiagnostik AG

Table of Contents

1. INTENDED USE	19
2. INTRODUCTION / CLINICAL RELEVANCE	19
3. MATERIAL SUPPLIED	21
4. MATERIAL REQUIRED BUT NOT SUPPLIED	21
5. PREPARATION AND STORAGE OF REAGENTS	22
6. SPECIMEN COLLECTION AND PREPARATION	23
7. ASSAY PROCEDURE	23
<i>PRINCIPLE OF THE TEST</i>	23
<i>SAMPLE PREPARATION PROCEDURE</i>	24
<i>TEST PROCEDURE</i>	25
8. RESULTS	26
9. LIMITATIONS	27
10. QUALITY CONTROL	28
<i>REFERENCE RANGE</i>	28
11. PERFORMANCE CHARACTERISTICS	28
<i>PRECISION AND REPRODUCIBILITY</i>	28
<i>SPIKING RECOVERY</i>	29
<i>DILUTION RECOVERY</i>	29
<i>ANALYTICAL SENSITIVITY</i>	29
<i>SPECIFICITY</i>	30
<i>CORRELATION WITH HPLC-MS</i>	30
12. PRECAUTIONS	30
13. TECHNICAL HINTS	31
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	31
15. REFERENCES	32

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of Tryptophan in EDTA plasma, serum and urine. For research use only. Not for use in diagnostic procedures.

2. INTRODUCTION / CLINICAL RELEVANCE

Tryptophan is an essential amino acid, which cannot be synthesized by humans and, therefore, must be obtained from diet. Tryptophan plays an important role in the synthesis of serotonin, as it is the only source for synthesizing serotonin (see figure 1).

Serotonin itself is released from the appendices of a widespread transmitter system with global modulating effects on different areas of life. The activity of this system depends not least on the actual tryptophan level. In some disease patterns a tryptophan deficiency is found, which is induced by an activation of the enzyme indoleamine 2,3-dioxygenase (IDO). The activation of IDO causes a shifting of the metabolism of tryptophan towards the kynurenine pathway and its degradation (see figure 1).


Figure 1: Synthesis of serotonin and melatonin *via* tryptophan (left), and shifting of the metabolism of tryptophan to the kynurenine-pathway *via* the enzyme indoleamine 2,3-dioxygenase (right).

As a consequence, the tryptophan level in plasma is reduced and it is not available for serotonin synthesis and the following melatonin production. As a result, fluctuations of mood, aggression, sleep and eating disorder, or depression can occur.

In recent years, an increasing number of disease patterns have been detected, in which either the tryptophan level is significantly lowered due to the activation of IDO or for other reasons, or which show an adverse ratio of tryptophan-to-kynurenine levels. In the former case, the tryptophan deficiency is associated with an increase in inflammation markers. Some indications are mentioned in the following:

- In case of **irritable bowel syndrome** (IBS) interferon-gamma activates the enzyme IDO causing a tryptophan deficiency.
- In **obese patients** inflammation markers, which show a chronic activation of the immune system, are increased as well. IDO is also activated here, and the plasma tryptophan levels are low.
- **Stress** activates IDO likewise and facilitates the degradation of tryptophan *via* the kynurenine pathway.

In addition, gastrointestinal dysfunctions like **fructose intolerance** are associated with difficulties in absorbing tryptophan from the gut.

In all these disease patterns fluctuation of mood, sleep and eating disorder or depression are known.

The symptoms of these diseases can be eased by administration of tryptophan. However, such a treatment requires the knowledge of the tryptophan level, which can be determined appropriately with the now available test.

3. MATERIAL SUPPLIED

Catalog No	Label	Kit Components	Quantity
K 7730MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 7730ST	STD	Standards, ready to use (0, 6, 20, 60, 200, 600 μ M)	6 x 200 μ l
K 7730KO1 K 7730KO2	CTRL 1 CTRL 2	Controls, ready to use (see specification for range)	2 x 200 μ l
K 7730WP	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 7730AK	AB	Anti-tryptophan antibody, lyophilized	2 x 1 vial
K 7730K	2.AB	POD antibody (concentrate, 200x)	1 x 90 μ l
K 7730CSP	2.ABDIL	Conjugate stabilizing buffer, ready to use	1 x 16 ml
K 7730RP	REABUF	Reaction buffer, ready to use	2 x 100 ml
K 7730DR	DER	Derivatization reagent	2 x 13.3 mg
K 7730LM	DMSO	Dimethylsulfoxide (DMSO)	2 ml
K 7730AP	ASYBUF	Assay buffer, ready to use	2 x 30 ml
K 7730TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 7730AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipets and 10-1000 μ l tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The wash buffer concentrate is stable at **2-8°C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month.**
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20°C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use. They can be re-frozen up to 3 times.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (13.3 mg) in 0.8 ml DMSO.** Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use.** When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.
- Dissolve the content of one vial of **anti-tryptophan antibody (AB) in 3 ml of diluted wash buffer.** The ELISA kit can be separated into two performances by providing two AB vials. Diluted anti-tryptophan antibody can be stored at **2-8°C for 4 weeks.**

- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL) (**e.g. 60 µl 2.AB + 12 ml 2.ABDIL, prepare only the required amount**). The undiluted POD antibody (2.AB) is stable at 2-8°C until the expiry date stated on the label. Diluted POD antibody can be stored at **2-8°C for 1 week**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

6. SPECIMEN COLLECTION AND PREPARATION

EDTA plasma, serum and urine

- Serum, EDTA-Plasma and urine are suited for this test system. The samples are stable for up to 6 hours at 2-8°C. For longer storage, blood- and urine samples must be frozen at -20°C. We recommend acidifying the urine samples.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- Samples with visible amounts of precipitates should be centrifuged.
- The EDTA plasma, serum and urine samples are diluted for derivatization (see sample preparation procedure).
- For sample preparation, a derivatization reagent (DER) for derivatization of tryptophan is added (details are given in the sample preparation procedure).

7. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for tryptophan derivatization. Afterwards, the treated samples and a polyclonal tryptophan antiserum are incubated in the wells of a microtiter plate coated with tryptophan derivative (tracer). During the incubation period, the target tryptophan in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The tryptophan in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-

bound antibody is inverse proportional to the tryptophan concentration in the sample.

During the second incubation step a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-tryptophan antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow color is inverse proportional to the tryptophan concentration in the sample; this means, high tryptophan concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. Tryptophan present in the patient samples is determined directly from this curve.

Sample preparation procedure

Dilute standards (STD), controls (CTRL) and samples (SAMPLE) by factor **1:41** as follows:

25 µl STD/CTRL/SAMPLE + 1000 µl REABUF (reaction buffer)

Please note: Samples from patients with tryptophan supplementation (e.g. in depletion studies) probably require further dilution, we recommend diluting these samples additionally 1:2 with REABUF.

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials) as follows:

- | | |
|----|---|
| 1. | Bring all reagents and samples to room temperature (15-30°C) and mix well. |
| 2. | Add 100 µl of diluted standards/ controls/ samples (STD/ CTRL/ SAMPLE) in the corresponding vial. |
| 3. | Add 25 µl of freshly prepared derivatization reagent (DER) into each vial (STD, CTRL, SAMPLE) and mix thoroughly by repeated inversion or several seconds on a vortex mixer. Incubate for 45 min at room temperature (15-30°C) on a shaker (180-240 rpm). |

4. Afterwards add **1000 µl of assay buffer (ASYBUF)** into each vial, mix well and incubate for **45 min at room temperature** (15-30°C) on a shaker (180-240 rpm).

2 x 50 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5. Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) in duplicate on a **protocol sheet**.
6. Take as many **microtiter strips (PLATE)** as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
7. Wash the microtiter strips **5x with 250 µl of diluted ELISA wash buffer before use**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
8. For the analysis in duplicate, take **2 x 50 µl** of derivatized **standards (STD)/ controls (CTRL)/ samples (SAMPLE)** out of the vials and add into the respective wells of the microtiter plate.
9. Add **50 µl of anti-tryptophan antibody (AB)** into each well. Cover the plate tightly.
10. Incubate overnight (**15-20 hours**) at **2-8°C**.
11. Aspirate the contents of each well. Wash each well **5 x** with **250 µl of diluted wash buffer**. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
12. Add **100 µl** of diluted **POD antibody (2. AB)** into each well.
13. Cover plate tightly and incubate for **1 hour at room temperature** (15-30°C) on a horizontal shaker (180-240 rpm).

14.	Aspirate the contents of each well. Wash each well 5 x with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
15.	Add 100 µl of TMB substrate (SUB) into each well.
16.	Incubate for 8-12 min at room temperature (15-30°C) in the dark*.
17.	Add 100 µl of stop solution (STOP) into each well, mix thoroughly.
18.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls, and samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results.**

Exception: The results from additionally diluted samples must be multiplied by this dilution factor.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.


3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

The concentration of controls and patient samples can be determined directly from the calibration curve in $\mu\text{mol/l}$. In the following, an example of a calibration curve is given; do not use it for the calculation of your results.

Example of calibration curve

9. LIMITATIONS

Samples with tryptophan concentrations above $600 \mu\text{mol/l}$ should be further diluted and re-assayed.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

Reference Range

Based on internal studies with plasma samples of apparently healthy persons (n=40) a mean value of 54 µmol/l was estimated. The standard deviation was 12 µmol/l.

EDTA plasma mean value ± 2x standard deviation: 54 ± 24 µmol/l

Normal range: 30 - 78 µmol/l

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-assay (n=12)

Sample	Tryptophan [µmol/l]	CV [%]
1	17.7	11.4
2	56.7	4.1

Inter-assay (n=5)

Sample	Tryptophan [µmol/l]	CV [%]
1	19.4	6.9
2	59.1	4.0

Spiking Recovery

One serum sample was spiked with different tryptophan concentrations and measured in this assay. The mean recovery rate for all concentrations was 100.8 % (n=5).

Spike [μmol/l]	Tryptophan expected [μmol/l]	Tryptophan measured [μmol/l]	Recovery [%]
0		44.8	
25	69.8	74.7	107.1
50	94.8	89.5	94.4

Dilution Recovery

One spiked serum sample was diluted and measured in this assay. The mean recovery was 99.8 % (n=5).

Dilution	Tryptophan expected [μmol/l]	Tryptophan measured [μmol/l]	Recovery [%]
original		44.8	
1+1	22.4	22.5	100.4
1+3	11.2	11.1	99.1

Analytical Sensitivity

The zero-standard was measured 20 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be 7 μmol/l.

Sample	Tryptophan mean value [OD]	2 x standard deviation (SD)	Detection limit [μmol/l]
Zero-standard	2.3	0.05	7


Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to tryptophan. The specificity is calculated in percent in relation to the tryptophan binding activity.

5-HTP (5-Hydroxy-Tryptophan)	< 0.5 %
L-Phenylalanine	< 0.1 %
L-Tyrosine	< 0.1 %

Correlation with HPLC-MS

10 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.98$.

12. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.

- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of sulfuric acid, which is a strong acid. Although diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control Samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure which is not coordinated with the producer may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *AM J Med Genet B Neuropsychiatr Genet.* 2009 Mar 5; 150B (2): 202-8.

Richard DM, Dawes MA, Mathias CW, Acheson A, Hill-Kapturczak N, Dougherty DM: L-Tryptophan: Basic metabolic functions, behavioural research and therapeutic indications. *International Journal of Tryptophan Research* 2009; 2 45-60.

Demisch K, Bauer J, Georgi K, Demisch L. Treatment of severe chronic insomnia with L-tryptophan: Results of a double blind cross-over study. *Pharmacopsychiatry.* 1987 Nov; 20 (6): 242-4.

Komer E, Bertha G, Flooh E, et al. Sleep inducing effect of L-tryptophan. *Er Neurol.* 1986 ; 25 (Suppl 2): 107-13.

Gendall KA, Joyce PR. Meal induced changes in tryptophan:LNAAs ratio: effects on craving and binge eating. *Eat Beh* 2000; 1 (1): 53-62.

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; 8 (3): 289-95.

Breum L, Rasmussen MH, Hilsted J, Femstrom JD. Twenty-four-hour plasma tryptophan concentrations and ratio are below normal in obese subjects and are not normalised by substantial weight reduction. *Am J Clin Nutr.* 2003; 77 (5): 1112-8.

Porter RJ, Mulder RT, Joyce PR, Miller AL, Kennedy M. Tryptophan hydrolase gene (TPH1) and peripheral tryptophan level in depression. *J Affect Disord.* 1;109(1-2):209-12.

Almeida-Montes LG, Valles-Sanches V, Moreno-Aguilar J et al. Relation of serum cholesterol, lipid, serotonin and tryptophan levels to severity of depression and suicide attempts. *J Psychiatry Neurosci.* 2000 Sep; 25 (4): 371-7.

Fitzgerald P et al: tryptophan catabolism in females with irritable bowel syndrome: relationship to interferon gamma, severity of symptoms and psychiatric co-morbidity; *Neurogastroenterol Motil* 2008; 20: 1291-7

Hideki Miura et al: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress*, Volume 11, Issue 3 2008, 198-209.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number