# TNF- $\alpha$ ELISA Kit

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
*For the in vitro determination of TNF- $\alpha$  in serum, plasma, stool and cell culture supernatant*

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---

<b>1. INTENDED USE</b>	<b>16</b>
<b>2. INTRODUCTION</b>	<b>16</b>
<b>3. TEST PRINCIPLE</b>	<b>16</b>
<b>4. MATERIAL SUPPLIED</b>	<b>17</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>17</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>18</b>
<b>7. PRECAUTIONS</b>	<b>19</b>
<b>8. SAMPLE PREPARATION</b>	<b>20</b>
<b>PLEASE NOTE</b>	<b>21</b>
<b>9. ASSAY PROCEDURE</b>	<b>21</b>
PROCEDURAL NOTES	21
TEST PROCEDURE	22
<b>10. RESULTS</b>	<b>23</b>
<b>11. LIMITATIONS</b>	<b>24</b>
<b>12. QUALITY CONTROL</b>	<b>24</b>
EXPECTED VALUES	24
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>24</b>
PRECISION AND REPRODUCIBILITY	24
SENSITIVITY	25
CROSS REACTIVITY	25
LINEARITY	25
<b>14. REFERENCES</b>	<b>25</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>26</b>

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**14. LITERATUR**

1. Cerami A, Beutler B (1988). "The history, properties, and biological effects of cachectin". *Biochemistry* 27 (20): 7575–7582.
2. Chen G, Goeddel DV (2002). "TNF-R1 signaling: a beautiful pathway". *Science* 296 (5573): 1634–5.
3. Locksley RM, Killeen N, Lenardo MJ (2001). "The TNF and TNF receptor superfamilies: integrating mammalian biology". *Cell* 104 (4): 487–501.

4. Swardfager W, Lanctôt K, Rothenburg L, Wong A, Cappell J, Herrmann N (2010). "A meta-analysis of cytokines in Alzheimer's disease". *Biol Psychiatry* 68 (10): 930-941.
5. Vilcek J, Lee TH (1991). "Tumor necrosis factor. New insights into the molecular mechanisms of its multiple actions". *J. Biol. Chem.* 266 (12): 7313–7316.
6. Wajant H, Pfizenmaier K, Scheurich P (2003). "Tumor necrosis factor signaling". *Cell Death Differ.* 10 (1): 45–65.

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Hersteller

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# TNF- $\alpha$ ELISA Kit

*For the in vitro determination of TNF- $\alpha$  in serum, plasma, stool  
and cell culture supernatant*


Valid from 23.11.2012

**K 9610**

96

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim  
Tel.: ++49 6251 70190-0  
Fax: ++ 49 6251 849430  
e.mail: [Info@immundiagnostik.com](mailto:Info@immundiagnostik.com)  
[www.immundiagnostik.com](http://www.immundiagnostik.com)

## 1. INTENDED USE

The *Immundiagnostik Assay* is intended for the quantitative determination of TNF- $\alpha$  in plasma, serum, stool and cell culture supernatant. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Tumor necrosis factor-alpha (TNF- $\alpha$ ) is a cytokine involved in systemic inflammation. The primary role of TNF- $\alpha$  is in the regulation of immune cells. TNF- $\alpha$  stimulates the acute phase reaction, induces apoptotic cell death, cellular proliferation and differentiation, inhibits tumor genesis and viral replication. Dysregulation of TNF- $\alpha$  production has been implicated in a variety of human diseases like cancer and Alzheimer.

TNF- $\alpha$  is secreted by macrophages, monocytes, neutrophils, T-cells as well as NK-cells (natural killer cells) following their stimulation by bacterial lipopolysaccharides. Human TNF- $\alpha$  is a non-glycosylated protein with a molecular weight of 17.5 kDa and a length of 157 amino acids. TNF- $\alpha$  shows a wide spectrum of biological activities. It causes cytolysis and/or cytostasis of many tumor cell lines *in vitro*. Within hours after injection, TNF- $\alpha$  leads to destruction of small blood vessels within malignant tumors. TNF- $\alpha$  enhances phagocytosis and cytotoxicity in neutrophilic granulocytes and modulates the expression of many other proteins.

Elevated TNF- $\alpha$  serum levels are found in patients suffering from Crohn's disease, ulcerating colitis or rheumatoid arthritis.

## 3. TEST PRINCIPLE

This Enzyme-Linked-Immunosorbent Assay (ELISA) allows the quantitative determination of TNF- $\alpha$ . In a first incubation step, TNF- $\alpha$  is bound to monoclonal mouse antibodies against TNF- $\alpha$ , which are immobilized on the surface of the microtiter plate. After a washing step, to remove all interfering substances, a second monoclonal anti-TNF- $\alpha$  antibody and a horseradish peroxidase labeled conjugate are added. The amount of the converted substrate by the peroxidase is directly proportional to the amount of bound TNF- $\alpha$  and can be determined photometrically at 450 nm.

#### 4. MATERIAL SUPPLIED

Cat. No.	Content	Kit Components	Quantity
K 9610MTP	PLATE	One holder with coated strips	12 x 8 wells
K 9610WP	WASHBUF	ELISA wash buffer concentrate, 10x	100 ml
K 9610A2	AB	2 <sup>nd</sup> Antibody (Mouse anti-hTNF- $\alpha$ , biotinylated)	1 x 150 $\mu$ l
K 9610K	CONJ	Conjugate (Streptavidin, Peroxidase- labeled)	1 x 200 $\mu$ l
K 9610ST	STDKONZ	Standard concentrate, lyophilized	3 x 1 vial
K 9610KO1	CTRL	Control, lyophilized	3 x 1 vial
K 9610KO2	CTRL	Control, lyophilized	3 x 1 vial
K 9610PV	STDBUF	Sample dilution buffer, ready-to-use	1 x 25 ml
K 9610TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9610AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm

\*Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2  $\mu$ m) with an electrical conductivity < 0.055  $\mu$ S/cm at 25°C ( $\geq$ 18.2 M $\Omega$  cm).

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that the reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to three times within the expiry date stated on the label.
- Reagents with a volume less than **100  $\mu$ l** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solution. The crystals must be re-dissolved in a water bath at 37°C before dilution. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The **AB** (2<sup>nd</sup> antibody) is diluted **1:101** in wash buffer (100  $\mu$ l **AB** + 10 ml wash buffer). Undiluted antibody is stable at 2-8°C until the expiry date given on the label. **Diluted antibody solution is not stable and cannot be stored.**
- The **CONJ** (conjugate) must be diluted **1:101** in wash buffer (100  $\mu$ l **CONJ** + 10 ml wash buffer). The undiluted conjugate is stable at **2-8°C** until the expiry date stated on the label. **Diluted conjugate is not stable and cannot be stored.**
- The lyophilized **CTRL** (controls) are stable at **2-8°C** until the expiry date stated on the label. The **CTRL** must be reconstituted with **500  $\mu$ l** of **STDBUF** (standard dilution buffer). Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted controls are not stable and cannot be stored.**
- The lyophilized **standard concentrate** (STDKONZ) must be reconstituted in **standard dilution buffer** (STDBUF). To prepare a calibration curve, dilute the reconstituted standard concentrate in **1:2** dilution steps. Reconstitution volume, dilution scheme and the resulting concentrations are described in the product specification data sheet.

The lyophilized standard concentrate (STDKONZ) is stable up to the date stated on the label when stored closed at **2–8°C**. **The reconstituted standard concentrate is not stable and cannot be stored.**

- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 8. SAMPLE PREPARATION

### Stool samples

#### *Extraction of the stool sample*

##### **1a. Stool Sample Application System (SAS)** (Cat. No.: K 6998SAS)

###### ***Stool sample tube – Instruction for use***

Please note that the dilution factor of the final stool suspension depends on the used amount of stool sample and the volume of the buffer.

###### **SAS with 0.75 ml Buffer:**

Applied amount of stool:	15 mg
Buffer Volume:	0.75 ml
Dilution Factor:	1:50

###### **Diluted wash buffer is used as extraction buffer.**

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For remarkably inhomogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) **Fill the empty sample tube with 0.75 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick exhibits notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off and leave 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.
- e) Allow sample to stand for app. 10 minutes until sediment has settled down. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the turquoise ring plus the dipstick. Discard cap and dipstick. Make sure, the sediment will not be dispersed again.

**1b. Sample preparation kit from Roche Diagnostics, Mannheim, Germany (Cat. No. 10 745 804 322)**

Alternatively, other stool sample preparation kits (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany) can be used. In the Roche sample preparation kit, 100 mg of stool sample are suspended in 5 ml of extraction buffer using a vibrator mixer (e.g. Vortex mixer). Centrifugation of the suspension is recommended.

**Dilution Factor (1a. or 1b.):** **1:50**

**PLEASE NOTE**

For TNF- $\alpha$  measurement in stool, please freeze the stool samples at -20°C immediately after collection (not later than two hours after collection). If samples have to be transported, please make sure to keep them frozen. The stability of the samples must be checked before using the method in routine.

**Serum, plasma**

Serum/Plasma can be used without further dilution. Samples must be stored at -20°C.

**9. ASSAY PROCEDURE***Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not to assemble wells of different microtiter plates for analysis, even if they are of the same batch.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. The Immundiagnostik AG can therefore not be held responsible for any damage.
- Carry out the assay with the actual manual delivered with the kit.

## Test procedure

Bring all reagents to room temperature (18-26°C) and mix well before use. Avoid direct sunlight during all incubation steps. Cover the microtiter plate during the different incubation steps. Carry out the tests in duplicate.

Wash the coated microtiter plate **5 x with 250  $\mu$ l** ELISA wash buffer before use. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess liquid.

1. Add **100  $\mu$ l** of **STD** (Standards), **SAMPLE** (Sample/s) or **CTRL** (Controls) into respective well.
2. Incubate for **2 hours** at room temperature (18-26°C) on a horizontal mixer.
3. Decant the contents of the plate and wash the cavities **5 x with 250  $\mu$ l** washing buffer solution. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess liquid.
4. Add **100  $\mu$ l** of diluted **AB** (2. antibody).
5. Incubate for **1 hour** at room temperature (18-26°C) on a horizontal mixer.
6. Decant the contents of the plate and wash the cavities **5 x with 250  $\mu$ l** washing buffer solution. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess liquid.
7. Add **100  $\mu$ l** of diluted **CONJ** (conjugate).
8. Incubate for **1 hour** at room temperature (18-26°C) on a horizontal mixer.
9. Decant the contents of the plate and wash the cavities **5 x with 250  $\mu$ l** washing buffer solution. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess liquid.
10. Add **100  $\mu$ l** of **SUB** (TMB substrate).
11. Incubate for **10-20 minutes** at room temperature (18-26°C).
12. Add **100  $\mu$ l STOP** (stop solution) and mix shortly.

13. Determine the **absorption** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

## 10. RESULTS

We recommend to use the 4-parameter-algorithm to calculate the results.

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

Alternatively, following algorithms can be used:

2. Point-to-point-calculation

We recommend for the optical density a linear ordinate and for the concentration a linear abscissa.

3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a manual control of the paired values should be made.

### Stool

To obtain the TNF- $\alpha$  concentration of a sample, multiply the estimated value by the dilution factor of **50**.

**Note:** If the sample is further diluted, this dilution must be considered by the calculation of the TNF- $\alpha$  concentration.

## 11. LIMITATIONS

Samples with TNF- $\alpha$  levels greater than the highest calibrator, should be diluted and re-assayed. This dilution must be considered by the calculation of the TNF- $\alpha$  concentration.

## 12. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### **Reference value for TNF- $\alpha$**

**Plasma (apparently healthy persons, n = 40):** < 20 pg/ml

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay**

The precision of the TNF- $\alpha$  ELISA test was calculated from 20 replicate determinations of one sample.

n= 20

Sample	TNF- $\alpha$ mean [pg/ml]	CV [%]
1	155	6.3

#### **Inter-Assay**

The reproducibility of the TNF- $\alpha$  ELISA test was calculated from data obtained by several persons measuring a normal range sample on different days.

n= 20

Sample	TNF- $\alpha$ mean [pg/ml]	CV [%]
1	142	8.2

### Sensitivity

The detection limit of TNF- $\alpha$  ELISA is calculated as  $B_0 + 2SD = 10$  pg/ml.

### Cross reactivity

GM-CSF (Granulocyte-macrophage colony-stimulating factor)	< 10%
M-CSF (Macrophage colony-stimulating factor)	< 10%
G-CSF (Granulocyte colony-stimulating factor)	< 10%

### Linearity

Linearity of the test was evaluated by dilution of TNF- $\alpha$ -containing material. Linearity was obtained in the range of 10 - 500 pg/ml.

## 14. REFERENCES

1. Cerami A, Beutler B (1988). "The history, properties, and biological effects of cachectin". *Biochemistry* 27 (20): 7575–7582.
2. Chen G, Goeddel DV (2002). "TNF-R1 signaling: a beautiful pathway". *Science* 296 (5573): 1634–5.
3. Locksley RM, Killeen N, Lenardo MJ (2001). "The TNF and TNF receptor superfamilies: integrating mammalian biology". *Cell* 104 (4): 487–501.
4. Swardfager W, Lanctôt K, Rothenburg L, Wong A, Cappell J, Herrmann N (2010). "A meta-analysis of cytokines in Alzheimer's disease". *Biol Psychiatry* 68 (10): 930-941.
5. Vilcek J, Lee TH (1991). "Tumor necrosis factor. New insights into the molecular mechanisms of its multiple actions". *J. Biol. Chem.* 266 (12): 7313–7316.

6. Wajant H, Pfizenmaier K, Scheurich P (2003). "Tumor necrosis factor signaling". *Cell Death Differ.* 10 (1): 45–65.

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number