Vitamin-D-Bindungsprotein ELISA Kit

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15. LITERATUR

Allgemeine Literatur

1. Bouillon, R., van Baelen, H. & de Moor, P., 1977. The measurement of the vitamin D-binding protein in human serum. *The Journal of clinical endocrinology and metabolism*, **45**(2), pp.225–31.
2. Fu, L. et al., 2009. Common genetic variants of the vitamin D binding protein (DBP) predict differences in response of serum 25-hydroxyvitamin D [25(OH)D] to vitamin D supplementation. *Clinical biochemistry*, **42**(10-11), pp.1174–7.
3. Haughton, M.A. & Mason, R.S., 1992. Immunonephelometric assay of vitamin D-binding protein. *Clinical chemistry*, **38**(9), pp.1796–801.
4. Heijboer, A.C. et al., 2012. Accuracy of 6 Routine 25-Hydroxyvitamin D Assays: Influence of Vitamin D Binding Protein Concentration. *Clinical Chemistry*, **58**(3), pp.543–548.
5. Jørgensen, C.S. et al., 2004. Gc globulin (vitamin D-binding protein) levels: an inhibition ELISA assay for determination of the total concentration of Gc globulin in plasma and serum. *Scandinavian journal of clinical and laboratory investigation*, **64**(2), pp.157–66.
6. Malik, S. et al., 2013. Common variants of the vitamin D binding protein gene and adverse health outcomes. *Critical reviews in clinical laboratory sciences*.
7. Schmidt-Gayk, H. et al., 1977. 25-hydroxy-vitamin-D in nephrotic syndrome. *Lancet*, **2**(8029), pp.105–8.
8. Thomas, L., 1982. Proteindiagnostik: Diagnose, Therapiekontrolle. 1st ed., Frankfurt am Main: Behringwerke, *Medizinische Information und Vertrieb*.

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9. Cauley, J.A. et al., 2011. Serum 25-hydroxyvitamin D and clinical fracture risk in a multiethnic cohort of women: the Women's Health Initiative (WHI). *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research*, **26**(10), pp.2378–88.
10. Comabella, M. et al., 2010. Cerebrospinal fluid chitinase 3-like 1 levels are associated with conversion to multiple sclerosis. *Brain : a journal of neurology*, **133**(Pt 4),

pp.1082–93.

11. Correale, J., Ysraelit, M.C. & Gaitán, M.I., 2010. Gender differences in 1,25 dihydroxyvitamin D3 immunomodulatory effects in multiple sclerosis patients and healthy subjects. *Journal of immunology (Baltimore, Md. : 1950)*, **185**(8), pp.4948–58.
12. Doorenbos, C.R.C. et al., 2012. Antiproteinuric treatment reduces urinary loss of vitamin D-binding protein but does not affect vitamin D status in patients with chronic kidney disease. *The Journal of steroid biochemistry and molecular biology*, **128**(1-2), pp.56–61.
13. Gressner, O. et al., 2007. Gc-globulin concentrations and C5 haplotype-tagging polymorphisms contribute to variations in serum activity of complement factor C5. *Clinical biochemistry*, **40**(11), pp.771–5.
14. Jeffery, L.E. et al., 2012. Availability of 25-hydroxyvitamin D(3) to APCs controls the balance between regulatory and inflammatory T cell responses. *Journal of immunology (Baltimore, Md. : 1950)*, **189**(11), pp.5155–64.
15. Jeng, L. et al., 2009. Alterations in vitamin D status and anti-microbial peptide levels in patients in the intensive care unit with sepsis. *Journal of translational medicine*, **7**, p.28.
16. Lee, S.-H. et al., 2011. Relationship between group-specific component protein and the development of asthma. *American journal of respiratory and critical care medicine*, **184**(5), pp.528–36.
17. Meguro, S. et al., 2011. Plasma 25-hydroxyvitamin d is independently associated with hemoglobin concentration in male subjects with type 2 diabetes mellitus. *International journal of endocrinology*, 2011, p.362981.
18. Mirković, K. et al., 2013. Urinary vitamin D binding protein: a potential novel marker of renal interstitial inflammation and fibrosis. *PloS one*, **8**(2), p.e55887.
19. Prytuła, A. et al., 2012. Urinary and dialysate losses of vitamin D-binding protein in children on chronic peritoneal dialysis. *Pediatric nephrology (Berlin, Germany)*, **27**(4), pp.643–9.
20. Turner, A.M. et al., 2013. Circulating DBP level and prognosis in operated lung cancer: an exploration of pathophysiology. *The European respiratory journal*, **41**(2), pp.410–6.
21. Wagner, D. et al., 2013. Randomized clinical trial of vitamin D3 doses on prostatic vitamin D metabolite levels and ki67 labeling in prostate cancer patients. *The Journal of clinical endocrinology and metabolism*, **98**(4), pp.1498–507.
22. Wang, X. et al., 2013. Vitamin D-binding protein levels in female patients with primary hyperparathyroidism. *Endocrine practice : official journal of the American*

College of Endocrinology and the American Association of Clinical Endocrinologists, 19(4), pp.609–13.

23. Wood, A.M. et al., 2011. Vitamin D-binding protein contributes to COPD by activation of alveolar macrophages. *Thorax, 66(3), pp.205–10.*

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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
Vitamin D binding protein ELISA kit

*For the in vitro determination of vitamin D binding protein
in serum, plasma and urine*

Valid from 25.06.2015

REF K 2314

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	17
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. STORAGE AND PREPARATION OF REAGENTS	18
6. STORAGE AND PREPARATION OF SAMPLES	19
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Test procedure</i>	20
8. RESULTS	21
9. LIMITATIONS	22
10. QUALITY CONTROL	22
<i>Reference range</i>	22
<i>Additional reference ranges</i>	23
11. PERFORMANCE CHARACTERISTICS	23
<i>Spiking Recovery</i>	23
<i>Dilution recovery</i>	24
<i>Precision and reproducibility</i>	24
<i>Analytical Sensitivity</i>	25
12. PRECAUTIONS	25
13. TECHNICAL HINTS	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	26
<i>General literature</i>	26
<i>Literature using the Immundiagnostik VDBP ELISA</i>	27

1. INTENDED USE

This Immundiagnostik assay is an enzyme immunoassay intended for the quantitative determination of free and not actin complex bound vitamin D binding protein in serum, plasma and urine. For *in vitro* diagnostic use only.

2. INTRODUCTION

Vitamin D binding protein (VDBP, also known as group-specific component / Gc protein) is a multifunctional serum protein which is formed in the liver. High estrogen levels, caused e. g. by pregnancy or hormonal birth control, stimulate its synthesis. VDBP can be found in plasma, ascites, liquor, urine and at the surface of various cell types.

In the blood, VDBP binds the bigger part of the circulating 25-OH vitamin D and brings it to the kidneys, where it is transformed into the hormone 1,25-(OH)₂ vitamin D.

In addition, VDBP binds monomeric actin at the ratio of 1:1. Actin is an intracellular protein which is available as monomer or filament. Massive tissue destruction or cell death cause the plasma level of actin to rise significantly, whereupon VDBP-actin complexes are formed which are removed quickly. Hence the VDBP level of trauma or sepsis patients decreases quickly, especially in patients with a high level of multi-organ failure.

Furthermore VDBP is a precursor of the immunomodulating protein Gc-MAF (Gc protein-derived macrophage activating factor) which increases the activity of macrophages against tumours.

Indications

- Prognostic factor for trauma patients
- Nephritic syndrome
- Tumours

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 2314	PLATE	Holder with precoated strips	12 x 8 wells
K 2314	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 2314	CONJ	POD antibody, (rabbit-anti-VDBP, peroxidase-labeled), pre-diluted	1 x 200 µl

Cat. No.	Label	Kit components	Quantity
K 2314	STD	Calibrators, lyophilized (60; 20; 6,6; 2,2; 0 ng/ml)	4 x 5 vials
K 2314	STDBUF	Standard dilution buffer	1 x 20 ml
K 2314	CTRL 1	Control, lyophilized	4 vials
K 2314	CTRL 2	Control, lyophilized	4 vials
K 2314	SAMPLEBUF	Dilution buffer, ready to use	2x 100 ml
K 2314	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 2314	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.

- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- The **lyophilized standards (STD)** and **controls (CTRL)** are stable at **2–8 °C** until the expiry date stated on the label. Before use, the standards and controls must be reconstituted with **500 µl** of **STDBUF** (standard dilution buffer). Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Reconstituted standards and controls are not stable and cannot be stored**.
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The CONJ is stable at **2–8 °C** until expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Serum and plasma samples

Dilute all plasma and serum samples 1:40 000 with SAMPLEBUF (sample dilution buffer). For example:

20 µl sample + **980 µl** SAMPLEBUF, mix well = **1:50 (dilution I)**

20 µl dilution I + **980 µl** SAMPLEBUF, mix well = **1:50 (dilution II)**

20 µl dilution II + **300 µl** SAMPLEBUF, mix well = **1:16 (dilution III)**

This results in a final dilution of 1:40 000.

For analysis, pipet **100 µl** of **dilution III** per well.

Other sample collectives should be diluted according to the expected VDBP concentration.

Urine

Urine samples have to be diluted 1:10 with SAMPLEBUF (sample dilution buffer). For example:

100 µl sample + **900 µl** SAMPLEBUF, mix well = **1:10**

Other sample collectives should be diluted according to the expected VDBP concentration.

7. ASSAY PROCEDURE

Principle of the test

This enzyme immunoassay is a sandwich assay for VDBP determination in serum, plasma and urine samples. The wells of the micro titer plate are coated with polyclonal anti-VDBP antibodies. In a first incubation step, the VDBP in the samples is bound to the coated polyclonal rabbit antibodies (in excess). To remove all unbound substances, a washing step is carried out. In a second incubation step, a polyclonal peroxidase-labeled rabbit anti-VDBP antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine. An acidic stopping solution is then added. The color converts to yellow. The intensity of the yellow color is directly proportional to the VDBP concentration in the sample. A dose response curve of the absorbance (at 450 nm) unit vs. concentration is generated.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Wash the precoated PLATE (microtiter plate) **5 times** with **250 µl** wash buffer per well. After the final washing step, remove residual buffer by tapping the plate on absorbent paper.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

Carry out the tests in duplicate.

1.	Add 100 µl STD (standard), CTRL (controls) and pre-diluted sample into the respective wells.
2.	Incubate for 1 hour shaking on a horizontal mixer at room temperature (15–30 °C).
3.	Discard the contents of the plate. Wash 5 times with 250 µl wash buffer per well. After the final washing step, remove residual buffer by tapping the plate on absorbent paper.

4.	Add 100 µl conjugate (peroxidase-labeled antibody).
5.	Incubate for 1 hour shaking on a horizontal mixer at room temperature (15–30 °C).
6.	Discard the contents of the plate. Wash 5 times with 250 µl wash buffer per well. After the final washing step, remove residual buffer by tapping the plate on absorbent paper.
7.	Add 100 µl SUB (TMB substrate).
8.	Incubate for 10–20 minutes at room temperature*.
9.	Add 100 µl of STOP (stop solution) and mix shortly.
10.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the

duplicate values should be evaluated manually.

Serum/plasma samples

The obtained VDBP levels of plasma and serum samples have to be multiplied with the dilution factor of 40 000.

Urine samples

The obtained VDBP levels of urine samples have to be multiplied with the dilution factor of 10.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

Analytical sensitivity × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Plasma / serum samples

200–550 mg/l

L. Thomas, 1982

You can find further reference ranges in the following publications in the references section: Bouillon (1977), Haughton (1992), Jorgensen (2004), Heijboer (2012).

We recommend each laboratory to establish its own reference range.

Urine samples

Doorenbos et al. analysed the urinary loss of VDBP in 24 h urine samples. In the healthy control group, they reported a value of 64 µg/24 h (23–111 µg/24 h).

Mirkovic et al. normalised the measured VDBP levels for albuminuria. Normalalbuminuria was reported as 0,44 mg/mg albumin (0,22–0,77 mg/mg albumin).

We recommend each laboratory to establish its own reference range.

Additional reference ranges

For some patient groups, other reference ranges of serum / plasma samples have been reported.

Pregnant women

Samples of pregnant women were measured to have a 30–80% higher reference range than the control groups (see publications cited for the serum reference range).

Liver diseases

According to Haughton et al., the reference range of patients with liver diseases is 35% lower than the one of healthy controls.

11. PERFORMANCE CHARACTERISTICS

Spiking Recovery

Two samples were spiked with VDBP calibrator and measured with this assay. (n = 2).

Sample	Unspiked Sample [ng/ml]	Spike [ng/ml]	VDBP expected [ng/ml]	VDBP measured [ng/ml]
A	2.2	5	7.2	7.7
	2.2	10	12.2	12.7
	2.2	20	22.2	25.3
B	6.7	2.5	9.2	8.7
	6.7	7.5	14.2	13.1
	6.7	15	21.7	22.5

Dilution recovery

Two patient samples were diluted with sample dilution buffer. The results are shown below (n = 2):

Sample	Dilution	VDBP expected [ng/ml]	VDBP measured [ng/ml]
A	1:5000	46.2	46.2
	1:10000	23.1	23.3
	1:20000	11.5	10.4
	1:40000	5.7	5.8
	1:80000	2.8	2.7
B	1:5000	38	38
	1:10000	19	20.2
	1:20000	9.5	8.7
	1:40000	4.7	4.3
	1:80000	2.3	2.4

Precision and reproducibility

Intra-Assay

The precision (intra-assay variation) of the Immundiagnostik VDBP ELISA test was calculated from 16 replicate determinations on each of one samples. (n = 16)

Sample	VDBP mean value [mg/dl]	CV [%]
1	24.2	5.0
2	42.9	3.2

Inter-Assay

The total precision (inter-assay variation) of the Immundiagnostik VDBP ELISA test was calculated from data on one sample obtained in 14 different assays by three technicians on two different lots of reagents over a period of three months.

Sample	VDBP mean value [mg/dl]	CV [%]
1	19.3	12.7

Analytical Sensitivity

The Zero-standard was measured 20 times. The detection limit was set as $B_0 + 2 SD$ and estimated to be 1,23 ng/ml.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

General literature

1. Bouillon, R., van Baelen, H. & de Moor, P., 1977. The measurement of the vitamin D-binding protein in human serum. *The Journal of clinical endocrinology and metabolism*, **45**(2), pp.225–31.
2. Fu, L. et al., 2009. Common genetic variants of the vitamin D binding protein (DBP) predict differences in response of serum 25-hydroxyvitamin D [25(OH)D] to vitamin D supplementation. *Clinical biochemistry*, **42**(10-11), pp.1174–7.
3. Haughton, M.A. & Mason, R.S., 1992. Immunonephelometric assay of vitamin D-binding protein. *Clinical chemistry*, **38**(9), pp.1796–801.
4. Heijboer, A.C. et al., 2012. Accuracy of 6 Routine 25-Hydroxyvitamin D Assays: Influence of Vitamin D Binding Protein Concentration. *Clinical Chemistry*, **58**(3), pp.543–548.
5. Jørgensen, C.S. et al., 2004. Gc globulin (vitamin D-binding protein) levels: an inhibition ELISA assay for determination of the total concentration of Gc globulin in plasma and serum. *Scandinavian journal of clinical and laboratory investigation*, **64**(2), pp.157–66.
6. Malik, S. et al., 2013. Common variants of the vitamin D binding protein gene and adverse health outcomes. *Critical reviews in clinical laboratory sciences*.
7. Schmidt-Gayk, H. et al., 1977. 25-hydroxy-vitamin-D in nephrotic syndrome. *Lancet*, **2**(8029), pp.105–8.

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Literature using the Immundiagnostik VDBP ELISA

9. Cauley, J.A. et al., 2011. Serum 25-hydroxyvitamin D and clinical fracture risk in a multiethnic cohort of women: the Women's Health Initiative (WHI). *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research*, **26**(10), pp.2378–88.
10. Comabella, M. et al., 2010. Cerebrospinal fluid chitinase 3-like 1 levels are associated with conversion to multiple sclerosis. *Brain : a journal of neurology*, **133**(Pt 4), pp.1082–93.
11. Correale, J., Ysraelit, M.C. & Gaitán, M.I., 2010. Gender differences in 1,25 dihydroxyvitamin D3 immunomodulatory effects in multiple sclerosis patients and healthy subjects. *Journal of immunology (Baltimore, Md. : 1950)*, **185**(8), pp.4948–58.
12. Doorenbos, C.R.C. et al., 2012. Antiproteinuric treatment reduces urinary loss of vitamin D-binding protein but does not affect vitamin D status in patients with chronic kidney disease. *The Journal of steroid biochemistry and molecular biology*, **128**(1-2), pp.56–61.
13. Gressner, O. et al., 2007. Gc-globulin concentrations and C5 haplotype-tagging polymorphisms contribute to variations in serum activity of complement factor C5. *Clinical biochemistry*, **40**(11), pp.771–5.
14. Jeffery, L.E. et al., 2012. Availability of 25-hydroxyvitamin D(3) to APCs controls the balance between regulatory and inflammatory T cell responses. *Journal of immunology (Baltimore, Md. : 1950)*, **189**(11), pp.5155–64.
15. Jeng, L. et al., 2009. Alterations in vitamin D status and anti-microbial peptide levels in patients in the intensive care unit with sepsis. *Journal of translational medicine*, **7**, p.28.
16. Lee, S.-H. et al., 2011. Relationship between group-specific component protein and the development of asthma. *American journal of respiratory and critical care medicine*, **184**(5), pp.528–36.
17. Meguro, S. et al., 2011. Plasma 25-hydroxyvitamin d is independently associated with hemoglobin concentration in male subjects with type 2 diabetes mellitus. *International journal of endocrinology*, 2011, p.362981.
18. Mirković, K. et al., 2013. Urinary vitamin D binding protein: a potential novel marker of renal interstitial inflammation and fibrosis. *PLoS one*, **8**(2), p.e55887.
19. Prytuła, A. et al., 2012. Urinary and dialysate losses of vitamin D-binding protein

- in children on chronic peritoneal dialysis. *Pediatric nephrology (Berlin, Germany)*, **27**(4), pp.643–9.
20. Turner, A.M. et al., 2013. Circulating DBP level and prognosis in operated lung cancer: an exploration of pathophysiology. *The European respiratory journal*, **41**(2), pp.410–6.
21. Wagner, D. et al., 2013. Randomized clinical trial of vitamin D3 doses on prostatic vitamin D metabolite levels and ki67 labeling in prostate cancer patients. *The Journal of clinical endocrinology and metabolism*, **98**(4), pp.1498–507.
22. Wang, X. et al., 2013. Vitamin D-binding protein levels in female patients with primary hyperparathyroidism. *Endocrine practice : official journal of the American College of Endocrinology and the American Association of Clinical Endocrinologists*, **19**(4), pp.609–13.
23. Wood, A.M. et al., 2011. Vitamin D-binding protein contributes to COPD by activation of alveolar macrophages. *Thorax*, **66**(3), pp.205–10.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by