WISP2 / CCN5 ELISA Kit

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Additional information about our products is available on our homepage

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Table of content	Page
1. INTENDED USE	17
2. INDICATIONS	17
3. INTRODUCTION	17
4. PRINCIPLE OF THE TEST	17
5. MATERIAL SUPPLIED	18
6. MATERIAL REQUIRED BUT NOT SUPPLIED	18
7. PREPARATION AND STORAGE OF REAGENTS	19
8. PREPARATION OF STANDARDS	20
9. SAMPLE PREPARATION	20
10. PRECAUTIONS	20
<i>PROCEDURAL NOTES</i>	20
<i>TEST PROCEDURE</i>	21
11. RESULTS	22
12. LIMITATIONS	22
13. QUALITY CONTROL	22
14. PERFORMANCE CHARACTERISTICS	23
<i>PRECISION AND REPRODUCIBILITY</i>	23
<i>SENSITIVITY</i>	24
<i>CROSS REACTIVITY</i>	24
15. REFERENCES	24
16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	25

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15. LITERATUR

Banerjee, S. u. a., 2008. CCN5/WISP-2 expression in breast adenocarcinoma is associated with less frequent progression of the disease and suppresses the invasive phenotypes of tumor cells. *Cancer Research*, 68(18), S.7606–7612.

Banerjee, S. u. a., 2003. WISP-2 Gene in Human Breast Cancer: Estrogen and Progesterone Inducible Expression and Regulation of Tumor Cell Proliferation. *Neoplasia (New York, N.Y.)*, 5(1), S.63–73.

Davies, S.R. u. a., 2010. Differential expression of the CCN family member WISP-1, WISP-2 and WISP-3 in human colorectal cancer and the prognostic implications. *International Journal of Oncology*, 36(5), S.1129–1136.

Dhar, G. u. a., 2007. Loss of WISP-2/CCN5 signaling in human pancreatic cancer: a potential mechanism for epithelial-mesenchymal-transition. *Cancer Letters*, 254(1), S.63–70.

- Fritah, A. u. a., 2008. Role of WISP-2/CCN5 in the Maintenance of a Differentiated and Noninvasive Phenotype in Human Breast Cancer Cells. *Molecular and Cellular Biology*, 28(3), S.1114–1123.
- Inadera, H. u. a., 2000. WISP-2 as a novel estrogen-responsive gene in human breast cancer cells. *Biochemical and Biophysical Research Communications*, 275(1), S.108–114.
- Jones, J.A. u. a., 2007. CCN5 expression in mammals. *Journal of Cell Communication and Signaling*, 1(2), S.127–143.
- Jun, J.-I. & Lau, L.F., 2011. Taking aim at the extracellular matrix: CCN proteins as emerging therapeutic targets. *Nature Reviews. Drug Discovery*, 10(12), S.945–963.
- Lake, A.C. u. a., 2003. CCN5 Is a Growth Arrest-Specific Gene That Regulates Smooth Muscle Cell Proliferation and Motility. *The American Journal of Pathology*, 162(1), S.219–231.
- Leask, A. & Abraham, D.J., 2006. All in the CCN family: essential matricellular signaling modulators emerge from the bunker. *Journal of cell science*, 119(Pt 23), S.4803–4810.
- Rittié, L. u. a., 2011. Spatial-temporal modulation of CCN proteins during wound healing in human skin in vivo. *Journal of Cell Communication and Signaling*, 5(1), S.69–80.
- Russo, J.W. & Castellot, John J., 2010. CCN5: biology and pathophysiology. *Journal of Cell Communication and Signaling*, 4(3), S.119–130.
- Sabbah, M. u. a., 2011. CCN5, a Novel Transcriptional Repressor of the Transforming Growth Factor β Signaling Pathway ∇ . *Molecular and Cellular Biology*, 31(7), S.1459–1469.
- Tanaka, I. u. a., 2005. Expression and regulation of WISP2 in rheumatoid arthritic synovium. *Biochemical and Biophysical Research Communications*, 334(4), S.973–978.
- Wiesman, K.C. u. a., 2010. CCN5, a secreted protein, localizes to the nucleus. *Journal of Cell Communication and Signaling*, 4(2), S.91–98.
- Yeger, H. & Perbal, B., 2007. The CCN family of genes: a perspective on CCN biology and therapeutic potential. *Journal of Cell Communication and Signaling*, 1(3), S.159–164.

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
WISP2/CCN5 ELISA Kit

For the determination of WISP2/CCN5 in serum, plasma and synovial liquid

Valid from 15.10.2012

K1015

1. INTENDED USE

The assay described here is intended for the quantitative determination of **WISP2/CCN5** in serum, plasma and synovial liquid.

2. INDICATIONS

Diseases of bone metabolism

3. INTRODUCTION

CCN (Cyr 61; CTGF; NOV) proteins are extracellular modulators of signal transduction between matrix associated proteins and molecules belong to the cell surface. WISP2/ CCN5 is participated in many physiological and pathophysiological processes. It can be measured in bone, muscle and dermis, also in the synovium of patients with rheumatoid arthritis (RA). At cancer e.g. breast, pancreas and colon cancer the expression level of WISP2/CCN5 is decreased when the tumor becomes more invasive. Additionally it could be shown (in vivo) that WISP2/CCN5 expression stopped the epithelial mesenchymal transition (EMT) and the migration of cancer cells.

4. PRINCIPLE OF THE TEST

This Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of **WISP2/CCN5** in serum, plasma and synovial liquid. The WISP2/CCN5 in the sample is bound to polyclonal antibodies (rabbit-anti-human WISP2/CCN5) coated on a microtiter plate. A wash step removes all unbound components of the sample. The detection of bound WISP2/CCN5 takes place with a biotinylated secondary antibody (monoclonal mouse-anti WISP2/CCN5). In a next step peroxidase-labeled Streptavidin is pipetted into the wells. The Peroxidase-Tetramethylbenzidine (TMB) reaction produces a blue color. After adding a stop solution, the color changes from blue to yellow. The intensity of the color reaction product is proportional to the concentration of the analyte in the samples, standards and controls. With the help of the standard curve the concentration of the sample can exactly be determined.

5. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K1015MTP	PLATE	One holder with coated strips	96 wells
K1015WP	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K1015K	CONJ	Conjugate	1 x 200 µl
K1015A2	2AB	2. Antibody	1 x 200 µl
K1015ST	STD	Standard, lyophilized	3 x 1 vial
K1015KO1	CTRL	Control, lyophilized	3 x 1 vial
K1015KO2	CTRL	Control, lyophilized	3 x 1 vial
K1015TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K1015SV	STDBUF	Standard Dilution buffer	1 x 100 ml
K1015PV	SAMPLEBUF	Sample Dilution buffer	1 x 100 ml
K1015AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

6. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra-pure water*
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm

*Immundiagnostik AG recommends the use of ultra-pure Water (water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity <0.055 µS/cm at 25°C (≥18.2 MΩ cm).

7. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 3 times within the expiry date stated on the label.
- The **ELISA WASHBUF** (wash buffer concentrate) has to be diluted with ultra-pure-water **1:10** before use (100 ml WASHBUF + 900 ml ultra-pure-water). Crystals could occur due to high salt concentration. The crystals must be re-suspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrate is stable at 2-8°C until the expiry date stated on the label. Diluted solutions can be stored at 2-8°C for 1 month.
- The **STD** (standard) and **CTRL** (controls) must be reconstituted with **500 µl ultra-pure-water**. Allow the vial content to dissolve for **2 minutes** and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standard and controls are not stable and cannot be used repeatedly.**
- The **2. AB** (secondary Antibody) must be diluted **1:100** in diluted wash buffer (e.g. 100 µl 2.AB + 990 µl wash buffer). Undiluted 2.AB is stable at 2-8°C until the expiry date stated on the label. **Diluted antibody solution is not stable and cannot be used repeatedly.**
- The **CONJ** (conjugate) must be diluted **1:100** in diluted wash buffer (e.g. 100 µl CONJ + 990 µl wash buffer). Undiluted CONJ is stable at 2-8°C until the expiry date stated on the label. **Diluted conjugate solution is not stable and cannot be used repeatedly.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2-8°C**.

8. PREPARATION OF STANDARDS

For the standard curve dilute the concentrate WISP2-STD (Standard, S7; concentration: 2400 ng/ml) with STDBUF (standard dilution buffer) **in 1:2 dilution steps** as described in the following example:

S7 **S7:** 2400 ng/ml

S6 = 250 µl S7 + 250 µl STDBUF **S6:** 1200 ng/ml

S5 = 250 µl S6 + 250 µl STDBUF **S5:** 600 ng/ml

S4 = 250 µl S5 + 250 µl STDBUF **S4:** 300 ng/ml

S3 = 250 µl S4 + 250 µl STDBUF **S3:** 150 ng/ml

S2 = 250 µl S3 + 250 µl STDBUF **S2:** 75 ng/ml

S1 **S1:** 0 ng/ml

SAMPLEBUF is used as zero standard **S1** (0 ng/ml).

9. SAMPLE PREPARATION

Serum and Plasma samples, synovial liquid

Samples are diluted **1:10** in **SAMPLEBUF** (sample dilution buffer), e.g. 30 µl Sample + 270 µl **SAMPLEBUF**.

Use **100 µl** of the final diluted sample per well.

The undiluted samples are stable up to 6 months when stored at -20°C and can be reanalyzed.

10. PRECAUTIONS

Procedural notes

- Do not mix different lot numbers of any kit component within the same assay.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- Carry out the assay with the actual manual delivered with the kit.

- The quality control guidelines should be observed.
- Carry out the determinations in the PLATE (micro titer plate) in duplicate.

TEST PROCEDURE

1. Wash the coated PLATE (microtiter plate) **5 x with 250 µl** ELISA wash buffer. After the last wash, tap the PLATE firmly on absorbent paper to remove excess liquid.
2. Add **100 µl of STD/CTRL/SAMPLE** (Standards/Controls/Sample) in duplicate into respective well.
3. Incubate for **2 hours**, shaking on a horizontal mixer, at room temperature (18-26°C).
4. Aspirate and wash the wells **5 x with 250 µl** diluted wash buffer. After the final washing step, tap firmly the microtiter plate on absorbent paper to remove excess liquid.
5. Add **100 µl** of diluted 2.AB in each well.
6. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature (18-26°C).
7. Aspirate and wash the wells **5 x with 250 µl** diluted wash buffer. After the final washing step, tap firmly the microtiter plate on absorbent paper to remove excess liquid.
8. Add **100 µl** of diluted CONJ (conjugate) in each well.
9. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature (18-26°C).
10. Aspirate and wash the wells **5 x with 250 µl** diluted wash buffer. After the final washing step, tap firmly the microtiter plate on absorbent paper to remove excess liquid.
11. Add **100 µl** of SUB (TMB substrate) in each well.
12. Incubate for **10-20 minutes** at room temperature (18-26°C).
13. Add **100 µl** of STOP (stop solution) in each well and mix shortly.
14. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

11. RESULTS

We recommend to use the 4-parameter-algorithm for evaluation of the results. Alternatively, the other two algorithms can be used.

1. **4-parameter-algorithm**

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. **Point-to-point-calculation**

We recommend for the optical density a linear ordinate and for the concentration a linear abscissa.

3. **Spline-algorithm**

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a manual control of the paired values should be made.

The estimated concentration must be multiplied by the dilution factor of the sample.

12. LIMITATIONS

Samples with a WISP2/CCN5 level greater than the highest calibrator, should be prepared in a higher dilution and re-assayed.

13. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected results

Normal range

The normal range was estimated based on sample values obtained from 200 apparently healthy persons.

WISP2/CCN5 in serum or plasma

Age 0-15	0 – 25 [ng/ml]
Age 15-20	0 – 50 [ng/ml]
Age 20-30	200 – 700 [ng/ml]
Age 30-40	400 – 900 [ng/ml]
Age > 50	800– 1800 [ng/ml]

Synovial liquid has a considerable higher WISP2 /CCN5 level.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) was calculated from 74 replicate measurements of two normal samples in the WISP2/CCN5 ELISA by one person.

Intra-Assay CV n= 74

Sample	WISP2 [ng/ml]	Intra-Assay CV [%]
1	1049	7
2	3517	4

The precision (inter-assay variation) was calculated from WISP2/CCN5 ELISA measurements of two normal samples on different days by two persons.

Inter-Assay CV n= 28

Sample	WISP2 [ng/ml]	Inter-Assay CV [%]
1	989	9
2	3554	7

Sensitivity

LoB Limit of Blank: 8,12 ng/ml

The limit of blank was determined as $B_0 + 1,645 SD_{\text{Blank}}$. Standard 0 (Sample buffer; sample without analyte) was measured 22 times in the WISP2 /CCN5 ELISA.

LoD Limit of Detection: 25,13 ng/ml

The lower limit of detection was determined as $LoB + 1,645 SD_{\text{low sample}}$. A sample with a low amount of analyte was measured 22 times in the WISP2/CCN5 ELISA.

Cross reactivity

There was insignificant cross-reactivity with the other proteins of the CCN-family (Cyr61; CTGF; NOV; WISP1; WISP3) and albumin.

Cross-Reactivity WISP2/CCN5 ELISA					
Albumin	CYR 61	CTGF	NOV	WISP1	WISP3
0.78%	0.01%	1.97%	0.16%	0.27%	3.19%

15. REFERENCES

- Banerjee, S. u. a., 2008. CCN5/WISP-2 expression in breast adenocarcinoma is associated with less frequent progression of the disease and suppresses the invasive phenotypes of tumor cells. *Cancer Research*, 68(18), S.7606–7612.
- Banerjee, S. u. a., 2003. WISP-2 Gene in Human Breast Cancer: Estrogen and Progesterone Inducible Expression and Regulation of Tumor Cell Proliferation. *Neoplasia (New York, N.Y.)*, 5(1), S.63–73.
- Davies, S.R. u. a., 2010. Differential expression of the CCN family member WISP-1, WISP-2 and WISP-3 in human colorectal cancer and the prognostic implications. *International Journal of Oncology*, 36(5), S.1129–1136.
- Dhar, G. u. a., 2007. Loss of WISP-2/CCN5 signaling in human pancreatic cancer: a potential mechanism for epithelial-mesenchymal-transition. *Cancer Letters*, 254(1), S.63–70.
- Fritah, A. u. a., 2008. Role of WISP-2/CCN5 in the Maintenance of a Differentiated and Noninvasive Phenotype in Human Breast Cancer Cells. *Molecular and Cellular Biology*, 28(3), S.1114–1123.

- Inadera, H. u. a., 2000. WISP-2 as a novel estrogen-responsive gene in human breast cancer cells. *Biochemical and Biophysical Research Communications*, 275(1), S.108–114.
- Jones, J.A. u. a., 2007. CCN5 expression in mammals. *Journal of Cell Communication and Signaling*, 1(2), S.127–143.
- Jun, J.-I. & Lau, L.F., 2011. Taking aim at the extracellular matrix: CCN proteins as emerging therapeutic targets. *Nature Reviews. Drug Discovery*, 10(12), S.945–963.
- Lake, A.C. u. a., 2003. CCN5 Is a Growth Arrest-Specific Gene That Regulates Smooth Muscle Cell Proliferation and Motility. *The American Journal of Pathology*, 162(1), S.219–231.
- Lake, A.C. & Castellot, John J., 2003. CCN5 modulates the antiproliferative effect of heparin and regulates cell motility in vascular smooth muscle cells. *Cell communication and signaling* [\[7\]](#): CCS, 1, S.5.
- Leask, A. & Abraham, D.J., 2006. All in the CCN family: essential matricellular signaling modulators emerge from the bunker. *Journal of cell science*, 119(Pt 23), S.4803–4810.
- Rittié, L. u. a., 2011. Spatial-temporal modulation of CCN proteins during wound healing in human skin in vivo. *Journal of Cell Communication and Signaling*, 5(1), S.69–80.
- Russo, J.W. & Castellot, John J., 2010. CCN5: biology and pathophysiology. *Journal of Cell Communication and Signaling*, 4(3), S.119–130.
- Sabbah, M. u. a., 2011. CCN5, a Novel Transcriptional Repressor of the Transforming Growth Factor β Signaling Pathway ∇ . *Molecular and Cellular Biology*, 31(7), S.1459–1469.
- Tanaka, I. u. a., 2005. Expression and regulation of WISP2 in rheumatoid arthritic synovium. *Biochemical and Biophysical Research Communications*, 334(4), S.973–978.
- Wiesman, K.C. u. a., 2010. CCN5, a secreted protein, localizes to the nucleus. *Journal of Cell Communication and Signaling*, 4(2), S.91–98.
- Yeger, H. & Perbal, B., 2007. The CCN family of genes: a perspective on CCN biology and therapeutic potential. *Journal of Cell Communication and Signaling*, 1(3), S.159–164.

16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide as bactericides. Sodium azide is toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number