# anti-Gliadin IgA ELISA

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
# anti-Gliadin sIgA ELISA

*For the in vitro determination of  
anti-Gliadin sIgA antibodies in stool*

Valid from 2015-05-07

**REF** K 9311

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>16</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>18</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>18</b>
<i>Sample storage</i>	18
<i>Extraction of the stool samples</i>	19
<i>Dilution of samples</i>	20
<b>7. ASSAY PROCEDURE</b>	<b>20</b>
<i>Principle of the test</i>	20
<i>Test procedure</i>	20
<b>8. RESULTS</b>	<b>22</b>
<b>9. LIMITATIONS</b>	<b>22</b>
<b>10. QUALITY CONTROL</b>	<b>22</b>
<i>Reference range</i>	22
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>23</b>
<i>Precision and reproducibility</i>	23
<i>Analytical Sensitivity</i>	23
<b>12. PRECAUTIONS</b>	<b>23</b>
<b>13. TECHNICAL HINTS</b>	<b>24</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>24</b>
<b>15. REFERENCES</b>	<b>25</b>

## 1. INTENDED USE

The assay is intended for the quantitative determination of anti-gliadin-sIgA antibodies in stool. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

The celiac/coeliac disease is caused by the gliadin fraction of wheat gluten and similar proteins of rye and barley. The disease is mainly manifested as chronic digestion insufficiency in children or young adults. In addition, patients with celiac disease have a greatly increased risk of developing malignant T-cell lymphoma of the small bowel, as T-cell lymphoma was found in 6 of 10 patients with coeliac sprue. Early diagnosis of celiac disease is important, because there is evidence that a gluten-free diet might help to prevent complications and malignancies (including intestinal lymphoma).

The measurement of anti-gliadin-sIgA-antibodies may be useful as a screening criterion before jejunal biopsy, the „Gold Standard“, and for the monitoring of the gluten-free diet treatment.

### Indication

- Celiac disease
- Food intolerance
- See also our proposal for stepwise diagnosis of gluten sensitivity on page 17.

### 3. MATERIAL SUPPLIED


Cat. No.	Label	Kit components	Quantity
K 9311	PLATE	Holder with precoated strips	12 x 8 wells
K 9311	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 9311	CONJ	Conjugate (Peroxidase-labeled), ready to use	1 x 15 ml
K 9311	CTRLNEG	Control negative, lyophilized	4 vials
K 9311	CTRLPOS	Control positive, lyophilized	4 vials
K 9311	CTRLCUTOFF	Cut-off control, lyophilized	4 vials
K 9311	SUB	TMB substrate, ready-to-use	1 x 15 ml
K 9311	STOP	ELISA stop solution, ready to use	1 x 15 ml
K 9311	IDK Extract®	Extraction buffer concentrate <i>IDK Extract® 2.5x</i>	1 x 100 ml
K 9311	DIL	Dilution buffer, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Horizontal microtiter plate shaker
- Calibrated precision pipettors and 5–1000 µl tips
- Foil to cover the microtiter plate
- Multi-channel pipets or repeater pipets
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- **Preparation of the extraction buffer:** The **extraction buffer concentrate IDK Extract®** must be diluted with ultra pure water **1:2.5** before use (100 ml *IDK Extract®* + 150 ml ultra pure water ), mix well. Crystals could occur due to high salt concentration in the stock solutions. Before dilution, the crystals must be redissolved at 37 °C in a water bath. The **IDK Extract®** is stable at **2–8 °C** until the expiry date stated on the label. Extraction buffer (1:2.5 diluted *IDK Extract®*) can be stored in a closed flask at **2–8 °C for three months**.
- The **lyophilized CTRLNEG, CTRLPOS and CTRLCUTOFF** (controls, negative, positive and cut-off) are stable at **2–8 °C** until the expiry date stated on the label. Reconstitution details are given in the data sheet. Diluted controls are not stable and can not be stored.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *Sample storage*

#### **Raw stool**

Raw stool samples can be stored for 4 weeks at -20 °C. Avoid repeated freezing and thawing.

## Stool suspensions

Stool extract can be stored for 7 days at 2–8°C or -20°C or for one day at room temperature (15–30°C). Avoid more than three freeze-thaw cycles.

### *Extraction of the stool samples*

**Diluted extraction buffer** is used as a sample extraction buffer. We recommend the following sample preparation:

## Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

### ***Stool sample tube – Instructions for use***

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

### ***SAS with 1.5 ml extraction buffer:***

Applied amount of stool: 15 mg

Buffer Volume: 1.5 ml

Dilution Factor: 1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty sample tube** with **1.5 ml** of ready-to-use extraction buffer before using it with the sample. Important: Allow the extraction buffer to reach room temperature.
- c) Unscrew the tube (orange part of cap) to open. Insert yellow dipstick into sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for app. 10 minutes improves the result.
- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.

- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

**Dilution I: 1:100**

### *Dilution of samples*

The suspension from the sample extraction (dilution I) is further diluted **1:50 with wash buffer**. For example:

- **20 µl** dilution I + **980 µl** wash buffer, mix well = **1:50 (dilution II)**  
This results in a final dilution of **1:5 000**.

For analysis, pipet **100 µl** of **dilution II** per well.

## **7. ASSAY PROCEDURE**

### *Principle of the test*

The antigen gliadin is immobilized on the microtiter plate. During the first incubation step, the human anti-gliadin-sIgA antibodies in the samples are bound by the immobilized antigen. After a washing step, the anti-gliadin-sIgA antibodies of the samples are quantitatively determined by addition of a peroxidase-labeled anti-sIgA antibody. Tetramethylbenzidine is used as a peroxidase substrate. The enzymatic reaction is stopped by an acidic stop solution. The absorbance of the color compound is determined photometrically at 450 nm. The measured absorbance is directly proportional to the amount of bound anti-gliadin-sIgA antibodies.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of STD /SAMPLE/CTRL (standards/sample/controls) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash the pre-coated microtiter plate <b>5 x with 250 µl wash buffer before use</b> . After the final washing step, the inverted microtiter plate should be tapped on absorbent paper.
2.	Add <b>100 µl</b> of <b>CTRLNEG</b> , <b>CTRLPOS</b> and <b>CTRLCUTOFF</b> (controls, negative, positive and cut-off) and diluted <b>SAMPLE</b> (patient samples).
3	Incubate for <b>1 hour</b> at room temperature <b>shaking</b> on a horizontal shaker.
4.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5.	Add <b>100 µl</b> of <b>conjugate</b> per well.
6.	Incubate for <b>1 hour</b> at room temperature <b>shaking</b> on a horizontal shaker.
7.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add <b>100 µl SUB</b> (TMB substrate).
9.	Incubate for <b>10–20 minutes</b> at room temperature*.
10.	Add <b>100 µl STOP</b> (ELISA stop solution) and mix shortly.
11.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of one sample exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

## 8. RESULTS

Samples with an optical density higher than the average optical density of the cut off control are positive.

$$\text{Cut off} = \text{OD}_{\text{cut off control}} = 100 \text{ U/l}$$

### Example

$$\text{OD}_{\text{patient sample}} = 0.685$$

$$\text{OD}_{\text{cut off control}} = 0.321 = 100 \text{ U/l}$$

$$\text{Concentration patient sample} = \frac{0.685 * 100 \text{ U/l}}{0.321} = 213.4 \text{ U/l}$$

**Attention:** Calculation is only valid for a sample dilution factor of 1:5000.

## 9. LIMITATIONS

The lower limit of the measurement range (see definition below) can be calculated as:

$$\text{LoB} \times \text{sample dilution factor to be used}$$

Samples with concentrations lower than the measurement range cannot be clearly quantified.

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

Based on an internal study with apparently healthy adults (n = 45) resulted in a reference value of < 100 U/l.

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay (n = 20)**

The precision (intra-assay variation) of the Immundiagnostik anti-gliadin-sIgA-antibody-ELISA test was calculated from 20 replicate determinations on each one of two samples.

Sample	anti-gliadin sIgA-Ab mean value [U/l]	CV [%]
1	527.7	2.9
2	348.9	3.9

#### **Inter-Assay (n = 12)**

The reproducibility of the results of two samples on different days was evaluated. The samples were measured by different technicians on different days using the anti-gliadin-sIgA-antibody-ELISA test.

Sample	anti-gliadin sIgA-Ab mean value [U/l]	CV [%]
1	425.3	14.9
2	237.6	8.7

### *Analytical Sensitivity*

The LoB (limit of blank) was evaluated according to the guideline CLSI EP17-A2 and resulted in 36.2 U/l.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions than sealed ones.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

### **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Control samples should be analyzed with each run.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Dieterich et al.: 1997; *Nat. Med.* **7** (3), 797
2. Rieken et al.: 1998; *Dtsch. med. Wschr.* **123**, 1454
3. Green et al.: 1998; *Clin. Persp. Gastroenterol.* November, **133**
4. Mothes, Th.: 1997; *Münsch. Med. Wschr.* **139**, 111

### Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by