

AmpliSens[®] *Chlamydia trachomatis*-EPh
PCR kit
Instruction Manual

AmpliSens[®]



Federal Budget Institute of
Science "Central Research
Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia

TABLE OF CONTENTS

1. INTENDED USE.....	3
2. PRINCIPLE OF PCR DETECTION	3
3. CONTENT	3
4. ADDITIONAL REQUIREMENTS	4
5. GENERAL PRECAUTIONS	5
6. SAMPLING AND HANDLING.....	5
7. WORKING CONDITIONS	6
8. PROTOCOL	6
9. DATA ANALYSIS	7
10. TROUBLESHOOTING	8
11. TRANSPORTATION	8
12. STABILITY AND STORAGE	9
13. SPECIFICATIONS	9
14. REFERENCES.....	9
15. QUALITY CONTROL	9
16. KEY TO SYMBOLS USED	10

1. INTENDED USE

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is an in vitro nucleic acid amplification test for qualitative detection of *Chlamydia trachomatis* DNA in the clinical material (urogenital and conjunctival swabs, urine sediment, prostate gland secretion, and synovial fluid) by using electrophoretic detection of the amplified products in agarose gel.

2. PRINCIPLE OF PCR DETECTION

Chlamydia trachomatis detection by the polymerase chain reaction (PCR) is based on the amplification of pathogen DNA specific region using specific *Chlamydia trachomatis* primers. After PCR, the amplified product is detected in agarose gel. AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is a qualitative test that contains the Internal Control (IC). It must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition. AmpliSens® *Chlamydia trachomatis*-EPh PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase by using a wax layer. Wax melts and the reaction mixture components mix only at 95 °C.

3. CONTENT

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is produced in 2 forms:

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit variant 100 R (0.5-ml tubes),

REF B1-100-R0,5-CE.

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit variant 100 R (0.2-ml tubes),

REF B1-100-R0,2-CE.

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit variant 100 R (0.2-ml tubes) in bulk¹,

REF B1-100-R0,2-CE-B.

¹ In bulk form contains unlabeled tubes. Tubes with identical reagent are packed in one bag with label.

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit variant 100 R includes:

Reagent	Description	Volume, ml	Quantity
PCR-mix-1-R <i>Chlamydia trachomatis</i> ready-to-use single-dose test tubes (<i>under wax</i>)	colorless clear liquid	0.005	110 tubes of 0.5 or 0.2 ml
PCR-mix-2 blue	blue clear liquid	1.2	1 tube
Mineral oil for PCR	colorless viscous liquid	4.0	1 vial
Positive Control DNA <i>Chlamydia trachomatis</i> (C₊ <i>Chlamydia trachomatis</i>)	colorless clear liquid	0.2	1 tube
DNA-buffer	colorless clear liquid	0.5	1 tube
Negative Control (C₋)*	colorless clear liquid	1.2	1 tube
Internal Control complex ICc**	colorless clear liquid	1.0	1 tube

* must be used in the extraction procedure as Negative Control of Extraction.

** add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture (see DNA-sorb-AM, **REF** K1-12-100-CE or DNA-sorb-B, **REF** K1-2-100-CE protocols).

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit variant 100 R is intended for 110 reactions (including controls).

4. ADDITIONAL REQUIREMENTS

- DNA extraction kit.
- Agarose gel detection kit.
- Disposable powder-free gloves.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol filters (up to 200 µl).
- Vortex mixer.
- Desktop microcentrifuge with a rotor for 2-ml reaction tubes (RCF max. 16,000 g).
- PCR box or biological cabinet.
- Vacuum aspirator with a flask for removing supernatant.
- Tube racks.
- 1.5-ml sterile polypropylene tubes.
- Refrigerator for 2–8 °C.
- Deep-freezer with temperature ≤ –16 °C.
- Waste bin for used tips.
- Permanent pen for labeling.

REF B1-100-R0,2-CE; **REF** B1-100-R0,5-CE; **REF** B1-100-R0,2-CE-B /

VER 18.08.09 – 14.10.13 /Page 4 of 11

- Thermostat with controlled temperature and capable of incubating at 25°C and 100 °C.
- Personal thermocyclers (for example, Gradient Palm Cycler (Corbett Research, Australia) or MaxyGene (Axygen Scientific, USA)).

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store and handle amplicons away from all other reagents.
- Thaw all components thoroughly at room temperature before starting detection.
- When thawed, mix the components and centrifuge briefly.
- Use protective gloves, laboratory coats, protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all samples and unused reagents in compliance with local authorities requirements.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills with 0.5% sodium hypochlorite solutions or other suitable disinfectant.
- Avoid contact with the skin, eyes, and mucous membranes. If skin, eyes and mucous membranes contact, immediately flush with water, seek medical attention
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional; it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area where you carried out the previous step.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING



Obtaining samples of biological materials for PCR-analysis, transportation, and storage is described in manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is intended for analysis of DNA extracted with DNA extraction kits from:

- Urogenital swabs;
- Urine sediment (use the first part of the stream);
- Prostate gland secretion;
- Synovial fluid;
- Conjunctival smears.

7. WORKING CONDITIONS

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. DNA extraction

It is recommended that the following nucleic acid extraction kits are used:

- DNA-sorb-AM, **REF** K1-12-100-CE.
- DNA-sorb-B (for the prostate gland secretion), **REF** K1-2-100-CE.



Extract DNA according to the manufacturer's protocol.

8.2. Preparing PCR

Total reaction volume is **25 µl**, the volume of DNA sample is **10 µl**.

1. Take the required quantity of tubes with **PCR-mix-1-R *Chlamydia trachomatis*** and wax for amplification of DNA extracted from clinical and control samples.
2. Add **10 µl of PCR-mix-2 blue** to the surface of the wax layer of each tube ensuring that it does not fall under the wax and mix with the reagents in the tube.
3. Add above 1 drop of **mineral oil for PCR** (~25 µl). This step can be omitted if a thermocycler with a constant-temperature cover is used.

8.2.2 Amplification

Use the prepared tubes for PCR. **Add 10 µl of DNA samples** obtained from clinical or control samples at the DNA extraction stage using tips with aerosol barrier above or under the oil layer.

Carry the **control amplification reactions**:

NCA -Add 10 µl of **DNA-buffer** to the tube for Negative Control of Amplification (NCA).

C+ -Add 10 µl of **Positive Control DNA *Chlamydia trachomatis*** to the tube for Positive Control of Amplification.

Run the following program in the thermocycler (see Table 1). When the temperature reaches 95°C (pause mode), insert tubes into the cells of the thermocycler and press the button to continue.

It is recommended to sediment drops from the walls of tubes by short centrifugation (1–3 s) before placing them into the thermocycler.

Table 1

Programming thermocyclers for *Chlamydia trachomatis* DNA amplification

	Thermocyclers with active temperature adjustment:						Thermocyclers with block temperature adjustment:		
	Terzik (DNA-Technology)			GeneAmp PCR System 2700 (Applied Biosystems), Gradient Palm Cycler (Corbett Research)			Uno II (Biometra), MiniCycler, PTC-100 (MJ Research)		
Step	Temperature	Time	Cycles	Temperature	Time	Cycles	Temperature	Time	Cycles
0	95 °C	pause		95 °C	pause		95 °C	pause	
1	95 °C	5 min	1	95 °C	5 min	1	95 °C	5 min	1
2	95 °C	10 s	42	95 °C	15 s	42	95 °C	1 min	42
	65 °C	10 s		65 °C	25 s		65 °C	1 min	
	72 °C	10 s		72 °C	25 s		72 °C	1 min	
3	72 °C	1 min	1	72 °C	1 min	1	72 °C	1 min	1
4	4 °C	storage		4 °C	storage		10 °C	storage	

Amplification in thermocyclers with block and active temperature adjustment continues for 2 h 30 min and 1 h 50 min, respectively.

After the reaction is completed, PCR tubes should be collected and transferred to the room for PCR product analysis.

Analysis of amplification products is performed by separation of DNA fragments in agarose gel.

The amplified samples can be stored for 16 h at room temperature, for 1 week at 2–8 °C, and for a long time at –16 °C (warm up samples to room temperature before electrophoretic run).

9. DATA ANALYSIS

It is recommended to use the following detection agarose kit:

- EPh variant 200, **REF** K5-200-CE.

Analysis of results is based on the presence or absence of specific bands of amplified DNA in agarose gel (1.7%). The length of specific amplified DNA fragments is as follows:

- *Chlamydia trachomatis*, 330 bp;
- Internal Control, 740 bp.



Put on a protective mask or use a glass filter while visualizing and photographing the gel.

9.1. Interpretation of results

Table 2

Results for controls

Control	Which step of test is controlled	Specific 330-bp bands in agarose gel	Specific 740-bp bands in agarose gel	Interpretation
C–	DNA extraction	No	Yes	OK
NCA	Amplification	No	No	OK
C+	Amplification	Yes	No	OK

- The sample is considered to be positive for *Chlamydia trachomatis* DNA if the 330-bp band is present in gel. The band of IC (740 bp) may be absent in samples with a high concentration of *Chlamydia trachomatis* DNA.
- The sample is considered to be negative for *Chlamydia trachomatis* DNA if the 330-bp band is absent whereas the 740-bp band is present in gel.

In addition to the specific bands, fuzzy bands corresponding to primer dimers may appear in lanes below the 100-bp level.

10. TROUBLESHOOTING

The results of analysis are not taken into account in the following cases:

- If the results obtained for the control samples results do not correspond to the results for controls listed in Table 2, the appropriate stage of the test should be repeated. Discard any reagents that may be suspect.
- If neither the 330-kb nor the 740-kb band is detected in lanes, the result of analysis of this sample is irrelevant. In this case, analysis of this sample should be repeated starting from the DNA extraction stage. This may be caused by clinical processing errors that led to the loss of DNA or inhibition of PCR.
- The appearance of nonspecific bands of different molecular weight in lanes may be caused by the lack of “hot start” or an inappropriate temperature regime in the thermocycler. In this case, the results of analysis are invalid.
- The appearance of the specific 330-bp band in lanes corresponding to negative controls (NCA, C–) suggests contamination of reagents or samples. In such cases, the results of analysis are considered to be invalid. Analysis of all samples must be repeated and measures to detect and eliminate the source of contamination must be taken.

11. TRANSPORTATION

AmpliSens® *Chlamydia trachomatis*-EPh PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the AmpliSens® *Chlamydia trachomatis*-EPh PCR kit are to be stored at 2–8 °C when not in use. They are stable until the expiration date on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

13. SPECIFICATIONS

13.1. Sensitivity

The analytical sensitivity of AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is not less than 5×10^3 genome equivalents per 1 ml of sample (GE/ml).



The claimed analytical characteristics of AmpliSens® *Chlamydia trachomatis*-EPh PCR kit are guaranteed only when additional reagent kits DNA-sorb-AM or DNA-sorb-B (for the prostate gland secret) and EPh, manufactured by FBIS CRIE, are used.

13.2. Specificity

The specificity of AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is ensured by selection of specific primers and stringent reaction conditions. It was confirmed in laboratory and clinical trials.








14. REFERENCES

1. Manual “Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology”, Moscow, 2008.

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of AmpliSens® *Chlamydia trachomatis*-EPh PCR kit is tested against predetermined specifications to ensure consistent product quality.

16. KEY TO SYMBOLS USED

REF	Catalogue number		Caution
LOT	Batch code		Sufficient for
RUO	Research use only		Expiration Date
VER	Version		Consult instructions for use
	Temperature limitation	NCA	Negative control of amplification
	Manufacturer	C-	Negative control of extraction
	Date of manufacture	C+	Positive control of amplification
FBIS CRIE	Federal Budget Institute of Science “Central Research Institute for Epidemiology”	ICc	Internal control complex

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
29.06.11 VV	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"
	Cover page	The phrase "For Professional Use Only" was added
	Content	The "Explanation of Symbols" section was renamed to "Key to Symbols Used"
	Key to Symbols Used	The explanation of symbols was corrected
12.07.11 VV	Content	New sections "Working conditions" and "Transportation" were added
	Stability and Storage	The information about the shelf life of open reagents was added
11.04.12 IvI	Title page, Key to symbols used	Symbol IVD <i>in vitro</i> diagnostic medical device was changed to RUO research use only
14.10.13 ME	Content	The form in bulk was added
	Footer	REF B1-100-R0,2-CE-B was added