



For Professional Use Only

**AmpliSens<sup>®</sup> *C.trachomatis* / *Ureaplasma* /  
*M.genitalium*-MULTIPRIME-FEP  
PCR kit  
Instruction Manual**

**AmpliSens<sup>®</sup>**



Ecoli s.r.o., Studenohorska 12  
841 03 Bratislava 47  
Slovak Republic  
Tel.: +421 2 6478 9336  
Fax: +421 2 6478 9040



Federal Budget Institute of  
Science "Central Research  
Institute for Epidemiology"  
3A Novogireevskaya Street  
Moscow 111123 Russia

## TABLE OF CONTENTS

1. INTENDED USE .....	3
2. PRINCIPLE OF PCR DETECTION .....	3
3. CONTENT .....	3
4. ADDITIONAL REQUIREMENTS .....	4
5. GENERAL PRECAUTIONS.....	5
6. SAMPLING AND HANDLING .....	6
7. WORKING CONDITIONS.....	6
8. PROTOCOL .....	6
9. DATA ANALYSIS .....	8
10. TROUBLESHOOTING.....	9
11. TRANSPORTATION.....	9
12. STABILITY AND STORAGE.....	9
13. SPECIFICATIONS.....	10
14. REFERENCES .....	10
15. QUALITY CONTROL.....	11
16. KEY TO SYMBOLS USED .....	12

## 1. INTENDED USE

**AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit is an *in vitro* nucleic acid amplification test for multiplex detection of DNA of *Chlamydia trachomatis*, *Ureaplasma* (*U.parvum* and *U.urealyticum*), and *Mycoplasma genitalium* in clinical materials (urogenital, rectal and pharyngeal swabs; conjunctival discharge; prostate gland secretion; and urine samples) by using end-point hybridization-fluorescence detection of amplified products.



The results of PCR analysis are taken into account in complex diagnostics of disease.

## 2. PRINCIPLE OF PCR DETECTION

*C.trachomatis / Ureaplasma / M.genitalium* detection by the multiplex polymerase chain reaction (PCR) is based on the amplification of pathogen genome specific region using specific *C.trachomatis / Ureaplasma / M.genitalium* primers. In **Fluorescent End-Point** PCR, the amplified product is detected by using fluorescent dyes. These dyes are linked to oligonucleotide probes which bind specifically to the amplified product during thermocycling. A multichannel rotor-type fluorometer is specially designed to detect fluorescence emission from the fluorophores in the reaction mixture after PCR. It allows detection of the accumulating product without re-opening the reaction tubes after the PCR run. **AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit is a qualitative test that contains the Internal Control (IC). It must be used in the extraction procedure in order to control the extraction of each individual sample and to identify possible reaction inhibition. **AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase by using a wax layer. Wax melts and reaction components mix only at 95 °C.

## 3. CONTENT

**AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit is produced in 2 forms:

AmpliSens<sup>®</sup> *C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP* PCR kit variant FEP (0.5 ml tubes),

**REF** B46-100-R0,5-FEP-CE.

AmpliSens<sup>®</sup> *C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP* PCR kit variant FEP (0.2 ml tubes),

**REF** B46-100-R0,2-FEP-CE.

**AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit variant FEP includes:

<b>Reagent</b>	<b>Description</b>	<b>Volume, ml</b>	<b>Quantity</b>
<b>PCR-mix-1-FL C.trachomatis / Ureaplasma / M.genitalium</b> ready-to-use single-dose test tubes ( <i>under wax</i> )	colorless clear liquid	0.01	110 tubes of 0.5 or 0.2 ml
<b>PCR-mix-2-FL-red</b>	red clear liquid	1.1	1 tube
<b>Mineral oil for PCR*</b>	colorless viscous liquid	4.0	1 dropper bottle
<b>PCR-mix-Background-red</b>	red clear liquid	0.6	1 tube
<b>Positive Control complex (C+)</b>	colorless clear liquid	0.2	1 tube
<b>DNA-buffer</b>	colorless clear liquid	0.5	1 tube
<b>Negative Control (C-)**</b>	colorless clear liquid	1.2	1 tube
<b>Internal Control-FL (IC)***</b>	colorless clear liquid	1.0	1 tube

\* must be used for thermocyclers without a constant-temperature lid.

\*\* must be used in the extraction procedure as Negative Control of Extraction.

\*\* add 10 µl of Internal Control-FL (IC) during the DNA extraction procedure directly to the sample/lysis mixture (DNA-sorb-AM, **REF** K1-12-100-CE).

**AmpliSens<sup>®</sup> C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit variant FEP is intended for 110 reactions, including controls.

#### 4. ADDITIONAL REQUIREMENTS

- DNA extraction kit.
- Transport medium.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol barriers (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with rotor for 2 ml reaction tubes.
- PCR box.
- Personal thermocyclers (for example, Terzik (DNA-Technology), Gradient Palm Cyclor (Corbett Research, Australia), MaxyGene (Axygen, USA), GeneAmp PCR System 2700 (Applied Biosystems, USA) or equivalent).

- Fluorometer (fluorescence detector; for example, ALA-1/4 (Biosan, Latvia) or equivalent).
- Disposable polypropylene microtubes for PCR (0.5- or 0.2-ml; for example, Axygen, USA).
- Refrigerator for 2–8 °C.
- Deep-freezer for  $\leq -16$  °C.
- Waste bin for used tips.

## 5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store and handle amplicons away from all other reagents.
- Thaw all components thoroughly at room temperature before starting detection.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats, protect eyes while samples and reagents handling. Thoroughly wash hands afterward.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all samples and unused reagents in compliance with local authorities requirements.
- Samples should be considered potentially infectious and handled in a biological cabinet in accordance with appropriate biosafety practices.
- Clean and disinfect all samples or reagent spills using a disinfectant such as 0.5 % sodium hypochlorite, or other suitable disinfectant.
- Avoid contact with the skin, eyes, and mucous membranes. If skin, eyes, and mucous membranes contact, immediately flush with water and seek medical attention.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional; it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

## 6. SAMPLING AND HANDLING



Obtaining samples of biological materials for PCR-analysis, transportation and storage is described in manufacturer's handbook [1]. It is recommended to read this handbook before starting work.

**AmpliSens® *C.trachomatis* / *Ureaplasma* / *M.genitalium*-MULTIPRIME-FEP** PCR kit is intended for analysis of DNA extracted by using DNA extraction kits from:

- urogenital swabs;
- rectal swabs;
- pharyngeal swabs;
- conjunctival discharge and prostate gland secretion;
- urine samples (use the first part of the stream).

## 7. WORKING CONDITIONS

**AmpliSens® *C.trachomatis* / *Ureaplasma* / *M.genitalium*-MULTIPRIME-FEP** PCR kit should be used at 18–25 °C.

## 8. PROTOCOL

### 8.1. DNA Extraction

It's recommended to use the following nucleic acid extraction kits:

- "DNA-sorb-AM", **REF** K1-12-100-CE.
- Other nucleic acid extraction kits recommended by Federal Budget Institute of Science "Central Research Institute for Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being (see manufacturer's **Guidelines**).



Extract DNA according to the manufacturer's instructions.

### 8.2. Preparing PCR

The total reaction volume is **30 µl**, the volume of DNA sample is **10 µl**.

#### 8.2.1 Preparing tubes for PCR

1. Prepare the required number of tubes with **PCR-mix-1-FL *C.trachomatis* / *Ureaplasma* / *M.genitalium*** for amplification of DNA from clinical and control samples.
2. Add **10 µl** of **PCR-mix-2-FL-red** to the surface of the wax layer in each tube ensuring that it does not fall under the wax and mix with **PCR-mix-1-FL *C.trachomatis* / *Ureaplasma* / *M.genitalium***.
3. Add above **1 drop of mineral oil for PCR (~ 25 µl)**.

4. Prepare 1 tube with **PCR-mix-1-FL *C.trachomatis* / *Ureaplasma* / *M.genitalium*** and mark them as **Background**. Add **20 µl** of **PCR-mix-Background-red** to the surface of the wax layer of each tube, ensuring that it does not fall under the wax and mix with **PCR-mix-1-FL *C.trachomatis* / *Ureaplasma* / *M.genitalium***. Add above **1** drop of **mineral oil for PCR**.



Use **mineral oil for PCR** if working with thermocyclers without a constant-temperature lid.



PCR-mix-Background-red is to be used if DNA was extracted using DNA-sorb-AM, **REF** K1-12-100-CE, or DNA-sorb-B, **REF** K1-2-100-CE. See the manufacturer's instruction if other nucleic acid extraction kits were used.

5. Using tips with aerosol barrier add **10 µl** of **DNA samples** obtained from clinical or control samples at the DNA extraction stage.
6. Carry out the control amplification reactions:
- NCA -Add **10 µl** of **DNA-buffer** to the tube labeled NCA (Negative Control of Amplification).
- C+ -Add **10 µl** of **Positive Control complex** to the tube labeled C+ (Positive Control of Amplification).



It is recommended to sediment drops from walls of tubes by short centrifugation (1–3 s) before placing them in the thermocycler.

## 8.2.2 Amplification

Run the following program on the thermocycler (see Table 1). When the temperature reaches 95 °C (pause regimen), insert tubes to cells of amplifier and press the button to continue.

Table 1

### Amplisens-1-FEP amplification program

Step	Terzik (DNA-Technology)			GeneAmp PCR System 2700 (Applied Biosystems)			Gradient Palm Cycler (Corbett Research), MaxyGene (Axygen)		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
0	95	pause		95	pause		95	pause	
1	95	5 min	1	95	5 min	1	95	5 min	1
2	95	2 s	35	95	20 s	20	95	2 s	24
	65	5 s		65	25 s		65	10 s	
	72	5 s		72	30 s		72	10 s	
3	95	2 s	9	95	20 s	24	95	2 s	20
	60	10 s		60	30 s		60	15 s	
	72	5 s		72	30 s		72	10 s	
4	95	2 s	1	95	20 s	1	95	2 s	1
	60	10 s		60	30 s		60	15 s	
5	10	storage		10	storage		10	storage	



Programming for other thermocyclers is described in manufacturer's guidelines [2]. It is recommended to read this handbook before starting work.

## 9. DATA ANALYSIS

Detection is conducted using a fluorescence detector.



Please read the fluorescence detector Operating Manual before using this kit.

Program the detector according to the manufacturer's manual and Important product information bulletin.

Accumulation of ***Chlamydia trachomatis* DNA** amplification product is detected in the **FAM/Green** channel, ***Ureaplasma spp. (U.parvum and U.urealyticum)* DNA** is detected in the **HEX/Yellow** channel, ***Mycoplasma genitalium* DNA** is detected in the **ROX/Orange** channel, and Internal Control-FL (IC) is detected in the **Cy5/Red** channel.

### Results interpretation

1. Principle of interpretation:

- *Chlamydia trachomatis* DNA is **detected** in a sample if its signal in the FAM channel is greater than the defined threshold value of the positive result.
- *Ureaplasma spp. (U.parvum and U.urealyticum)* DNA is **detected** in a sample if its signal in the HEX channel is greater than the defined threshold value.
- *Mycoplasma genitalium* DNA is **detected** in a sample if its signal in the ROX channel is greater than the defined threshold value of the positive result.
- *Chlamydia trachomatis* DNA, *Ureaplasma spp.* DNA and *Mycoplasma genitalium* DNA are **not detected** in a sample if their signals in FAM and HEX channels are less than the defined threshold value of the negative result while the signal in the Cy5 channel is greater than the defined threshold value.
- The result is **invalid** if the signal of a sample in FAM, HEX, ROX and Cy5 channels is less than the defined threshold values for these channels.
- The result is **equivocal** if the signal of a sample in the FAM and/or ROX channel(s) is greater than the defined threshold value of the negative result but less than the threshold value of the positive result (the signal is between thresholds).



If the result is invalid or equivocal, PCR should be repeated once again.



2. Result of the analysis is considered reliable only if the results obtained for Positive and Negative Controls of Amplification as well as Negative Control of Extraction are correct (Table 2).

Table 2

**Results for controls**

Control	Stage for control	Result of automatic interpretation				Interpretation
		FAM channel (samples)	HEX channel (samples)	ROX channel (samples)	Cy5 channel (IC)	
<b>C-</b>	DNA extraction	< threshold of negative result	< threshold	< threshold of negative result	> threshold	“-“ or “OK”
<b>NCA</b>	Amplification	< threshold of negative result	< threshold	< threshold of negative result	< threshold	“nd”
<b>C+</b>	Amplification	> threshold of positive result	> threshold	> threshold of positive result	> threshold	“+” or “OK”

**10. TROUBLESHOOTING**

Results of analysis are not being registered in the following cases:

- If the signal of C+ in the FAM, HEX, and ROX channels is less than the threshold of the positive result, run PCR and detection for all samples in which the signal in either channel was less than the threshold of the positive result.
- If the signal of the C– and/or NCA in the FAM, HEX, and ROX channels is greater than the threshold of the positive result, run PCR test starting from the extraction stage for all samples in which the signal in either channel was greater than the threshold of the positive result.

If you have any further questions or if encounter problems, please contact our Authorized representative in the European Community.

**11. TRANSPORTATION**

**AmpliSens® C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit should be transported at 2–8 °C for no longer than 5 days.

**12. STABILITY AND STORAGE**

All components of the **AmpliSens® C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit are to be stored at 2–8 °C when not in use. All components of the **AmpliSens® C.trachomatis / Ureaplasma / M.genitalium-MULTIPRIME-FEP** PCR kit are to be stable until the expiration date. The shelf life of reagents before and after the first use is the same, unless otherwise stated.



PCR-mix-1-FL *C.trachomatis* / *Ureaplasma* / *M.genitalium* is to be kept away from light.

## 13. SPECIFICATIONS

### 13.1. Sensitivity

Clinical material	Nucleic acid extraction kit	PCR kit	Microorganism	Sensitivity, GE/ml <sup>1</sup>
Cervical urethral swabs <sup>2</sup> or	DNA-sorb-AM	PCR kit variant FEP	<i>Chlamydia trachomatis</i>	5x10 <sup>2</sup>
			<i>Ureaplasma</i> spp.	10 <sup>3</sup>
			<i>Mycoplasma genitalium</i>	10 <sup>3</sup>
Urine <sup>3</sup>	DNA-sorb-AM	PCR kit variant FEP	<i>Chlamydia trachomatis</i>	10 <sup>3</sup>
			<i>Ureaplasma</i> spp.	2x10 <sup>3</sup>
			<i>Mycoplasma genitalium</i>	2x10 <sup>3</sup>



Analytical Sensitivity of each microorganism doesn't change even at high concentrations of two other microorganisms.

### 13.2. Specificity

Specificity of **AmpliSens<sup>®</sup> *C.trachomatis* / *Ureaplasma* / *M.genitalium*-MULTIPRIME-FEP** PCR kit is ensured by selection of specific primers and probes, as well as the selection of strict reaction conditions. The primers and probes have been checked for possible homologies to all sequences published in gene banks by sequence comparison analysis. Specificity of **AmpliSens<sup>®</sup> *C.trachomatis* / *Ureaplasma* / *M.genitalium*-MULTIPRIME-FEP** PCR kit was confirmed in laboratory clinical trials.

Nonspecific responses were absent while testing human DNA samples and DNA samples of microorganisms: *Gardnerella vaginalis*, *Lactobacillus* spp., *Escherichia coli*, *Staphylococcus* spp., *Streptococcus* spp., *Candida albicans*, *Ureaplasma urealyticum*, *Ureaplasma parvum*, *Mycoplasma hominis*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Neisseria* spp., *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, *Treponema pallidum*, *Toxoplasma gondii*, HSV 1 and 2, CMV, HPV.

## 14. REFERENCES

- Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics", developed by Federal State Institute of Science "Central Research Institute of Epidemiology" of Federal Service for Surveillance on Consumers' Rights

<sup>1</sup> The quantity of genome equivalents of microorganism per 1 ml of the sample from transport medium.

<sup>2</sup> Cervical and urethral swabs are to be placed into Transport Medium for Swabs (**REF** 956-CE, 987-CE) or Transport Medium with Mucoytic (**REF** 952-CE, 953-CE).

<sup>3</sup> Pretreatment is needed.














Protection and Human Well-Being, Moscow, 2008.

2. Guidelines “End-Point PCR Detection of STIs and Other Reproductive Tract Infections”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology” of Federal Service for Surveillance on Consumers’ Rights Protection and Human Well-Being, Moscow.

## 15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of **AmpliSens® *C.trachomatis* / *Ureaplasma* / *M.genitalium*-MULTIPRIME-FEP** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

## 16. KEY TO SYMBOLS USED

	Catalogue number		Sufficient for
	Batch code		Expiration Date
	<i>In vitro</i> diagnostic medical device		Consult instructions for use
	Version		Keep away from sunlight
	Temperature limitation	<b>NCA</b>	Negative control of amplification
	Manufacturer	<b>C-</b>	Negative control of extraction
	Date of manufacture	<b>C+</b>	Positive control of amplification
	Authorised representative in the European Community	<b>IC</b>	Internal control
	Caution		

### List of Changes Made in the Instruction Manual

<b>VER</b>	<b>Location of changes</b>	<b>Essence of changes</b>
23.06.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"