# Pregnenolonsulfat LC-MS/MS Kit

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Hersteller

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

# Pregnenolone sulfate LC-MS/MS Kit

*For the determination of pregnenolone sulfate  
in EDTA plasma and serum*

Valid from 13.11.2013

REF

KM2000

IVD

**Immundiagnostik AG**, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	18
7. CHROMATOGRAPHIC CONDITIONS	19
8. MS/MS METHOD	20
9. CALCULATION	21
10. EXAMPLES OF CHROMATOGRAMS	21
11. QUALITY CONTROL	22
12. PERFORMANCE CHARACTERISTICS	23
13. PRECAUTIONS	23
14. TECHNICAL HINTS	23
15. DISPOSAL	23
16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24
17. REFERENCES	24

## 1. INTENDED USE

The described pregnenolone sulfate LC-MS/MS-Application is intended for the quantitative determination of pregnenolone sulfate in serum and EDTA-plasma. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Steroids and neurosteroids are synthesized from cholesterol, the basic building block of all steroid hormones. While the classic steroids are produced in tissues like adrenal gland, gonads or placenta, the synthesis of neurosteroids takes place in the central and peripheral nervous system. In the 1980s, Baulieu et al. demonstrated that several steroids, such as dehydroepiandrosterone and pregnenolone, are synthesized *de novo* in the central nervous system and present in higher concentrations in the brain than in the blood. These steroids are universally referred as neurosteroids and are differently distributed in the brain and the pituitary gland.

Neurosteroids are formed from cholesterol via pregnenolone and progesterone (see Fig.1). Pregnenolone and progesterone are precursors not only of the synthesis of gluco- and mineral corticoids in the adrenal gland and of the sexual hormones in the gonads or placenta, but also of neurosteroids in the central nervous system, i.e. pregnenolone is an endogenous, naturally produced steroid that is a precursor of many body hormones.

Numerous studies have demonstrated that neurosteroids interact with neurotransmitter receptors, influence the brain excitability, thereby functioning as potent allosteric modulators of several neurotransmitter receptors: the endogenous neurosteroid allopregnanolone, a positive modulator of the  $\gamma$ -aminobutyric acid type A (GABA<sub>A</sub>) receptor complex, activates the GABA<sub>A</sub> receptor functions lowered by various stresses, and has anxiolytic and anti-stress effects. In contrast, negative modulators, like dehydroepiandrosterone sulfate and pregnenolone sulfate, inhibit almost completely the GABA<sub>A</sub>-receptor.

Due to their lipophilicity, most of the free steroids can easily pass the blood-brain barrier. A number of steroids are conjugated as sulfate or fatty acid esters and occur in concentrations much higher than these of the free steroids. In contrast to the free steroids, they can not pass through the blood-brain barrier and represent the active forms involved in the control of various physiological and pathophysiological processes.

Neurosteroids affect a number of processes in the central nervous system, which have a powerful effect on our thinking and feeling. Furthermore, neurosteroids influence the social and sexual behavior. At the same time, they are promising pharmaceutical targets for important indications like epilepsy, anxiety disorders and dementia.

For these reasons, the determination of neurosteroids, especially of pregnenolone and pregnenilone sulfate, is expected to be helpful for understanding the roles of neurosteroids, for diagnosis of psychic or mental disorders and for the development of new steroidal therapeutic agents.

**Indications**

- Cognitive ability
- Weakness of memory
- Mood swings like depression


Fig. 1: Biosynthesis of neuroactive steroids, Weiss M &amp; Hess S, 2000

### 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KM2000LA	MOPHA A	Mobile phase A	1000 ml
KM2000LB	MOPHA B	Mobile phase B	1000 ml
KM2000KA	CAL	Calibrator (lyophilized, see label for concentration)	5 vials (each 300 µl)
KM2000KO	CTRL 1 CTRL 2	Control 1 and 2 (lyophilized, see product specification for concentration)	5 vials (each 300 µl)
KM2000FR	PREC	Precipitation solution (-20°C, contains internal standard)	50 ml
KM2000DIL	DIL	Dilution solution (see sample preparation)	120 ml
KM2000RE	RECSOL	Reconstitution solution	5 ml

As a first step for the application of the Immundiagnostik AG pregnenolone sulfate LC-MS/MS-Kit, a tuning is necessary to estimate the optimal LC-MS/MS-settings as well as to assess the sufficiency of the sensitivity.

The UPLC separation column (KM2000RP) can be ordered separately from Immundiagnostik AG.

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml Eppendorf reaction tubes
- Precision pipettors and disposable tips to deliver 10–1000 µl
- Cooling centrifuge at 4°C for 1.5 ml Eppendorf reaction tubes and capable of 10000 g
- Vortex-Mixer
- LC-MS/MS-Equipment
- RP-C18 column, e.g. Acquity BEH 18, 1,7 µm (2,1 x 50 mm)

### 5. PREPARATION AND STORAGE OF REAGENTS

The controls (CTRL 1 and CTRL 2), calibrator (CAL) and precipitation solution (PREC) are stable at -20°C until the expiry date. All other test reagents are stable until the expiry date when stored at 2–8°C.

### Calibrator and controls

Dissolve **calibrator** (CAL) and **controls** (CTRL1, CTRL2) in 300 µl of reconstitution solution (RECSOL).

## 6. SAMPLE PREPARATION

Serum and EDTA-plasma samples are suited for the assay.

Control samples should be analyzed with each run.

1.	Prior to use in the assay allow all samples and reagents, except for the precipitation solution (PREC), to come to room temperature (18–26 °C). Mix well samples and reagents before use
2.	Pipette into 1.5 ml Eppendorf reaction tubes: <b>200 µl</b> Sample, Calibrator (CAL) or Control (CTRL 1, CTRL 2) + <b>400 µl</b> -20 °C cold precipitation solution (PREC incl. internal standard), mix well
3.	Leave the samples for <b>15 minutes at 2–8 °C</b> and then centrifuge for <b>15 minutes at 4 °C</b> and 10000 g
4.	Dilute the supernatant with DIL <b>1:50</b> (e.g. 10 µl supernatant + 490 µl DIL) and transfer into LC-MS/MS suitable vials

## 7. CHROMATOGRAPHIC CONDITIONS

### *UPLC-Method*

<b>Column material:</b>	e.g. Acquity BEH C18; 1.7 µm
<b>Column dimension:</b>	2.1 x 50 mm
<b>Flow rate:</b>	0.5 ml/min
<b>Column temperature:</b>	45 °C
<b>Injection volume:</b>	50 µl
<b>Running time:</b>	3.5 min

<b>Gradient:</b>	0 min	70% A	30% B
	0.5 min	70% A	30% B
	1.0 min	50% A	50% B
	1.5 min	0% A	100% B
	3.5 min	0% A	100% B

It is recommended that a guard column/filter is used to extend column's life.

After the analysis, the separation column should be washed with ca. 20 ml of 50 % methanol. The column can be stored in 50 % methanol.

## 8. MS/MS METHOD

The MS/MS method listed here is an example for a Waters Quattro Premier XE tandem mass spectrometer.

<b>Mode:</b>	MRM
<b>Polarity:</b>	ESI-
<b>Capillary (kV):</b>	3
<b>Cone (V):</b>	var.
<b>Extractor (V):</b>	-3
<b>RF Lens (V):</b>	0
<b>Source Temperature (°C):</b>	130
<b>Desolvation Temperature (°C):</b>	450
<b>Cone Gas Flow (L/Hr):</b>	50
<b>Desolvation Gas Flow (L/Hr):</b>	950
<b>Collision Gas Flow (mL/Min):</b>	0.25

### *MRM transitions (m/z)*

#### **Pregnenolone sulfate**

395.02 > 96.5	Cone voltage: 50	Collisions energy: 35
395.02 > 79.6	Cone voltage: 50	Collisions energy: 90

#### **Internal Standard pregnenolone sulfate-d6:**

399,1 > 96,6	Cone voltage: 50	Collisions energy: 30
399,1 > 79,6	Cone voltage: 50	Collisions energy: 30

## 9. CALCULATION

$$\frac{\text{peak intensity sample} \times \text{calibrator concentration}}{\text{peak intensity internal standard of the sample}} \times F = \text{sample concentration}$$


$$F = \frac{\text{peak intensity internal standard of calibrator}}{\text{peak intensity calibrator}}$$

**Note:** The peak area can be used as an alternative to the peak intensity for analysis.


## 10. EXAMPLES OF CHROMATOGRAMS

*Examples of chromatograms for pregnenolone sulfate*


**Blank:**


**Calibrator:**

## Internal standard

## Sample:

## 11. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Normal range*

#### **EDTA plasma and serum**

11.8–115 ng/ml

We recommend each laboratory to establish its own norm concentration range.

## 12. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-Assay (n = 10)**

Sample	Pregnenolon sulfate [ng/ml]	CV [%]
1	17.5	7.6

#### **Inter-Assay (n = 6)**

Sample	Pregnenolon sulfate [ng/ml]	CV [%]
1	17.2	7.5
2	150.7	4.8

## 13. PRECAUTIONS

- For in vitro diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

## 14. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.

## 15. DISPOSAL

Mobile phases (MOPHA A, MOPHA B) dilution solution (DIL) and precipitation reagent (PREC) must be disposed as non-halogenated solvents.

## 16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

## 17. REFERENCES

1. Hu, Z Y, Bourreau, E, Jung-Testas, I, Robel, P, Baulieu, E E: Neurosteroids: oligodendrocyte mitochondria convert cholesterol to pregnenolone, *Proc Natl Acad Sci USA* **84** (1987), 8215-8219
2. Rupprecht, R, Holsboer, F: Neuroactive steroids: mechanism of action and neuro-psychological perspectives, *Trends Neurosci* **22** (1999), 410-416
3. Mehta, A K, Ticku, M K: An update on GABAA receptors, *Brain Res Brain Res Rev* **29** (1999), 196-217
4. Lambert, J J, Belelli, D, Hill-Venning, C, Peters, J A: Neurosteroids and GABAA receptor function, *Trends Pharmacol Sci* **16** (1995), 295-303
5. Park-Chung, M, Malayev, A, Purdy, R H, Gibbs, T T, Farb, D H: Sulfated and unsulfated steroids modulate  $\gamma$ -aminobutyric acidA receptor function through distinct sites, *Brain Res* **830** (1999), 72-87
6. Monaghan, E P, Navalta, L A, Shum, L, Ashbrook, D W, Lee, D A: Initial human experience with ganaxolone, a neuroactive steroid with antiepileptic activity, *Epilepsia* **38** (1997), 1026-1031
7. Weiss M & Hess S: Neuroaktive Steroide: Wirkungen und Risiken, *Pharmazie in unserer Zeit* **29. Jahrg** ( 2000) 6: 372-376
8. Jääntti SE, Tammimäki A, Raattamaa H, Piepponen P, Kostianen R, Ketola RA: De-termination of steroids and their intact glucuronide conjugates in mouse brain

by capillary liquid chromatography-tandem mass spectrometry. *Anal Chem.* **82** (2010) Apr 15; (8):3168-75

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by