# Vitamin B<sub>6</sub> HPLC Kit

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
# Vitamin B<sub>6</sub> HPLC Kit

For the determination of Vitamin B<sub>6</sub> in plasma and serum

Valid from 07.11.2013

**REF** KC2100

CAL
CTRL 1
CTRL 2

**Immundiagnostik AG**  
Stubenwald-Allee 8a  
64625 Bensheim, Germany

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>15</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>4. MATERIAL SUPPLIED</b>	<b>16</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>17</b>
<b>7. PRECAUTIONS</b>	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>18</b>
<b>9. ASSAY PROCEDURE</b>	<b>18</b>
<i>Procedural notes</i>	18
<i>Sample and standard preparation</i>	19
<i>Chromatographic conditions</i>	19
<b>10. TREATMENT OF THE COLUMN</b>	<b>20</b>
<b>11. RESULTS</b>	<b>20</b>
<i>Calculation</i>	20
<i>Typical chromatogram</i>	20
<b>12. QUALITY CONTROL</b>	<b>21</b>
<i>Expected values</i>	21
<i>Controls</i>	21
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>21</b>
<i>Precision and reproducibility</i>	21
<i>Linearity</i>	21
<i>Detection limit</i>	21
<i>Recovery</i>	21
<b>14. DISPOSAL</b>	<b>22</b>
<b>15. TROUBLESHOOTING</b>	<b>22</b>
<b>16. REFERENCES</b>	<b>23</b>
<b>17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>24</b>

## 1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of vitamin B<sub>6</sub> in plasma and serum.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin B<sub>6</sub> is the term for pyridoxin, pyridoxal and pyridoxamin and the appropriate phosphate products. All forms can be transformed into the active form pyridoxal-5-phosphate (here described as "vitamin B6").

Vitamin B6 functions as a coenzyme and is essential for more than 50 reactions in the protein metabolism thereby synthesizing, transforming or degrading amino acids. Vitamin B6 supports the resorption of amino acids and their transport into the cells. Furtheron vitamin B6 contributes to the synthesis of neuro transmitters and amine products (histamin).

Due to the fact that vitamin B6 contributes to a variety of different reactions lack of vitamin B6 results in various clinical pictures.

### Applications:

- determination of vitamin B<sub>6</sub> status
- homocysteinaemia
- skin diseases
- movement disorders
- Anaemia, depression

## 3. PRINCIPLE OF THE TEST

The first step in the determination of Vitamin B6 includes the sample preparation with additional derivatisation. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is used for derivatisation (20 min at 60°C) thereby transforming the vitamin B6 into a fluorescent product. The sample is cooled, centrifuged and injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a „reversed phase“ column; one run lasts 10 minutes. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak heights.

## Summary

The application of vitamin B<sub>6</sub> for HPLC makes it possible to determine the vitamin in an easy, fast and precise method. The kits includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As many other parameters the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test series the one-point calibration is much more economic than 6-point calibration for immuno assays).

## 4. MATERIAL SUPPLIED

Cat.No	Label	Kit Components	Quantity
KC2100LM	MOPHA	Mobile phase	1 x 1000 ml
KC2100KA	CAL	Calibrator (500 µl lyophilized; concentration, see product specification)	5 vials
KC2100FR	PREC	Precipitating reagent	1 x 5 ml
KC2100RE	RECSOL	Reconstitution solution	1 x 10 ml
KC2100DL	DER	Derivatisation solution	3 x 8,5 ml
KC2100KO	CTRL 1 CTRL 2	Control 1 and 2 (250 µl lyophilized; concentration, see product specification)	2 x 3 vials

HPLC column (KC2100RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C18-column
- Heatable shaker or water bath
- Vortex mixer

## 6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **CAL** (calibrator) in **500 µl RECSOL** (reconstitution solution). The concentration of vitamin B<sub>6</sub> might have minor changes from lot to lot.
- Reconstitute **CTRL 1, CTRL 2** (control 1 and 2) in **250 µl RECSOL** (reconstitution solution).
- **All other test reagents are stable at 2-8 °C**, the **CAL** (calibrator) and **CTRL 1, CTRL 2** (control 1 and 2) **at -20 °C**, up to the date of expiry stated on the label.

## 7. PRECAUTIONS

- Human materials used in kit components were tested for hepatitis-B (HBsAg), HCV, HIV-1 and HIV-2 antibodies and found to be negative. However, for safety reasons, all kit components should be treated as potentially infectious.
- The **PREC** (precipitating reagent) contains acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.
- As the **DER** (derivatisation solution) contains KCN, it should be pipetted under an extractor hood.

## 8. SPECIMEN COLLECTION AND PREPARATION

EDTA-plasma and serum are suitable for use in the assay.

Vitamin B6 is light- and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately.

The samples are stable in the dark at 2-8°C for 1 week. For longer storage samples should be frozen at -20°C.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

## 10. TREATMENT OF THE COLUMN


After the analysis the column should be flushed with 30 ml aqua bidest (1 ml/min) and stored in 50% methanol in aqua bidest (ca. 30 ml, flow 0.5 ml/min). Before use, the system should be equilibrated with ca. 30 ml Eluent.

## 11. RESULTS

### *Calculation*

$$\text{Concentration sample} = \frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height calibrator}}$$

### *Typical chromatogram*

## 12. QUALITY CONTROL

### *Expected values*

**Vitamin B<sub>6</sub> normal values:**                      **4.3 – 17.9 ng/ml**                      **(n=90)**

We recommend that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### *Controls*

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

<b>Intra-Assay CV:</b>	8.3 % (6.6 ng/ml)	[n = 6]
	6.0 % (15.1 ng/ml)	[n = 6]
<b>Inter-Assay CV:</b>	5.9 % (10.5 ng/ml)	[n = 13]
	5.9 % (27.2 ng/ml)	[n = 13]

### *Linearity*

up to 250 ng/ml

### *Detection limit*

0.8 ng/ml

### *Recovery*

98.3 %

## 14. DISPOSAL

The DER (derivatisation solution) must be oxidized with hydrogen peroxide, the pH value adjusted to 6-8, and disposed as aqueous salt solution. The MOPHA (mobile phase) and the PREC (precipitation solution) must be neutralized with NaOH to neutral pH and disposed as salt solution.

**Important:** Reaction will produce heat, be careful! Please refer to the appropriate national guidelines.

## 15. TROUBLESHOOTING

Problem	Possible reasons	Solution
No signal	No or defect connection to evaluation system.	Check signal cord and connection.
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating	Injector dirty	Clean injector
peaks	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min

<b>Problem</b>	<b>Possible reasons</b>	<b>Solution</b>
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not calm	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

## 16. REFERENCES

Ambrosch A. et al. (2000). Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic Med* 18; 185-192.

Dierkes J. et al. (2001). Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis* 158; 161-164.

Dierkes et al. (2001). Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *J Renal Nut* 11; 67-72.

## 17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

### Used Symbols:

Store at

Catalog Number

In Vitro Diagnostic Device

No. of tests

Manufacturer

Use by

Chargennummer

Controls

Calibrator

**Immundiagnostik AG**

Stubenwald-Allee 8a  
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)