# 25-OH-Vitamin D<sub>3</sub>/D<sub>2</sub> LC-MS/MS Kit

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## 19. LITERATUR

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D<sub>3</sub> in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone-

rone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009).

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

*In-Vitro*-Diagnostikum

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Hersteller

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
# 25-OH Vitamin D<sub>3</sub>/D<sub>2</sub> LC-MS/MS Kit

*For the determination of 25-OH vitamin D<sub>3</sub> and  
25-OH vitamin D<sub>2</sub> in plasma and serum*

Valid from 2015-02-04

**REF** KM5000

**IVD** **CE**

**Immundiagnostik AG**

# Table of Contents

<b>1. INTENDED USE</b>	<b>19</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>19</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>19</b>
<b>4. MATERIAL SUPPLIED</b>	<b>20</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>23</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>23</b>
<b>7. SPECIMEN COLLECTION AND PREPARATION</b>	<b>24</b>
<b>8. ASSAY PROCEDURE</b>	<b>24</b>
<i>Sample preparation 1</i>	24
<i>Sample preparation 2</i>	25
<i>Sample preparation 3</i>	26
<b>9. CHROMATOGRAPHIC CONDITIONS</b>	<b>27</b>
<b>10. MS/MS METHOD (LISTED AS AN EXAMPLE FOR A WATERS QUATTRO PREMIER XE TANDEM MASS SPECTROMETER)</b>	<b>27</b>
<i>MRM transitions (m/z)</i>	28
<b>11. RESULTS</b>	<b>28</b>
<b>12. LIMITATIONS</b>	<b>28</b>
<b>13. QUALITY CONTROL</b>	<b>29</b>
<i>Reference ranges for 25(OH) vitamin D<sub>3</sub></i>	29
<b>14. PERFORMANCE CHARACTERISTICS</b>	<b>30</b>
<i>Precision and reproducibility</i>	30
<i>Sensitivity</i>	31
<b>15. DISPOSAL</b>	<b>31</b>
<b>16. PRECAUTIONS</b>	<b>31</b>
<b>17. TECHNICAL HINTS</b>	<b>31</b>
<b>18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>31</b>
<b>19. REFERENCES</b>	<b>32</b>

## 1. INTENDED USE

The reversed phase HPLC application is intended for the quantitative determination of 25-OH vitamin D<sub>3</sub> and 25-OH vitamin D<sub>2</sub> in serum and plasma. This assay is designed for in vitro diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

D-vitamins and calciferols arise from provitamins by the UV radiation of sunlight catalysed splitting of the B ring of the steran backbone. The two most important D vitamins are vitamin D<sub>3</sub> and vitamin D<sub>2</sub>. Contrary to vitamin D<sub>2</sub>, which has to be added via food, vitamin D<sub>3</sub> can be produced in the liver.

In the skin formed or together with vitamin D<sub>2</sub> by food ingested vitamin D<sub>3</sub> is bound to a vitamin D binding protein in the plasma, transported into the liver and hydroxylated in position 25 to form 25-OH-D. More than 95% of 25-OH-D is 25-OH-D<sub>3</sub>. 25-OH-D<sub>2</sub> is only detectable in patients with medication of vitamin D<sub>2</sub>.

## 3. PRINCIPLE OF THE TEST

For the determination of 25-OH vitamin D<sub>3</sub> and 25-OH vitamin D<sub>2</sub> samples, we developed three sample extraction methods depending on your LC-MS/MS-system (sample preparation 1, 2, 3).

### Sample preparation 1:

Suitable for a LC-MS/MS-system with high sensitivity.

### Sample preparation 2 and 3:

Suitable for all LC-MS/MS-systems. The test procedure is a SPE via 1 ml cartridges for test procedure 2 and for test procedure 3 a SPE via 96-well-plate.

## 4. MATERIAL SUPPLIED

### For sample preparation 1:

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by *protein precipitation*

Cat. No	Content	Kit Components	Quantity
KM5000LA	MOPHA A	Mobile Phase A	1000 ml
KM5000LB	MOPHA B	Mobile Phase B	1000 ml
KM5000KA	CAL	Calibrator (lyoph. 1,2 ml; concentration is given on the label)	5 vials
KM5000KO	CTRL1 CTRL2	Control 1 and 2 (lyoph. 0,6 ml; for concentration see product specification)	2 x 5 vials each
KM5000RE	RECSOL	Reconstitution solution	10 ml
KM5000AC	ACTSOL	Activation reagent	2,5 ml
KM5000FR	PREC	Precipitation reagent (-20 °C)	45 ml
KM5000SO	SOL A	Solution A	2 x 50 ml
KM5000IS	INT STD	Internal Standard	1,1 ml

For the preparation of the 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> application, tuning is necessary to find the ideal LC-MS/MS settings. The tuning solution (KM5000TU, 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> and INT STD as pure substances), the column (KM5000RP) and all other individual components can be ordered separately. Please ask for our single component price list.

**For sample preparation 2:**

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by purification via cartridges (Waters Oasis HLB 30 mg, 1cc)

Cat. No	Content	Kit Components	Quantity
KM5000LA	MOPHA A	Mobile Phase A	1000 ml
KM5000LB	MOPHA B	Mobile Phase B	1000 ml
KM5000KA	CAL	Calibrator (lyoph. 1,2 ml; concentration is given on the label)	5 vials
KM5000KO	CTRL1 CTRL2	Control 1 und 2 (lyoph. 0,6 ml; for concentration see product specification)	2x 5 vials each
KM5000RE	RECSOL	Reconstitution solution	10 ml
KM5000AC	ACTSOL	Activation reagent	2,5 ml
KM5000FR2	PREC	Precipitation reagent (-20 °C)	55 ml
KM5000SO2	SOL A	Solution A	25 ml
KM5000WL2	WASHSOL	Washing solution	160 ml
KM5000EL2	ELUSOL	Elution solution (contains acetonitril)	210 ml
KM5000IS	INT STD	Internal Standard	1,1 ml

For the preparation of the 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> application, tuning is necessary to find the ideal LC-MS/MS settings. The tuning solution (KM5000TU, 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> and INT STD as pure substances), the column (KM5000RP), Oasis® cartridges HLB 30 mg, 1cc (KM5000CK) and all other individual components can be ordered separately. Please ask for our single component price list.

**For sample preparation 3:**

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by *purification via 96-well plate (Oasis HLB  $\mu$ Elution-plate)*

**Attention:** number of determinations is limited to 96.

Cat. No	Content	Kit Components	Quantity
KM5000LA	MOPHA A	Mobile Phase A	1000 ml
KM5000LB	MOPHA B	Mobile Phase B	1000 ml
KM5000KA	CAL	Calibrator (lyoph. 1,2 ml; concentration is given on the label)	5 vials
KM5000KO	CTRL1 CTRL2	Control 1 und 2 (lyoph. 0,6 ml; for concentration see product specification)	2x 5 vials each
KM5000RE	RECSOL	Reconstitution solution	10 ml
KM5000AC	ACTSOL	Activation reagent	2,5 ml
KM5000FR3	PREC	Precipitation reagent (-20 °C)	20 ml
KM5000SO3	SOL A	Solution A	10 ml
KM5000WL3	WASHSOL	Washing Solution	20 ml
KM5000EL3	ELUSOL	Elution solution (contains acetonitril)	10 ml
KM5000IS	INT STD	Internal Standard	1,1 ml

For the preparation of the 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> application, tuning is necessary to find the ideal LC-MS/MS settings. The tuning solution (KM5000TU, 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> and INT STD as pure substances), the column (KM5000RP), Oasis® 96-well- $\mu$ Elution plate HLB (KM5000PL) and all other individual components can be ordered separately. Please ask for our single component price list.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Centrifuge (4 °C)
- Vortex mixer
- 1.5 ml reaction tubes (Eppendorf)
- Various pipettes
- LC-MS/MS equipment
- Glass tubes, LC-MS/MS suitable
- 500 ml graduated cylinder, LC-MS/MS suitable
- UPLC column, e. g. Acquity BEH C18, 1.7 µm

### For sample preparation 2 you also need:

- Methanol p. a.
- Ultra pure water\*
- Oasis HLB cartridges 30 µm, 1cc (KM5000CK)
- Sample evaporation unit

### For sample preparation 3 you also need:

- Methanol p. a.
- Ultra pure water\*
- Oasis HLB 96-well µElution plate (KM5000PL)
- Sample evaporation unit for µElution plate

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 6. PREPARATION AND STORAGE OF REAGENTS

Internal Standard (**INT STD**) and precipitation reagent (**PREC**) are stable at **-20 °C** until the expiry date stated on the label.

### Mobile phase and solution A

Before use, a 0,1 % activation reagent (ACTSOL) must be added to the mobile phases (MOPHA A, MOPHA B) and solution A (SOL A):

e. g. 500 ml MOPHA + 500 µl ACTSOL

The prepared solutions can be used within 2 weeks. For this reason, it is recommended to prepare only the desired amount necessary for each assay.

**WARNING:** The activation reagent (ACTSOL) must be added under the fume hood. All vials to be used must be absolutely clean, detergent-free and preferably made of a LC-MS/MS suitable glass.

## Calibrators and controls

- The **lyophilized calibrator** (CAL), plasma with a specific concentration of 25-OH vitamin D<sub>3</sub> and 25-OH vitamin D<sub>2</sub>, must be stored at **-20 °C** until use. The content of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> changes slightly from lot to lot, the exact concentration is stated on the label. Before use, the calibrator has to be **reconstituted in 1.2 ml reconstitution solution** (RECSOL). Reconstituted calibrator **cannot be stored**.
- The **lyophilized controls** (CTRL1, CTRL2) must be stored at **-20 °C** until use. Before use, they have to be **reconstituted** in each **0.6 ml of RECSOL**. Reconstituted controls **cannot be stored**.

All other test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 7. SPECIMEN COLLECTION AND PREPARATION

Serum and plasma could be used. The samples should be stored at -20 °C until testing. Samples are stable for at least four weeks when stored at -20 °C.

## 8. ASSAY PROCEDURE

### Sample preparation 1

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by *protein precipitation*.

1.	Pipet <b>200 µl</b> of calibrator ( <b>CAL</b> ), controls ( <b>CTRL1, CTRL2</b> ) or sample in 1.5 ml reaction tubes
2.	Add <b>10 µl</b> of internal standard ( <b>INT STD</b> ) and mix
3.	Add <b>400 µl</b> of precipitation reagent ( <b>PREC</b> ) (ice-cold)
4.	Vortex for <b>1 min</b>
5.	Incubate <b>10 min</b> at <b>4 °C</b>
6.	Centrifuge for <b>10 min</b> at 10 000 rpm at <b>4 °C</b>
7.	Dilute supernatant <b>1:10</b> with <b>activated solution A</b> (SOL A) (e. g. 100 µl supernatant + 900 µl SOL A)
8.	Inject <b>50 µl</b> into the LC-MS/MS system

### Sample preparation 2

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by purification via cartridges (Waters Oasis HLB 30 mg, 1cc).

1.	Equilibrate cartridges with 1 ml methanol and 1 ml ultra pure water, respectively
2.	Pipet <b>500 µl</b> of calibrator ( <b>CAL</b> ), controls ( <b>CTRL1</b> , <b>CTRL2</b> ) or <b>sample</b> in 1.5 ml reaction tubes
3.	Add <b>10 µl</b> of internal standard ( <b>INT STD</b> ) and mix
4.	Add <b>500 µl</b> of precipitation reagent ( <b>PREC</b> ) (ice-cold)
5.	Vortex for <b>1 min</b>
6.	Incubate <b>10 min</b> at <b>4 °C</b>
7.	Centrifuge for <b>10 min</b> at 10 000 rpm at <b>4 °C</b>
8.	Pipet <b>500 µl</b> of <b>supernatant</b> in cartridges
9.	Wash <b>3 x</b> with <b>500 µl</b> washing solution ( <b>WASHSOL</b> )
10.	Elute <b>2 x</b> with <b>1 ml</b> elution solution ( <b>ELUSOL</b> )
11.	Evaporate eluent
12.	Pipet <b>200 µl</b> of <b>activated solution A</b> (SOL A) and vortex for 1 min
13.	Inject <b>50 µl</b> into the LC-MS/MS system

### Sample preparation 3

Determination of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub> by *purification via 96-well-plate (Oasis HLB μElution plate)*

1.	Equilibrate μElution plate with 0.2 ml methanol and 0.2 ml ultra pure water, respectively 0.2 ml water, respectively
2.	Pipet <b>200 μl</b> of calibrator ( <b>CAL</b> ), controls ( <b>CTRL1, CTRL2</b> ) or <b>sample</b> in 1.5 ml reaction tubes
3.	Add <b>10 μl</b> of internal standard ( <b>INT STD</b> ) and mix
4.	Add <b>200 μl</b> of precipitation reagent ( <b>PREC</b> ) (ice-cold)
5.	Vortex for <b>1 min</b>
6.	Incubate <b>10 min</b> at <b>4 °C</b>
7.	Centrifuge for <b>10 min</b> at 10 000 rpm at <b>4 °C</b>
8.	Pipet <b>200 μl</b> of <b>supernatant</b> in μElution wells
9.	Wash <b>3 x</b> with <b>200 μl</b> washing solution ( <b>WASHSOL</b> )
10.	Elute with 100 μl elution solution ( <b>ELUSOL</b> )
11.	Add <b>100 μl</b> of <b>activated solution A</b> (SOL A) to the eluent
12.	Inject <b>50 μl</b> into the LC-MS/MS system

## 9. CHROMATOGRAPHIC CONDITIONS

Column material: e.g. Acquity BEH C18; 1.7 µm

Column dimension: 2.1 x 50 mm

Flow: 0.3 ml/min

Column temperature: 45 °C

Inject volume: 50 µl

Run time: 5 min

Gradient:

Time	% A	% B
0 min	100	0
2 min	0	100
2,6 min	0	100
2,9 min	100	0
4,0 min	100	0

After analysis, the separation column can be stored in the mobile phase, tightly closed.

## 10. MS/MS METHOD (LISTED AS AN EXAMPLE FOR A WATERS QUATTRO PREMIER XE TANDEM MASS SPECTROMETER)

Mode: MRM

Polarity: ESI+

Capillary (kV): 3

Cone (V): var.

Extracor (V): 4

RF Lens (V): 0

Source temperature (°C): 130

Desolvation Temperatur (°C): 450

Cone Gas Flow (l/h): 50

Desolvation Gas Flow (l/h): 950

Collision Gas Flow (ml/min): 0,15

*MRM transitions (m/z)***25-OH vitamin D<sub>3</sub>**

401.5 > 383.2	Cone voltage: 22	Collision energy: 10
401.5 > 159.0	Cone voltage: 22	Collision energy: 28

**25-OH vitamin D<sub>2</sub>**

413.4 > 395.3	Cone voltage: 20	Collision energy: 10
413.4 > 83.1	Cone voltage: 20	Collision energy: 25

**Internal standard (25-OH-vitamin-D<sub>3</sub>-d6)**

407.5 > 389.4	Cone voltage: 22	Collision energy: 10
407.5 > 159.1	Cone voltage: 22	Collision energy: 28

**11. RESULTS****Calculation**

$$\frac{\text{Peak height sample} \times \text{Calibrator concentration}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration of the sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

**Tip:** Alternatively, the peak area instead of the peak height can be used for quantification.

**Conversion factors:**

*for 25-OH vitamin D<sub>3</sub>*

1 ng/ml = 2,5 nmol/l

1 nmol/l = 0,4 ng/ml

*for 25-OH vitamin D<sub>2</sub>*

1 ng/ml = 2,42 nmol/l

1 nmol/l = 0,412 ng/ml

**12. LIMITATIONS**

This test cannot be used with whole blood samples.

### 13. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

#### *Reference ranges for 25(OH) vitamin D<sub>3</sub>*

##### **Information from ASBMR 2011**

Deficiency (seriously deficient)	< 20 ng/ml	or	< 50 nmol/l
Insufficiency (deficient)	20–29 ng/ml	or	50–74 nmol/l
Sufficiency (adequately supplied)	> 30 ng/ml	or	> 75 nmol/l

Society of Osteologists SACHSEN E. V.

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#### **Reference intervals for 25(OH)-vitamin D<sub>3</sub> (ng/ml), males and females\***

Age	n	2.5%	97.5%
0 to < 3 months	131	5	42
3 to < 6 months	135	9	60
6 months to < 1 year	147	18	58
1 to < 3 years	394	15	54
3 to < 10 years	619	14	46
10 to < 13 years	286	11	50
13 to < 15 years	275	10	44
15 to < 18 years	390	8	45
> 18 years	421	8	56

\*All seasons

(Soldin et al., 2009)

We recommend each laboratory to establish its own reference range.

#### **Note**

The vitamin D production in the skin is high variable and depends on the season and daily time, degree of latitude, age, sun protection etc. The normal ranges depend on

the method used (e.g. vitamin D release from the vitamin D binding protein, VDBP) and serve only as orientation.

### Literature references

The following literature references for the vitamin D reference values can be found in the references on page 28: Grant et al., Soldin et al., Visser et al., Wicherts et al.

## 14. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### 25-OH vitamin D<sub>3</sub>

##### *Intra-Assay (n = 6)*

Sample	25-OH vitamin D <sub>3</sub> [nmol/l]	CV [%]
1	84.4	2.9

##### *Inter-Assay (n = 18)*

Sample	25-OH vitamin D <sub>3</sub> [nmol/l]	CV [%]
1	83.4	6.6

#### 25-OH vitamin D<sub>2</sub>

##### *Intra-Assay (n = 6)*

Sample	25-OH vitamin D <sub>2</sub> [nmol/l]	CV [%]
1	43.9	7.5

##### *Inter-Assay (n = 18)*

Sample	25-OH vitamin D <sub>2</sub> [nmol/l]	CV [%]
1	48.2	12.1

### *Sensitivity*

Detection limit of 25-OH vitamin D<sub>3</sub>: 0,64 nmol/l

Detection limit of 25-OH vitamin D<sub>2</sub>: 5,4 nmol/l

It should be noted that the determination of the detection limit depends not only on the application method but also on the instrument.

## **15. DISPOSAL**

The mobile phases (MOPHA A, MOPHA B), solution A (SOL A), precipitation reagent (PREC), activation reagent (ACTSOL), washing solution (WASHSOL) and elution reagent (ELUSOL) must be disposed as non-halogenated solvent. Please refer to the appropriate national guidelines.

## **16. PRECAUTIONS**

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.

## **17. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date stated on kit label.
- The assay should always be performed according the enclosed manual.

## **18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.

- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 19. REFERENCES

1. Chapuy, M. C. et al. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *The Journal of clinical endocrinology and metabolism* **81**, 1129–33 (1996).
2. Scharla, S. H. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporosis international* **8** Suppl 2, S7–12 (1998).
3. Grant, W. B. & Holick, M. F. Benefits and requirements of vitamin D for optimal health: a review. *Alternative medicine review : a journal of clinical therapeutic* **10**, 94–111 (2005).
4. Visser, M., Deeg, D. J. H., Puts, M. T. E., Seidell, J. C. & Lips, P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *The American journal of clinical nutrition* **84**, 616–22; quiz 671–2 (2006).
5. Wicherts, I. S. et al. Vitamin D status predicts physical performance and its decline in older persons. *The Journal of clinical endocrinology and metabolism* **92**, 2058–65 (2007).
6. Wielders, J. P. M. & Wijnberg, F. A. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clinical chemistry* **55**, 1584–5 (2009).
7. Soldin, O. P., Sharma, H., Husted, L. & Soldin, S. J. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clinical biochemistry* **42**, 823–7 (2009)

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for &lt;n&gt; tests

Lot number

Use by