Vitamin K₁ HPLC Kit

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Hersteller

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
Vitamin K₁ HPLC Kit

For the determination of Vitamin K₁ in plasma and serum

Valid from 03.01.2013

REF

KC 2400

CAL
INT STD
CTRL 1
CTRL 2
STD

Immundiagnostik AG

IVD

Content

1. INTENDED USE	17
2. INTRODUCTION	17
3. PRINCIPLE OF THE TEST	17
4. MATERIAL SUPPLIED	18
5. MATERIAL REQUIRED BUT NOT SUPPLIED	19
6. PREPARATION AND STORAGE OF REAGENTS	19
7. PRECAUTIONS	20
8. SPECIMEN COLLECTION AND PREPARATION	20
9. ASSAY PROCEDURE	21
Procedural notes	21
Sample and standard preparation	21
Chromatographic conditions	22
10. TREATMENT OF THE COLUMN	22
11. RESULTS	22
Calculation	22
Typical chromatogram	23
12. LIMITATIONS	23
13. QUALITY CONTROL	23
Reference value	23
Controls	23
14. PERFORMANCE CHARACTERISTICS	24
Precision and reproducibility	24
Linearity	24
Detection limit	24
15. DISPOSAL	24
16. TROUBLESHOOTING	25
17. REFERENCES	26
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26

1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of Vitamin K₁ in plasma and serum. This Assay is designed for *in vitro* diagnostic use only.

2. INTRODUCTION

Vitamin K₁ is a derivative of 2-methyl-1, 4-naphtochinone and is found in green plants. Vitamin K₁ is insoluble in water and readily soluble in ether, n-hexane and chloroform.

It's an essential co-factor in the posttranslational carboxylation reaction of glutamic acid residues (GLU) to γ -carboxyglutamic acid residues (GLA) in a number of blood clotting factors and also in some other proteins e.g. osteocalcin. The adjacent carboxyl groups of the GLA-residues provide the Vitamin K dependent proteins with characteristic calcium- and phospholipid-binding properties that are essential for their activation and function. A decrease in Vitamin K₁ is reported in osteoporotic patients. Other clinical symptoms of Vitamin K₁ deficiency are clotting disorders, which manifest themselves as bleeding in the skin, in the mucous membranes, in muscles and in internal organs. Symptoms of deficiencies normally appear within few days. Vitamin K₁ is rapidly metabolized and only minor amounts are stored in the organism.

Applications:

- Determination of Vitamin K₁ status
- Vitamin K deficiency induced by:
 - obstructive liver disease
 - obstructive icterus
 - malabsorption due to celiac disease,
 - pancreatitis, diarrhea, antibiotic abuse
- Blood clotting disorder
- Bone metabolism disorders
- Haemorrhagic disorders of newborns

3. PRINCIPLE OF THE TEST

The newly developed Vitamin K HPLC-application is the first commercially available kit. After a solid phase extraction on SPE-cartridges, serum or plasma samples are precipitated. The supernatant is then extracted with an organic solvent and evaporated. After resuspension the sample is measured in an isocratic HPLC-system. A post-column reduction reactor reduces Vitamin K and enables the measurement of Vitamin K with a fluorescence detector. An internal standard is added before the solid phase extraction step to ensure the high quality of the measurement.

Summary:

The application of Vitamin K₁ for HPLC makes it possible to determine the Vitamin in an easy, fast and precise way. The kits includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns. As with many other parameters the advantage of HPLC measurements are the simultaneous handling of many analytes in a single test. The „complete HPLC-system“ enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test runs the one-point calibration is much more economic than 6-point calibration for immuno assays).

4. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
KC2400LM	MOPHA	Mobile phase	3 x 1000 ml
KC2400KA	CAL	Calibrator, lyophilized	8 vials
KC2400IL	STD	Isopropanolic standard	10 ml
KC2400IS	INT STD	Internal Standard	1 ml
KC2400FR	PREC	Precipitating reagent	200 ml
KC2400EX	EXTSOL	Extraction solution	400 ml
KC2400ZI	ZINC	Zinc (Important: store under argon)	20 g
KC2400ZN	ZUB	Accessories for post-column reduction-reactor	5 pieces
KC2400KO	CTRL 1 CTRL 2	Control 1 and 2; 1.1 ml lyophilized	2 x 3 vials

HPLC column (KC2400RP), post-column reduction-reactor (KC2400NR), SPE cartridges (KC2400CK) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Glass tubes for centrifugation, V-bottom (10 ml)
- SPE cartridges, C18 (KC 2400CK)
- MERCK-Cartridge Holder Manu-Kart (ordering no. KC2400RK)
- post-column reduction-reactor (KC2400NR)
- Evaporator
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C18-column – Superspher RP 18,4 µm, 125 x 4,6 mm (KC2400RP)
- Vortex mixer

6. PREPARATION AND STORAGE OF REAGENTS

Preparation of the calibrator and controls


- Reconstitute **CAL** (calibrator) in aqua bidest. The volume is given on the label. One vial is for a single use only; discard the material, which have not been used. The content of Vitamin K₁ might have minor changes from lot to lot.
- Reconstitute **CTRL 1** and **CTRL 2** (control 1 and 2) in 1.1 ml aqua bidest.
- All other test reagents are stable at 2-8 °C, up to the date of expiry stated on the label.

Preparation of the post-column reduction reactor

For getting well detectable peaks it is necessary to exchange the zinc particles in the post-column reduction reactor each day. The zinc particles are used up after 12 h by oxidation. Filling the column is very easy and takes just 10 min of time.

For assembling the column see the following figure.

1. Cap nut
2. Stainless steel inlet
3. PTFE seal
4. Stainless steel sieve (grey)
5. Glass fiber sieve (3 pieces, white)
6. Stainless steel sieve (grey)
7. column tube

1. Close one side of the column according to the figure above.
2. Fill in the zinc-particles with a funnel while knocking the column slightly on the table, so that the packing will not show any cavities.
3. Close the upper side of the column.

The post-column reduction reactor should be mounted in the HPLC-system as described by the following picture:

Note: The performance of the reactor must be checked (see Sample and standard preparation, Point 10).

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The supplied reagents contain different dangerous chemical reagents like methanol (mobile phase), ethanol (precipitating reagent), n-hexan (extraction solution), isopropanol and acid (mobile phase). Although diluted, it still must be handled with care and should only be handled with gloves, eye protection, and appropriate protective clothing. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Plasma and serum can both be used for analysis. The samples must be cooled immediately.

The samples are stable at 2-8°C for 1 week. For longer storage samples should be frozen at -20°C.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

Sample and standard preparation

1. The cartridge is rinsed with 3 ml methanol and 3 ml aqua bidist.
2. Add 10 µl INT STD (internal standard) to 1 ml patient sample, CAL (calibrator) or CTRL 1 and CTRL 2 (control 1 and 2) and pipette the mixture on the SPE cartridge. Let it soak through by vacuum.
3. Collect the break-through in a glass vial.
4. Add 2 ml PREC (precipitation reagent), vortex for 1 min and centrifuge for 10 min at 3.500 x g.
5. Pipette the supernatant in a fresh glass vial and add 4 ml of EXTSOL (extraction solution).
6. Vortex for 2 min and centrifuge for 5 min at 3.500 x g.
7. Pipette the upper phase in a new glass vial and evaporate to dryness.
8. The dried sample is stable for 8 days at 4-8°C
9. Connect the post-column reduction reactor in the HPLC-system, as described above and wait for equilibration (30-45 min).
10. Check the performance of the reactor by the injection of 100 µl of STD (isopropanolic standard) and determine the signal to noise ratio , which should be greater than 25 .
11. Resuspend the dried sample in 150 µl MOPHA (mobile phase) and inject 100 µl in the HPLC-system.

Chromatographic conditions

For determination of the retention time, inject **100 µl** of the **STD** (isopropanolic standard) in the HPLC-system.

Column material :	Superspher 100 RP 18; 4 µm
Column dimension:	125 mm x 4.6 mm
Flow rate:	0.9 - 1.2 ml/min
Fluorescence detector:	ex.: 248 nm em.: 418 nm
Temperature:	30°C
Running time:	approx. 20 min

Immundiagnostik recommends to use a guard-column to enlarge lifetime of the column.

10. TREATMENT OF THE COLUMN

After analysis the column should be flushed with 30 ml aqua bidest (1,0 ml/min) and stored in 50% methanol in aqua bidest (approx. 30 ml, flow 0.5 ml/min). Before use, the system should be equilibrated with ca. 50 ml MOPHA (mobile phase).

Important: Do not re-circulate the MOPHA (mobile phase) in this test system.


11. RESULTS

Calculation

$$\frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

Typical chromatogram

12. LIMITATIONS

We recommend not to measure hemolytic and lipaemic patient samples.

13. QUALITY CONTROL

Reference value

0.22 - 2.28 ng/ml (Mean value = 1.29 ng/ml) (n = 19)

We recommend that each laboratory should develop its own reference range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV: 4.5 % (1.7 ng/ml) [n = 6]

Inter-Assay CV: 5.6 % (1.7 ng/ml) [n = 6]

Linearity

up to 25 ng/ml

Detection limit

0.15 ng/ml

15. DISPOSAL

MOPHA (mobile phase), **STD** (isopropanolic standard), **INTSTD** (internal standard,) **EXT-SOL** (extraction solution) and **PREC** (precipitating reagent) must be disposed as non-halogenated solvent. Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reasons	Solution
No signal	No or defect connection to evaluation system.	Check signal cord and connection.
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

17. REFERENCES

- Saupe J. (1992). Vitamin K und Indikationen zur Bestimmung im Plasma. *Klin. Lab.* 38: 51-54
- Szulc P, Arlot M., Chapuy M.C., Duboef F., Meunier P.J., Delmas P.D. (1994). Serum undercarboxylated osteocalcin correlates with hip bone mineral density in elderly woman. *J Bone Miner Res* 9:1591-95
- Hodges S.J., Akesson K., Vergnaud P., Obrant K., Delmas P.D. (1993). Circulating levels of vitamin K1 and K2 decreased in elderly women with hip fracture. *J Bone Miner Res* 8:1241-45
- Haroon Y, Bacon DS, Sadowski JA. (1986) Liquid-chromatographic determination of vitamin K1 in plasma, with fluorometric detection. *Clin Chem.* Oct;32(10):1925-9.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.

Used Symbols:

Store at

Catalog Number

In Vitro Diagnostic Device

No. of tests

Manufacturer

Use by

Lot number

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com