Zink-Protoporphyrin/ Protoporphyrin HPLC Kit

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TABLE OF CONTENTS

1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	16
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
CHROMATOGRAPHIC CONDITIONS	17
10. TREATMENT OF THE COLUMN	18
11. RESULTS	18
CALCULATION	18
TYPICAL CHROMATOGRAM	18
12. QUALITY CONTROL	19
EXPECTED VALUES	19
CONTROLS	19
13. PERFORMANCE CHARACTERISTICS	19
PRECISION AND REPRODUCIBILITY	19
14. DISPOSAL	19
15. TROUBLESHOOTING	20
16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	21

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
Manual

Zinc-Protoporphyrin/ Protoporphyrin HPLC Kit

*For the determination of Zinc-Protoporphyrin/Protoporphyrin in
EDTA whole blood*

Valid from 2015-01-09

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of Zinc-Protoporphyrin and Protoporphyrin in EDTA whole blood. This Assay is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Zinc protoporphyrin (ZPP) is a metabolite formed in erythrocytes during hemoglobin synthesis. By iron deficiency during erythropoiesis, instead of incorporating a ferrous ion to form the hem precursor, zinc becomes an alternative metal ion in the protoporphyrin complex. As a result, zinc protoporphyrin is produced instead of hemoglobin. As a consequence, in cases of iron deficiency anemia the zinc protoporphyrin concentration in the erythrocytes is elevated.

Indications

- Lead poisoning
- Iron deficiency
- Sickle cell anemia
- Sideroblastic anemia
- Anemia of chronic disease
- Vanadium exposure

3. PRINCIPLE OF THE TEST

The first step in the determination of Zinc-Protoporphyrin and Protoporphyrin is precipitation of the higher molecular components. After their removal by centrifugation, the supernatant is injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a reversed phase column. One run lasts 5 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas by the external standard method.

Parallel measurements of the hemoglobin concentration are required, because the Zinc-Protoporphyrin concentration is related to the hemoglobin concentration.

Summary

This HPLC technique provides an easy, fast and precise method for quantitative determination of Zinc-Protoporphyrin and Protoporphyrin. Except the column, the kit contains all reagents necessary for sample preparation and separation in ready-to-use form.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
KC2700LM	MOPHA	Mobile phase	1000 ml
KC2700KA	CAL	Calibrator (lyophilized, 500 µl)	5 vials
KC2700IS	INT STD	Internal Standard (lyophilized, 2 ml)	4 x 1 vials
KC2700RE	RECSOL	Reconstitution solution	7 ml
KC2700FR	PREC	Precipitation reagent	150 ml
KC2700KO	CTRL 1 CTRL 2	Control 1 and 2 (250 µl lyophilized; concentration, see specification)	2 x 3 vials

HPLC column (KC2700RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C₁₈-column
- Vortex Mixer

6. PREPARATION AND STORAGE OF REAGENTS

- The **calibrator** (CAL) is stable at -20 °C until the expiry date stated on the label. Before use, the calibrator has to be **reconstituted**: use 500 µl of the included reconstitution solution (RECSOL) per vial. Unused leftovers are not stable and have to be discarded. The concentration of Zinc-Protoporphyrin might slightly change in between different lots, the exact amount is stated on the label.
- The **controls** (CTRL1, CTRL2) are stable at -20 °C until the expiry date stated on the label. Before use, the controls have to be **reconstituted**: use 250 µl of the included reconstitution solution (RECSOL) per vial. Unused leftovers are not stable and have to be discarded.
- The **internal standard** (INT STD) is stable at -20 °C until the expiry date stated on the label. Before use, the internal standard has to be **reconstituted**: use 2 ml of the included precipitation solution (PREC) per vial. Unused leftovers are not stable and have to be discarded.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at 2–8 °C.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA whole blood is suited for this analysis. Before analysis, the erythrocytes are lysed by freezing and thawing in order to release Zinc-Protoporphyrin and Protoporphyrin.

Zinc-Protoporphyrin as well as Protoporphyrin is quite light sensitive. Samples should be transported light-protected.

10. TREATMENT OF THE COLUMN

After analysis, the column should be flushed with 30 ml aqua bidest. (1 ml/min) and stored in 50% methanol in aqua bidest (approx. 30 ml, flow 0.5 ml/min). Before use, the system should be equilibrated with ca. 30 ml eluent.


11. RESULTS

Calculation

$$\frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

Typical chromatogram

12. QUALITY CONTROL

Expected values

Zinc-Protoporphyrin: < 40 µmol/mol haem

We recommend each laboratory to establish its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV:

Zinc-Protoporphyrin	3.03 % (349 nmol/l)	[n = 10]
Protoporphyrin	5.54 % (198 nmol/l)	[n = 10]

Inter-Assay CV:

Zinc-Protoporphyrin	9.28 % (310 nmol/l)	[n = 8]
Protoporphyrin	9.57 % (79 nmol/l)	[n = 7]

14. DISPOSAL

The MOPHA (mobile phase) and PREC (precipitation reagent) must be disposed as non-halogenated solvents.

Please refer to the appropriate national guidelines.

15. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Avoid contact with skin or mucous membranes.
- For *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		