



For Professional Use Only

AmpliSens[®] *hRSV-FEP*

PCR kit

Instruction Manual

AmpliSens[®]



Ecoli s.r.o., Studenohorska 12
841 03 Bratislava 47
Slovak Republic
Tel.: +421 2 6478 9336
Fax: +421 2 6478 9040



Federal Budget Institute of
Science "Central Research
Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia

TABLE OF CONTENTS

1. INTENDED USE	3
2. PRINCIPLE OF PCR DETECTION	3
3. CONTENT	3
4. ADDITIONAL REQUIREMENTS	4
5. GENERAL PRECAUTIONS.....	5
6. SAMPLING AND HANDLING	5
7. WORKING CONDITIONS.....	7
8. PROTOCOL	7
9. DATA ANALYSIS	8
10. TROUBLESHOOTING.....	9
11. TRANSPORTATION.....	9
12. STABILITY AND STORAGE.....	9
14. REFERENCES	10
15. QUALITY CONTROL.....	10
16. KEY TO SYMBOLS USED	11

1. INTENDED USE

AmpliSens® hRSV-FEP PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of *Human Respiratory Syncytial Virus* RNA in the clinical materials (nasal and oropharyngeal swabs, nasopharyngeal or tracheal sputum or aspirate, and autopsy material) by using end-point hybridization-fluorescence detection of amplified products.



The results of PCR analysis are taken into account in complex diagnostics of disease.

2. PRINCIPLE OF PCR DETECTION

Human Respiratory Syncytial Virus detection by the polymerase chain reaction (PCR) is based on the amplification of the pathogen genome specific region using special *Human Respiratory Syncytial Virus* primers. In **Fluorescent End-Point** PCR, the amplified product is detected by using fluorescent dyes. These dyes are linked to oligonucleotide probes that bind specifically to the amplified product during thermocycling. A multichannel rotor-type fluorometer is specially designed to detect fluorescence emission from the fluorophores in a reaction mixture after PCR. It allows detection of the accumulating product without re-opening the reaction tubes after the PCR run. **AmpliSens® hRSV-FEP** PCR kit is a qualitative test that contains the Internal Control (IC). It must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition. **AmpliSens® hRSV-FEP** PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase by using a wax layer. Wax melts and reaction components mix only at 95 °C.

3. CONTENT

AmpliSens® hRSV-FEP PCR kit is produced in 2 forms:

AmpliSens® hRSV-FEP PCR kit (tubes of 0.2 ml), **REF** V37-50-R0,2-FEP-CE.

AmpliSens® hRSV-FEP PCR kit (tubes of 0.5 ml), **REF** V37-50-R0,5-FEP-CE.

AmpliSens® hRSV-FEP PCR kit includes:

Reagent	Description	Volume (ml)	Amount
PCR-mix-1-FEP/FRT hRSV ready-to-use single-dose test tubes (<i>under wax</i>)	colorless clear liquid	0.008	55 tubes of 0.5 or 0.2 ml
PCR-mix-2-FL	colorless clear liquid	0.77	1 tube
PCR-mix-Background	colorless clear liquid	0.5	1 tube
Mineral oil for PCR	colorless viscous liquid	4.0	1 dropper bottle
Positive Control cDNA hRSV-Flu (C+hRSV-Flu)	colorless clear liquid	0.1	1 tube
Positive Control STI-88 (CS+)	colorless clear liquid	0.1	1 tube
TE-buffer	colorless clear liquid	0.5	1 tube
Negative Control (C-)*	colorless clear liquid	1.2	1 tube
Internal Control STI-rec (IC)**	colorless clear liquid	0.12	5 tubes

* must be used in the extraction procedure as Negative Control of Extraction.

** add 10 µl of Internal Control STI-rec (IC) during the RNA extraction procedure directly to the sample/lysis mixture (see RIBO-sorb **REF** K2-1-Et-50-CE, RIBO-prep, **REF** K2-9-Et-50-CE protocols).

AmpliSens® hRSV-FEP PCR kit is intended for 55 reactions, including controls.

4. ADDITIONAL REQUIREMENTS

- RNA extraction kit.
- Reverse transcription kit.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol barriers (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with a rotor for 2-ml reaction tubes.
- PCR box.
- Personal thermocyclers (for example, Gradient Palm Cyclor (Corbett Research, Australia), GeneAmp PCR System 2700 (Applied Biosystems, USA), MyCycler (Bio-Rad, USA), Uno-2 (Biometra, Germany), or equivalent).
- Fluorometer (for example, ALA-1/4 (Biosan, Latvia), or equivalent).
- Disposable polypropylene microtubes for PCR (0.5- or 0.2-ml; for example, Axygen, USA).
- Refrigerator for 2–8 °C.
- Deep-freezer for ≤ –16 °C.

- Waste bin for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store and handle amplicons away from all other reagents.
- Thaw all components thoroughly at room temperature before starting detection.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats, protect eyes while samples and reagents handling. Thoroughly wash hands afterward.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all samples and unused reagents in compliance with local authorities' requirements.
- Samples should be considered potentially infectious and handled in a biological cabinet in accordance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid contact with the skin, eyes, and mucosa. If skin, eyes, or mucosa contact, immediately flush with water and seek medical attention.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING



Obtaining samples of biological materials for PCR-analysis, transportation and storage is described in manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

AmpliSens® hRSV-FEP PCR kit is intended for the analysis of RNA extracted with RNA extraction kits from nasal and oropharyngeal swabs, nasopharyngeal or tracheal sputum or aspirate, and autopsy material.

Nasal swabs are taken using a probe with a dry cotton swab. Make the patient blow his nose if it is filled with mucus. Insert the probe gently along the external nasal wall by 2–3 cm

till the inferior nasal concha. Then move the probe slightly lower, insert in the inferior nasal meatus under the inferior nasal concha, rotate, and pull it back along the external nasal wall.

When material is obtained, place the working area of the probe with a cotton swab into a sterile disposable tube with 500 µl of transport medium for storage and transportation of respiratory swabs. Break off the terminal part of the probe or cut it off to allow tight closing of the tube cup. Close the tube with the solution and the working area of the probe.

Oropharyngeal swabs are taken using a probe with a dry cotton swab. Obtain swabs by rotating the probe over the surface of tonsils, palatine arches, and the posterior wall of the pharynx after gargling of the oral cavity with water.

When material is obtained, place the working area of the probe with a cotton swab into a sterile disposable tube with 500 µl of transport medium for storage and transportation of respiratory swabs. Break off the terminal part of the probe or cut it off to allow tight closing of the tube cup. Close tube with the solution and the working area of the probe.



It is recommended to combine nasal and throat swabs. To do this, working parts of probes are placed into one tube with 500 µl of medium for storage and transportation of respiratory swabs and studied as one sample.

The sputum is collected into a sterile disposable container after preliminary gargling of the oral cavity with water. Nasopharyngeal or tracheal aspirates are taken by the conventional method and placed in a sterile disposable container.

Sectional material is placed in a sterile disposable container. Material is to be frozen after sampling or analyzed within 1 h.



Only one freeze–thaw cycle is allowed.

Pretreatment of material

Any operation with studied material is carried out in compliance with local authority's requirements.

All pretreatment manipulations are carried out using adjustable pipettes with disposable tips with aerosol barriers and disposable polypropylene tubes (1.5 ml or 10.0 ml). Disposable plastic dishes (tubes and tips) are to be discarded into a special container with a suitable disinfectant and utilized in compliance with local authority's requirements.

Swabs are used without preliminary preparation.

Sputum must be treated with Mucolysin reagent **REF** 180-CE according to the Mucolysin instruction manual. 100 µl of the pretreated sputum is used for RNA extraction.

Section material is homogenized with a sterile porcelain mortar and pestle, with subsequent preparation of a 10 % suspension in sterile saline or phosphate buffer. Transfer the suspension to 1.5-ml tube and centrifuge for 5 min at 10,000 rpm. 100 µl of the supernatant is used for RNA extraction. To repeat the test later, freeze the remained

suspension at $\leq -16\text{ }^{\circ}\text{C}$.

Additional reagents manufactured by Federal Budget Institute of Science “Central Research Institute for Epidemiology” are required:

1. Transport Medium for Storage and Transportation of Respiratory Swabs (unaliquoted)

REF 959-CE is used for sampling and storage of nasal and oropharyngeal swabs.

2. Mucolysin reagent **REF** 180-CE is used for pretreatment of sputum and aspirates.

7. WORKING CONDITIONS

AmpliSens® hRSV-FEP PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. RNA extraction

It is recommended to use the following nucleic acid extraction kits:

- RIBO-sorb, **REF** K2-1-Et-50-CE.

- RIBO-prep, **REF** K2-9-Et-50-CE.



Extract RNA according to the manufacturer's instructions.



Add **10 µl** of **Internal Control STI-rec** to each tube.

8.2. Reverse transcription

It is recommended to use the following kit for complementary DNA (cDNA) synthesis from RNA:

- REVERTA-L, **REF** K3-4-50-CE.



Carry out the reverse transcription according to the manufacturer's instructions.

8.3. Preparing the PCR

The total reaction volume is **25 µl**, the volume of cDNA sample is **10 µl**.

8.3.1. Preparing tubes for PCR

1. Prepare the required number of tubes with **PCR-mix-1-FEP/FRT hRSV** and wax for amplification of cDNA from clinical and control samples.
2. Add **7 µl** of **PCR-mix-2-FL** to the surface of the wax layer of each tube ensuring that it does not fall under the wax and mix with **PCR-mix-1-FEP/FRT hRSV**.
3. Add above **1** drop of **mineral oil for PCR** (about **25 µl**).
4. Prepare 2 tubes with **PCR-mix-1-FEP/FRT hRSV** and mark them as **Background**. Add **17 µl** of **PCR-mix-Background** to the surface of the wax layer of each tube ensuring that it does not fall under the wax and mix with **PCR-mix-1-FEP/FRT hRSV**. Add above **1** drop of **mineral oil for PCR**.

5. Using tips with aerosol barrier, add **10 µl** of **cDNA samples** obtained in RNA reverse transcription reaction.
6. Carry out the control amplification reactions:
 - NCA** - Add **10 µl** of **TE-buffer** to the tube labeled NCA (Negative Control of Amplification).
 - C+_{hRSV-Flu}** - Add **10 µl** of **Positive Control cDNA hRSV-Flu** to the tube labeled C+_{hRSV-Flu} (Positive Control of Amplification).
 - CS+** - Add **10 µl** of **Positive Control STI-88** to the tube labeled CS+.

8.3.2. Amplification

Run the following program in the thermocycler (see Table 1). When the temperature reaches 95 °C (pause mode), insert tubes into cells of the thermocycler and press the button to continue.

It is recommended to sediment drops from walls of tubes by short centrifugation (1–3 s) before placing them in the thermocycler.

Table 1

Programming thermocyclers at cDNA amplification of *Human Respiratory Syncytial virus*

	Thermocyclers with active temperature adjustment:						Thermocyclers with block temperature adjustment: Biometra		
	Terzik (DNA-Technology)			GeneAmp PCR System 2700 (Applied Biosystems), Gradient Palm Cycler (Corbett Research), MyCycler (Bio-Rad)					
Step	Tempe- rature	Time	Cycles	Tempe- rature	Time	Cycles	Tempe- rature	Time	Cycles
0	95 °C	pause		95 °C	pause		95 °C	pause	
1	95 °C	5 min	1	95 °C	5 min	1	95 °C	5 min	1
2	95 °C	10 s	42	95 °C	10 s	42	95 °C	25 s	42
	54 °C	20 s		54 °C	25 s		54 °C	40 s	
	72 °C	10 s		72 °C	25 s		72 °C	25 s	
3	72 °C	1 min	1	72 °C	1 min	1	72 °C	1 min	1
4	10 °C	storage		10 °C	storage		10 °C	storage	

9. DATA ANALYSIS

Detection is conducted with an ALA-1/4 fluorescence detector.



Please read Ala-1 Operating Manual before use of this kit.

Program the detector according to the manufacturer's manual and Guidelines.

9.1. Results interpretation

1. When the analysis is complete the results are automatically shown in the table as follows:

pos – positive result;

neg – negative result;

eq – equivocal result (signal is in grey zone);

nd – invalid result (specific signal and IC signal are absent in the sample).

2. The result of analysis is considered reliable only if the results obtained for Positive and

Negative Controls of Amplification and Negative Control of Extraction are correct (Table 2).

Table 2

Results for controls

Control	Stage for control	Result of automatic interpretation		Interpretation
		FAM channel (IC)	HEX channel (samples)	
C–	RNA extraction	+	RSv - neg	OK
NCA	Amplification	-	RSv - nd	OK
C+_{hRSV-Flu}	Amplification	-	RSv - pos	OK
CS+	Amplification	+	RSv - neg	OK

10. TROUBLESHOOTING

Results of analysis are not taken into account in the following cases:

1. PCR and detection of the samples with the **nd** result (except for NCA) are to be repeated. If the same result is obtained, analysis should be repeated starting from the RNA extraction stage. The **nd** result for NCA is normal.
2. PCR and detection of the samples with the **eq** result are to be repeated. If the same result is obtained, the samples are considered to be **positive**.
3. If no positive signal is detected for Positive Controls of Amplification, this may indicate that the temperature profile programming of the used instrument was incorrect, the configuration of the PCR reaction was incorrect, the storage conditions for kit components did not comply with the manufacturer's instruction, or that the reagents kit has expired. Check programming of the used Instrument, storage conditions, and the expiration date of the reagents and then repeat PCR analysis.
4. If a positive signal is detected for C– in the HEX channel or for NCA in both channels, this indicates contamination of reagents or samples. In such cases, the results of analysis must be considered as **invalid**. Analysis should be repeated and measures to detect and eliminate the source of contamination are to be taken.

If you have any further questions or if you encounter problems, please contact our Authorized Representative in the European Community.

11. TRANSPORTATION

AmpliSens® hRSV-FEP PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the **AmpliSens® hRSV-FEP** PCR kit are to be stored at 2-8 °C when not in use. All components of the **AmpliSens® hRSV-FEP** PCR kit are stable until labeled expiration date. The shelf life of reagents before and after the first use is the same, unless otherwise stated.



PCR-mix-1-FEP/FRT *hRSV* is to be stored away from light.

13. SPECIFICATIONS

13.1. Sensitivity

Analytical Sensitivity of **AmpliSens® hRSV-FEP** PCR kit is not less than 1×10^3 genome equivalents per 1 ml of sample (GE/ml).



The claimed analytical features of **AmpliSens® hRSV-FEP** PCR kit are guaranteed only when additional reagents kits RIBO-sorb, RIBO-prep and REVERTA-L (manufactured by Federal Budget Institute of Science “Central Research Institute for Epidemiology”) are used.

13.2. Specificity

The analytical specificity of **AmpliSens® hRSV-FEP** PCR kit is assured by selection of specific primers and probes as well as stringent reaction conditions. The primers and probes were checked for possible homologies to all sequences published in gene banks by sequence comparison analysis. The clinical specificity of **AmpliSens® hRSV-FEP** PCR kit was confirmed in laboratory clinical trials.














14. REFERENCES

1. Handbook “Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics”, developed by Federal State Institute of Science “Central Research Institute for Epidemiology” of Federal Service for Surveillance on Consumers’ Rights Protection and Human Well-Being, Moscow, 2008.

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of **AmpliSens® hRSV-FEP** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

16. KEY TO SYMBOLS USED

	Catalogue number		Sufficient for
	Batch code		Expiration Date
	<i>In vitro</i> diagnostic medical device		Consult instructions for use
	Version		Keep away from sunlight
	Temperature limitation	NCA	Negative control of amplification
	Manufacturer	C-	Negative control of extraction
	Date of manufacture	C+_{hRSV-Flu}, CS+	Positive control of amplification
	Authorised representative in the European Community	IC	Internal control
	Caution		

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
13.12.10	Cover page	The phrase “For Professional Use Only” was added
	Intended use	The phrase “The results of PCR analysis are taken into account in complex diagnostics of disease” was added
	Content	New sections “Working Conditions” and “Transportation” were added
		The “Explanation of Symbols” section was renamed to “Key to Symbols Used”
	Stability and Storage	The information about the shelf life of reagents before and after the first use was added
		Information that PCR-mix-1-FEP/FRT <i>hRSV</i> is to be kept away from light was added
	Key to Symbols Used	The explanation of symbols was corrected
	Text	Positive Control STI was changed to Positive Control STI-88
		Positive Control cDNA <i>hRSV</i> -Flu (C+) was changed to Positive Control cDNA <i>hRSV</i> -Flu (C+ _{<i>hRSV</i>-Flu})
22.06.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science “Central Research Institute for Epidemiology”