



For Professional Use Only

MC-EPh PCR kit Instruction Manual

AmpliSens[®]



Federal Budget Institute of Science
“Central Research Institute for
Epidemiology”
3A Novogireevskaya Street
Moscow 111123 Russia

TABLE OF CONTENTS

1. INTENDED USE 3
2. PRINCIPLE OF PCR DETECTION 3
3. CONTENT 3
4. ADDITIONAL REQUIREMENTS 3
5. GENERAL PRECAUTIONS..... 4
6. SAMPLING AND HANDLING 5
7. WORKING CONDITIONS..... 6
8. PROTOCOL 6
9. DATA ANALYSIS 7
10. TROUBLESHOOTING..... 8
11. TRANSPORTATION..... 8
12. STABILITY AND STORAGE..... 8
13. SPECIFICATIONS..... 9
14. REFERENCES 9
15. QUALITY CONTROL..... 9
16. KEY TO SYMBOLS USED 10

1. INTENDED USE

MC-EPh PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of *Mycoplasma* microorganisms DNA in the biological material by using electrophoretic detection of the amplified products in agarose gel.



This kit is used only for scientific research.

2. PRINCIPLE OF PCR DETECTION

Mycoplasma detection by the polymerase chain reaction (PCR) is based on the amplification of pathogen DNA specific region using special *Mycoplasma* primers. After PCR the amplified product is detected in agarose gel. **MC-EPh PCR kit** uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase by using a wax layer. Wax melts and reaction components mix only at 95 °C.

3. CONTENT

MC-EPh PCR kit is produced in 2 forms:

MC-EPh PCR kit variant 50 R (tubes 0.2 ml) **REF** B72-50-R0,2-CE;

MC-EPh PCR kit variant 50 R (tubes 0.5 ml) **REF** B72-50-R0,5-CE.

MC-EPh PCR kit variant 50 R includes:

Reagent	Description	Volume (ml)	Quantity
PCR-mix-1-R MYC-COM ready-to-use single-dose test tubes (<i>under wax</i>)	colorless clear liquid	0.005	55 tubes of 0.5 or 0.2 ml
PCR-mix-2 blue	blue clear liquid	0.6	1 tube
Mineral oil for PCR	colorless viscous liquid	2.0	1 dropper bottle
Positive Control DNA <i>Mycoplasma hominis</i> (C+_{M.hominis})	colorless clear liquid	0.2	1 tube
DNA-buffer	colorless clear liquid	0.5	1 tube
Negative Control (C-)*	colorless clear liquid	1.2	1 tube

* must be used in the extraction procedure as Negative Control of Extraction.

MC-EPh PCR kit variant 50 R is intended for 55 reactions, including controls.

4. ADDITIONAL REQUIREMENTS

- DNA extraction kit.

- Agarose gel detection kit.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol barriers (up to 200 µl).
- PCR box.
- Tube racks.
- Personal thermocyclers (for example, Terzik (DNA-Technology, Russia)).
- Refrigerator with temperature 2-8 °C.
- Deep-freezer with temperature not more than minus 16 °C.
- Waste bin for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store and handle amplicons away from all other reagents.
- Thaw all components thoroughly at room temperature before starting detection.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats, protect eyes while samples and reagents handling. Thoroughly wash hands afterward.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all samples and unused reagents in compliance with local authorities requirements.
- Samples should be considered potentially infectious and handled in a biological cabinet in accordance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills using a disinfectant such as 0.5% sodium hypochlorite, or other suitable disinfectant.
- Avoid contact with the skin, eyes and mucosa. If skin, eyes and mucosa contact immediately flush with water, seek medical attention.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional, it should begin in the Extraction Area,

move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.

- Observe the rules and instructions to prevent semination of the environment objects.



Some components of this kit contain Sodium Azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING

MC-EPh PCR kit is intended for analysis of DNA extracted with DNA extraction kits from the biological material:

- Nasal swabs, conjunctive swabs, outflows.
- Joint synovial fluid.
- Dotter substance, embryo allantoic fluid.
- Parenchymatous organs, trachea, air sacs.
- Whole blood.
- Medicinal serum, cell cultures (1×10^6 cells for 1 ml).

6.1. Material sampling

Blood with 6 % EDTA (20:1).

The parts of tissues and organs – about 1x1x1 (cm) or less.

Materials are to be analyzed the next day after sampling if they are stored at the temperature 2-8 °C.



Materials can be stored at the temperature not more than minus 16 °C for one month.

6.2. Preparation of samples

Whole blood samples with EDTA, medicinal serum and cell cultures are used for DNA extraction without treatment.

Other fluid samples (1.5-ml). Centrifuge at 10000 rpm for 5 min. Carefully remove and discard the supernatant using a tip with aerosol barrier leaving about 200 µl of the liquid on the sediment. If the pellet is too small, add more material (1.5 ml) and repeat centrifugation. Suspend the pellet in supernatant and use 100 µl for DNA extraction.

Parenchymatous organs, trachea, air sacs should be ground in a sterile porcelain mortar or glass homogenizer. Add equal volume of saline and thoroughly homogenize. Incubate at room temperature for 30 min. Transfer 100 µl of the upper part of prepared suspension to a sterile tube for DNA extraction.

7. WORKING CONDITIONS

MC-EPh PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. DNA extraction

It is recommended to use the following nucleic acid extraction kits:

- "DNA-sorb-B", **REF** K1-2-50-CE



Carry out the amplification according to the manufacturer instruction.
The volume of clinical sample is 100 µl.



Add 100 µl of Negative Control (C-) into tube used as Negative Control of Extraction.



Add 90 µl of Negative Control (C-) and 10 µl of Positive Control DNA *Mycoplasma hominis* (C+) into tube used as Positive Control of Extraction.

8.2. Preparing the PCR

Total reaction volume - **25 µl**, volume of cDNA sample - **10 µl**.

8.2.1. Detection of *Mycoplasma* DNA

1. Prepare the required number of PCR tubes with **PCR-mix-1-R MYC-COM** and wax for amplification of DNA from clinical and control samples.
2. Add **10 µl** of **PCR-mix-2 blue** to the surface of the wax layer of each tube ensuring that it does not fall under the wax and mix with **PCR-mix-1-R MC**.
3. Add above **1 drop** of **mineral oil for PCR** (about **25 µl**).
4. Using tips with aerosol barrier add **10 µl** of **DNA samples** obtained from clinical or control samples.
5. Carry out the control amplification reactions:

NCA -Add **10 µl** of **DNA-buffer** to the tube labeled NCA (Negative Control of Amplification).

C+ -Add **10 µl** of **Positive Control DNA *Mycoplasma hominis*** to the tube labeled C+.

8.2.2. Amplification of *Mycoplasma* DNA

Run the following program on the thermocycler (see Table 1). When the temperature reaches 95 °C (pause regimen), insert tubes into the cells of amplifier and press the button to continue.

It is recommended to precipitate drops from walls of tubes by short vortexing (1–3 sec) before their insertion in a thermocycler.

Programming thermocyclers for *Mycoplasma* DNA amplification

Step	Thermocyclers with active temperature adjustment						Thermocyclers with block temperature adjustment		
	Terzik (DNA-Technology), Omn-E (Hybaid)			GeneAmp PCR System 2700 (Applied Biosystems), Palm Cycler (Corbett Research)			Ampli-3 (Biocom), Biometra, MiniCycler, PTC-100 (MJ Research)		
	Temperature, °C	Time	Cycles	Tempera- ture, °C	Time	Cycles	Tempera- ture, °C	Time	Cycles
0	95	pause		95	pause		95	pause	
1	95	5 min	1	95	5 min	1	95	5 min	1
2	95	10 s	41	95	10 s	41	95	1 min	41
	61	10 s		61	25 s		61	1 min	
	72	10 s		72	25 s		72	1 min	
3	72	1 min	1	72	1 min	1	72	1 min	1
4	10	storage		4	storage		10	storage	

After the reaction is finished PCR tubes must be collected and sent to the room for PCR products analysis.

Analysis of amplification products is performed by separation of DNA fragments in agarose gel.

The amplified samples can be stored for 16 h at room temperature, for 1 week at 2-8 °C (be sure to warm the samples to room temperature before running electrophoresis).

9. DATA ANALYSIS

It's recommended to use the following agarose kit for electrophoretic detection:

- "EPh" variant 200, **REF** K5-200-CE.

Analysis of results is based on the presence or absence of specific bands of amplified DNA in agarose gel (1.7 %). The lengths of specific amplified DNA fragments are:

- *Mycoplasma* – 509 bp (base pairs)



Put the protective mask or use the glass barrier while watching and photographing the gel.

Start analysis from results for controls (see Table 2).

Results for controls

Control	Controlled step	Specific bands in the agarose gel 509 bp	Interpretation
C-	DNA extraction	No	OK
PCE	DNA extraction	Yes	OK
NCA	Amplification	No	OK
C+	Amplification	Yes	OK

- The sample is considered to be positive if the band of 509 bp is present in agarose gel.
- The sample is considered to be negative if the band of 509 bp is absent.
- Besides specific bands, the fuzzy bands of primer dimers can be seen in lanes, they are situated lower than the level of 100 bp.

10. TROUBLESHOOTING

Analysis results are not taken into account in the following cases:

1. If the results of control samples do not correspond to the listed above (Table 2), then the tests should be repeated.
2. If the specific 509-bp band is absent in lanes corresponding to positive control (PCE, C+), it may be mistake of preparation of reagents, amplification or program error of thermocycler.
3. If nonspecific bands appear in lanes at different levels, it may be caused by the lack of “hot start” or a false temperature regime in thermocycler.
4. If the specific 509-bp band appears in lanes corresponding to negative control (NCA, C-), it means that reagents or samples are contaminated. In such cases, analysis results must be considered as irrelevant. Test analysis should be repeated and measures for detecting contamination source must be taken.

11. TRANSPORTATION

MC-EPh PCR kit should be transported at 2–8 °C for no longer than 5 days. Once received, the PCR kit should be dekitted according to the indicated storage conditions.

12. STABILITY AND STORAGE

All components of **MC-EPh** PCR kit are to be stored are to be stored at 2–8 °C when not in use.. All components of the **MC-EPh** PCR kit are stable until labeled expiration date. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

13. SPECIFICATIONS

13.1. Sensitivity

Analytical Sensitivity of **MC-EPh** PCR kit is no less than 1×10^3 copies per 1 ml of a sample (cop/ml).



The claimed analytical features of MC-EPh PCR kit are guaranteed only when additional kits of reagents, “DNA-sorb-B” and “EPh” (manufactured by Federal State Budget of Science Central Research Institute for Epidemiology) are used.

13.2. Specificity

Specificity of **MC-EPh** PCR kit is ensured by selection of specific primers and strict reaction conditions as well as laboratory and clinical trials.














14. REFERENCES

1. García, M., M. W. Jackwood, S. Levisohn, and S. H. Kleven. Detection of *Mycoplasma gallisepticum*, *M. synoviae*, and *M. iowae* by multi-species polymerase chain reaction and restriction fragment length polymorphism. Avian Dis 39: 606–616.1995.

15. QUALITY CONTROL

In compliance with Federal Budget Institution of Science Central Research Institute for Epidemiology ISO 13485 – certified Total Quality Management System, each lot of **MC-EPh** PCR kit is tested against predetermined specifications to ensure consistent product quality.

16. KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Sufficient for
	Research use only		Expiration Date
	Version		Consult instructions for use
	Temperature limitation		Keep away from sunlight
	Upper limit of temperature	NCA	Negative control of amplification
	Manufacturer	C-	Negative control of extraction
	Date of manufacture	C+	Positive control of amplification
FBIS CRIE	Federal Budget Institute of Science "Central Research Institute for Epidemiology"	IC	Internal control

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
25.01.12 LA	Cover page	The phrase "For Professional Use Only" was added
	Content	New sections "Working Conditions" and "Transportation" were added
		The "Explanation of Symbols" section was renamed to "Key to Symbols Used"
	Stability and Storage	The information about the shelf life of reagents before and after the first use was added
	Key to Symbols Used	The explanation of symbols was corrected
Cover page, text	The name of Institution was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"	
11.02.14 ChA	Front page, Key to symbols used	Symbol IVD was changed to RUO
13.10.14 PM	Content	"Positive Control DNA <i>Mycoplasma hominis</i> (C+)" was changed to "Positive Control DNA <i>Mycoplasma hominis</i> (C+ _{<i>M.hominis</i>})"