---

# CERTEST

## *E. coli* O111

---

ONE STEP *E. coli* O111  
CARD TEST

---

CERTEST BIOTEC S.L.

---

# CERTEST *E. coli* O111

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## ENGLISH

### INTENDED USE

CerTest *E. coli* O111 one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Escherichia coli* O111 in stool samples and *E. coli* O111 suspected colonies in stool culture.

CerTest *E. coli* O111 card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of *Escherichia coli* O111 infection and it could be used to identify of suspected isolates of *E. coli* O111 from selective media (stool culture).

### INTRODUCTION

Enterohaemorrhagic *Escherichia coli* (EHEC) infection presents with a wide spectrum of clinical manifestations, including asymptomatic carriage, non-bloody diarrhoea, haemorrhagic colitis (HC), the haemolytic-uremic syndrome (HUS), and thrombotic thrombocytopenic purpura (TTP).

Not only is *Escherichia coli* O157:H7 an important agent for haemorrhagic colitis, it is also one of the leading causes of bacterial diarrhoea. Transmission of *Escherichia coli* (EHEC) is primarily food-borne. Undercooked meat is the most common culprit, dairy products and secondary person-to-person spread is also important.

O157:H7 is epidemiologically the most frequent serotype; however the improvement of microbiological detection methods has reported an increasing frequency of other serotypes such as O145:H2, O111:H2, O26:H11, O103:H2, and O121:H19 (non-O157 STEC).

The organism produces at least two Shiga-like toxins. These toxins are thought to have direct pathogenic significance in *Escherichia coli* O157:H7 infection. This infection is usually diagnosed from a positive stool culture, from the presence of Shiga toxins, or both. Timely collection (within 7 days of illness onset) of a stool sample for culture is imperative for a high recovery rate.

### TEST PRINCIPLE

CerTest *E. coli* O111 is based on the principle of a qualitative immunochromatographic assay for the determination of *Escherichia coli* O111 in stool samples and *E. coli* O111 suspected colonies in stool culture.

The strip consists of a nitrocellulose membrane pre-coated with antibodies on the test line (T), in the results window, against *Escherichia coli* O111 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (antibodies anti-*Escherichia coli* O111) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Escherichia coli* O111 antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Escherichia coli* O111 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Escherichia coli* O111 antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Escherichia coli* O111 antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *E. coli* O111 test. Do not use any other commercial kit component.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


## SPECIMEN COLLECTION AND PREPARATION

### Stool samples

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

### Suspected *E. coli* O111 colonies in stool culture


The Sorbitol-MacConkey Agar (SMAC) is the method of choice for the isolation of non-O157 STEC (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation in SMAC the non-O157 STEC colonies will be small and pink.

1. Examine Sorbitol MacConkey plates after 1 day incubation. Select non-O157 STEC typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected non-O157 STEC colonies (2) and add them to the collection tube.

Non-O157 STEC colonies in SMAC

2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

## MATERIALS

### MATERIALS PROVIDED

- CerTest *E. coli* O111 card tests
- Instructions for use
- Collection tubes with diluent

### MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer
- Sorbitol MacConkey Agar (SMAC)

## TEST PROCEDURE


Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the collection tube in order to assure good sample dispersion.
2. Remove the CerTest *E. coli* O111 card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>E. coli</i> O111	Interpretation of the results
1.	- GREEN	There is no <i>Escherichia coli</i> O111 presence. No infection caused by <i>Escherichia coli</i> O111.
2.	+ GREEN-RED	There is <i>Escherichia coli</i> O111 presence. Infection caused by <i>Escherichia coli</i> O111.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

**NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

**INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

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## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

## LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. CerTest *E. coli* O111 should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of *Escherichia coli* O111 in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical method or by PCR) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, an *Escherichia coli* O111 determination should be carried out on a sample from an enrichment culture.

## EXPECTED VALUES

*Escherichia coli* O157:H7 causes 73,000 illnesses in the United States annually. That means 8,598 cases, 17% require hospitalizations, 4% hemolytic uremic syndrome cases and 0.5% deaths. The main transmissions routes are: 52% foodborne, 21% unknown, 14% person to person, 31% waterborne, 3% animal contact and 0.3% laboratory related.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity (detection limit)

Detection limit value of CerTest *E. coli* O111 is  $2.37 \times 10^4$  CFU/mL.

### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (CerTest *E. coli* O111, CerTest) and these results were confirmed using culture (positive results). The results were as follows:

IC test: CerTest <i>E. coli</i> O111		Culture		
		+	-	Total
+		5	0	5
-		0	60	60
Total		5	60	65

	Sensitivity	Specificity	PPV	NPV
CerTest <i>E. coli</i> O111 vs culture	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Escherichia coli* O111 using CerTest *E. coli* O111.

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### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *E. coli* O111; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter coli</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Salmonella enteritidis</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i>	<i>Salmonella paratyphi</i>	<i>Shigella sonnei</i>
<i>E. coli</i> O26	<i>Salmonella typhi</i>	<i>Staphylococcus aureus</i>
<i>E. coli</i> O157	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica</i>
<i>Helicobacter pylori</i>	<i>Shigella boydii</i>	

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