CERTEST
Rotavirus+
Adenovirus

ONE STEP
Rotavirus+Adenovirus
COMBO CARD TEST

CERTEST BIOTEC S.L.

For Information Purposes Only

CERTEST Rotavirus+Adenovirus

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ENGLISH

INTENDED USE

CerTest *Rotavirus+Adenovirus* one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Rotavirus* and *Adenovirus* in stool samples.

CerTest *Rotavirus+Adenovirus* combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Rotavirus* and/or *Adenovirus* infection.

INTRODUCTION

Rotavirus, *Adenovirus* and *Astrovirus* are most common and major causes of severe gastroenteritis in infants and young children. Pattern also observed in adults. They are transmitted by faecal-oral contact. The main symptoms of viral gastroenteritis are watery diarrhoea and vomiting. The affected person may also have headache, fever and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness (*Rotavirus* 3 days, *Adenovirus* 5-8 days and *Astrovirus* 3 days).

TEST PRINCIPLE

CerTest *Rotavirus+Adenovirus* is based on the principle of a qualitative immunochromatographic assay for the determination of *Rotavirus* and *Adenovirus* in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Rotavirus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Rotavirus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Adenovirus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Adenovirus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is *Rotavirus* positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Rotavirus* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is *Adenovirus* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*Adenovirus* monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Rotavirus* antibodies present on the membrane of strip A (test line) and the anti-*Adenovirus* antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *Rotavirus* and *Adenovirus* presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Rotavirus* and anti-*Adenovirus* antibodies present on the membranes (test lines) will not capture the antigen-red -coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

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
- Components provided in the kit are approved for use in the CerTest *Rotavirus+Adenovirus* combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

MATERIALS

MATERIALS PROVIDED

- CerTest *Rotavirus+Adenovirus* combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Rotavirus+Adenovirus* combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5) and 4 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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

Rotavirus Strip A-procedure


Adenovirus Strip B-procedure

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

INVALID RESULTS

	A (Rotavirus)	B (Adenovirus)	Interpretation of the results
1.	- GREEN	- GREEN	There is no Rotavirus or Adenovirus presence. No viral infection caused by Rotavirus or Adenovirus.
2.	+ GREEN-RED	+ GREEN-RED	There is Rotavirus and Adenovirus presence. Viral infection caused by Rotavirus and Adenovirus.
3.	+ GREEN-RED	- GREEN	There is Rotavirus presence. Viral infection caused by Rotavirus.
4.	- GREEN	+ GREEN-RED	There is Adenovirus presence. Viral infection caused by Adenovirus.
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

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QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. After one week of infection, the number of viruses in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
5. CerTest *Rotavirus+Adenovirus* should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
6. Positive results determine the presence of *Rotavirus* and/or *Adenovirus* in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
7. A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit values. If the symptoms or situation still persist, a *Rotavirus* and/or *Adenovirus* determination should be carried out with another technique (for example microscopy).

EXPECTED VALUES

Acute gastroenteritis is a common disorder in young children; moreover, the associated dehydration is a leading cause of admission to hospital in industrialized countries. Acute diarrhea is a major health problem throughout the world and main source of mortality in developing countries. Enteric viruses have been recognized as the most significant etiological agents of the disease and yet, four categories of viruses are being considered clinically relevant: group A *Rotavirus* (family *Reoviridae*), *Norovirus* (family *Caliciviridae*), *Adenovirus* and *Astrovirus*. Several studies proved co-infections in infants 46% with acute watery diarrhea.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (CerTest *Rotavirus+Adenovirus*, CerTest) and these results were compared with a commercial available test (Ridascreen®*Rotavirus* ELISA Test, r-Biopharm) for strip A and the results for strip B were confirmed using PCR. The results were as follows:

CerTest <i>Rotavirus+Adenovirus</i> (<i>Rotavirus</i>)	Ridascreen® <i>Rotavirus</i> ELISA Test		
		+	-
+	18	1	19
-	0	43	43
Total	18	44	62

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Rotavirus+Adenovirus</i> (<i>Rotavirus</i>) vs Ridascreen® <i>Rotavirus</i> ELISA Test	>99%	98%	>94%	>99%

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The results showed a high sensitivity and specificity to detect *Rotavirus* and *Adenovirus* using CerTest *Rotavirus+Adenovirus*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Rotavirus+Adenovirus*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Adenovirus</i> (Strip A)	<i>Escherichia coli</i> O157:H7	<i>Salmonella typhimurium</i>
<i>Astrovirus</i>	<i>Giardia lamblia</i>	<i>Salmonella typhi</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Shigella boydii</i>
<i>Campylobacter jejuni</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Clostridium difficile</i>	<i>Norovirus</i>	<i>Shigella flexneri</i>
<i>Cryptosporidium parvum</i>	<i>Rotavirus</i> (Strip B)	<i>Shigella sonnei</i>
<i>Enterovirus</i>	<i>Salmonella enteritidis</i>	<i>Staphylococcus aureus</i>
<i>Entamoeba histolytica</i>	<i>Salmonella paratyphi</i>	<i>Yersinia enterocolitica</i>

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